GIMA

Gima S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy
Italia: tel. 199 400 401 - fax 199 400 403 Export: tel. +39 02
953854209/221/225 fax +39 02 95380056 gima@gimaitaly.
com - export@gimaitaly.com www.gimaitaly.com

PROFESSIONAL MEDICAL PRODUCTS

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The contrast sensitivity chart was developed, based on a definition of contrast ($C = \text{background luminance minus letter luminance, divided by background luminance}$) and of contrast sensitivity, which is equal to the reciprocal of the contrast ($S = 1/C$). The chart has 7 lines, each of which has three triplets. Each triplet has three letters with the same contrast.

The chart, therefore, shows 21 different contrasts, progressively decreasing from top to bottom and from left to right.

As the contrasts progress across the values of the triplet, 21 different degrees of contrast sensitivity are established. The variation in contrast between one triplet and another is always constant, at all points of the table and, consequently, the modifications in sensitivity are also always constant. Each triplet can be identified by a number above it, shown on the back of the chart.

This number represents the contrast sensitivity required to see that triplet and corresponds to the exponent given on the base of 10 to obtain the value of that particular sensitivity.

The test should be carried out in good lighting conditions, preferably under natural lighting, at a distance varying from 2 to 3 meters, depending on the visual acuity of the patient being tested.

The separate triplets are to be considered as seen, if at least two out of the three letters that make them up are perceived.

The eyes are, firstly, examined separately, beginning from the stronger one, and then together. The date, distance of reading and identification number of the last triplet perceived should be registered, for tests regarding OD, OO and OS. The difference in contrast sensitivity of the two eyes examined separately and the contrast sensitivity of the two eyes examined together varies, with results between 0.1 and 0.3 approximately. The test may be repeated at regular intervals, following treatment in order to evaluate the results or, independently from this, to monitor whether the cataract condition is progressive, stationary or regressive over time.

Over-frequent testing is unnecessary, however the interval between tests should always be decided by the doctor, in relation to the requirements of each single patient.

With every three triplets, contrast sensitivity is doubled or is halved. It doubles, moving from top left towards bottom right and, vice versa, it is halved, going from bottom right towards top left, because the contrast is halved or doubled respectively. If we express sensitivity using its logarithm, it increases or decreases every three triplets by 0.3, as each triplet differs from the one before and the one after for a value of 0.1, expressed using the sensitivity logarithm. Overall, it is a full version of the chart with all the triplets that go to make it up, as well as with the logarithmic values of sensitivity relative to each triplet.

Suggested method for testing contrast sensitivity using the attached chart.

The test should be carried out under good and constant lighting, at the best distance chosen by the patient, which will vary from 2 to 3 meters. It is advisable to administer the test first to the strongest eye, then to the other and finally to both eyes together, registering the date of the test, distance of reading and identification number of the last triplet perceived (the identification number is that found above