

Chapter11 Memory Function

The sphygmomanometer is designed to store the blood pressure, pulse rate values and the date and time when measured, which are up to 100 groups. If there have been stored 100 groups, the earliest results will be deleted when saving the 101 group of measurement results.

11.1 Review Memory Values

1. In the main interface (interface when starting-up), press **【MEMORY】** button to review the latest measured values in big-font with the serial number from 1 to 100.

2. Press **【UP】** / **【DOWN】** button to circularly switch the former measurement values.

*The right figure shows that there is no measurement result can be displayed.

3. Press **【LIST】** button to display the data list interface.

4. Press **【TREND】** button to display trend interface.

End to display the measurement values:

Press **【EXIT】** button to return to the main interface or hold **【ON/OFF】** button to turn the device off.

11.2 Delete Memory Values

Users can delete all memory values of a user instead of separately delete one memory value.

1. Press **【MENU】** button to enter the system menu, select **【DELETE DATA】** option to enter its interface, select the user whose data to be deleted, after confirming again, all measurement results of the selected user will be deleted.

2. Finish Operation

Select **【EXIT】** to return to the previous menu, or hold **【ON/OFF】** button to turn the device off.

Chapter12 SpO₂ Measurement Function

(This chapter is only suitable for European Union market)

Precautions during SpO₂ Measurement:

 **Note** 

⊗Make sure the nail covers the light. The probe cable should be on the backside of the hand. Improper probe placement or improper contact with the test site will influence the measurement.

⊗SpO₂ value always displays in the fixed place.

⊗The test site shall not use external coloring agent (such as nail polish, dyestuff or color skin care products, etc.), otherwise it will affect the measurement.

⊗The fingers which are too cold or too thin may affect the measure accuracy, please insert the thicker finger such as thumb or middle finger deeply enough into the probe.

⊗The SpO₂ probe is suitable for children and adults (not suitable for infants and newborns). The device may not applicable for all patients. If you are unable to achieve stable readings, stop using it.

⊗Data averaging and signal processing will delay the SpO₂ displaying and data values transmitting. The update time of measurement data is less than 30 seconds, when signal attenuation, weak perfusion or other interference appears, it will result in time increasing of dynamic data averaging, which depends on the PR value.

⊗The PLETH waveforms are not normalized, which is used as the indicator for signal incompleteness. So the accuracy of the measured values may decrease when the waveform does not tend to smooth and stable. When the waveform tends to smooth and stable, the reading is optimal value, and the waveform at the moment is the most standard one..

⊗The temperature for the contact surface of the device with the body is less than 41℃, and this temperature value is measured by a temperature measuring device.

⊗The device does not provide over-limit alarm function, so it is inapplicable for using in places where need such function.

⊗The SpO₂ probe has been calibrated before leaving factory. It does not need to be calibrated during maintenance.

⊗SpO₂ probe is calibrated to show functional oxygen saturation.

⊗The SpO₂ probe and photoelectric receiving tube should be arranged in a way with the subject's arteriole in a position there between. Make sure the optical path is free from any optical obstacles like rubberized fabric, to avoid inaccurate measurement.

⊗As the measure is taken on the basis of arteriole pulse, substantial pulsating blood flow of subject is required. For a subject with weak pulse due to shock, low ambient/body temperature, major bleeding, or use of vascular contracting drug, the SpO₂ waveform (PLETH) will decrease. In this case, the measurement will be more sensitive to interference.

⊗Percentage modulation of infrared signal as the indication of pulsating signal strength, patient simulator has been used to verify its accuracy under conditions of low perfusion. SpO₂ and PR values are different due to low signal conditions, compare them with the known SpO₂ and PR values of input signal.

⊗The claim for SpO₂ accuracy should be supported by clinical research measurements covering the entire spectrum. By artificially inducing to different stable oxygen levels, make it in the range of 70 % ~ 100 % of SpO₂. Use secondary standard SpO₂ measuring equipment for comparison to collect SpO₂ values together with the tested product, compose paired data groups for accuracy analysis.

⊗Clinical report records data of 12 healthy volunteers, including 6 females and 6males. Volunteers' age ranges from 21 to 29. Skin color is distributed from dark to light, including 3 dark black skins, 2 medium black skins, 5 light skins, 2 white skins.

⊗When using the device, please keep it away from the instruments that can generate strong electric or magnetic field. Use the device in an inappropriate environment may cause interference to the surrounding radio equipment or affect its working.

⊗If necessary, please log on our company's official website to download the list of SpO₂ probes and extension cords that can be used in conjunction with this device.

 **Warning** 

⊗Check if the cable of SpO₂ probe is in normal condition before measuring. After unplugging the SpO₂ probe cable from the socket, "SpO₂%" and "bnp" on the screen will disappear.

⊗Do not use the SpO₂ probe once the package or the probe is found damaged. Instead, you shall return it to the vendor.

⊗The supplied SpO₂ probe is only suitable for use with this device. This device can only use the SpO₂ probe described in this manual. It is the responsibility of the operator to check the compatibility of the device and the SpO₂ probe (and extension cable) before use. Incompatible accessories may result in device performance decreasing or cause injury to the patient.

⊗SpO₂ probe is a medical product that can be used repeatedly.

⊗The measured value may be normal seemingly for the testee who has anemia or dysfunctional hemoglobin (such as carboxyhaemoglobin (COHb), methaemoglobin (MetHb) and sulphaemoglobin (SuHb)), but the testee may appear hypoxia, it is recommended to perform further assessment according the clinical situations and symptoms.

⊗Pulse oxygen has only reference significance for anemia and toxic hypoxia, as some patients with severe anemia still show better pulse oxygen measurements.

⊗Measurement accuracy can be affected by the interference of electrosurgical equipment.

⊗Do not install the SpO₂ probe on an extremity with arterial catheter or receiving intravenous injection.

⊗Do not perform SpO₂ measuring and NIBP measuring on same limb at one time, because obstruction of blood flow during NIBP measuring may adversely affect the reading of SpO₂ value.

⊗Excessive movement (active or passive) of the subject or severe activity can affect the measurement accuracy.

⊗Excessive ambient light may affect the measuring results, such as surgical light (especially xenon light sources),

bilirubin lamp, fluorescent lamp, infrared heater and direct sunlight, etc. In order to prevent interference from ambient light, make sure to place the probe properly and cover the probe with opaque material.

⊗The measured value may be inaccurate during defibrillation and in a short period after defibrillation, as the SpO₂ probe not have defibrillation-proof function.

⊗The person who is allergic to silicone, PVC, TPU, TPE or ABS can not use this device.

⊗For some special patients, it should be a more prudent inspecting in the measurement part. The probe can not be clipped on the edema and tender tissue.

⊗Do not stare at the luminescent component directly when the device is turned on (infrared light is invisible), even if for maintenance purpose, or it may have bad influence to the eyes.

⊗Uncomfortable or painful feeling may appear if using the SpO₂ probe ceaselessly, especially for the microcirculation barrier patients. It is recommended that the measurement should not be taken at the same position for over 2 hours. Continuous, long measurements may increase the risk of unwanted changes in the skin characteristics, such as exceptionally sensitive, reddish, blistering or oppressive necrosis, especially for neonates or the patients with perfusion disorder and change or immature skin form. It should be paid special attention to check the placement position of the probe according to the skin quality change, correct optical alignment and attachment method. Check the attachment position periodically and change it when the skin quality decreases. A more frequent check may be required due to the difference of patient's state.

⊗Some models of functional tester or patient simulator can measure the accuracy of the device that reproduces the calibration curve, but it can not be used to evaluate the accuracy of this device.

⊗Please refer to related medical literature for detailed clinical restrictions and contraindications,

⊗This device is not used for treatment purpose.

⊗Do not use the SpO₂ probe during MRI and CT scanning, as the induced current may cause burns.

⊗When the device is ON, if power interrupt for more than 30s, the SpO₂ probe needs no operation after the power is restored, after the device is turned on, ensure that the SpO₂ probe can be used normally.


⊗The probe can be used before/after sport, but not recommended to use during exercising.

⊗Patient simulator has been used to verify the pulse rate accuracy, it is stated as the root-mean-square difference between the PR measurement value and the value set by simulator.

Chapter13 SpO₂ Measurement Method

(This chapter is only suitable for European Union market)

1) Attach the SpO₂ probe to the appropriate site of the patient finger as the following figure.



Place SpO₂ probe

2) Plug the connector of the SpO₂ probe cable into the USB socket in the lower right of the device. The main interface will switch to SpO₂ interface. This operation brings no affection to other functions.If no finger inserted into the probe or finger placement is improper, the SpO₂ interface displays finger-out prompt; when the probe cable is unplugged from the USB socket, the SpO₂ interface will display probe-off prompt, a few seconds later, it automatically return to the main interface.

3) The pulse sound can be set to on or off during measuring the SpO₂. The operation steps are: Press **【MENU】** to enter the system menu, select **【POROMPT SETUP】** to enter the setting interface, and select **【SpO₂ PROMPT】** to enter the SpO₂ prompt interface, in which you can set the pulse sound to on or off. Note: Please make sure the prompt switch is set to ON state before turning on the pulse sound. (set the prompt switch to on or off in the **【POROMPT SETUP】** interface)

 **Note** 

SpO₂ display range: 0 % ~ 100 %, PR display range: 30 bpm (beats/min) ~ 250 bpm (beats/min)

A SpO₂ probe fault may exist in the following situations, please stop using it and contact our after-sales service:

⊗After reconnecting the probe cable with the USB port multiple times, the device main interface does not show any response;

⊗When probe cable is connected normally and the finger properly placed in the probe, the SpO₂ interface still displays finger-out or probe-off prompt;

⊗The SpO₂ interface prompts probe error;

Measurement Limitation

During operation, the accuracy of SpO₂ readings can be affected by:

- High-frequency electromagnetic interference such as interference from electrosurgical apparatus connected to the system.
- Intravenous dyestuff.
- Excessive patient movement.
- External light.
- Improper SpO₂ probe installation or incorrect contact position of the patient.
- Temperature of SpO₂ probe(optimal temperature range: 28℃ ~ 40℃).
- Place the SpO₂ probe on an limb that has a blood pressure cuff, arterial catheter, or intravascular line.
- Concentrations of dysfunctional hemoglobin, such as carboxyhemoglobin(COHb) and methemoglobin(MetHb).
- SpO₂ is too low, Bad circular perfusion of the part being measured.
- Intravascular staining agents (such as indocyanine green or methylene blue), skin pigmentation.
- It is required to use SpO₂ probe which is provided by our company, contact with our sale department when necessary.

Chapter14 Error Message

Error message will be displayed in the screen if there is something wrong when measuring. The causes and solutions are shown as follows:

Error Message	Causes	Solutions
Self-test failure System failure	Function abnormal	Please contact us
Loose cuff	Cuff is not connected correctly.	Correctly connect cuff (refer to Chapter 10)
Air leakage	Cuff plug fall off	Make sure the cuff plug is securely inserted in the windpipe (refer to Chapter 10)
Air pressure error	Air pressure error	Refer to the troubleshooting
Weak signal	The pulse signal is too weak or the cuff is loose.	Correctly connect cuff (refer to Chapter 10)
Overpressure	Cuff is blocked or squeezed	Correctly connect cuff (refer to Chapter 10)
Excessive movement Over range Saturated signal	The signal extent is too big owing to the arm or body moving or other reasons when measuring	Keep arm, body still, measure again
Time out	It takes too much time	

Chapter15 Troubleshooting

Abnormal Phenomenons	Causes	Solutions
BP measurement values too high or too low.	Cuff is not connected correctly.	Correctly connect cuff (refer to Chapter 10)
	Talk or move arms when measuring	Keep quiet and restart a measurement

	The turnout clothing presses the arm	Take off the clothing which presses the arm, and restart a measurement
No pressure	Cuff leakage	Buy a new cuff
	The cuff windpipe is not correctly connected with cuff	Correctly connect
	Cuff is not inflated	Stop using the device and contact us
Cuff deflates in short time	Loose cuff	Correctly apply cuff
It can not carry on measurement when press the measurement button		Switch on the power once again and restart a measurement
Power off suddenly when inflating	No use for a long time, the power of batteries can be exhausted owing to the changed temperature	Replace all four batteries with new ones.
Hold the on/off button but can not start the device	Power of batteries can be exhausted	Replace all four batteries with new ones.
	The battery polarities is reversed	Check the battery installation for proper placement of the battery polarities.
Cuff inflation start before press the measurement button or never stop inflating when measuring		Pull out the cuff to deflate. Stop using the device and contact us.
Cuff never deflation		Pull out the cuff to deflate. Stop using the device and contact us.
Air pressure error	No deflation or deflation error or inflation without stop	Pull out the cuff to deflate. Stop using the device and contact us.
	Others	Keep arm, body still, measure again.
No press value displayed or the value unchanged or change erratically when cuff inflated		Pull out the cuff to deflate. Stop using the device and contact us.
Other phenomenon		Switch on the power once again and restart an operation.
		Replace the batteries. If no, please contact us.

Chapter16 Keys and Symbols

Your device may not contain all the following symbols.

Signal	Description	Signal	Description
	Caution: read instructions (warnings) carefully		Follow instructions for use
SIS	Systolic pressure	DIA	Diastolic pressure
INFO	Information	PR	Pulse rate (bpm)
IP22	It means this pulse oximeter is protected against harmful effects of dripping water when tilted at 15°.	CEM	Electromagnetic compatibility
ADU	Adult	P/N	Material code of manufacturer
	Type BF applied parts		Recyclable
	Open the Prompt sound indication		Close the Prompt sound indication
	Lot number		Expiration date
	This way up		Fragile, handle with care
	Keep in a cool, dry place		Atmospheric pressure limit
	Temperature limit		Humidity limit
	Manufacturer		Manufacture Date
	UDI		Unique device identifier
	Batteries Power		Pulse rate (bpm)
	1.No NIBP data to review 2.An indicator of signal inadequacy		1.No Pulse rate 2.An indicator of signal inadequacy
	WEEE disposal		Medical Device compliant with Directive 93/42/EEC
	Serial number		Class II applied
	Authorized representative in the European community		Interface for connecting cuff
	MR Unsafe,can not be used in MRI		Socket for power adapter
	No SpO ₂ Alarms		Artery indicator label
	Medical device		Imported by
	Keep away from sunlight		Product code

Chapter17 Maintenance, Cleaning and Keeping

*Please do obey the precautions and correct operating methods in this user manual. Otherwise, we will not responsible for any fault.

 **Warning** 

Remove the batteries before cleaning. The accessories and main unit must be separated for cleaning. Maintenance is not allowed during device using.

Do not squeeze the rubber tube on the cuff.

 **Caution** 

- High pressure disinfection to the device and accessories is not allowed.
- Do not let water or cleaning agent flow into the socket to avoid device damage.
- Do not soak the device and accessories in liquid.
- If any damage or deterioration of the device and accessories is found, please do not use it.

Maintenance:

- Clean the device and accessories regularly. It is recommended to clean them every one month.
- Before cleaning the device, remove batteries and disconnect it from the AC power. The accessories and main unit must be separated for cleaning. Do not maintain or repair the device during use.
- When cleaning the device, dip a clean cloth in isopropyl alcohol (70%), wring it out fully, and wipe the main unit, cuff and cuff windpipe separately for about 3 minutes, then use the other clean cloth moistened with distilled water, wring it out fully and respectively wipe the main unit, cuff and cuff windpipe for about 2 minutes. Repeat above 5 times until there is no obvious residual cleaning agent. Avoid isopropyl alcohol or water entering the main unit during cleaning. After cleaning, place the product in a dry and ventilated place to dry.
- Visual inspection to make sure the product is thoroughly cleaned. If any residue exists, please repeat the entire process described above.
- The device should be inspected and calibrated periodically (or obey the requirements of the hospital). It is available to inspect in the state specified inspection institution or by professional personal, or you can contact our company.Long press "USER" button in main interface for 5s to enter the calibration interface.

 **Advice** 

- Do not use gasoline, volatile oil, diluent, etc. to wipe the device.
- Do not clean or wet the cuff.

Storage:

 **Advice** 

- Do not expose the device in direct sunlight for long time, otherwise the display screen maybe damaged.
- The basic performance and safety of the device are not affected by the dust or cotton wool in home environment, while the device shall not be placed where with high temperature, humidity or dusty.
- Aged cuff may result in inaccurate measurement, please replace the cuff periodically according to the user manual.
- To avoid device damage, keep the device out the reach of children and pets.
- Avoid the device close to extreme high temperature such as fireplace, otherwise the device performance may be affected.
- Do not store the device with chemical medicine or corrosive gas.
- Do not place the device where there is water.
- Do not place the device where with slope, vibration or impact
- Take the batteries out if the device is not to be used for three months or longer.

Chapter18 NIBP Specification

Name	Electronic Sphygmomanometer		
Display mode	2.8" color LCD Display		
The degree of protection against ingress of water	IP22		
Measurement method	Oscillometric method		
Working mode	Automatic		
Operation mode	Continuous operation		
Cuff pressure range	0-297 mmHg(0.39.6 kPa)		
Overpressure protection	adult mode	295±5 mmHg(39.33±0.67 kPa)	
	pediatric mode	240±5 mmHg(32±0.67 kPa)	
	neonatal mode	145±5 mmHg(19.33±0.67 kPa)	
Measurement Range	pressure	adult	SYS: 30~270 mmHg(4~36 kPa) DIA: 10~220 mmHg(1.3~29.3 kPa)
		pediatric	SYS: 30~235 mmHg(4~31.3 kPa) DIA: 10~195 mmHg(1.3~26 kPa)
		neonatal	SYS:30~135 mmHg(4~18 kPa) DIA: 10~100 mmHg(1.3~13.3 kPa)
		Pulse: 40~240 bpm	
Inflation	adult	160±5 mmHg(21.33±0.67 kPa)	
	pediatric	120±5 mmHg(16±0.67 kPa)	
	neonatal	70±5 mmHg(9.33±0.67 kPa)	
Resolution	Pressure: 1mmHg(0.1 kPa)		
	Pulse: 1 bpm		
Accuracy	Static pressure: ±3 mmHg(±0.4 kPa)		
	Pulse: ±5 bpm or ±5% select larger		
Error	The BP value measured by the device is equivalent with the measurement value of Stethoscopy, perform clinical verification in accordance with the requirements in ISO 81060-2: 2013, whose error meets the followings: Maximum mean error: ±5 mmHg Maximum Standard deviation: 8 mmHg		
Operating temperature/ humidity	+5℃~40℃ 15 %RH~85 %RH(Non-condensing)		
Transport	Transport by general vehicle or according to the order contract, avoid pounded, shake and splash by rain and snow in transportation.		
Storage	Temperature: -20℃~+55℃; Relative humidity: ≤95 %; No corrosive gas and drafty.		
Atmospheric pressure	700 hPa~1060 hPa		
Power supply	4 "AA" alkaline batteries, AC Adapter(AC, 100 V-240 V, optional)		
Rated current	≤600 mA		
Battery life	When the temperature is 23℃, limb circumference is 270 mm, the measured blood pressure is normal, 4 "AA" alkaline batteries can be used about 300 times.		
Dimensions	130(L)*110(W)*80(H) mm		
Unit Weight	300 gram(without batteries)		
Safety classification	Class II equipment (power supplied by power adapter)/Internally powered equipment (power supplied by batteries) Type BF applied part		
Service life	The service life of the device is five years or 10000 times of BP measurement.		
Date of manufacturer	See the label		
Accessories	Standard Configure: Adult Cuff: limb circumference 22-32 cm (upper arm center) four "AA" alkaline batteries,User Manual Optional Configure: four "AA" alkaline batteries AC Adapter Input: voltage: AC 100 V~240 V frequency: 50 Hz/60 H Rated current: AC 150 mA Output: DC 5.0 V±0.2 V 1.0 A Power Adapter cable Cuff the range of limb circumference is 6-11 cm (middle part of upper arm) the range of limb circumference is 10-19 cm (middle part of upper arm) the range of limb circumference is 18-26 cm (middle part of upper arm) the range of limb circumference is 22-30 cm (middle part of upper arm) the range of limb circumference is 22-43 cm (middle part of upper arm) the range of limb circumference is 32-43 cm (middle part of upper arm)		