**Test Monofase Cocaina Strisce Urina**  
**One Step Cocaine Test Strip (Urine)**  
**1 Etape Test de Cocaïne sur Bandelette (Urine)**  
**Einstufen Kokain Teststreifen (Urin)**  
**Prueba de Cocaína en Un Solo Paso en Tira (Orina)**  
**Tira para Teste de Cocaína Em Um Só Passo (Urina)**  
**Τεστ Κοκαΐνης Ταινίας Ενός Σταδίου (Ούρων)**  
**فحص أحادي الطور كوكاينا شرائح بول**

**USO PROFESSIONALE**  
**PROFESSIONAL USE**

**MANUALE D'USO**  
**OPERATOR'S MANUAL**  
**MANUEL D'UTILIZATION**  
**BEDIENUNGSANLEITUNG**  
**MANUAL DE USO**  
**MANUAL DE USO**  
**Εγχειρίδιο χρήσης**  
**دليل للإرشادات**

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### One Step Cocaine Test Strip (Urine)

A rapid, one step test for the qualitative detection of Cocaine metabolites in human urine. For professional IVD (In Vitro Diagnostics) use only.

#### INTENDED USE

The COC One Step Cocaine Test Strip (Urine) is a rapid chromatographic immunoassay for the qualitative detection of Cocaine metabolite, Benzoylcegonine, in human urine at a cut-off concentration of 300 ng/mL. This test will detect other related compounds, please refer to the Analytical Specificity table in this package insert. **This assay provides only a qualitative, preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.**

#### PRINCIPLE

The COC One Step Cocaine Test Strip (Urine) is an immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody. During testing, a urine specimen migrates upward by capillary action. A drug-positive urine specimen will not generate a colored line in the test line region because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

#### REAGENTS

The test contains mouse monoclonal anti-Benzoylcegonine antibody-coupled particles and Benzoylcegonine-protein conjugate. A goat antibody is employed in the control line system.

#### PRECAUTIONS

- For professional IVD (In Vitro Diagnostics) use only. Do not use after the expiration date.
- The test should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test should be discarded according to local regulations.

#### STORAGE AND STABILITY

Store as packaged in the sealed pouch either at room temperature or refrigerated (2-30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

#### SPECIMEN COLLECTION AND PREPARATION

##### Urine Assay

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear supernatant for testing.

##### Specimen Storage

Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

## MATERIALS

### Materials Provided

- Test strips
- Package insert


### Materials Required But Not Provided

- Specimen collection container
- Timer

## DIRECTIONS FOR USE

Allow the test, urine specimen, and/or controls to equilibrate to room temperature (15-30°C) prior to testing.

1. Bring the pouch to room temperature before opening it. Remove the test strip from the sealed pouch and use it as soon as possible.
2. With arrows pointing toward the urine specimen, immerse the test strip vertically in the urine specimen for at least 10-15 seconds. Do not pass the maximum line (MAX) on the test strip when immersing it. See the illustration below.
3. Place the test strip on a non-absorbent flat surface, start the timer and wait for the colored line(s) to appear. Read results at 5 minutes. Do not interpret the result after 10 minutes.

## INTERPRETATION OF RESULTS

(Please refer to the illustration above)

**NEGATIVE:**\* Two lines appear. One colored line should be in the control line region (C), and another apparent colored line should be in the test line region (T). This negative result indicates that the Benzoylcegonine concentration is below the detectable level (300 ng/mL).

\*NOTE: The shade of color in the test line region (T) may vary, but it should be considered negative whenever there is even a faint colored line.

**POSITIVE:** One colored line appears in the control line region (C). No line appears in the test line region. This positive result indicates that the Benzoylcegonine concentration is above the detectable level (300 ng/mL).

**INVALID:** Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

## QUALITY CONTROL

A procedural control is included in the test. A colored line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

**LIMITATION**

1. The COC One Step Cocaine Test Strip (Urine) provides only a qualitative, preliminary analytical result. A secondary quantitative analytical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.<sup>3,4</sup>
2. It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
3. Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
4. A positive result does not indicate level of intoxication, administration route or concentration in urine.
5. A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
6. Test does not distinguish between drugs of abuse and certain medications.
7. Certain foods and food supplements can give a positive result.

**PERFORMANCE CHARACTERISTICS**

**Analytical Sensitivity**

A drug-free urine pool was spiked with Benzoylcegonine at the following concentrations: 0 ng/mL, 150 ng/mL, 225 ng/mL, 300 ng/mL, 375 ng/mL and 450 ng/mL. The result demonstrates >99% accuracy at 50% above and 50% below the cut-off concentration. The data are summarized below:

Benzoylcegonine Concentration (ng/mL)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0	30	30	0
150	-50%	30	30	0
225	-25%	30	30	0
300	Cut-off	30	4	26
375	-25%	30	0	30
450	-50%	30	0	30

**Analytical Specificity**

The following table lists compounds that are positively detected in urine by the COC One Step Cocaine Test Strip (Urine) at 5 minutes.

Compound	Concentration (ng/mL)
Benzoylcegonine	300
Cocaine HCl	780
Cocaethylene	12,500
Ecgonine HCl	32,000

**Cross-Reactivity**

A study was conducted to determine the cross-reactivity of the test with compounds in either drug- negative urine or Benzoylcegonine positive urine. The following compounds show no interference when tested with the COC One Step Cocaine Test Strip (Urine) at a concentration of 100 µg/mL.

**Non Cross-Reacting Compounds**

Acetaminophen	Diazepam	Methadone	Prednisone
Acetophenetidin	Diclofenac	Methoxyphenamine	Procaine
N-Acetylprocainamide	Diflunisal	(±)3,4-Methylenedioxy-amphetamine	Promazine
Acetylsalicylic acid	Digoxin	(±)-3,4-Methylenedioxy-methamphetamine	Promethazine
Aminopyrine	Diphenhydramine	Morphine-3-β -D glucuronide	D,L-Propranolol
Amitypyline	Doxylamine	Morphine Sulfate	D-Propoxyphene
Amobarbital	Ecgonine methylester	Nalidixic acid	D-Pseudoephedrine
Amoxicillin	(-)-α -Ephedrine	Naloxone	Quinidine
Ampicillin	Erythromycin		Quinine
L-Ascorbic acid	β-Estradiol		Ranitidine
D,L-Amphetamine sulfate	Estrone-3-sulfate		Salicylic acid

Apomorphine	Ethyl-p-aminobenzoate	Naltrexone	Secobarbital
Aspartame	Fenoprofen	Naproxen	Serotonin
Atropine	Furosemide	Niacinamide	Sulfamethazine
Benzilic acid	Gentisic acid	Nifedipine	Sulindac
Benzoic acid	Hemoglobin	Norcodein	Temazepam
Benzphetamine	Hydralazine	Norethindrone	Tetracycline
Bilirubin	Hydrochlorothiazide	D-Norpropoxyphene	Tetrahydrocortisone,
(±)-Brompheniramine	Hydrocodone	Noscapine	3-acetate
Caffeine	Hydrocortisone	D,L-Octopamine	Tetrahydrocortisone
Cannabidiol	O-Hydroxyhippuric acid	Oxalic acid	3-(β-D glucuronide)
Cannabinol	p-Hydroxy-	Oxazepam	Tetrahydrozoline
Chloralhydrate	Methamphetamine	Oxolinic acid	Thebaine
Chloramphenicol	3-Hydroxytyramine	Oxycodone	Thiamine
Chlordiazepoxide	Ibuprofen	Oxymetazoline	Thioridazine
Chlorothiazide	Imipramine	Papaverine	D,L-Tyrosine
(±)-Chlorpheniramine	Iproniazid	Penicillin-G	Tolbutamide
Chlorpromazine	(±)-Isoproterenol	Pentobarbital	Triamterene
Chlorquine	Isoxsuprine	Perphenazine	Trifluoperazine
Cholesterol	Ketamine	Phencyclidine	Trimethoprim
Clomipramine	Ketoprofen	Phenelzine	Trimipramine
Clonidine	Labetalol	Phenobarbital	Tryptamine
Codeine	Levorphanol	Phentermine	D,L-Tryptophan
Cortisone	Loperamide	L-Phenylephrine	Tyramine
(-)-Cotinine	Maprotiline	β-Phenylethylamine	Uric acid
Creatinine	Meperidine	Phenylpropanolamine	Verapamil
Deoxycorticosterone	Meprobamate	Prednisolone	Zomepirac
Dextromethorphan			

## BIBLIOGRAPHY

1. Stewart DI, T Inoba, M Ducassen, W Kalow. Clin. Pharmacol. Ther. 1979; 25: 264
2. Ambre J.J. Anal. Toxicol. 1985; 9: 241
3. Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982: 488
4. Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986

## Index of Symbols

	Attention, see instruction for use		Test per kit		Manufacturer
	For <i>in vitro</i> diagnostic use only		Use by		Do not reuse
	Store between 2-30° C		Lot Number		Catalog <b>24540</b>
	Keep away from humidity		Protect against direct sunlight		Read the user instructions

 GIMA Spa  
 Via Marconi, 1 - 20060  
 Gessate (MI) - Italia  
 Imported product

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Acetaminophen	Diazepam	Methadone	Prednisone
Acetophenetidin	Diclofenac	Methoxyphenamine	Procaine
N-Acetylprocainamide	Diflunisal	(±)3,4-Methylenedioxy-	Promazine
Acetylsalicylic acid	Digoxin	amphetamine	Promethazine
Aminopyrine	Diphenhydramine	(±)-3,4-Methylenedioxy-	D,L-Propranolol
Amitriptyline	Doxylamine	methamphetamine	D-Propoxyphene
Amobarbital	Ecgonine methylester	Morphine-3-β -D	D-Pseudoephedrine
Amoxicillin	(-)-α -Ephedrine	glucuronide	Quinidine

Ampicillin	Erythromycin	Morphine Sulfate	Quinine
L-Ascorbic acid	β-Estradiol	Nalidixic acid	Ranitidine
D,L-Amphetamine sulfate	Estrone-3-sulfate	Naloxone	Salicylic acid
Apomorphine	Ethyl-p-aminobenzoate	Naltrexone	Secobarbital
Aspartame	Fenopropfen	Naproxen	Serotonin
Atropine	Furosemide	Niacinamide	Sulfamethazine
Benzilic acid	Gentisic acid	Nifedipine	Sulindac
Benzoic acid	Hemoglobin	Norcodein	Temazepam
Benzphetamine	Hydralazine	Norethindrone	Tetracycline
Bilirubin	Hydrochlorothiazide	D-Norpropoxyphene	Tetrahydrocortisone, 3-acetate
(±)-Brompheniramine	Hydrocodone	Noscapine	Tetrahydrocortisone
Caffeine	Hydrocortisone	D,L-Octopamine	3-(β-D glucuronide)
Cannabidiol	O-Hydroxyhippuric acid	Oxalic acid	Tetrahydrozoline
Cannabinol	p-Hydroxy-	Oxazepam	Thebaine
Chloralhydrate	Methamphetamine	Oxolinic acid	Thiamine
Chloramphenicol	3-Hydroxytyramine	Oxycodone	Thioridazine
Chlordiazepoxide	Ibuprofen	Oxymetazoline	D,L-Tyrosine
Chlorothiazide	Imipramine	Papaverine	Tolbutamide
(±)-Chlorpheniramine	Iproniazid	Penicillin-G	Triamterene
Chlorpromazine	(±)-Isoproterenol	Pentobarbital	Trifluoperazine
Chlorquine	Isoxsuprine	Perphenazine	Trimethoprim
Cholesterol	Ketamine	Phencyclidine	Trimipramine
Clomipramine	Ketoprofen	Phenelzine	Tryptamine
Clonidine	Labetalol	Phenobarbital	D,L-Tryptophan
Codeine	Levorphanol	Phentermine	Tyramine
Cortisone	Loperamide	L-Phenylephrine	Uric acid
(-)-Cotinine	Maprotiline	β-Phenylethylamine	Verapamil
Creatinine	Meperidine	Phenylpropanolamine	Zomepirac
Deoxycorticosterone	Meprobamate	Prednisolone	
Dextromethorphan			

## BIBLIOGRAFIA

1. Stewart DI, Inoba T, Ducassen M, Kalow W. Clin. Pharmacol. Ther. 1979; 25: 264
2. Ambre J. J. Anal. Toxicol. 1985; 9: 241
3. Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982: 488
4. Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986

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Acetaminophen	Diazepam	Methadone	Prednisone
Acetophenetidin	Diclofenac	Methoxyphenamine	Procaine
N-Acetylprocainamide	Diflunisal	(±)3,4-Methylenedioxy- amphetamine	Promazine
Acetylsalicylic acid	Digoxin	(±)-3,4-Methylenedioxy- methamphetamine	Promethazine
Aminopyrine	Diphenhydramine	Morphine-3-β -D glucuronide	D,L-Propranolol
Amitriptyline	Doxylamine	Morphine Sulfate	D-Propoxyphene
Amobarbital	Ecgonine methylester	Nalidixic acid	D-Pseudoephedrine
Amoxicillin	(-)-☼ -Ephedrine	Naloxone	Quinidine
Ampicillin	Erythromycin	Naltrexone	Quinine
L-Ascorbic acid	β-Estradiol	Naproxen	Ranitidine
D,L-Amphetamine sulfate	Estrone-3-sulfate	Niacinamide	Salicylic acid
Apomorphine	Ethyl-p-aminobenzoate	Nifedipine	Secobarbital
Aspartame	Fenoprofen	Norcodein	Serotonin
Atropine	Furosemide	Norethindrone	Sulfamethazine
Benzilic acid	Gentisic acid	D-Norpropoxyphene	Sulindac
Benzoic acid	Hemoglobin	Noscapine	Temazepam
Benzphetamine	Hydralazine	D,L-Octopamine	Tetracycline
Bilirubin	Hydrochlorothiazide	Oxalic acid	Tetrahydrocortisone, 3-acetate
(±)-Brompheniramine	Hydrocodone	Oxazepam	Tetrahydrocortisone 3-(β-D glucuronide)
Caffeine	Hydrocortisone	Oxolinic acid	Tetrahydrozoline
Cannabidiol	O-Hydroxyhippuric acid	Oxycodone	Thebaine
Cannabinol	p-Hydroxy- Methamphetamine	Oxymetazoline	Thiamine
Chloralhydrate	3-Hydroxytyramine	Papaverine	Thioridazine
Chloramphenicol	Ibuprofen	Penicillin-G	D,L-Tyrosine
Chlordiazepoxide	Imipramine	Pentobarbital	Tolbutamide
Chlorothiazide	Iproniazid	Perphenazine	Triamterene
(±)-Chlorpheniramine	(±)-Isoproterenol	Phencyclidine	Trifluoperazine
Chlorpromazine	Isoxsuprine	Phenelzine	Trimethoprim
Chlorquine	Ketamine	Phenobarbital	Trimipramine
Cholesterol	Ketoprofen	Phentermine	Tryptamine
Clomipramine	Labetalol	L-Phenylephrine	D,L-Tryptophan
Clonidine	Levorphanol	β-Phenylethylamine	Tyramine
Codeine	Loperamide	Phenylpropanolamine	Uric acid
Cortisone	Maprotiline	Prednisolone	Verapamil
(-)Cotinine	Meperidine		Zomepirac
Creatinine	Meprobamate		
Deoxycorticosterone			
Dextromethorphan			

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