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*National Institutes of Health
The Practical guide. Identification, Evaluation and Treatment of Overweight and Obesity in Adults (October 2000)

Professional TAPE to measure abdominal circumference

The measuring tape is a device used to measure the distribution of body fat (subcutaneous, visceral and on the torso).

ABS body with patented ergonomic handle, measuring tape made of special plastic material: nondeformable, latex free and double colored. Double colored tape shows risk indications for men and women on two different sides. This device can be used to measure any body circumference through the push button and flexible lock that allows an easy reading of data.

The excellent flexibility of the thermoplastic arm allows the tape to be used in the pediatric field as a cephalometer to measure the circumference of children's head.


This device can also be used in sport activities to accurately measure the arm circumference and to check the hypertrophy of athlete. In the medical field it can assess the loss and gain of lean body mass after traumas or surgical procedures.

**ABDOMINAL CIRCUMFERENCE
RISK OF METABOLIC COMPLICATIONS**

RISK LEVEL	MEN	WOMEN
LEVEL 1 ALERT ZONE 1	≥ 94 cm	≥ 80 cm
LEVEL 2 ACTION ZONE 2	≥ 102 cm	≥ 88 cm

1 WHO Tech Rep Ser 2000;894(3):i-xii, 1-253. Mod da 1
2 Han, TS et al. - BMJ 1995; 311:1401-1405 Mod da 2

**TAPE POSITION
FOR A PRECISE MEASURE
OF ABDOMINAL CIRCUMFERENCE***

1. Patient in straight position.
2. Detect the upper level of hip and tip of right iliac crest.
3. Place the tape horizontally around abdomen at the superior border of right iliac crest.
4. Before taking the measure make sure that tape is tight on the body, the tape must not squeeze the skin and be parallel to the floor.
5. Take the measure after a normal expiration.

*National Institutes of Health
The Practical guide. Identification, Evaluation and Treatment of Overweight and Obesity in Adults (October 2000)

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