www.tecniwork.it

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NUMERO VERDE
800 90 40 55
GRATUITO

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TecniWork®

PODOSCOPIO AP500

A LUCE POLARIZZATA LED
WITH LED POLARISED LIGHT

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Tecniwork has conceived and implemented this podoscope thanks to its long-lasting experience in the field of podiatry. This podoscope is originated from the request of professionals of the field who needed to work with a practical and at the same time efficient and resistant instrument for their job, avoiding the inconveniences of bulky and obsolete appliances.

Tecniwork has then reached its target, developing this innovative device, thanks to the research of new materials and to the application of the most advanced technology.

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Company quality system

Tecniwork has a certified company quality system according to the UNI EN ISO 9001:2008 norm which guarantees the effectiveness of the company organization system, which undergoes controls by an accredited external body auditor, every year.

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GB CONTENTS

Chapter 1 - GENERAL DESCRIPTION p. 14

- 1.1 General Identification data p. 14
- 1.2 General characteristics p. 15
- 1.3 Technical features and voltage p. 15
- 1.4 Classification and reference norms p. 16
- 1.5 Electrical features and voltage p. 17
- 1.6 Components identification Electrical schema p. 18
- 1.7 General safety warning p. 18

Chapter 2 - INSTALLATION p. 18

- 2.1 Transport and stocking p. 18
- 2.2 Conditions for installation p. 19
- 2.3 Electric connection p. 21

Chapter 3 - INSTRUCTIONS FOR USE p. 21

- 3.1 Use of the podoscope p. 21

Chapter 4 - MAINTENANCE p. 21

- 4.1 Ordinary maintenance p. 21
- 4.2 Dismantling p. 22
- 4.3 Technical assistance service p. 22

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Chapter 1 – GENERAL DESCRIPTION

1.1 GENERAL IDENTIFICATION DATA

This instruction handbook is considered as an essential document of podoscope AP 500. Therefore they must not be separated in case of cessation to third parties, in accordance with EC Directive 93/42/EC and following norms.

This handbook gives the operator all the relevant information to carry out the following operations correctly:

- INSTALLATION
- USAGE
- MAINTENANCE

Instructions for a correct installation are contained in the par. 2.2 CONDITIONS FOR THE INSTALLATION.

FOLLOW THE INSTRUCTIONS FOR USE

READ CAREFULLY THIS MANUAL FOR THE INSTALLATION
THE USAGE AND THE MAINTENANCE OF THE PODOSCOPE

COMPANY NAME

COMMERCIAL TRADEMARK: TECNIWORK

TYPOLOGY: PODOSCOPE

MODEL: AP500

MANUFACTURER: TECNIWORK S.p.A. - Via R. Benini, 8 - 50013 Campi Bisenzio - FI - Italy

IDENTIFICATION AND LABELS - SYMBOLS EXPLICATION

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ON EACH DEVICE YOU WILL FIND A SIMILAR PLATE WITH THE MANUFACTURER'S IDENTIFICATION.

PROTECTION AGAINST INDIRECT CONTACTS:
DEVICE OF CLASS II

DEVICE WITH APPLIED
PART TYPE B

FOLLOW THE INSTRUCTIONS FOR USE

ATTENTION: DANGER

1.2 DESCRIPTION AND GENERAL CHARACTERISTICS

The particular system of light diffusion of the PODOSCOPE AP500 enables you to study the patient's foot plant and to mark the parts with the heaviest pressure.

Its structure is very simple: the methacrylate structure is an excellent light conductor and it is also endowed with an inner light source: in this way you can discover variations of colour intensity on the foot plant, corresponding to higher or lower plant pressure.

The base is of ABS.

A separable net cable with plug, a jack, and a led strip are its electric components.


1.3 TECHNICAL FEATURES AND VOLTAGE

MATERIAL	Transparent methacrylate
TENSION	220-240 V single-phase
FREQUENCY	50/60 Hz
POWER	7,5 VA
ILLUMINATION	Led Strip Light
WAVELENGTH	520 nm
RISK CLASS EN 62471	Exempt
NR. OF INCLUDED PLUGS	Nr. 1
2 FUSES	5 x 20 mm: 2 x T 0,5 A - 250 V

CONDITIONS OF USE
temperature: from 10°C up to +40°C
umidity: from da 20% up to 85%

STOCKING CONDITIONS
temperature: from 0°C up to 60°C
umidity: from da 10% up to 90%

1.4 CLASSIFICATION AND REFERENCE NORMS

The podoscope AP500 allows to evaluate the foot print, supplying secondary indirect information on the heel alignment and on the toes situation. It allows to evaluate the signs and the symptoms of a particular condition of the foot, in order to formulate a judgement on the health conditons and establish a therapy or a chiropody treatment.

The device is designed and realized according to the law prescriptions concerning electromedical equipment, and it's therefore indicated for medical purposes and used in clinics, consultation rooms and sports centres by operators with knowledge of biomechanics, like doctors in sport medicine, orthopedic specialists and physiotherapists.

REFERENCE HARMONIZED NORMS:

- CEI EN 60601-1 "Safety of electromedical devices".
- CEI EN 60601-1-2 "Collateral harmonized norm for electromedical devices
Electromagnetic compatibility".

CLASSIFICATION ACCORDING TO SAFETY NORM EN 60601-1:

- Portable device
- Class II
- Type B
- Ordinary protection against humidity
- Do never use with inflammable anaesthetics
- Voltage: 220-240 V ~ 50/60Hz.

The device highlights areas of the foot plant with heaviest pressure thanks to different intensity of light. In order to visualize them the patient just needs to step on the podoscope barefoot. The particular features of this device grant adequate protection and isolation against direct or indirect electric shocks and against electric or mechanic surcharge.

The thermic power implied in this podoscope avoid dangerous overheating for materials and components.

Electric isolation and well-balanced distances between the surfaces avoid the formation of electric arches. In addition to this, the podoscope has a solid mechanical structure and an excellent stability on the ground. It has no sharp edges, therefore there exists no danger of mechanical nature both for operators and patients. Each part of the inner circuits are assembled and fixed in such a solid way so as to stand up to the often adverse transport conditions.

As to indirect shocks, the device is classified in Class II.

As to current leakage it can be classified in class B.

No essential performance is defined for this device.

1.5 ELECTRICAL FEATURES OF THE PODOSCOPE

The conventional led power supply functions with single-phase current (220-240 V - 50/60 Hz); the mark is ILC NanoLED.


The power cord - H05VV-F - 2 mt (6.56 ft) long, section 2 x 0.75 bears the mark IMQ, has a two-pole plug and it is connected to the podoscope through a connector.

The internal connectors of the power cord satisfy the norm CEI 20-20 (type H05V).


1.6 COMPONENTS IDENTIFICATION

PODOSCOPE AP 500

- | | |
|--|-------------------------|
| 1- Transparent methacrylate surface | 4- Power cord jack |
| 2- Basement in ABS (acrylonitrile butadiene styrene with wiring) | 5- On/Off switch |
| 3- Mirror | 6- Fuse holder |
| | 7- Identification plate |

ELECTRIC SCHEMA

1.7 GENERAL SAFETY WARNING

As far as safety is concerned, this podoscope has undergone all the necessary checks and inspections as foreseen. In order to respect the safety prerequisites during the whole life of your device, we suggest that you always take the following precautions:

- A) Always unplug the podoscope during maintenance operations
- B) Do not move the podoscope while using it
- C) Each time before the usage, check the functioning of the cable and of its electric parts.
- D) Operators must not intervene or take any initiatives on the podoscope which are not foreseen in this handbook.

ATTENTION: DO NOT CHARGE THE PODOSCOPE WITH LOADS OVER 135 KG (297 LB) IN ORDER TO AVOID ANY RISK OF BREAKING

ATTENTION: NO CHANGE OF THE DEVICE IS ALLOWED.

ATTENTION: NOT COMPLYING WITH WARNINGS WILL RELIEVE TECNIWORK COMPANY FROM ANY KIND OF RESPONSIBILITY

Chapter 2 – INSTALLATION

2.1 TRASPORT AND STOCKING

Before shipping the podoscope AP500 is covered by a film and packed in a cardboard with protective panels inside. The box contains the assembled podoscope and its cable (which must be obviously connected to the apparatus before using it); the caoutchouc ref. AP504 is not included.

In case of a longer storage, please store the whole package clean and dry.

2.2 CONDITIONS FOR INSTALLATION


- 1- Unpack the podoscope completely.
- 2- Disposal of packing elements (plastic foams, polystyrene, polyethylene, cardboard) must follow your local waste disposal guidelines.
- 3- Verify that the podoscope has not been damaged during transport. In case, please contact the manufacturer immediately.

IMPORTANT: Please take notice of the fact that light stripes and/or microscopic luminescent dots on the transparent surface of the podoscope are due to the particular nature of its material and do not impair the equipment functionality.

RIGHT POSITIONING AND SAFETY INSTALLATION OF THE PODOSCOPE

Position the podoscope so as to leave enough room around it to be able to safely operate and to connect and disconnect easily the device to the power supply network.

We suggest to attach a handle on the wall as indicated in the scheme below to make easier for patients to go on and down the podoscope.

ELECTROMAGNETIC COMPATIBILITY (EMC)

The podoscope needs special precautions concerning electromagnetic compatibility (EMC) and must be installed and operated according to the EMC information contained in this manual.

In order to avoid possible risks of electromagnetic interferences, do not use radio frequency portable or mobile devices near the podoscope.

Generally the podoscope should not be used near or together other equipments; in the case this would not be possible, it's necessary to observe its operativity to verify the normal functioning.

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSIONS

The equipment AP 500 is intended for use in the electromagnetic environment specified below.
The customer or the user of the equipment AP 500 should assure that it is used in such an environment.

EMISSIONS TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT – GUIDANCE
RF emissions CISPR 11	Group1	The equipment “AP 500” must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
RF emissions CISPR 11	Class B	The “AP 500” is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Compliant	
Voltage fluctuations / Flicker emissions IEC 61000-3-3	Compliant	

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY

The equipment AP 500 is intended for use in the electromagnetic environment specified below.
The customer or the user of the equipment AP 500 should assure that it is used in such an environment.

IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT – GUIDANCE
Electrostatic discharge (ESD) IEC 61000-4-2	±2 kV contact ±4 kV contact ±6 kV contact ±2 kV air ±4 kV air ±8 kV air	±2 kV contact ±4 kV contact ±6 kV contact ±2 kV air ±4 kV air ±8 kV air	Floor should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at less 30 %.
Electrical fast transient/burst IEC 61000-4-5	±2 kV for power supply lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for cycle 70 % U_T (30 % dip in U_T) for cycle <5 % U_T (>95 % dip in U_T) for 5 s	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for cycle 70 % U_T (30 % dip in U_T) for cycle <5 % U_T (>95 % dip in U_T) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of AP 500 requires continued operation during power mains interruptions, it is recommended that the AP 500 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic field should be at levels characteristic of a typical location in a typical commercial or hospital environment.

2.3 ELECTRIC CONNECTION

Connect your podoscope verifying that the indications on the plate correspond to the local electric voltage, then plug it into the socket which must be supplied with a surcharge safety device.

The use of accessories, transducers and cables, except the ones sold by the producer as spare parts, can increase the emissions and decrease the EM (electromagnetic) device immunity.

Chapter 3 – INSTRUCTIONS FOR USE

3.1 USE OF THE PODOSCOPE

Plug the podoscope in to light it up, then switch on (5) to position I.

Let the patient step on the transparent basis of the podoscope and make sure that the feet are right in the centre of it. Then proceed with a sight examination on the mirror.

If have not been supplied with article nr. AP504 (caoutchouc) the examination can be performed only if the patient is barefoot, otherwise lay the caoutchouc on the transparent surface. This material is also useful to check the validity of arch support treatments.

WARNING: DON'T LET STEP ON PATIENTS WITH INJURED FOOT SKIN

At the end of the examination, press (5) again to put it back to the 0 position and turn it off.

IMPORTANT: PLEASE HELP ELDERLY PATIENTS STEP ON THE PODOSCOPE AND REMAIN STEADY.

Chapter 4 – MAINTENANCE

4.1 ORDINARY MAINTENANCE

The ordinary maintenance described in this paragraph can be carried out directly by the operator.

- 1) Disinfect the surface of the podoscope after each use with disinfectants indicated for this material: **never use alcohol nor solvents for this purpose.**
- 2) Clean the podoscope at least once a week, always after having unplugged it. Use clean clothes, slightly moistened with water and detergent for glasses. **Never use alcohol nor solvents for this purpose.**
- 3) The frequency of the checks on internal wiring may vary according to external conditions (i.e. presence of dust) and repeated usage.
- 4) As a general rule, we recommend a sight inspection every six months or at least within one year after the purchase of the podoscope, to ascertain whether there is no wear and tear in the cable or in the electric connections.

If you have to replace one of the fuses, open the fuse holder placed in the position indicated in the picture on page 17. Then get out the blown fuse with a suitable instrument and replace it with a T 0,5 A – 5x20 250V fuse. If you have to replace the power cord, the new one must have the same characteristics and certification of the original one. In the case of troubles or needs to replace spare parts, always contact the producer.

WARNING: ALL MAINTENANCE SHOULD BE DONE ONLY AFTER DISCONNECTING THE DEVICE FROM THE ELECTRICAL MAINS.

LED REPLACEMENT

In case of troubles with LED, please contact the Technical Assistance Service. Do not make direct operations, as Warranty shall become void.

WARNING: ALL THE OPERATIONS TO REPLACE LED SHOULD BE MADE ONLY BY THE AUTHORIZED TECHNICAL ASSISTANCE

4.2 DISMANTLING

CAUTION! THE WRONG DISPOSAL OF ELECTRIC AND ELECTRONIC EQUIPMENT MAY INVOLVE SANCTIONS.

The symbol on the plate of the device indicates the separated collection of electric and electronic equipment. At the end of life of the device, do not dispose it as mixed solid municipal waste, but dispose it referring to a specific collection centre located in your area or returning it to the distributor, when buying a new device of the same type to be used with the same functions. This procedure of separated collection of electric and electronic devices is carried out forecasting a European environmental policy aiming at safeguarding, protecting and improving environment quality, as well as avoiding potential effects on human health due to the presence of hazardous substances in such equipment or to an improper use of the same or of parts of the same.

4.3 TECHNICAL ASSISTANCE SERVICE

For any inquiry, please call:

FAX: +39 055 8991.801

TEL: +39 055 8991.71

Chapter 5 – WARRANTY

The podoscope AP500 has a 12-month guarantee from the date of shipping; however, the manufacturer cannot be held responsible in case of wrong installation or maintenance, negligence, misuse or disregard of the norms and warnings contained in this handbook or in case of unduly operations by unauthorised assistance personnel.

This warranty gives you the right to prompt replacement of damaged components, parts or material which may result faulty.

Custom duties, import/export costs and charges covering the supply of spare parts to our factory are not included in the guarantee.

The guarantee will expire after 12 months form the shipping date: therefore any substitution or assistance will be charged according to the current prices.

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in accordance with the European Directive 93/42/CEE and subsequent modifications and integrations

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declares under its responsibility that the product

Marca:	TECNIWORK
Modello:	AP500
S/N:	

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It is in compliance with the prescription of the communitarian directive 93/42/CEE,

recepita in Italia con D.L. n° 46/97 e D.Lgs 37/10,

accepted in Italy with the D.L. n° 4697 and DLgs 37/10

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and subsequent amendments and additions

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
it has been classified in class I and has been declared in conformity with Medical Device Directive according to Annex VII

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Moreover, it declares that the device has been planned and built in conformity with the technical norms:

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