ITALIANO

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ENGLISH

ANTI-DECUBITUS LINE IN SILICONE HOLLOW FIBRE

WARNINGS

All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

Before carrying out any operation, visually check the integrity of the device. If the product shows any abnormalities, immediately contact the point of sale.

Before carrying out any operation, visually check the integrity of the device. If the product shows any abnormalities, immediately contact the point of sale.

Erregi S.a.s. is not liable in case of improper use of the device.

DESCRIPTION

Cushions for the prevention of decubitus in 100% cotton, with padding made from silicone hollow fibre enclosed in a cotton cover with Velcro fastening system. The high resilience and breathability allow an effective action in the prevention of the onset of pressure ulcers.

INDICATIONS

Decubitus prevention for bedridden patients in order to decrease pressure peaks, reduce shear and friction forces and promote air circulation.

PRECAUTIONS FOR USE

Medical device indicated for the prevention of decubitus, to be used following a clinical evaluation carried out by a specialist doctor. Place the cushion on the wheelchair seat, making sure that the Velcro fasteners, where present, are positioned on the backrest side. Then secure the cushion to a rigid wheelchair structure using the Velcro fasteners. It is therefore recommended to regularly monitor the affected area, to move the patient frequently, and to pay particular attention to hygiene.

For this purpose, it is advisable to consult the protocols and guidelines for the prevention of pressure ulcers. Always wear personal protective clothing when using the device. Non-fireproof device.

CONTRAINDICATIONS

No contraindications were found except in cases of ascertained sensitivity to the component materials.

MATERIAL

Padding: Silicone hollow fibre

Cover: 100% cotton

WASHING INSTRUCTIONS

The manufacturer is not liable for any changes in the technical features of the aid in the event of failure to observe the washing methods

STORAGE

Store away from sources of light, heat and humidity.

DISINFECTION

Use medical surgical products.

DISPOSAL

Dispose of the device in compliance with the current regulations on environmental protection and recycling.

	Product code		Medical Device
	Lot number		Manufacturer
	Medical Device compliant with Regulation (EU) 2017/745		Consult instructions for use
	Keep away from sunlight		Keep in a cool, dry place
	Date of manufacture		Not made with natural rubber latex
	Unique Device Identifier		

DO NOT remove the label from the Medical Device, as it allows tracing the supplied product (Manufacturer, UDI, batch number, product code).

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

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