MASCHERA IN SILICONE BESMED BESMED SILICONE MASK

Manuale d'uso - User manual

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Model:

SM-86503 - Transparent Silicone Mask, Size 0, Infant

SM-86513 - Transparent Silicone Mask, Size 1, Infant Large

SM-86523 - Transparent Silicone Mask, Size 2, Child

SM-86131 - Silicone Mask, Size 3, Child Large

SM-86141 - Silicone Mask, Size 4, Adult

SM-86151 - Silicone Mask, Size 5, Adult Large

Product Description

Besmed Silicone mask is designed for connection to manual resuscitators for use as an adjunct to artificial respiration and cardiopulmonary resuscitation, automatic resuscitators, and ventilators. Six different sizes are available, from infant to large adult. Anti-slip design for easy holding and an adjustable soft-cuff for promoting better seal and patient comfort.

Intended Use

Besmed Silicone mask is intended to be used with manual resuscitators, automatic resuscitators, and ventilators. It provides a seal such that positive pressure from the positive pressure source is directed into the patient's nose and mouth.

Besmed Silicone mask Technical Specifications

Product Name	Silicone Mask	
Manufacturer	Besmed Health Business Corp.	
Reference No.	SM-86503, SM-86513, SM-86523	SM-86131, SM-86141, SM-86151
Representative Photo		
Patient category	Adult/Child/Infant	

Material		Silicone		
Color		Transparent		
Component type		Accessory		
Device type		Reusable		
Patient Usage type		Single patient reuse		
Patient Weight		Size 0	<2 kg	
		Size 1	2 ~10 kg	
		Size 2	10 ~25 kg	
		Size 3	25 ~60 kg	
		Size 4	60 ~100 kg	
		Size 5	>100 kg	
Specifications	Engagement of 22 mm latching sockets to cones	Did not disconnected at (50 ± 5) N for 10 s		
	Leakage from 22 mm latching sockets	No signs of leakage at (8 ± 0.5) kPa		
	Drop procedure for 22 mm latching sockets	No signs of damage		
Package		Non-sterile packages		
Connection diameter		22 mm		
Operation temperature		25±5 °C, Room Temperature		
Storage temperature		-20°C ~ 60 °C		
Relative humidity		30 - 70% RH		
Shelf life		5 years		
Attestation		Test Item	Guideline	Result
			Conical connectors	ISO 5356-1
		Cytotoxicity	EN ISO 10993-5:2009	PASS
		Sensitization	EN ISO 10993-10:2023	PASS
		Skin Irritation	EN ISO 10993-23:2021	PASS
		Usability	IEC 62366-1	PASS

Intended users

Medical personnel received the training of respiratory care.

Indication

Besmed Silicone mask is used with manual resuscitators ventilators, with an appropriate perimeter of Besmed Silicone mask reduces the pressure sensitivity and enhances facial contour with comfort.

Instruction for use

1. Check for clean airflow passage prior to use.
2. Connect to breathing circuit or resuscitation with appropriate fitting.
3. Set the gas to the proper flow and check for gas flow through the device.
4. Position the elastic strap behind the neck.
5. Make sure all connections are compatible and fit firmly
6. Gently pull the ends of the strap until the mask is secure.

Contraindications

Besmed Silicone mask is possibly not suitable for the use to the patients with recent facial, oral, skull surgery, or trauma, and patients who have claustrophobia, vomiting, or active GI bleeding. Inability to maintain an adequate seal between the mask and the patient's face, which may result in inadequate oxygenation.

Risks/Side Effects

1. Inappropriate size of mask and strap tension may cause discomfort.
2. Mask related complication includes claustrophobia, ulceration, facial rash, or itching.
3. Cramped hands and tired muscles and limits the user's ability to perform other tasks especially in obese and patients with a larger head

Cleaning and Sterilization**1. Cleaning**

1.1 Hand wash the mask thoroughly to remove adhering substances from each component while soaking in a commercially available, anionic detergent such as ANIOSYME DD1 for 10 mins. Use 1:200 ratio of dilution in tepid water for each mask keeping the temperature 35 °C.

1.2 Rinse the mask with 5 liters of water and proceed to the sterilization step.

2. Sterilization

2.1 Sterilize the mask, SM-86503, SM-86513, and SM-86523 with (Autoclave max. 134°C), for 20 mins.

2.2 Sterilize the mask, SM-86131, SM-86141, and SM-86151 with (Autoclave max. 121°C), for 20 mins.

2.3 Reprocess limit is 40 cycles max.

2.4 Air dry out of direct sunlight for 25 min. Make sure the masks are dry before use. If not air-dried, please extend the air-drying time.

2.5 Inspect all components; if any tears or any damage, please replace them.

Warning

- Read this instruction before use.
- Do not use if the package is open or damage or contaminated.
- Any skin redness, itching or any abnormal situation happens, stop using this mask.
- Do not use any alcohol, bleaching chemical or cleaning solution that contains conditioners or lubricant.
- Any serious incident that has occurred in relation to the device should be reported to Besmed or its Authorized Representative/competent authorities and the Regulatory Authorities of the country.

Disposal

Dispose of in accordance with local regulations

LEGEND OF SYMBOLS

	Medical Device compliant with Directive 93/42/EEC		Product code
	Imported by		Manufacturer
	Please read instructions carefully		Date of manufacture
	Lot number		Disposable device, do not re-use
	Expiration date (see box / package)		Keep in a cool, dry place
	Keep away from sunlight		Authorized representative in the European community
	Not made with natural rubber latex		No-DEHP formulation

 A symbol showing a thermometer with a vertical line and a horizontal line at the top, indicating a temperature limit.	Temperature limit	 A triangular warning symbol with the text 'NON STERILE' inside.	Non-sterile
 A square symbol with the letters 'MD' inside, indicating a medical device.	Medical Device		

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.