Abbott

NycoCard™

U-ALBUMIN

EN	5-8	EL	21-24	HU	37-40	PL	53-56
CS	9-12	ES	25-28	IT	41-44	PT	57-60
DA	13-16	FI	29-32	NL	45-48	RO	61-64
DE	17-20	FR	33-36	NO	49-52	RU	65-68

REF 1116819

1116930 Rev. A 2019/02

SYMBOLS/ SYMBOLY/ SYMBOLER/ SYMBOLE/ ΣΥΜΒΟΛΑΙ/ SÍMBOLOS/ SYMBOLIT/ SYMBOLES/
JELMAGYARÁZAT/ SIMBOLI/ SYMBOLEN/ SYMBOLER/ SYMBOLE/ SÍMBOLO/ SIMBOLURI/
СИМВОЛЫ

Conformity to the European directive 98/79/EC on *in vitro* diagnostic medical devices
Splnění evropské směrnice 98/79/ES o zdravotnických prostředcích pro diagnostiku *in vitro*
Overensstemmelse med EU-direktivet 98/79/EC om *in vitro* diagnostisk udstyr til medicinsk brug
Konformität mit der europäischen Richtlinie 98/79/EC über medizinische Geräte zur *In-vitro*-Diagnostik
Συμμορφώνεται με την ευρωπαϊκή οδηγία 98/79/EK περί ιατροτεχνολογικών βοηθημάτων που χρησιμοποιούνται στη διάγνωση *in vitro*
Conformidad con la Directiva europea 98/79/CE sobre dispositivos médicos para el diagnóstico *in vitro*
Yhdenmukainen Euroopan direktiivin 98/79/EC kanssa koskien *in vitro* diagnostisia lääkinällisiä laitteita
Conformité avec la directive européenne 98/79/CE sur les appareils médicaux de diagnostic *in vitro*
Az Európai Parlament és a Tanács 98/79/EK irányelve az *in vitro* diagnosztikai orvostechnikai eszközökről
Conformità con la direttiva europea 98/79/CE sui dispositivi medicali per la diagnostica *in vitro*
Voldoet aan de Europese richtlijn 98/79/EG betreffende medische hulpmiddelen voor *in-vitro*diagnostiek
Overensstemmelse med EU-direktiv 98/79/EC for *in vitro* diagnostisk medisinsk udstyr
Zgodność z dyrektywą 98/79/WE w sprawie wyrobów medycznych używanych do diagnostyki *in vitro*
Conformidade com a Diretiva Europeia 98/79/CE de dispositivos médicos de diagnóstico *in vitro*
Conformitate cu directiva europeană 98/79/CE privind dispozitivele medicale pentru diagnostic *in vitro*
Знак соответствия европейской директиве 98/79/EC по медицинским устройствам, предназначенным для диагностики *in vitro*

Lot number
Číslo šarže
Lot-nummer
Chargen-Nr.
Αριθμός παρτίδας
Número de lote
Eränumero
Numéro de lot
Gyártási tételszám
Numero di lotto
Partijnummer
Lotnummer
Numer partii
Número de lote
Număr de lot
номер партии

Caution
Varovani
Forsigtig
Achtung
Προσοχή
Advertencia
Varoitus
Faire attention
Figyelmeztetés
Attenzione
Opgeet
Advarsel
Ostrzeżenie
Atenção
Avertizare
Предупреждение

In Vitro diagnostic medical device
Diagnostické lékařské zařízení pro testy *in vitro*
In vitro diagnostisk medicinsk udstyr
Medizinisches Gerät zur *In-vitro*-Diagnostik
Ιατροτεχνολογικό βοήθημα που χρησιμοποιείται στη διάγνωση *In Vitro*
Dispositivo médico para el diagnóstico *in vitro*
In vitro diagnostinen lääkinällinen laite
Dispositif médical de diagnostic *in vitro*
In vitro diagnosztikai orvosi eszköz
Dispositivo medico per la diagnostica *in vitro*
In-vitro-diagnostisch medisch apparaat
Medisinsk udstyr til *in vitro* diagnostikk
Wyrob medyczny do diagnostyki *in vitro*
Dispositivo médico de diagnóstico *in vitro*
Dispozitiv medical pentru diagnostic *in vitro*
Медицинское устройство для диагностики *in vitro*

Catalogue number
Katalogové číslo
Katalognummer
Bestellnummer
Αριθμός καταλόγου
Número de catálogo
Luettelonumero
Numéro de référence du catalogue
Katalógusszám
Numero di catalogo
Catalogusnummer
Katalognummer
Numer katalogowy
Número do catálogo
Număr de catalog
Номер по каталогу

2°C 8°C

Storage temperature 2-8°C (36-46°F)
Skladovací teplota 2-8 °C (36-46 °F)
Opbevaringstemperatur 2-8°C (36-46°F)
Lagertemperatur 2-8°C (36-46°F)
Θερμοκρασία φύλαξης 2-8°C (36-46°F)
Temperatura de almacenamiento 2-8°C (36-46°F)
Säilytyslämpötila 2-8°C (36-46°F)
Température de conservation 2-8°C (36-46°F)
Tárolási hőmérséklet 2-8 °C (36-46 °F)
Temperatura di conservazione 2-8°C (36-46°F)
Bewaartemperatuur 2-8°C (36-46°F)
Lagringstemperatur 2-8°C (36-46°F)
Temperatura przechowywania 2-8°C (36-46°F)
Temperatura de conservação 2-8°C (36-46°F)
Temperatură de depozitare 2-8°C (36-46°F)
Температура хранения 2-8 °C (36-46 °F)

Consult the instruction for use
Viz návod k použití
Se brugervejledningen
Bitte die Gebrauchsanleitung beachten
Συμβουλευτείτε τις οδηγίες χρήσης
Consulte las instrucciones de uso
Katso lisätietoja käyttöohjeesta
Consulter le mode d'emploi
További információt a használati útmutatóban talál
Consultare le istruzioni per l'uso
Raadpleeg de gebruiksaanwijzing
Se bruksanvisningen
Sprawdź w instrukcji użytkowania
Consulte as instruções de utilização
Consultați instrucțiunile de utilizare
См. инструкцию по использованию

Contents sufficient for 24 tests
Obsah postačuje pro 24 testů.
Indholdet er tilstrækkeligt til 24 test
Inhalt ausreichend für 24 Tests
Περιεχόμενα επαρκή για 24 δοκιμές
Contenido suficiente para 24 pruebas
Sisältö riittää 24 testimäärälle
Contenu suffisant pour 24 dosages
24 vizsgálatra elegendő tartalom
Contenuto sufficiente per 24 test
Inhoud voldoende voor 24 testen
Innholdet rekker til 24 tester
Zawiera ilość wystarczającą na 24 badań
Conteúdo suficiente para 24 testes
Conține o cantitate suficientă pentru 24 teste
Содержимого достаточно для проведения 24 анализов

Expiry date (year-month-day)
Datum expirace (rok-měsíc-den)
Udløbsdato (år-måned-dag)
Verfallsdatum (Jahr-Monat-Tag)
Ημερομηνία λήξης (έτος-μήνας-ημέρα)
Fecha de caducidad (año-mes-día)
Viimeinen käyttöpäivämäärä
(vuosi-kuukausi-päivä)
Date de péremption (année-mois-jour)
Lejárati dátum (év-hónap-nap)
Data di scadenza (anno-mese-giorno)
Uiterste gebruiksdatum (jaar-maand-dag)
Utløpsdato (år-måned-dag)
Data przydatności (rok-miesiąc-dzień)
Prazo de validade (ano-mês-dia)
Data expirării (an-lună-ziuă)
Срок годности (год-месяц-день)

Manufacturer
Výrobce
Producent
Hersteller
Κατασκευαστής
Fabricante
Valmistaja
Fabricant
Gyártó
Produttore
Fabrikant
Produsent
Producent
Fabricante
Producător
Изготовитель

Date of Manufacture (year-month-day)
Datum výroby (rok-měsíc-den)
Produktionsdato (år-måned-dag)
Herstellungsdatum (Jahr-Monat-Tag)
Ημερομηνία κατασκευής
(έτος-μήνας-ημέρα)
Fecha de fabricación (año-mes-día)
Valmistuspäivämäärä (vuosi-kuukausi-päivä)
Date de fabrication (année-mois-jour)
Gyártás dátuma (év-hónap-nap)
Data di fabbricazione (anno-mese-giorno)
Productiedatum (jaar-maand-dag)
Produktsjonsdato (år-måned-dag)
Data produkcji (rok-miesiąc-dzień)
Data de fabrico (ano-mês-dia)
Data fabricației (an-lună-ziuă)
Дата изготовления (год-месяц-день)

For professional near-patient testing and laboratory use. For use with the NycoCard™ Reader II.

⚠ The NycoCard Reader II is calibrated for use with the NycoCard Tests only.

PRODUCT DESCRIPTION

Intended use

NycoCard™ U-Albumin is a rapid *in vitro* test for measurement of low albumin concentrations in human urine.

Test principle

NycoCard U-Albumin is a solid phase, sandwich-format, immunometric assay. The test device contains a membrane coated with immobilised albumin specific monoclonal antibodies. When the diluted sample is applied to the test device, the sample flows through the membrane, and immobilised antibodies on the membrane capture the albumin molecules. Albumin trapped on the membrane will bind the gold-antibody conjugate then added, in a sandwich-type reaction. Unbound conjugate is removed from the membrane by the washing solution. The paper layer underneath the membrane absorbs excess liquid. Due to the bound gold particles the membrane appears red with colour intensity proportional to the albumin concentration of the sample. The colour intensity is measured quantitatively by using the colour densitometer NycoCard Reader II.

Kit contents, 24 test kit

TD/Test Device	1 x 24 units
Plastic device containing a membrane coated with monoclonal anti-albumin anti- bodies.	
R1/Dilution Liquid	1 x 24 x 1.0 mL
Phosphate buffer (pH 5.6), organic solvent (<10%) and a small amount of a yellow pigment.	
R2/Conjugate	1 x 2.0 mL
Borate buffered solution containing monoclonal anti-albumin antibodies labelled with ultra-small gold particles.	
R3/Washing Solution	1 x 2.0 mL
Phosphate buffered NaCl solution (pH 7.4).	

Material required

(not supplied with the kit)

- Pipette (50 µL) and pipette tips for dilution of sample and application of diluted sample, R2 and R3
- NycoCard Reader II for measurement of the test result

Warnings and precautions

- For *in vitro* diagnostic use.
- R1, R2 and R3 contain sodium azide (<0.1%), a toxic agent. The organic solvent in the R1/dilution liquid can be irritating to the skin, throat and eyes.

TEST CHARACTERISTICS

Analytical specificity

Monoclonal antibodies specific to human albumin are used in the test. No other human urine components have been found to cross-react with albumin in the NycoCard U-Albumin Test System.

Standardization

NycoCard U-Albumin is calibrated with internal urine standards. These standards are assayed against ERM®-DA470 reference preparation.

Measuring range

Measuring range albumin: 5-200 mg/L

Measuring interval: 1 mg/L

Precision

In controlled laboratory testing, a precision of <10% expressed by the coefficient of variation (CV) is usually obtained.

Limitations of the test

No interference is observed from glucose (50 mmol/L), creatinine (60 mmol/L), nitrite (10 mmol/L), NaCl (500 mmol/L), acetone (2.4 g/L), IgA (0.5 g/L), IgG (0.5 g/L), hemoglobin (0.05 g/L), myoglobin (0.5 g/L), β_2 -microglobulin (250 mg/L), bilirubin (50 mg/L), urea (200 g/L) or pH.

Whole blood HSA (human serum albumin) in the sample will give falsely elevated results at hemoglobin concentrations >0.5 mg/L in the urine sample.

STABILITY AND STORAGE

Unopened kits

The expiry date of the kit applies to storage at 2-8°C in original container. Exposure of kit components to temperatures above 25°C should be avoided. Do not freeze.

Opened kits

TD/Test devices should be stored refrigerated (2-8°C) and can be used directly from the refrigerator without equilibration to room temperature. The test devices are stable for 8 weeks when stored at room temperature during the working day. The sealing foil should be removed immediately before use.

R1/Dilution liquid can be stored refrigerated or at room temperature (2-25°C). The dilution liquid must be equilibrated to room temperature before use.

R2/Conjugate should be stored refrigerated (2-8°C) and can be used directly from the refrigerator without equilibration to room

temperature. The conjugate is stable for 8 weeks when stored at room temperature (15-25°C).

R3/Washing solution can be stored refrigerated or at room temperature (2-25°C). The washing solution can be used directly from the refrigerator without equilibration to room temperature.

Urine samples can be stored refrigerated (2-8°C) up to 14 days without significant change in measured albumin concentration. Frozen urine samples can be analyzed with the NycoCard U-Albumin Test System. The frozen samples are stable for 12 weeks at -20°C and can be frozen and thawed once only.

Urine samples, diluted with R1 can be stored for 14 days before analysis with NycoCard U-Albumin, refrigerated or at room temperature (2-25°C).

TEST PROCEDURE

Important procedural notes!

- Do not interchange components from different kit lots.
- Do not touch the test device membrane with the pipette tip.
- Change pipette tip between each pipetting step.
- Bring the test tubes with R1/dilution liquid to room temperature (15-25°C) before use.

Sample material

Human urine.

Pre-treatment of the urine sample or control prior to dilution is not required. Samples that appear turbid after dilution should be centrifuged or filtered before application to the test device. For screening of microalbuminuria a random urine sample as the first morning urine can be used^{1,2,3}.

Internal quality control

NycoCard U-Albumin Control should be used to confirm the efficacy of the reagents and correct performance of the test. The measured values should be within the acceptance limits stated on the vial label.

Test procedure (illustrations on page 69)

1 Sample dilution

Add 50 µL urine sample or C/control to a test tube with R1/dilution liquid. Mix well.

Note! R1/dilution liquid must have room temperature (15-25°C) when used.

2 Application of sample

Apply 50 µL diluted urine or diluted control to a TD/test device. Allow the diluted sample to soak completely into the membrane (approximately 50 sec.).

Note! Avoid air bubbles.

3 Application of R2/conjugate

Apply 50 µL of R2/conjugate to the test device. Allow the conjugate to soak completely into the membrane (approximately 50 sec.).

Note! Avoid air bubbles.

4 Application of R3/washing solution

Immediately apply 50 µL of R3/washing solution to the test device. Allow the washing solution to soak completely into the membrane (approximately 50 sec.).

Note! Avoid air bubbles.

5 Read the results

Read the results within 5 minutes using the NycoCard Reader II. Follow the instructions given in the NycoCard Reader II User Manual.

INTERPRETATION OF RESULTS

Interpret the NycoCard U-Albumin Test results with careful consideration to the patient's medical history, clinical examinations and other laboratory results. If the test result is questionable or if clinical signs and symptoms appear inconsistent with the test result, re-test the sample or confirm the result using another method. Analyse control materials frequently to verify the performance of the NycoCard Reader II Test System.

Reference range

<20 mg/L

Albumin is normally excreted in the urine at the rate 5-20 µg/min (up to 30 mg/24 hours), giving a urine albumin concentration up to 20 mg/L with a normal urine volume.

Increased values

The condition microalbuminuria is present when the urine albumin excretion is persistently elevated to 30-300 mg/24 hours, giving a urine albumin concentration of 20-200 mg/L with a normal urine volume⁴.

The urine albumin excretion rate (UAER) is given as µg/min when timed urine is used:

$$\frac{C_{\text{albumin}} \times V_{\text{urine}}}{\text{Sampling time, min}} = \mu\text{g/min}$$

C_{albumin} = albumin concentration, mg/L
 V_{urine} = total urine volume, mL

TROUBLE SHOOTING

Problem	Possible cause	Corrective action
The colour of the test device membrane is not homogeneous (white spots).	An air bubble has been formed during application of the sample.	Retest the sample and avoid air bubbles.
Unexpected low results.	Incorrect pipetting volume.	Check the pipette and pipette tips. Retest the sample.
	The R1/dilution liquid has not been equilibrated to room temperature (15-25°C) before use.	Bring R1/dilution liquid to room temperature and retest the sample.
Unexpected high results.	Incorrect pipetting volume.	Check the pipette and pipette tips. Retest the sample.
	A turbid diluted sample will give significantly increased flow time (time before the sample is absorbed into the membrane).	Clarify the diluted sample by centrifugation or filtration. Retest the sample.
	A reagent has been applied to the test device before the preceding reagent has soaked completely into the test device membrane.	Retest the sample. Make sure that all reagents soak completely into the test device before applying the following reagent.

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© 2019 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.