CATETERE DI ASPIRAZIONE STERILE MONOUSO
STERILE SUCTION CATHETER FOR SINGLE USE
CATHÉTER D'ASPIRATION STÉRILE À USAGE UNIQUE
CATÉTER DE SUCCIÓN ESTÉRIL DE UN SOLO USO
CATETER DE SUÇÇÃO ESTÉRIL PARA UTILIZAÇÃO ÚNICA
STERILER ABSAUGKATHETER FÜR DEN EINMALGEBRAUCH
ΑΠΟΣΤΕΙΡΩΜΕΝΟΣ ΚΑΘΗΤΗΡΑΣ ΑΝΑΡΡΟΦΗΣΗΣ ΜΙΑΣ ΧΡΗΣΗΣ
STERYLNY CEWNIK SSĄCY DO JEDNORAZOWEGO UŻYTKU
STERILNÍ SACÍ KATÉTR NA JEDNO POUŽITÍ
STERIL SUGKATETER FÖR ENGÅNGSBRUK
STERILNI SUKCIJSKI KATETER ZA ENKRATNO UPORABO
STERIILI, KERTAKÄYTTÖINEN IMUKATETRI
STERILNÝ ODSÁVACÍ KATÉTER NA JEDNO POUŽITIE
CATETER DE ASPIRAȚIE STERIL DE UNICĂ FOLOSINȚĂ
STERIELE ZUIGKATHETER VOOR EENMALIG GEBRUIK
STERILNI USISNI KATETER ZA JEDNOKRATNU UPORABU
STERIL SUGKATETER FÖR ENGÅNGSBRUK
STERILT SUGEKATETER TIL ENGANGSBRUG
СТЕРИЛЕН СМУКАТЕЛЕН КАТЕТЪР ЗА ЕДНОКРАТНА УПОТРЕБА
STERILUS VIENKARTINIO NAUDOJIMO ATSIURBIMO KATETERIS
STERILS SŪKŠANAS KATETRS VIENREIZĒJAI LIETOŠANAI
STERIILNE IMEMISKATEETER ÜHEKORDSEKS KASUTAMISEKS

Gima 22110 - 22112 - 22113 - 22114 - 22115 - 22116 - 22117 - 22118 - 22120 - 22122 - 22123 - 22124 - 22125

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ENGLISH

All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

Product Name: Sterile Suction Catheter for single use

Indications: It is necessary to use sputum suction to help the patient clean the respiratory tract, so as to prevent the patient from mistakenly sucking sputum, resulting in choking and suffocation.

Intend purpose: It is necessary to use sputum suction to help the patient clean the respiratory tract, so as to prevent the patient from mistakenly sucking sputum, resulting in choking and suffocation.

Contraindications: Banning the use of base of skull fracture patients.

Caution:

A. When using the product, it should strictly implement the requirements of aseptic operation specifications and relevant regulations, and can be used by trained doctors, nursing staff. In the course of use, the operator or patient should be treated in a timely manner according to medical regulations.

B. If you find any possible quality doubt during use, please stop using it immediately and notify our company quickly.

C. Sterilized by EO gas. Product is sterile.

D. For single use. Discard after use.

E. Please checking the individual package before using it.

F. Don't use if package is opened or damaged.

G. Sterile and validity for 5 year. Don't use if overdue.

H. Discard after use. After use, the products shall be collected, transported, stored and disposed to ensure that the products are harmless to human health and the environment. Or the product can be disposed in accordance with the local medical regulations for waste disposal.

Patient group: The suction catheter is suitable for patients with respiratory tract obstruction who cannot spit autonomously.

Intended users: It can be used by trained doctors, nursing staff.

Clinical benefits: Suction catheter can aspirate sputum from the patient's respiratory tract, indirect improvement of body function, alleviate symptoms.

Any residual risk and any undesirable side-effects: none
install information: none

Performance characteristics: The product is provided in a sterile state and sterilized with ethylene oxide. Single-use

Physical performance:

a. The tube, the head end and the hole are free of impurities, the surface is smooth, the top is round and without pitting point, the inner cavity is clean and smooth, if there is a print line, it should be clear and complete, not easy to fall off.

b. Security of construction: The adapter of the suction catheter and the tube shall be able to withstand a tensile force of 5N/15N/20N.

c. Shaft performance: When the adapter of the suction catheter is connected to the vacuum source, block the patient end, set the vacuum degree of the vacuum control device (if any) to the maximum, apply a negative pressure of 40 kPa at 23 °C ± 2 °C for 15 seconds, and the tube body shall not be flattened.

d. Vacuum control device performance: The residual vacuum of the suction catheter shall not exceed 0.33 kPa. Chemical performance: The residue of ethylene oxide in each product shall be no more than 10ug/g.

Biological performance: The product should be sterile.

Structural composition: The product is composed of tube and adapter.

Method:

a. Before sputum suction, check whether the connection of each part of the suction catheter is perfect and whether it can operate normally; Whether the suction catheter is unobstructed.

b. Open the small package and take the suction catheter out.

c. Connect the adapter of the suction catheter with the vacuum negative pressure system.

d. Try to aspirate with normal saline to check whether the catheter is unobstructed. Adjust suction pressure according to the patient's condition.

e. Tilt the patient's head to one side and lean

back slightly. Open the mouth of the patient with a tongue depressor, suction catheter is inserted into the throat from the buccal portion of the mouth, and the suction catheter is inserted into the trachea when the patient inhales. If suction operation is difficult from the oral, it can be inserted through the nasal cavity; A tracheotomy patient with a tracheal intubation or tracheal cannula can be inserted by a tracheal intubation or a tracheal cannula for Sputum suction.

Storage: The product should be stored in a well-ventilated room temperature.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

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