GIMA

PROFESSIONAL MEDICAL PRODUCTS

LUCE DI WORTH WORTH LIGHT TEST TEST DE WORTH WORTH-TEST TEST DE WORTH TESTE DE WORTH ΤΕΣΤ ΦΩΤΟΣ ΤΗΣ WORTH WORTH فحص ضوء وورث

Manuale utente - User manual - Notice d'utilisation - Gebrauchsanweisung
Manual de uso - Guia para utilização - Οδηγίες χρήσης -

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All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

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GIMA 31289

BASCH MEDICAL PVT. LTD.
Plot No.152, Jaggi Garden, Baldev Nagar,
Ambala City - 134003 - INDIA
MADE IN INDIA

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Purpose: Worth light Test is used to assess a patient's flat Fusional ability.

Indications: Flat Fusional Testing is Indicated any time stereopsis falls between (50) Sec of Arc, on those patients with suspected strabismus and on preschool Children. The Worth Test should also be used when evaluating cases of reduced monocular visual acuity that does not improve with the Pinhole Test.

Procedure

1. With the best refraction correction worn by the patient, place the Anaglyph glasses over the patient's correction, with red filter over the right eye.
2. In a slightly dimmed room, turn on the flashlight and hold the Worth Test with the red dot orientated up at approximately 16" from the patient and slightly below the LOS.
3. Conduct the following Monocular Check first:
 - a. Cover the right eye; ask how many dots the patient sees. They should report 3 green.
 - b. Then Cover the left eye and ask how many dots they see now. They should report 2 red.
4. Next, Conduct the binocular test:
 - a. With both eyes uncovered ask a third time how many dots they see. If normal flat fusion is present they should report 4.

Pediatric Note: This test can reliably be conducted on preschoolers as young as 2 if they are allowed to point to the dots.

5. Abnormal responses:
 - a. If the patient reports only 2 red dots under binocular conditions, this indicates that they are suppressing the left eye.
 - b. If the patient reports seeing 3 green dots under binocular conditions, then they are suppressing the right eye.
 - c. If they report 5 dots, they are diplopic. The type of diplopia can then be determined by asking which side are the green dots.
- If the green dots are located on the right, the patient has an eso deviation; to the left, an exo deviation.
 - a. If the green dots are reported above or below the red dots then a vertical deviation exists. A report of the green dots above the red dots would be seen with a right hyper deviation.
 - b. With a report of 6 or more dots, one should question the patient's reliability.
6. The binocular view is repeated at 5 and 10 feet, then repeat all these distances under greatly reduced light as both suppression and ocular deviation can be different under varying lighting conditions.

Battery Description	
Battery Type Required for this Instrument	No of Batteries Required
Alkaline Battery : Size "AA"	3

Disposal: The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

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