
CUCITRICE CUTANEA MONOUSO DISPOSABLE SKIN STAPLER AGRAFEUSE CUTANÉE JETABLE GRAPADORA DE PIEL DESECHABLE GRAMPEADOR DE PELE DESCARTÁVEL EINWEG-HAUTHEFTER

ITALIANO

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ENGLISH

Our special designed and manufactured disposable skin stapler is the specific instrument for surgical skin suture. With novel structure and application security, our skin stapler has the advantages of quick speed of sewing up during the operation, little histological reaction, good wound matching and healing, small surgical scar, and no pain for removing, etc. Our disposable skin stapler enables you not to only work with higher efficiency, but also avoid cross-infection effectively.

1. Main Structure

- The product is composed of skin stapler and staple. Each skin stapler has some quantities of staples inside.(Please check the detail refer to the product label). The skin stapler is composed of the component of body, handle, cartridge, staple push device and staple turn back device.
- The raw material of the staple is the stainless steel for medical usage.

2. Working Principle

For the aim of skin suture, after the handle was pushed by outside force, the staples will be shaped after extrusion under the mechanical transmission of the each part.

3. Capability Index

- The suture force of the skin stapler should be no more than 55N
- The metal components of the skin stapler have the good function of anti-eroded.
- he skin stapler is sterilized with ethylene oxide. The leftover of the ethylene oxide should be less than 10 mg/kg before the products are released from factory.

4. Application

Disposable skin stapler is used to do the wound suture of skin during the surgical operation.

5. How-to-use

- Remove the Tyvek paper and take out the skin stapler from blister package.
- With a pair of tweezers lift the skin of both sides of the wound, bring them together and hold it.
- Make sure that the front end of the stapler is well over the wound and the staple is vertical to the longitudinal axis of the incision. Gently put the stapler onto the incision till the fore-arrow levels with the middle line of the incision. Do not press hard on the wound using the stapler as it would make the removal of the staples difficult and painful.
- Caution : To clasp the handle and do not loosen it until the staple is out, remove the stapler and get ready for the next step.

See Images 1-2-3.

6. Limit of use

The use of this medical devices is exclusively for qualified medical personnel.

7. Precautions of use

- The skin stapler is sterilized with ethylene oxide.
- Sterile medical device, do not use in case of damaged single blister packaging or if the expiry date is over.
- Make sure the product can not be polluted when it was taken out of blister package.
- The product is limited to disposable use. Please destroy/dispose the finished product according to the related regulation.
- The skin stapler is only to be used for wound suture of at least 5 mm depth. If the depth of the wound suture is less than 5mm, the skin stapler is forbidden to be used for suture.

- If the hypodermic part under the incision is too thick, it is advised to first do an hypodermic suture and then a skin suture.
- The advised distance between each staple is of 5 to 10 mm. Please use the specific staple remover to remove the staples after 7-10 days. Or please refer to the doctor's prescribing.
- Please use the skin stapler within the period of validity.
- Contraindications**
If the wound depth is less than 5mm, the use of stapler for skin closure is contraindicated. The use of skin stapler is contraindicated in patients with known sensitivities or allergies to the metals contained in 317L stainless steel, i.e. chromium, nickel, iron, molybdenum, copper. Do not be used for internal tissue closure. This product is not designed, sold, or intended for use except as indicated.
- Storage conditions**
The skin stapler must be stored in a cool dry place, with a humidity level ≤ 80% with no erosive gases. To be manipulated with care. Do not store other products on top of skin staplers' cartons.
- 10. Period of Validity**
5 Years since sterile date. Please see expiry date of each product on the packaging.
- 11. Service after Sold**
If you have any quality complains, please contact with suppliers or manufacturer.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

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Ningbo Advan Electrical Co., Ltd.
Holding Corp. GmbH (Europe)
Gima S.p.A. Via Marconi, 1