TRUSILK®

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ENGLISH

DESCRIPTION

Trusilk is a nonabsorbable, sterile, surgical suture, composed of an organic protein, called fibroin. This protein is derived from the domesticated species Bombyx mori (B mori) of the family Bombycidae. Natural Silk contains approximately 30% extraneous material or gum which is removed by specially developed "Degumming" process. This process is essential for a compact braid while ensuring that the filaments retain their natural body and elasticity. Trusilk dyed with black Hematoxylin and coated with wax/ silicon. Trusilk meets all the requirements as per United States Pharmacopeia for Nonabsorbable Surgical suture.

INTENDED PURPOSE

Trusilk is indicated for use in soft tissue approximation and or ligation, including use in cardiovascular, ophthalmic and neurological tissues.

ACTION

Trusilk elicits an initial inflammatory reaction in tissues which is followed by gradual encapsulation of the suture by fibrous connective tissue. While Trusilk is not absorbed, progressive degradation of the proteinous silk fibres in vivo may result in gradual loss of the sutures over a period of time.

CONTRA INDICATIONS

Use of this suture is contra indicated in patients with known sensitivities or allergies to Silk. Due to the gradual loss of tensile strength which may occur over prolonged periods in vivo, Trusilk sutures should not be used where permanent retention of tensile strength is required as in fixation of vascular prostheses.

WARNINGS

Do not reuse. Discard open, unused sutures. Prolonged contact of this suture or any other suture with salt solutions such as those found in urinary and biliary tracts, may result in calculus formation. Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing Trusilk for wound closures, as risk of wound dehiscence may vary with the site of application and the suture material used. As any foreign material in the presence of bacterial contamination may enhance bacterial infectivity, acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

PRECAUTIONS

In handling this suture or any other suture material care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. As with any suture material, adequate knot security requires the accepted surgical techniques of flat and square ties, with additional throws as warranted by surgical circumstances and the experience of the surgeon. Discard used needles in 'Sharps' containers.

ADVERSE REACTIONS

Adverse reactions associated with the use of the device, Trusilk include wound dehiscence, gradual loss of all tensile strength over time, allergic response in patients who are known to be sensitive to silk, calculus formation in urinary & biliary tracts when prolonged contact with salt solutions such as urine and bile juice, enhanced bacterial infectivity, acute inflammatory tissue reaction, pain, edema and erythema at the wound site.

SUPPLY

Trusilk is available in various U.S.P. sizes, The suture is supplied sterile in pre-cut lengths, both non- needed and attached to various needle type, shape and length, which are packed in a printed box quantity as indicated on the box label. Non sterile reels are packed in 25 metres, six reels packed in a box.

STORAGE

Recommended storage conditions: Below 30°C away from moisture and di-

rect sunlight. Do not use after expiry.

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**Healthium Medtech Limited**
**MED DEVICES LIFESCIENCES B.V.**