MONOGLYDE®

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ENGLISH

DESCRIPTION

Monoglyde (Poliglecaprone25), is a monofilament synthetic absorbable suture, composed of a Polymer (Glycolide co – caprolactone). Monoglyde sutures are available undyed in beige color and in violet color, dyed with D & C violet. Poliglecaprone 25 copolymer has been found to be nonantigenic, nonpyrogenic and elicits only a slight tissue reaction during absorption. MONOGLYDE sutures comply with USP requirements, except for diameter.

INTENDED PURPOSE

Monoglyde sutures are indicated for use in general soft tissue approximation, but not for use in cardiovascular or neurological tissues, microsurgery or ophthalmic surgery.

ACTIONS

Monoglyde suture is a monofilament which elicits a minimal acute inflammatory reaction in tissues and in-growth of fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of Monoglyde sutures occurs by means of hydrolysis. Absorption begins as loss of tensile strength followed by loss of mass. Implantation studies in rats indicate that Monoglyde suture retains approximately 45% of its original tensile strength 7 days post implantation and approximately 15% of its original tensile strength at 14 days post implantation. All of the original tensile strength is lost by 21 days post implantation. The absolute strength remaining 14 days post implantation meets or exceeds that historically observed with plain and chromic gut sutures. Absorption of MONOGLYDE absorbable synthetic suture is essentially complete approximately 90 days.

CONTRA INDICATIONS

Monoglyde sutures being absorbable, should not be used where extended approximation of tissue under stress is required, such as in fascia.

WARNINGS

Users should be familiar with the surgical procedures and techniques involving absorbable sutures before employing Monoglyde suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Physicians should consider the in vivo performance when selecting a suture for use in patients. The use of this suture may be inappropriate in elderly, malnourished

or debilitated patients or in patients suffering from conditions which may delay wound healing. Do not reuse. Discard opened packages and unused sutures. As with any foreign body, prolonged contact of any suture with salt solutions such as those found in the urinary or biliary tracts may act transiently as a foreign body. Acceptable surgical practice should be followed for the management of contaminated or infected wounds. As this is an absorbable suture material the use of supplemental sutures should be considered by the surgeon in the closure of the sites which may undergo expansion, stretching or distention or which may require additional support.

PRECAUTIONS

Skin sutures which must remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated. Sub cuticular sutures should be placed as deeply as possible to minimize the erythema and induration normally associated with absorption. Under some circumstances, notably orthopedic procedures, immobilization of joints by external support may be employed at the discretion of the surgeon. Consideration should be taken in the use of absorbable sutures in tissue with poor blood supply as suture extrusion and delayed absorption may occur. In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. Monoglyde suture knots must be properly placed to be secure. Adequate knot security requires the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstances and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments. Avoid prolonged exposure to elevated temperature.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharp" containers.

ADVERSE REACTIONS

Adverse effects associated with the use of synthetic absorbable sutures include wound dehiscence, failure to provide adequate wound support in closure of the sites where expansion, stretching, or distention occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients or in patients suffering from conditions which may delay wound healing, infection, localized irritation when skin sutures are left in place for more than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, and transitory local irritation at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood-borne pathogens.

SUPPLY

Monoglyde is available in various U.S.P. sizes, The suture is supplied sterile in pre-cut lengths, both non-needled and attached to various needle type, shape and length. which are packed in a printed box quantity as indicated on the box label.

STORAGE

Recommended storage conditions: Below 30°C away from moisture and direct sunlight. Do not use after expiry.

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**Healthium Medtech Limited**
**MED DEVICES LIFESCIENCES B.V.**