Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

TRUSYNTH®

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ENGLISH

DESCRIPTION

Trusynth is an absorbable, sterile, surgical suture, composed of a co-polymer, made from 90% Glycolide and 10% Lactide, Polyglactin 910. Trusynth is colored violet with D & C Violet No. 2, conforming to U.S. Code of Federal regulations. Trusynth is available as undyed also. For added lubrication and smoothness, Trusynth is coated with unique combination of Polyglactin 370 and Calcium Stearate. Trusynth meets all the requirements, established by the United States Pharmacopeia for Absorbable Surgical suture (Synthetic) except for diameter.

INTENDED PURPOSE

Trusynth is indicated for use in general soft tissue approximation and or Ligation, including use in ophthalmic procedures but not for use in cardio vascular and neurological tissues.

ACTION

Trusynth suture gets absorbed in the site of implementation, owing to the hydrolysis, where the polymer degrades to Glycolic acid and Lactic acid. It is slowly absorbed and metabolized in the body. Absorption begins with the loss of tensile strength followed by loss of mass. Approximately 70% of the original tensile strength is retained up to 14 days after surgery. Approximately 15% of the original tensile strength is retained up to 28 days after surgery. Absorption of Trusynth is minimal upto 10% in two weeks, 25% of the mass gets absorbed in four weeks, upto 75% in two months and absorption is essentially complete Approximately 80 days.

CONTRA INDICATION

Trusynth suture being absorbable, should not be used where extended approximation of tissue is required.

WARNING

Do not resue, Discard open unused sutures. Prolonged contact of this suture or any other suture with salt solutions such as those found in urinary and binary tracts, may result in calculus formation. As an absorbable suture, it may act transiently as a foreign body. As any foreign material in the presence of bacterial contamination may enhance bacterial infectivity, acceptable surgical practice must be followed with respect to storage with respect to drainage and closure of infected or contaminated wounds. USE of Trusynth suture may be inappropriate in patient suffering from conditions which may be delay wound healing.

PRECAUTIONS

As this is an absorbable suture, the surgeon in the closure of sites undergoing expansion, stretching or distention, which may require additional support, should consider the use of supplemental non absorbable suture. Skin suture, which remain in place for more than 7 days may cause localized irritation and should be snipped off or removed. In handling Trusynth suture or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. As with any suture material, adequate knot security requires the accepted surgical technique of flat and square ties with additional throws as indicated by surgical circumstance and experience of the surgeon. Discard used needle in sharps containers.

ADVERSE REACTIONS

Adverse effects associated with the use of Trusynth include allergic response in certain patients, transient local irritation at the wound site, transient inflammatory foreign body response, erythema and induration during the absorption process of subcuticular sutures.

SUPPLY

Trusynth is available in various U.S.P. sizes, The suture is supplied sterile

in pre-cut lengths, both non- needed and attached to various needle type, shape and length. which are packed in a printed box quantity as indicated on the box label.

STORAGE

Recommended storage conditions: Below 30° C away from moisture and direct sunlight. Do not use after expiry.

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Healthium Medtech Limited
MED DEVICES LIFESCIENCES B.V.