Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

TRULON®

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ENGLISH

DESCRIPTION

Trulon is a nonabsorbable, sterile, surgical suture, composed of long chain-aliphatic polymers, Nyon 6 and Nylon 6.6. The monofilament Polyamide sutures attached needles are available in black and blue, colors dyed with FDA approved dyes. Polyamide 6.6 is formed by Polycondensation of hexamethylene diamine and adipic acid. Polyamide 6 is formed by polymerisation of epsilon caprolactum. Blue Nylon sutures is dyed with a blue pigment, colour Index No.69800 and Black Nylon suture is dyed with Haematin HCK, colour Index No. 75290. Trulon meets all the requirements, established by the United States Pharmacopeia for Nonabsorbable Surgical Suture.

INTENDED PURPOSE

Trulon is indicated for use in soft tissue approximation and or ligation, including use in cardiovascular, ophthalmic and neurological tissues.

ACTION

Trulon suture elicits a minimal acute inflammatory reaction in tissue which is followed by gradual encapsulation of suture by fibrous connective tissue. While it is not absorbed, progressive hydrolysis of Nylon in vivo may result in gradual loss of tensile strength over a period of time.

CONTRA INDICATIONS

Due to gradual loss of tensile strength which may occur over prolonged periods in vivo Trulon suture should not be used where permanent retention of tensile strength is required as in fixation of intraocular lenses or synthetic vascular grafts.

WARNINGS

Do not reuse. Discard open, unused sutures. Prolonged contact of this suture or any other suture with salt solutions such as those found in urinary and biliary tracts, may result in calculus formation. Users should be familiar with surgical procedures and techniques involving non absorbable sutures before employing Trulon for wound closures, as risk of wound dehiscence may vary with the site of application and the suture material used. As any foreign material in the presence of bacterial contamination may enhance bacterial infectivity; acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

PRECAUTIONS

In handling this suture or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. As with any suture material, adequate knot security requires the accepted surgical techniques of flat and square ties, with additional throws as warranted by surgical circumstances and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting with Trulon being a monofilament. Discard used needles in ‘Sharps’ containers.

ADVERSE REACTIONS

Adverse effects associated with the use of Trulon include wound dehiscence, gradual loss of tensile strength over time, minimal acute inflammatory tissue reaction, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile juice etc. occurs, enhanced bacterial infectivity, pain, edema and erythema at wound site.

SUPPLY

Trulon is available in various U.S.P. sizes, The suture is supplied sterile in pre-cut lengths, both non- needed and attached to various needle type, shape and length. which are packed in a printed box quantity as indicated on the box label.

STORAGE

Recommended storage conditions: Below 30° C away from moisture and direct sunlight. Do not use after expiry.

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