TRULENE®

FILO DA SUTURA CHIRURGICA NON RIASSORBIBILE USP (MONOFILAMENTO POLIPROPILENE)

NONABSORBABLE SURGICAL SUTURE U.S.P (MONOFILAMENT POLYPROPYLENE)

SUTURE CHIRURGICALE NON RÉSORBABLE USP (POLYPROPYLENE MONOFILAMENT)

NICHT-RESORBIERBARES CHIRURGISCHES NAHTMATERIAL U.S.P (MONOFILES POLYPROPYLEN)

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DESCRIPTION

Trulene is a nonabsorbable, sterile, surgical suture, composed of Polypropylene, a synthetic linear Polyolefin. The molecular formula is (C3H6) n. Trulene is colored blue with Phthalocyanine Copper blue, Color Index No.74160, and U.S.F.D.A. approved dye. Trulene is monofilament and it is uncoated. Trulene meets all the requirements as per United States Pharmacopeia for nonabsorbable surgical suture.

INTENDED PURPOSE

Trulene is indicated for use in soft tissue approximation and or ligation, including use in cardiovascular, ophthalmic and neurologic tissues.

ACTION

TRULENE being a monofilament suture is remarkably smooth and easy to handle. It does not have interstices in its structure. It elicits an initial inflammatory reaction in tissues which is followed by gradual encapsulation of the suture by fibrous connective tissue. Trulene is neither absorbed, nor it is subjected to degradation or weakening by the action of tissue enzymes. Due to its relative biological inertness, it is recommended for use where the least possible suture reaction is desired. Because of its lack of adherence to tissue, Trulene is effective as a pull out suture. Trulene, being a monofilament has been successfully employed in surgical wounds which could subsequently become infected or contaminated where it can minimize later sinus formation and suture extrusion.

CONTRA INDICATIONS

None known.

WARNINGS

Do not resterilise. Discard open, unused sutures. Prolonged contact of this suture or any other suture with salt solution such as those found in urinary and biliary tracts may result in calculus formation. Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing Trulene for wound closures, as risk of wound dehiscence may vary with the site of application and the suture material used. As any foreign material in the presence of bacterial contamination may enhance bacterial infectivity, acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

PRECAUTIONS

In handling this suture or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holder. As with any suture material, adequate knot security requires the accepted surgical techniques of flat and square ties, with additional throws as warranted by surgical circumstances and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting with Trulene being a monofilament. Discard used needles in 'Sharps' containers.

ADVERSE REACTIONS

Adverse effects associated with the use of Trulene include minimal initial inflammatory tissue reaction and transient local irritation at the wound site, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile juice etc. Like all foreign bodies, TRULENE may potentiate an existing infection.

SUPPLY

TRULENE is available in various U.S.P. sizes, The suture is supplied sterile in pre-cut lengths, both non- needed and attached to various needle type, shape and length. which are packed in a printed box quantity as indicated on the box label.

STORAGE

Recommended storage conditions: Below 30°C away from moisture and direct sunlight. Do not use after expiry.

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Healthium Medtech Limited
EC REP MED DEVICES LIFESCIENCES B.V.