Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

TRUGLYDE®

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ENGLISH

DESCRIPTION

Truglyde Suture is an absorbable, sterile, synthetic surgical suture, composed of homopolymers of Glycolide(100%). Truglyde is colored violet with D & C Violet No.2, conforming to U.S.Code of Federal regulations. Truglyde is available undyed also. For added lubrication and smoothness, Truglyde is coated with unique combination of polycaprolactone and calcium stearate. Truglyde meets all the requirements as per United States Pharmacopeia for Absorbable Surgical suture (Synthetic).

INTENDED PURPOSE

Truglyde is indicated for use in soft tissue approximation and or ligation, including use in plastic and ophthalmic procedures but not for use in cardiovascular and neurological tissues.

ACTION

Truglyde suture gets absorbed in the site of implantation, owing to the hydrolysis, where the polymer degrades to Glycolic acid. It is slowly absorbed and metabolized in the body. Absorption begins with loss of tensile strength followed by loss of mass. Approximately 70% of the original tensile strength is retained up to 14 days after surgery. Approximately 15% of the original tensile strength is retained up to 28 days after surgery. Absorption of Truglyde suture is minimal up to 10% in two weeks. 25% of the mass gets absorbed in four weeks, unto 75% in two months and absorption is essentially complete Approximately 70 days.

CONTRA INDICATION

Truglyde suture being absorbable, should not be used where extended approximation of tissue is required.

WARNING

Do not reuse. Discard open unused sutures. Prolonged contact of this suture or any other suture with salt solutions such as those found in urinary and biliary tracts, may result in calculus formation. As an absorbable suture, it may act transiently as a foreign body. As any foreign material in the presence of bacterial contamination may enhance bacterial infectivity, acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds. The use of this suture may be inappropriate in patients suffering from conditions which may delay wound healing.Discard used needle in sharps container.

PRECAUTIONS

As this is an absorbable suture, the surgeon in the closure of sites undergoing expansion, stretching or distention, which may require additional support, should consider the use of supplemental non absorbable suture. Under some circumstances, notably orthopaedic Procedure, immobilisation by external support may be employed at the discretion of the surgeon. Skin suture, which remain in place for more than 7 days may cause localized irritation and should be snipped off or removed. Discard used needle in sharps containers.

ADVERSE REACTIONS

Adverse effects associated with the use of Truglyde include allergic response in certain patients, transient local irritation at the wound site, transient inflammatory foreign body response, erythema and induration during the absorption process of subcuticular sutures.

SUPPLY

Truglyde is available in various U.S.P. sizes, The suture is supplied sterile in pre-cut lengths, both non- needed and attached to various needle type, shape and length. which are packed in a printed box quantity as indicated on the box label.

STORAGE

Recommended storage conditions: Below 30°C away from moisture and direct sunlight. Do not use after expiry.

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