GIMA

PROFESSIONAL MEDICAL PRODUCTS

Gima S.p.A.

Via Marconi, 1 - 20060 Gessate (MI) Italy

gima@gimaitaly.com - export@gimaitaly.com

www.gimaitaly.com

TRUBOND™

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ENGLISH

DESCRIPTION

Trubond is a nonabsorbable, sterile, surgical suture, composed affine filaments of Polyester, Poly Ethylene Tere Phthalate. The fine polyesterfibres are braided to produce firm suture that remains soft and pliable. Trubond is colored green with D&C Green No.5, USFDA. Approved dye. Trubond is available undyed also. For added lubrication and smoothness, Trubond is coated with Silicon. Trubond meets all the requirements, established by the Unitedci States Pharmacopeia far Nonabsorbable Surgical suture.

INTENDED PURPOSE

Trubond is indicatoci far use in soft tissue approximation and or ligation, including use in cardiovascular, ophthalmic and neurological tissues.

ACTION

Trubond being a braided suture, enables far secured knots. Trubond is soft and pliable. It is very easy to handle. It elicits an initial inflammatory reaction in tissues which is followed by gradual encapsulation of the suture by fibrous connective tissue. Trubond is not absorbed, nor is it subjected to degradation or weakening bythe action of tissue enzymes. Due to its total biological inertness, it is recommended far use where the least possible suture reaction is desired. Because of its lack of adherence to tissue, Trubond is effective as a pull out suture.

CONTRA INDICATIONS

NoneKnown.

WARNINGS

Do not reuse. Discard open, unused sutures. Prolonged contact of this suture or anyothersuture with salt solutions such as those found in urinary and biliary tracts, may result in calculus formation. Users should be familiarwith surgical procedures and techniques involving no absorbable sutures before employing Trubond far wound closures, as risk of wound dehiscence may vary with the site of application and the suture material used. As any foreign material in the presence of bacterial contamination may enhance bacterial infectivity acceptable surgical practice must be followed with respect to drainage and closure of infected orcontaminated wounds.

PRECAUTIONS

In handling Trubond suture or any other suture material care should be taken to avoid damage from handling. Avo id crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. As with any suture material, adequate knot security requires the accepted surgical techniques of flat and square ties, with additional throws as warranted by surgical circumstances and the experience of the surgeon. Discard used needles in 'Sharps' containers.

ADVERSE REACTIONS

Adverse effects associateci with the use of Trubond include minimal initial inflammatory tissue reaction and transient local irritation at the wound site, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile juice etc. Trubond may patentiate an existing infection.

SUPPLY

Trubond is available in various U.S.P. sizes, The suture is supplied sterile in pre-cut lengths, both non- needed and attached tovarious needle type, shape and length. which are packed in a printed box quantity as indicatoci on the box la bel.

STORAGE

Recommended storage conditions: Below 30° C awayfrom moisture and directsunlight. Do not use afterexpiry.

DIPOSAL

Discard used sutures and needles contaminateci with blood in the container meant infection waste. Unused expired pouches should be incineratoci or Disposal should be done as per loca l regulation.

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MC REF **MED DEVICES LIFESCIENCES B.V.**
**Healthium Medtech Private Limited**