TRUBARB™

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ENGLISH

DESCRIPTION

The Polydioxanone Knotless Tissue-Closure Device is comprised of barbed dyed Polydioxanone (PDO) suture material, armed with a surgical needle on one end and end stopper at the other end. The Polydioxanone Knotless Tissue-Closure Device barbs are oriented in one direction to allow for tissue approximation without the need to tie surgical knots. Polydioxanone Knotless Tissue-Closure Device is consisting of dyed (violet) polyester, poly(p-dioxanone). The empirical molecular formula of which is (C4H6O3)X. The pigment for the violet dye is D&C Violet No. 2. Polydioxanone has been found to be nonantigenic, nonpyrogenic and to elicit only a slight tissue reaction during absorption. While the formation of barbs in the Polydioxanone Knotless Tissue-Closure Device reduces the tensile strength relative to non-barbed suture material of the same size, tying of knots in non-barbed suture materials also reduce their effective strength. For this reason, the strength of the Polydioxanone Knotless Tissue-Closure Device can be compared to USP knot strength of non-barbed sutures. Additionally, USP designations for diameter are used to describe the Polydioxanone Knotless Tissue-Closure Device prior to barbing, except for minor variation in suture diameter with a maximum overage of 0.14mm.

Tensile Strength/Size Equivalency

Polydioxanone Knotless Tissue-Closure Device size	Pre-Barbing Suture Size	Equivalent Non-Barbed Size(USP)/ Tensile Strength(Kgf)
1	2	1/5.08
0	1	0/3.90
2-0	0	2-0/2.68
3-0	2-0	3-0/1.77
4-0	3-0	4-0/0.95
5-0	4-0	5-0/0.68

ACTIONS

Two important characteristics describe the in vivo performance of absorbable sutures: first, tensile strength retention and second, the absorption rate (loss of mass). Polydioxanone Knotless Tissue-Closure Device has been formulated to minimize the variability of these characteristics and to provide wound support through the critical wound healing period and an extended healing period.

The results of implantation studies in animals using Polydioxanone Knotless Tissue-Closure Device indicate that

Days Implantation	Approximate % Original Strength Remaining
14 days	75%
28 days	65%
42 days	55%

Absorption of Polydioxanone Knotless Tissue-Closure Device is essentially complete between 180 to 220 days.

INTENDED PURPOSE

Polydioxanone Knotless Tissue-Closure Device is indicated for use in soft tissue approximation where use of Polydioxanone absorbable suture is appropriate.

CONTRAINDICATIONS

Polydioxanone Knotless Tissue-Closure Device is not to be used where prolonged (beyond six weeks) approximation of tissues under stress is required and is not to be used in conjunction with or for fixation of prosthetic devices (e.g. heart valves or synthetic grafts) that are non-absorbable in nature.

WARNINGS

Do not resterilize. Discard opened, unused Polydioxanone Knotless Tissue-Closure Device and associated surgical needles. Users should be familiar with surgical procedure and techniques involving absorbable sutures before employing Polydioxanone Knotless Tissue-Closure Device for wound Closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Physicians should consider the in vivo performance (under ACTIONS section) when selecting a suture for use in patients. The use of this suture may be inappropriate in elderly, malnourished or debilitated patients, or in patients suffering from conditions which may delaywound healing. The safety and effectiveness of Polydioxanone Knotless Tissue-Closure Device has not been established for use in fascial closures (including abdominal wall, thoracic and extremity fascial closures), gastrointestinal anastomoses, cardiovascular tissue, neural tissue, osseous tissue, tendinous tissue, ophthalmic surgery, or for use in microsurgery, therefore this product should not be used for these purposes.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation. As an absorbable suture, Polydioxanone Knotless Tissue-Closure Device may act transiently as a foreign body. Acceptable surgical practice should be followed for the management of contaminated or infected wounds. As this is an absorbable suture material, the use of supplemental

nonabsorbable sutures should be considered by the surgeon in the Closure of the sites which may undergo expansion, stretching or distention or which may require additional support.

PRECAUTIONS

Polydioxanone Knotless Tissue-Closure Device contains unidirectionally oriented barbs to anchor tissues and does not require knots to approximate opposing edges of a wound. Tying knots on the barbed section of the material will damage the barbs and potentially reduce the suture tensile strength and barb effectiveness. Additionally, when completing placement, an additional backstitch or bite of tissue lateral to the end of the incision is required to lock the device in place. Avoid contacting the Polydioxanone Knotless Tissue-Closure Device and associated needles with other materials (e.g. surgical gauze, drapes, etc.) in the surgical field to prevent ensnaring on the barbs. If the barbs catch, carefully pull the material in the opposite direction of the needle to disengage it from the barbs.

Care should be taken to avoid damage when handling. Avoid crushing or crimping the suture material with surgical instruments, such as needle holders and forceps. Do not pull the Polydioxanone Knotless Tissue-Closure Device out of the package by the needles as this can cause the barbs to catch on one another. Do not attempt to remove memory in the polymer. by running fingers down the suture material as this can damage the barbs.

Infections, erythema, foreign body reactions, transient inflammatory reactions and in rare instances wound dehiscence are typical or foreseeable risks associated with any suture and hence are also potential complications associated with Polydioxanone Knotless Tissue-Closure Device.

When using Polydioxanone Knotless Tissue-Closure Device subcutaneously, the device should be placed as deeply as possible in order to minimize erythema and induration normally associated with absorption. Acceptable surgical practice should be followed with respect to drainage and Closure of infected wounds.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharps" containers.

ADVERSE REACTIONS

Adverse effects associated with the use of this device may include, wound dehiscence, failure to provide adequate wound support in closure of the site where expansion, stretching or distension occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients or in patients suffering from conditions which may delay wound healing, infection, minimal acute inflammatory tissue reaction, localized irritation when skin sutures are left in place for greater than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs and transitory local infection at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens. Due to prolonged suture absorption, some irritation and bleeding may occur.

APPLICATION

Polydioxanone Knotless Tissue-Closure Device is designed to be used in continuous suture patterns and intended to be used without anchoring knots to begin or terminate the device line Use as required per surgical procedure.

Secure the fixation end stopper portion to robust tissue by taking a bite in the designated tissue, passing the needle and pulling tautly around the anchoring tissue.

After Polydioxanone Knotless Tissue-Closure Device is anchored into the tissue, grasp the strand and approximate tissue to the desired tension. Bites, or passes through the tissue, can be taken in a continuous suturing technique to approximate the wound. Care should be taken to utilize the Polydioxanone Knotless Tissue-Closure Device on the barbed segments only. Do not attempt to approximate wounds using the non-barbed segment near the

needle, as the barbs are required for successful wound approximation with the Polydioxanone Knotless Tissue-Closure Device. To complete the closure in subcuticular closure, take at least one pass in the reverse direction. Then pass the needle across the incision and take a split thickness bite perpendicular to the incision and exit the skin.

To complete the closure in all other tissue layers, take at least two passes in the reverse direction to complete the closure. Then gently pull on the free end of the device and cut flush with the surface of the tissue.

STERILITY

Polydioxanone Knotless Tissue-Closure Device is sterilized by ethylene oxide gas. Do not resterilize. Do not use if package is opened or damaged. Discard opened unused sutures. Do not use after expiration date.

HOW SUPPLIED

Polydioxanone Knotless Tissue-Closure Device is available sterile, in various sizes in unidirectional barb configurations attached with various needle type, shape and length which are packed in a printed box, quantity as indicated on the label.

STORAGE

Recommended storage conditions: Below 30° C away from moisture and direct sunlight. Do not use after expiry.

DISPOSAL

Discard used sutures and needles contaminated with blood in the container meant for infection waste. Unused expired pouches should be Incinerated or Disposal should be done as per local regulation.

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	Healthium Medtech Private Limited Plot No.1600, R-6 West, Sri City (SEZ), Chervi Village, Irrugulam Post, Satyavedu Mandal, Chittoor District, Andhra Pradesh- 517588, India Email : sales@healthiummedtech.com Mfg. Lic. No.: 42/CT/AP/2012/S/R
---	--

	IT - Codice prodotto GB - Product code FR - Code produit DE - Erzeugniscodex ES - Código producto PT - Código produto GR - Κωδικός προϊόντος
	IT - Dispositivo monouso, non riutilizzare GB - Disposable device, do not re-use FR - Ne pas réutiliser DE - Für einmaligen Gebrauch, nicht wiederverwenden ES - Dispositivo monouso, no reutilizable PT - Dispositivo descartável, não reutilizar GR - Προϊόν μιας χρήσεως. Μην το χρησιμοποιείται εκ νέου
	IT - Numero di lotto GB - Lot Number FR - Numéro de lot DE - Chargennummer ES - Número de lote PT - Número de lote GR - Αριθμός παρτίδας
	IT - Data di fabbricazione GB - Date of Manufacturing FR - Date de fabrication DE - Herstellungsdatum ES - Fecha de fabricación PT - Data de fabrico GR - Ημερομηνία παραγωγής
	IT - Data di scadenza GB - Expiration date FR - Date d'échéance DE - Ablaufdatum ES - Fecha de Caducidad PT - Data de validade GR - Ημερομηνία λήξεως
	IT - Sterilizzato con ossido di etilene GB - Sterilized using ethylene oxide FR - Stérilisé à l'oxyde d'éthylène DE - Sterilisiert mit Ethylenoxid ES - Esterilizado con óxido de etileno PT - Esterilizado com óxido de etileno GR - Αποστειρωμένο με αιθυλενοξείδιο
	IT - Rappresentante autorizzato nella Comunità europea GB - Authorized representative in the European community FR - Représentant autorisé dans la Communauté européenne DE - Autorisierter Vertreter in der EG ES - Representante autorizado en la Comunidad Europea PT - Representante autorizado na União Europeia GR - Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Ένωση
	IT - Fabbrikante GB - Manufacturer FR - Fabricant DE - Hersteller ES - Fabricante PT - Fabricante GR - Παραγωγός
	IT - Non ri-sterilizzare GB - Do not resterilize FR - Ne pas restériliser DE - Nicht erneut sterilisieren ES - No reesterilizar PT - Não reesterilize GR - Μην αποστειρώνετε
	IT - Conservare in luogo fresco ed asciutto GB - Keep in a cool, dry place FR - À conserver dans un endroit frais et sec ES - Conservar en un lugar fresco y seco DE - An einem kühlen und trockenen Ort lagern PT - Armazenar em local fresco e seco GR - Διατηρείται σε όρσοερό και στεγνό περιβάλλον
	IT - Conservare al riparo dalla luce solare GB - Keep away from sunlight FR - À conserver à l'abri de la lumière du soleil DE - Vor Sonneneinstrahlung geschützt lagern ES - Conservar al amparo de la luz solar PT - Guardar ao abrigo da luz solar GR - Κρατήστε το μακριά από ηλιακή ακτινοβολία
	IT - Limite superiore di temperatura GB - Upper limit of temperature FR - Limites supérieure de température DE - Obergrenze der Temperatur ES - Limitaciones superiorde temperatura PT - Limitação superior de temperatura GR - Άνωτερο όριο θερμοκρασίας
	IT - Non utilizzare se l'imballaggio è danneggiato GB - Don't use if package is damaged FR - Ne pas utiliser si le colis est endommagé DE - Nicht verwenden, wenn das Paket beschädigt ist ES - No usar si el paquete está dañado PT - Não use se o pacote estiver danificado GR - Μην το χρησιμοποιείτε αν η συσκευασία είναι κατεστραμμένη
	IT - Leggere le istruzioni per l'uso GB - Consult instructions for use FR - Consulter les instructions d'utilisation DE - Gebrauchsanweisung beachten ES - Consultar las instrucciones de uso PT - Consulte as instruções de uso GR - Διαβάστε προσεκτικά τις οδηγίες χρήσης
	IT - Dispositivo medico conforme alla Direttiva 93/42/CEE GB - Medical Device complies with Directive 93/42/EEC FR - Dispositif médical conforme à la directive 93/42 / CEE DE - Medizinprodukt gemäß Richtlinie 93/42/CEE ES - Dispositivo médico según a la Directiva 93/42 / CEE PT - Dispositivo médico em conformidade com a Diretiva 93/42/CEE GR - Ιατρική συσκευή σύμφωνα με την οδηγία 93/42 / CEE
	IT - Identificatore univoco del dispositivo GB - Unique device identifier FR - Identifiant unique de l'appareil ES - Identificador de dispositivo único DE - Eindeutige Gerätekenntung PT - Identificador exclusivo do dispositivo GR - Μοναδικό αναγνωριστικό συσκευής
	IT - Attenzione: Leggere e seguire attentamente le istruzioni (avvertenze) per l'uso GB - Caution: read instructions (warnings) carefully FR - Attention: lisez attentivement les instructions (avertissements) DE - Achtung: Anweisungen (Warnings) sorgfältig lesen ES - Precaución: lea las instrucciones (advertencias) cuidadosamente PT - Cuidado: leia as instruções (aviso) cuidadosamente GR - Προσοχή: διαβάστε προσεκτικά τις οδηγίες (ενστάσεις)

REF BSN0316U (GIMA 22890)
BSN0326U (GIMA 22891)
BSN0315U (GIMA 22892)
BSN0325U (GIMA 22893)
BSN0615U (GIMA 22894)
BSN0624U (GIMA 22895)
BSN1100U (GIMA 22896)

EC REF MED DEVICES LIFESCIENCES B.V.
Kraijenhoffstraat 137 A, 1018RG Amsterdam, Netherlands
Email: info@meddevices.net Phone: +31-202254558

Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

Healthium Medtech Private Limited
Plot No.1600, R-6 West, Sri City (SEZ),
Chervi Village, Irugulam Post, Satyavedu Mandal, Chittoor District,
Andhra Pradesh- 517588, India
Email : sales@healthiummedtech.com
Mfg. Lic. No.: 42/CT/AP/2012/S/R