TRULENE MESH™

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ENGLISH

DESCRIPTION

Trulene mesh is constructed of knitted filaments of extruded Polypropylene, identical in composition with that used in Trulene polypropylene suture. Non absorbable surgical suture U.S.P. Trulene Mesh is available undyed in a variety of square and rectangular sizes for surgical procedures. This material, when used as sutures has been reported to be nonreactive and to retain its strength for an indefinite period in clinical use.

Trulene mesh is knitted by a process which interlinks each fiber junction and which provides for elasticity in both directions. This construction permits the mesh to be cut into any desired shape or size without unraveling. The fiber junctions are not subject to the same work – fatigue exhibited by more rigid metallic meshes. This bi-directional elastic properly allows adaptation to various stresses encountered in the body.

ACTIONS

Trulene mesh is a non absorbable mesh used to span and reinforce traumatic or surgical wounds to provide extended support during and following wound healing. This mesh elicits minimum inflammatory reaction, which is transient and is followed by the deposition of a thin fibrous layer of tissue which can grow through the interstices of the mesh, thus incorporating the mesh into the adjacent tissue. The mesh remains soft, pliable and the wound healing is not noticeably impaired. The mesh is not absorbed nor is subject to degradation or weakening by the action of tissue enzymes.

INTENDED PURPOSE

Trulene Mesh used for the repair of hernia. Types of Hernia: Inguinal Hernia and Ventral Hernia.

CONTRAINDICATIONS

When this mesh is used in infants or children with future growth potential and pregnant women, the surgeon should be aware that this product will not stretch significantly as the patient grows. Trulene Mesh in contaminated wounds should be used with the understanding that subsequent infection may require removal of the material.

WARNING

Trulene mesh is supplied as a sterile product. Re sterilization is not recommended.

RISK OF REUSE

Reprocessing or sterilization of this device may create the risk of contamination and patient infection. Do not reuse, reprocess or resterilize the Trulene mesh.

PRECAUTIONS

Users should be familiar with surgical procedures and techniques involving the use of non-absorbable meshes before employing Trulene mesh in wound repair or reinforcement and brief the patients about contraindications and precautions. A minimum of 6.5 mm inside the edge mesh extend beyond the suture line.

ADVERSE REACTIONS

Potential adverse reactions with Trulene mesh implantation are those typically associated with surgical implantable materials which including seroma formation, inflammation, adhesion formation, fistula formation, extrusion and potentiation of infection.

INSTRUCTIONS FOR USE

Mesh is implanted according to currently accepted surgical mesh procedures.

Some surgeons prefer to suture an uncut section of mesh that is considerably larger than the defect into position over the wound. The opposite sides are then sutured to assure proper closure under correct tension. When the margin sutures have all been placed, the extra mesh is trimmed away.

Adequate mesh fixation is required to minimize post-operative complications and recurrence. The fixation technique method and products used should follow the current standard of care. Careful attention to fixation and spacing will help prevent excessive tension or disruption between the mesh materials and connective tissue.

When fixating with non-absorbable sutures or other mechanical fixation devices, a safe from the edge of the mesh of not less than 6.5mm (¼") must be maintained. 6.5 mm to 12.5mm (¼" to ½") should be left between fixation points.

SUPPLY

Trulene mesh is available in various sizes which are packed in a printed box quantity as indicated on the box label.

STORAGE

Recommended storage conditions: Below 30°C away from moisture and direct sunlight. Do not use after expiry.

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