English Instructions for Use

Suction Filters and Sets.

Intended use: The hydrophobic filter protects surgical suction vacuum pump equipment from contamination.

Cautions

CAUTION: A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided, may result in harm to the patient.

Any use of the medical device requires full understanding and strict observation of all portions of these Instructions for Use. The medical device may only be used for the purpose specified under "Intended Use". Observe all CAUTION statements throughout this manual and all statements on medical device labels. The manufacturer disclaims any liability for patient injury caused by noncompliance with these statements.

Following installation of the filter, check that all connections are secure.

The medical device IFU is not available for each individual filter. Only one copy of the Instructions for Use is included in the distribution package and should therefore be kept in an accessible location for users.

Medical device is only to be used by healthcare professionals.

User has the responsibility to report any serious incident that has occurred in relation to the medical device, to the manufacturer,

EC REP and the competent authority of the Member State in which the user and/or the patient is established.

To avoid contamination and soiling, the product should remain packaged until ready to be used.

Do not use the product if the packaging and/or product is damaged.

If product was kept outside the stated storage temperature, product functionality cannot be guaranteed. Storage Temperature: +5°C to 40°C.

Store the product at room temperature, in dry place, free of dust and not exposed to direct sunlight.

Handle filter with care. Abusive handling can damage the filter.

For Your Safety and that of Your Patients, follow the Instruction for Use of the machine!

Before installation, check that all system components are free of obstructions and of foreign bodies. Otherwise, use is limited, or faulty operation is possible.

Inappropriate connection to other devices can result in a patient hazard.

Medical device products have been designed, tested and manufactured exclusively for disposable use and for a period of use not exceeding 24 hours following the first opening of the primary packaging.

Do not use after expiry date stated on primary packaging.

Manufacturing date

Expiry date

Monitor filter performance during the use. If any changes to filter performance or appearance are observed, change the filter to new one.

Always keep spare filter ready to exchange with the one in use, if any of the below are observed:

- Increase in resistance;
- Filter is visibly soiled or damaged or if it rattles;
- Excessive condensation;
- Associated alarm.

This product is for single patient use only.

The products must not be reused, cleaned or sterilised / not to be sterilised after use in order to re-use. Exceptions from this rule will be stated within the Instruction For Use for the particular products code if applicable. Reuse, cleaning or sterilization may lead to failure of the products, risk infection and injury to the patient.

This product is latex and PHT free.

Batch code symbol.

Prescription only.

Following use, the products must be disposed off in accordance with the local hospital hygiene and waste disposal regulations for hazardous infectious clinical waste.

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EC REP GVS S.p.A. - www.gvs.com
Via Roma, 50 - 40069 Zola Predosa ITALY

REF														
	Product Code; Code du produit; Artikelnummer; Codice prodotto; Código de producto; Productcode; Código de produto; Produktkod; Produktkode; Produktkode; Tuotekoodi; Kod produktu; Код на произукта; Šifra proizvoda; Kód výrobku; Toole kood; Κυβικός προιοντος; Termékkód; Cód Tárge; Izstrādājuma kods; Gaminio kodas; Kodići tal-prodotti; Cod produs; Kód výrobku; Kóda izdelka.	2000/53XXXX Suction Filter	2200/02XXXX Suction Filter	2200/02DIKXXXX Suction Set	2200/06XXXX Suction Filter	2200/16XXXX Suction Filter	2200/21XXXX Suction Filter	2200/26XXXX Suction Filter	2200/35XXXX Suction Filter	2200/36XXXX Suction Filter	2200/55XXXX Suction Filter	2200/70XXXX Suction Filter	2200/902XXXX 90mm Suction Filter	2200/911XXXX 90mm Suction Filter
Product description; Description du produit; Produktbeschreibung; Descrizione prodotto; Descripcón de producto; Productomschrijving; Descrição do produto; Produktbeskrivning; Produktbeskrivelse; Produktbeskrivelse; Tuutteen kuvaus; Opis produktu; Описание на продукта; Opis proizvoda; Popis výrobku; Tootekirjeldus; Περιγραφή προϊόντος; Termékleírás; Tuaris ar Tháirge; Izstrādājuma apraksts; Gaminio aprašymas; Deskripcijoni tal-prodott; Descrereera produsului; Popis výrobku; Opis izdelka. Suction Filters and Sets; Filtres et ensembles d'aspiration; Saugfilter und Sets; Filtri e set di aspirazi-one; Filtrae y conjuntos de succión; Aanzuigfilters en sets; Filtros e conjuntos de sucção; Sugfliter och uppsättningar; Sugeflitr og -set; Sugefiltr og sett; Imusuodattimet ja -sarjat; Filtry i zestawy ssące; Всмукателни филтри и комплекти; Uisäni filtri i pribor; Saci filtry a sady; Imemisfiltrid ja komplektid; Φίλτρα και σετ αναρρόφησης; Szívósűrök és -kész-letek; Foirme agus Scagairí Súcháin; Iesüklšanas filtri un komplekti; Siurbimo filtrai ir rinkiniai; Filtri u Sett-ijet tal-Gbid tal-Arja; Filtre și seturi de aspirare; Sacie filtre a súpravy; Sesalni filtri in kompleti.	Inlet/Outlet Connectors: 8mm-11mm hose barbed both sides, machine side indicated by markings on outer ring. Flow direction indicated by an arrow on the yellow overmould ring of the filter. Approximate dimensions: 54mm diameter x 51mm height. Weight: 11g (approx.). Unidirectional filter.	Inlet/Outlet Connectors: OD 8mm (approx.) hose barbed and ID 4mm (approx.) both sides. Approx. dimensions: 65mm diameter x 51.5mm length. Flow direction indicated by "IN" marking on filter face. Weight 17gm (approx.). Unidirectional filter.	Inlet/Outlet Connectors of 2200/02: OD 8mm (approx.) hose barbed and ID 4mm (approx.) both sides with luer connector. Approx. dimensions of 2200/02: 65mm diameter x 51.5mm length. Approx. dimensions of the tube: 30mm length, 9mm OD x 6mm ID. Luer connector: OD 5mm / ID 3.6mm from the tubing side and OD 7.7mm / ID 4.3mm on the opposite side. Approx. length: 23.8mm. Flow direction indicated by "IN" marking on filter face. Weight of filter: 17gm (approx.). Unidirectional filter.	Inlet/Outlet Connectors: 5 - 9.5mm hose barbed - both sides. Approx. dimensions: 65mm diameter x 65.8mm height. Flow direction indicated by "IN" marking on filter face. Weight 17gm (approx.). Unidirectional filter.	Inlet/Outlet Connectors: 11mm to 15mm hose barbed both sides. Approx. dimensions: 64.7mm diameter x 55.7mm height. Flow direction indicated by "IN" marking on filter face. Weight 17gm (approx.). Unidirectional filter.	Inlet/Outlet Connectors: IN 11mm -15mm hose barbed, OUT 15mm Male ISO. Approx. dimensions: 65mm diameter x 56mm length. Flow direction indicated by "IN" marking on filter face. Weight 17gm (approx.). Unidirectional filter.	Inlet/Outlet Connectors: IN - 8mm hose barbed, OUT - 15mm Male ISO. Approx. dimensions: 65mm diameter x 56mm length. Flow direction indicated by "IN" marking on filter face. Weight 17gm (approx.). Unidirectional filter.	Inlet/Outlet Connectors: IN - 8mm hose barbed, OUT - 7mm extended hose taper. Approx. dimensions: 65mm diameter x 52mm length. Flow direction indicated by "IN" marking on filter face. Weight 17gm (approx.). Unidirectional filter.	Inlet/Outlet Connectors: IN - 7mm hose plain, OUT - 8mm hose barbed. Approx. dimensions: 65mm diameter x 52mm length. Flow direction indicated by "IN" marking on filter face. Weight 17gm (approx.). Unidirectional filter.	Inlet/Outlet Connectors: IN - 8.6mm hose plain, OUT - 12mm hose straight. Approx. dimensions: 65mm diameter x 42.2mm height. Flow direction indicated by "IN" marking on filter face. Weight: 17gm (approx.). Unidirectional filter.	Inlet/Outlet Connectors: IN - 8.1mm conical, OUT - 12mm hose barbed. Approx. dimensions: 65mm diameter x 60.3mm height. Flow direction indicated by "IN" marking on filter face. Weight: 17gm (approx.). Unidirectional filter.	Inlet/Outlet Connectors: 11mm to 15mm hose barbed both sides. Approx. dimensions: 92mm diameter x 61mm height. Flow direction indicated by "IN" marking on filter face. Weight: 47gm (approx.). Unidirectional filter.	Inlet/Outlet Connectors: IN - 11mm hose plain, OUT - 11-15mm hose barbed. Approx. dimensions: 92mm diameter x 61.8mm height. Flow direction indicated by "IN" marking on filter face. Weight: 47gm (approx.). Unidirectional filter.	
 Product supplied non sterile; Produit livré non stérile; Das Produkt wird nicht steril geliefert; Prodotto non sterile; Producto suministrado sin esterilizar; Produkt niet-steriel aangeleverd; Produto fornecido não-esteril; Produkten levereras ej steril; Produkt leveret ikke-sterilt; Produkt leveret ikke-sterilt; Tuote toimittaan ei-steriilinä; Produkt dostarczany w stanie niejałowym; Προϊόντος se доставя нестерилно; Proizvod isporucen nesterilan; Výrobek je dodáván v nesterilním stavu; Toode tamitakse mittesterilisena; Το προϊόν παρέχεται μη στείρο; A termék nem sterilen szállítva; Soláthraithear an táirge neamhsterilúil; Izstrādājums piegādāts nesterils; Gaminys tiekiamas nesterilus; Prodott forunt miux stërlil; Produs furnizat nesteril; Výrobok sa dodáva v nesterilnom stave; Izdelek je dobavljen nesteriliziran.	Yes; Oui; Ja; Ji; Si; Si; Ja; Sim; Ja; Ja; Ja; Joo; Tak; DA; Da; Anò; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Anò; Da.	Yes; Oui; Ja; Ji; Si; Si; Ja; Sim; Ja; Ja; Ja; Joo; Tak; DA; Da; Anò; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Anò; Da.	Yes; Oui; Ja; Ji; Si; Si; Ja; Sim; Ja; Ja; Joo; Tak; DA; Da; Anò; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Anò; Da.	Yes; Oui; Ja; Ji; Si; Si; Ja; Sim; Ja; Ja; Ja; Joo; Tak; DA; Da; Anò; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Anò; Da.	Yes; Oui; Ja; Ji; Si; Si; Ja; Sim; Ja; Ja; Ja; Joo; Tak; DA; Da; Anò; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Anò; Da.	Yes; Oui; Ja; Ji; Si; Si; Ja; Sim; Ja; Ja; Ja; Joo; Tak; DA; Da; Anò; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Anò; Da.	Yes; Oui; Ja; Ji; Si; Si; Ja; Sim; Ja; Ja; Ja; Joo; Tak; DA; Da; Anò; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Anò; Da.	Yes; Oui; Ja; Ji; Si; Si; Ja; Sim; Ja; Ja; Ja; Joo; Tak; DA; Da; Anò; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Anò; Da.	Yes; Oui; Ja; Ji; Si; Si; Ja; Sim; Ja; Ja; Ja; Joo; Tak; DA; Da; Anò; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Anò; Da.	Yes; Oui; Ja; Ji; Si; Si; Ja; Sim; Ja; Ja; Ja; Joo; Tak; DA; Da; Anò; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Anò; Da.	Yes; Oui; Ja; Ji; Si; Si; Ja; Sim; Ja; Ja; Ja; Joo; Tak; DA; Da; Anò; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Anò; Da.	Yes; Oui; Ja; Ji; Si; Si; Ja; Sim; Ja; Ja; Ja; Joo; Tak; DA; Da; Anò; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Anò; Da.	Yes; Oui; Ja; Ji; Si; Si; Ja; Sim; Ja; Ja; Ja; Joo; Tak; DA; Da; Anò; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Anò; Da.	Yes; Oui; Ja; Ji; Si; Si; Ja; Sim; Ja; Ja; Ja; Joo; Tak; DA; Da; Anò; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Anò; Da.
Ambient Conditions: Conditions ambiantes : Umgebungsbedingungen : Condizioni ambientali : Condiciones ambientales : Omgevingscondities : Condições ambientais : Omgivningsförhållanden : Omgivelsetsbetingelser : Krav til omgivelsene : Ympäristöolosuhteet : Warunki otoczenia : Условия на околната среда; Uyjetil okoline; Okolini podmnky; Keskkonnatingimused; Συνθήκες περιβάλλοντος; Környezeti feltételek; Coonniollacha Comhthimpeallacha; Apkārtējie apstākļi; Aplinkos sąlygos; Kundizacijonijet; Ambjentali; Conditiì ambientale; Oskoltè podmienky; Okolje.	Temperature during storage; Température de stockage; Temperatur bei Lagerung; Temperatura durante l'immagazzinamento; Temperatura durante el almacenamiento; Temperatuur tijdens opslag; Temperatura em armazenamento; Temperatur under förvaring; Temperatur under opbevaring; Temperatur under oppbevaring; Lämpötila varastoinnin aikana; Temperatura podczas składowania; Температурна по архене на съхранение; Temperatura tijekom skladištenja; Teplota během skladování; Temperatuur ladustamise ajal; Θερμότητα κατά την αποθήκευση; Hőmérséklet tárolás közben; Teocht le linn na stórála; Temperatura uzglabāsanas laikā; Temperatūra sandėliavimo metu; Temperatura waqt il-hażna; Temperatura în timpul depozitării; Teplota počas składowania; Temperatura med skladiščenjem.	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)
Performance Data: Caractéristiques ventilatoires : Leistungsdaten : Dati di rendimento : Características de rendimiento : Karakteristieke waarden : Dados de desempenho : Prestandauppfigter : Ydelse : Kapasitetsdata : Suorituskykytiedot : Parametry sprężu; Даниа за ефективността; Podaci o izvedbi; Údaje o výkonu; Toimivusandmed; Δεδομένα απόδοσης; Teljesítmény adat; Sonrai Feidhmíochta; Veikispējas dati; Pasirodymo data; Dejlta tai-Prestazjonji; Date de performanță; Údaje o výkonnosti; Podatki o uspešnosti.	Air Flow; Flux d													