DIATERMO MB 200 T

ELECTROSURGICAL UNIT

USER MANUAL

SUMMARY

IMPORTANT	3
INTRODUCTION	4
STANDARD AND OPTIONAL COMPOSITION	4
GENERAL DESCRIPTION	6
INTENDED USE	7
TARGET USER	7
TARGET PATIENT POPULATION	7
ELECTROPHYSICAL PRINCIPLES	8
OPERATIVE TECHNICS	11
MONOPOLAR CUT	11
MONOPOLAR COAGULATION	11
BIPOLAR COAGULATION	12
CONTRAINDICATIONS	13
SAFETY	14
GENERAL	16
INSTALLATION	19
PATIENT SAFETY	20
CORRECT PATIENT POSITIONING	20
HIGH-FREQUENCY ELECTROSURGERY IN LAPAROSCOPY	22
PUTTING INTO SERVICE	24
CONNECTION AND USE OF ACCESSORIES	25
CONNECTORS AND CONTROLS	26
LABEL ON REAR PANEL	26
TECHNICAL DATA	26
MEANING OF GRAPHIC SYMBOLS	27
BOX LABEL	28
FRONT PANEL	29
OPERATING MODES	30
CONTROL AND POWER ON	30
PROGRAMS	31
SURGERY	32
MONOPOLAR SECTION	33
BIPOLAR SECTION	35

AUTOSTART AND AUTOSTOP	35
DISPENSING SCREEN.....	37
FLASH TIMER SECTION	37
PROGRAMS SECTION	38
NEUTRAL ELECTRODE CONTROL	39
LIGHTS	40
SETTINGS.....	41
UPDATE	42
BACK PANEL	43
POWER SUPPLY MODULE OF THE EQUIPMENT AND VOLTAGE SELECTOR	43
POWER SWITCH.....	43
TECHNICAL CHARACTERISTICS	44
HARDWARE REQUIREMENTS	45
MAINTENANCE.....	46
GENERAL.....	46
CLEANING THE ENCLOSURE.....	46
CLEANING AND STERILIZATION OF ACCESSORIES	46
GUIDE TO TROUBLESHOOTING	47
REPAIRS	47
FUSE REPLACEMENT	47
PRE-USE UNIT CHECK.....	48
SAFETY AND MEASUREMENT FUNCTION CHECKS	48
DIAGRAMS	49

IMPORTANT

This instruction manual represents an essential part of the HF electrosurgery unit and must be kept available to the operating staff at all times. It is of paramount importance to carefully read and fully understand all instructions and guidelines before attempting any use of an active electrode.

It is imperative to diligently follow all warnings and safety instructions. Please ensure that this documentation is included with the device in case it is transferred to another team.

In the event of technical assistance being required, please contact LED SpA.

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Code	Description	DIATERMO MB 200 T
310-110-05	BIPOLAR - Bipolar Forceps 11,5cm TIP0.5mm	○
310-112-05	BIPOLAR - Bipolar Forceps Curved 11,5cm TIP0.5mm	○
310-140-10	BIPOLAR - Bipolar Forceps 20cm TIP 1mm	○
310-140-20	BIPOLAR - Bipolar Forceps 20cm TIP 2mm	○
310-142-10	BIPOLAR - Bipolar Forceps Curved 20cm TIP 1mm	○
310-142-20	BIPOLAR - Bipolar Forceps Curved 20cm TIP 2mm	○
310-180-10	BIPOLAR - Bipolar Forceps Angled 20cm TIP 1mm	○
310-180-20	BIPOLAR - Bipolar Forceps Angled 20cm TIP 2mm	○
310-182-10	BIPOLAR - Bipolar Forceps Angled Curved 20cm TIP 1mm	○
310-185-10	BIPOLAR - Bipolar Forceps Angled Curved 20cm TIP 1mm	○
310-510	BIPOLAR - Bipolar electrode 20cm – direct	○
310-550	BIPOLAR - Bipolar electrode 20cm – curved	○
310-590	BIPOLAR - Bipolar electrode 20cm – curved 2	○
330-134-20	MONOPOLAR - Monopolar Forceps 20cm TIP2mm	○
330-160	MONOPOLAR - Monopolar Surgical Scissors 18cm	○
500500.L1	ELECTRODE - Straight thin wire electrode (5pcs) 5cm	○
500500.L1/L	ELECTRODE - Straight thin wire electrode (5pcs) 10cm	○
500500.L10	ELECTRODE - Bent ball electrode ø 3mm (5pcs) 5cm	○
500500.L10/L	ELECTRODE - Bent ball electrode ø 3mm (5pcs) 10cm	○
500500.L11	Needles for micro-surgery (10Pcs)	○
500500.L2	ELECTRODE - Bent thin wire electrode (5pcs) 5cm	○
500500.L2/L	ELECTRODE - Bent thin wire electrode (5pcs) 10cm	○
500500.L3	ELECTRODE - Loop electrode ø 4mm (5pcs) 5cm	○
500500.L3/L	ELECTRODE - Loop electrode ø 4mm (5pcs) 10cm	○
500500.L4	ELECTRODE - Loop electrode ø 8mm (5pcs) 5cm	○
500500.L4/L	ELECTRODE - Loop electrode ø 8mm (5pcs) 10cm	○
500500.L5	ELECTRODE - Bent hook electrode (5pcs) 5cm	○
500500.L5/L	ELECTRODE - Bent hook electrode (5pcs) 10cm	○
500500.L6	ELECTRODE - Bent thick wire electrode (5pcs) 5cm	○
500500.L6/L	ELECTRODE - Bent thick wire electrode (5pcs) 10cm	○
500500.L7	ELECTRODE - Drop electrode (L7) (5pcs) 5 cm	○
500500.L7/L	ELECTRODE - Drop electrode (L7) (5pcs) 10cm	○
500500.L8	ELECTRODE - Noose electrode (L8) (5pcs) 5 cm	○
500500.L8/L	ELECTRODE - Noose electrode (L8) (5pcs) 10cm	○
500500.L9	ELECTRODE - Straight ball electrode ø 3mm (5pcs) 5cm	○
500500.L9/L	ELECTRODE - Straight ball electrode ø 3mm (5pcs) 10cm	○
6429A	NEUTRAL - Steel Neutral Electrode 24x16cm	○
755VL	Disposable handle with finger switches (F4797)	○
F7520	Electrode cleaning sponge 47x50mm	○
F7915	Conductive rubber neutral electrode without cable	○
F7930	Conductive rubber split neutral electrode without cable	○
TR003	Trolley 3 Shelves	○
TR003W	Trolley 3 Shelves wide	○
TR004	Trolley 4 Shelves	○
TR005	Trolley 5 Shelves	○
TR005W	Trolley 5 Shelves wide	○

■/ PZ = STANDARD ○= OPTIONAL

GENERAL DESCRIPTION

The high-frequency electrosurgical device is capable of providing current for monopolar cutting, coagulated cutting with different levels of coagulation, both in monopolar and bipolar modes. In bipolar coagulation mode, it is possible to activate a tissue impedance detection system with automatic start and automatic stop once coagulation is achieved (AUTOSTART - AUTOSTOP functions). Additionally, a specific function can be used to perform vascular and venous synthesis and coagulation by sealing vessels with radiofrequency.

A total of ten different operating modes and power levels can be stored and recalled for easy use. You have the option to use single-plate neutral reference electrodes or electrodes with a divided conductive zone to monitor patient impedance stability during surgery.

Unit control is through the display. The power input is located on the rear panel. The units are equipped with automatic control systems that monitor internal parameters and signal any damage or errors. The operational parameters used are constantly stored so that the last selected parameters are recalled each time the unit is turned on or the operating method is changed.

Volume sound adjustment is customizable, allowing each operator to choose the most suitable level for working environmental conditions.

Units can be controlled through handpieces with or without buttons, thanks to a double pedal command. It is possible to connect a bipolar forceps to the unit for bipolar functions.

INTENDED USE

Medical device intended for temporary use for surgical operations in which cutting and/or coagulation of soft tissues is required, with a monopolar and/or bipolar technique, for survey minor and/or major in open and/or intra-operative percutaneous and/or endoscopic and/or laparoscopic.

The **DIATERMO MB 200 T** equipment is conceived for being used in the following sectors:

Description	DIATERMO
	MB 200T
Electrosurgical unit code	GMA10100.T40
Casualty Surgery	●
Dermatology	●
Dental	-
Endoscopy	●
Gastroenterology	●
General Surgery	●
Gynecology	●
First Aid	●
Neurosurgery	●
Orthopedics	●
Otorhinolaryngology	●
Pediatric Surgery	●
Plastic Surgery	●
Pneumology	●
Thorax Surgery	●
Trans Urethral Resection (TUR)	●
Urology	●
Vascular Surgery	●

● = Usable
- = Not Usable

INTENDED USER

Device for professional use. Use of the equipment is restricted to medical personnel with medical degrees specializing in high frequency electrosurgery.

INTENDED PATIENT POPULATION

Adults - Men and women (≥18 years), excluding patients in the contraindications section.

ELECTROPHYSICAL PRINCIPLES

In electrosurgical interventions the traditional use of surgical blade is substituted by an electrosurgical needle that allows for fast and effective cut and coagulation of the targeted tissue.

The electrosurgical needle operates on the principle of converting electrical energy into heat and consists of the following components:

- A radiofrequency sinusoidal oscillator (0.4 - 4MHz).
- A wave packet generator with a packet repetition frequency of 15 – 30 kHz.
- A mixer for transferring the waveform to the power amplification block, either for cutting, coagulation, or a signal obtained from an appropriate combination of the two.
- A power amplifier block capable of supplying the required power in terms of current and transmitting the amplified signal to the electrodes through a transformer.
- A safety circuit for the return electrode, designed to detect any cable interruptions and deactivate the radiofrequency delivery.
- A specially shaped active electrode (handpiece).
- A return electrode (neutral) that completes the circuit through the patient.

Electric current flowing through biological tissue usually can cause:

1. *Joule Effect*
2. *Faradic Effect*
3. *Electrolytic Effect*

1. *Joule Effect*

In biological tissue, when passed through by the electric current delivered by the electrosurgical scalpel, heating (Joule effect) is produced, which is dependent on tissue-specific electrical resistance, current density, and application time and can result in various cellular transformations.

$$Q = I^2 \times R \times T$$

The influence of the thermal effect (Joule effect) is realized by:

- **Current Intensity and output power**
- **Modulation level**
Parameters that can be interpreted from the waveform of the high-frequency current produced by the generator.
- **Electrode shape**
Pointed or rounded as required, it is very small in size; therefore, the current density on the tip surface [$A \cdot m^{-2}$] is very high. Thin-section electrodes create a 'high current density, and high temperature, promoting cutting action. Those with a large surface area create a lower current density, and lower temperature, realizing a coagulation effect.

- **State of active electrode**

Thermal effects can be related to the resistance of the human body to which the contact resistance of the electrode must be added. It is essential to keep the active electrodes perfectly clean in order not to have a reduction in the effects.

- **Characteristics of the tissue**

The resistive characteristics change according to the biological tissues.

Biological Tissue (range from 0,3 to 1 MHz)	Metals
Blood $0,16 \times 10^3 \Omega$	Silver $0,16 \times 10^{-5} \Omega$
Muscle, kidney, heart $0,2 \times 10^3 \Omega$	Branch $0,17 \times 10^{-5} \Omega$
Liver, spleen $0,3 \times 10^3 \Omega$	Gold $0,22 \times 10^{-5} \Omega$
Brain $0,7 \times 10^3 \Omega$	Aluminium $0,29 \times 10^{-5} \Omega$
Lung $1,0 \times 10^3 \Omega$	
Fat $3,3 \times 10^3 \Omega$	

(Example of specific resistances of organic and metallic materials)

Based on the temperature achieved and according to the pulse forms used, different techniques of using radiofrequency current on the human body are recognized as follows:

- **Coagulation**

Temperatures of 60 to 70 °C in the area around the active electrode cause slow heating of the intracellular fluid, the water contained in the cell evaporates, and a clotting action is achieved that stops bleeding.

- **Cut**

Temperatures above 100 °C in the area surrounding the active electrode result in the vaporization of the intracellular fluid and explosion of the cell. The vapor present around the electrode triggers an intercellular chain reaction in the direction in which the active electrode is handled, also transmitting the vaporization energy to the immediately surrounding tissues.

Electrotomy is, therefore, not mechanical resection. If the temperature reaches 500 °C, tissue carbonization occurs with a cauterizing action.

- **Mixed currents**

These are obtained by combining the effects of coagulation and electrotomy. A reduction in bleeding occurs during a cutting procedure, or as a cut that develops a substantial layer of eschar. The high frequencies used by the electrosurgical scalpel, however, do not allow the electromagnetic field to penetrate matter and cause the current to flow through the conductor more on the outermost surface, decreasing exponentially and becoming negligible in the center of the conductor's cross-section. This effect, called the 'skin effect,' results in a decrease in the useful cross-sectional area for the passage of a current, and an increase in the electrical resistance of the material, and becomes a major problem in the neutral electrode. In fact, in this electrode the current density is very high (KA/m^2) at the edge, where excessive temperature rise due to the 'Joule effect' causes burns to the patient. Therefore, it is no accident that burns to the patient, which has occurred

in surgical procedures, have the shape of the edge of the neutral electrode. To reduce the risk of burns, it is necessary to dose the delivered power ($I^2 \cdot t$) appropriately and follow the rules for applying the neutral electrode to the patient (see SAFETY chapter).

2. Faradic Effect

Pulsed electric current causes neuro-muscular stimulation, originating from the stimulation of the physiological process of ion exchange, which is responsible for the transmission of stimuli that cause muscle spasms and cardiac phenomena of extrasystole and ventricular fibrillation.

The effect of these stimuli is known as the faradic effect and is expressed by:

$$R = I / \sqrt{F}$$

The physiological system of stimulus transmission follows a limiting curve in which pulsed or low-frequency currents generate a pacing pulse. With the high-frequency alternating current (above 200 kHz), which is used in electrosurgery, there are no neuromuscular reactions (the change of polarity is so fast that it does not affect the patient at the level of neuro-muscular reactions), let alone electrolyte damage to the body.

For this reason, all high-frequency generating equipment for surgical use (electrosurgery) works on base frequencies above 300 kHz so as not to introduce electrical stimulation.

3. Electrolytic Effect

The use of high-frequency currents reduces the electrolytic effect (ionic separation) in tissues due to the very short unidirectional conduction period of the current.

OPERATIVE TECHNICS

MONOPOLAR CUT

Monopolar cutting is a technique used to dissect biological tissue by directing a high-frequency, high-density electrical current from the active electrode's tip. This current generates intense molecular heat within the cells, causing them to rupture. As the electrode moves through the tissue, it shears and destroys cells sequentially. This movement prevents the lateral spread of heat within the tissue, confining the destruction to a single line of cells.

For precise cutting with minimal thermal impact and limited hemostasis, a pure sine wave current without modulation is preferred. This allows for precise control and safe use without harming adjacent bone. However, some level of modulation in the current is desirable to facilitate effective coagulation during the cutting process, making electrosurgery an advantageous technique.

The following rules help the operator achieve a good cut:

- Keep the tissue moist but not wet.
- Maintain the electrode perpendicular to the tissue.
- Activate the output circuit before making contact with the tissue.
- Keep the electrode tip clean (for this purpose, you can use optional electrode-cleaning sponges with code F7520).
- Allow the tissue to cool before making a new cut.

When the output power level is adequate, you can expect to achieve:

- No resistance when moving the electrode through the tissue.
- No change in the color of the cut surfaces.
- No residual tissue fibers on the electrode.

MONOPOLAR COAGULATION

When an increase in temperature occurs, caused by the heat generated through the Joule effect in the tissue, thermal coagulation takes place. This process involves the partial solidification of organic liquids and the precipitation of colloidal substances. Specifically, within the blood, fibrin forms, which, as it solidifies, can obstruct blood vessels.

To achieve coagulation using an electrosurgical scalpel, it's essential to provide the active electrode with an intermittent current. This prevents excessive heat generation, which could lead to cell explosion and tissue cutting, allowing for controlled heating instead. This controlled heating causes the water within the cells to escape without destroying them. However, even with intermittent current, if the current intensity is too high, it can still result in a cutting effect.

Active electrodes well-suited for coagulation purposes include those with ball, plate, or lance shapes, used laterally.

Coagulation can be achieved through two different methods:

- **Coagulation by Desiccation**

This is achieved by supplying the electrode with low voltages to prevent the generation of sparks (ensuring that the action obtained is pure coagulation without any cutting effect). The electrode is placed in direct contact with the tissue, and the amount of heat generated on contact dries it out. Typically, coagulated cellular surfaces act as an insulating layer, preventing heat from subsequent current applications from penetrating too deeply. The current normally used for coagulation is modulated. Depending on the modulation percentage, you get a balance between precision in cutting, effectiveness in hemostasis, and tissue destruction. Greater modulation of the current leads to a more jagged cut, greater depth of tissue destruction, but more effective coagulation.

The following rules help the operator achieve good coagulation:

- Select a ball electrode or a thick wire electrode.
- Locate the bleeding vessel after excess blood has been dried from the area.
- Gently touch the bleeding vessel before activating the electrode.
- Cease electrode activation as soon as the tissue whitens to avoid damaging it.
- Keep the electrode tip clean (you can use optional electrode-cleaning sponges with code F7520).

- **Coagulation with Anatomical Forceps by Clamping**

The most commonly used coagulation technique involves blocking the blood flow by applying clamping pressure at the end of the forceps. After clamping the tissue portion or the blood vessel where coagulation is needed, the active electrode is brought into contact with the proximal metal part of the forceps. The high-frequency activation must occur after this contact (forceps - active electrode) to avoid the Faradic effect (initiation of an electrical discharge that uses air as a conductor), which could cause electrical shock, burns to the operator, etc.

BIPOLAR COAGULATION

In contrast to the monopolar technique, the bipolar technique focuses the high-frequency current on a very small portion of tissue. This method employs bipolar forceps with various sizes and shapes on their distal ends, serving as both the active and neutral electrodes. By clamping the tissue to be operated on between these forceps, the high-frequency current flows from one end to the other, using the tissue itself as an electrical bridge.

Bipolar coagulation, achieved using this technique, effectively controls bleeding in small blood vessels within body tissues situated between the two clamp tips. When the current density is reduced, it primarily dries the cell surface without deep penetration, resulting in coagulation.

The bipolar technique is considered safer due to the predictable and consistent direction of the high-frequency current. It eliminates the uncertainties and potential errors associated with unknown current paths, and it requires much lower power levels compared to the monopolar technique. Consequently, it is commonly used in delicate surgical procedures.

It is crucial to maintain the cleanliness of the distal ends of the forceps during surgery as they tend to accumulate coagulated tissue, which can impede current flow and cause tissue adhesion. While the use of a neutral electrode, mandatory in the monopolar technique, is not necessary in bipolar procedures, it is often allowed for practical reasons during the initial preparation phase.

CONTRAINDICATIONS

The use of electrosurgery is contraindicated in patients:

- pacemaker carriers
- with stimulation electrodes
- with metallic prostheses
- with serious blood pressure imbalances
- with serious diseases of the nervous system
- with serious kidney failure
- in state of pregnancy.

In the context of electrical surgery, burns due to high frequency are the main causes of burns caused to the patient, but they are not the only ones involved. One can also get necroses by compression, allergic reactions to disinfectants, gas sparks or flammable liquids.

Some of the causes of burns are to be attributed to:

- insufficient medical equip training about all modalities to avoid or reduce the risks of burns by using HF electrosurgical units
- use of disinfectants with high alcohol content
- incorrect position of the patient during the electrosurgical operation
- contact between active electrode and the skin
- contact with liquid
- long application of HF currents
- incorrect application of the patient-plate.

To avoid or reduce the risks associated with the use of high frequency electrosurgery, it is necessary to respect the rules and safety measures illustrated in the following chapter.

SAFETY

WARNING: Electrosurgery carries inherent risks, and improper usage of any component within the electrosurgical system can result in severe burns to the patient. It is absolutely crucial that you meticulously read and comprehensively understand all instructions before attempting to utilize an active electrode. Neither the manufacturer nor any retailers can be held liable for any harm or damage, whether direct or indirect, caused to individuals or equipment due to the incorrect use of the device and its accompanying accessories.

The accessories provided with this unit are designed to be compatible specifically with this unit and may not work with other electrosurgical units. Prior to connecting any additional accessories to this unit, the user must confirm that these accessories possess insulation characteristics that align with those of this unit (please refer to the "TECHNICAL CHARACTERISTICS" section for details). It is strongly allowed to assess the packaging integrity of sterile accessories before their initial use.

ATTENTION


- **DO NOT USE** on patients with electronic implants such as cardiac pacemakers without consulting a qualified professional (e.g., a cardiologist). There is a potential risk of interference with the functioning of the electronic implant or damage to the implant itself.
- **DO NOT USE** in the presence of flammable anaesthetics or oxidizing gases (such as nitrous oxide (N₂O) and oxygen) or near volatile solvents (such as ether or alcohol) as explosions may occur.
- **DO NOT PLACE** instruments near or in contact with flammable materials (such as gauze or surgical drapes). Activated or heated instruments can cause fires.
- When not in use, store instruments in a clean, dry, and highly visible area away from direct patient contact. Inadvertent contact with the patient can result in burns.
- **INSPECT** instruments and cables for damage before each use, especially the insulation of laparoscopic/endoscopic instruments. This inspection can be carried out visually under magnification or with a high-voltage insulation testing device. Insulation failures can lead to burns or other injuries to the patient or the operator.
- The surface of the active electrode may remain sufficiently hot to cause burns even after RF current is deactivated.
- Due to concerns about the potential carcinogenic and infectious properties of electrocautery byproducts (such as tissue smoke plumes and aerosols), protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures.
- Only connect adapters and accessories to the electrosurgical unit when the power is **OFF**. Failure to do so may result in patient or operating room personnel injury or electric shocks.
- If the device is powered with argon, warnings regarding gas embolisms must be included.
- If the instrument is reusable, a warning should be included that visual inspection alone may not be sufficient to ensure intact insulation.
- **DO NOT ACTIVATE** the instrument when it is not in contact with the target tissue, as this could cause injuries due to capacitive coupling with other surgical equipment.
- **ASPIRATE** fluids from the area before activating the instrument. Conductive fluids (e.g., blood or saline) in direct contact with or in proximity to an active electrode can carry electrical current or heat away from the target tissues, potentially causing unintended patient burns.
- **DO NOT USE** with hybrid systems, i.e., a combination of metal and plastic, when using monopolar active components. This can result in burns at alternative sites due to capacitive coupling. Use only all-metal or all-plastic systems.
- Before increasing the intensity, verify the adhesion of the neutral electrode and its connections. Apparent low power or device malfunction at normal operating settings may indicate improper neutral electrode application or poor contact in its connections.

- This unit has a CQM system; please note that the loss of secure contact between the neutral electrode and the patient will not trigger an alarm unless a compatible monitoring neutral electrode (split neutral electrode) is used.
- **CAUTION:** Set the intensity to the lowest level necessary to achieve the desired effect.
- **CAUTION:** Keep the active electrodes clean. Accumulated eschar may reduce the tool's effectiveness. Do not activate the instrument during cleaning. Operating room personnel may be injured.
- Any serious incidents related to the device must be reported to LED SpA (via Selciatella n.40, 04011 Aprilia (LT) ITALY) and the competent authority:
Ministero della salute – Direzione generale dei dispositivi medici e del servizio farmaceutico
Viale Giorgio Ribotta, 5 – Roma
E-mail: segr.dgfdm@sanita.it
Tel.: +39 06 5994 3199 / +39 06 5994 3207

GENERAL


The following precautions are crucial for minimizing the risk of inadvertent burns:

- Ensure the secure and complete attachment of the neutral electrode to the patient's body, preferably at the extremities and as close to the surgical site as possible. Avoid connecting the neutral electrode to bony protrusions, prosthetic devices, areas with scar tissue, regions susceptible to fluid accumulation, or areas with a thick layer of subcutaneous fat. The application site should be free from hair, dry, and clean. Avoid using alcohol for skin cleaning. Except for veterinary medicine applications, refrain from using electrode gel.

- When using single-use neutral electrodes, always adhere to the provided expiry dates.
- For multi-use electrodes, ensure that the fixation systems in place guarantee stability during use.
- When applying the neutral electrode, avoid a transverse path and instead favour a vertical or diagonal path, especially when using a bipartite neutral electrode. This helps distribute current evenly across the surface of the neutral electrode and reduces the risk of patient burns.
- If you encounter difficulty in correctly applying the neutral electrode, consider using the bipolar technique instead of the monopolar approach, if feasible.
- To prevent the patient from coming into contact with earthed metallic parts or components with significant grounding capacity (such as an operating table or supports), it is allowed to use an antistatic drape for this purpose.

- To avoid skin-to-skin contact (e.g., between the arm and trunk, between the legs, or on the breasts), it's allowed to insert dry gauze. Additionally, ensure that body areas prone to profuse sweating are kept dry.

1. Active Electrode – 2. Neutral Electrode

3. Dry Gauze – 4. Antistatic Drape

- When using both an electrosurgical scalpel and a physiological monitoring device on the same patient, you must place all monitoring electrodes as far away from the surgical electrodes as possible. Needle monitoring electrodes are discouraged. In any case, it is allowed to use monitoring systems that incorporate high-frequency current-limiting devices.
- Position surgical electrode cables in a manner that prevents contact with the patient or other conductive materials. Active electrodes that are not in use should be isolated from the patient.
- Consider utilizing bipolar techniques when operating on body parts with a relatively small cross-sectional area. This helps prevent unintended coagulation.
- Set the output power level to the lowest effective setting for the intended purpose, minimizing the risk of excessive tissue damage.
- If the electrosurgical unit exhibits an obvious low output level or operates incorrectly, even when set up for normal power delivery, this could indicate issues with the application of the neutral electrode or poor contact in the neutral electrode connections. Therefore, it is essential to verify the proper placement and connections of the neutral electrode before considering higher power settings.
- Avoid the use of flammable anesthetics or oxidizing gases, such as nitrous oxide (N₂O) and oxygen, especially in chest or head operations, unless they can be safely aspirated. Whenever possible, opt for non-flammable substances for cleaning and disinfection purposes. If flammable substances are used for cleaning, disinfection, or as solvents for adhesives, they should be allowed to completely evaporate before using the electrosurgical unit. There is a risk of flammable solutions accumulating under the patient or in cavities like the umbilicus and vagina. Any fluid in these areas should be removed before using the device. It's important to consider the presence of endogenous gases as well.

- Be aware that certain materials, such as cotton wool or gauze impregnated with oxygen, may ignite due to sparks produced by the appliance under normal conditions. Take necessary precautions to prevent such incidents.
- Patients with pacemakers or pacing electrodes are at risk of interference with their pacemaker's functionality or potential pacemaker damage when exposed to electrosurgical equipment. If any uncertainty arises, consult the cardiology department.
- Electrosurgical equipment emits high-frequency energy radiation that can affect other medical devices, unrelated electronics, telecommunications systems, and navigation systems. To prevent interference, you must maintain a minimum distance of at least 1.5 meters between the electrosurgical equipment and other devices.
- Regularly inspect accessories, with special attention to electrode cables and any endoscopy accessories, to ensure there is no damaged insulation or other defects that could compromise their safety or effectiveness.
- To connect accessories compatible with the equipment's characteristics, you must compare the insulation characteristics of the accessories (information provided by the manufacturers) with the specifications of the supplied unit (as outlined in the Technical Characteristics section). This step ensures proper compatibility and safe operation.
- **Caution:** Equipment failure could lead to an inadvertent increase in power output.
Note: Stimulation of the patient's muscles or nerves may be caused by low-frequency currents resulting from electrical sparking between the electrodes and the patient's tissue. If neuromuscular stimulation occurs during surgery, take the following steps:
 1. Pause the surgery immediately.
 2. Thoroughly inspect all connections to the generator to identify any potential issues or loose connections.
 3. If the problem persists and cannot be resolved through connection checks, it is imperative to have the generator inspected by qualified personnel for necessary maintenance and troubleshooting.

INSTALLATION

- Electrical safety is guaranteed only when the equipment is correctly connected to a reliable power supply network with proper grounding, in compliance with current safety standards. It is essential to ensure this fundamental safety requirement, and if there are any doubts, seek a comprehensive inspection of the system by qualified personnel. The manufacturer cannot be held responsible for potential damage caused by the absence of an efficient earth connection in the installation. Operating the equipment without a protective earth connection is strictly prohibited.
- Before connecting the equipment, verify that the provided voltage, as indicated on the rear panel, matches the voltage available in your mains power supply.
- In the event of any incompatibility between the power socket and the equipment's power cable, only replace it with a suitable type. It is not allowed to use adapters, multiple connections, or extension cables. If their use becomes necessary, it is mandatory to employ single or multiple adapters that comply with current safety standards.
- Protect the equipment from exposure to outdoor elements such as rain and direct sunlight. The apparatus must be shielded to prevent the infiltration of liquids.
- Do not keep the equipment plugged in unnecessarily. Turn it off when it is not in use to conserve energy and ensure safe operation.
- This equipment is not designed for use in explosive environments. Avoid using it in such environments where there may be a risk of ignition or explosion.
- The equipment should be used only for its intended purpose. Any other use should be considered improper and potentially dangerous. The manufacturer cannot be held responsible for any damage resulting from improper, incorrect, or unreasonable use.
- Modifying or attempting to modify the equipment is dangerous and should not be done. Altering the characteristics of the equipment can lead to unsafe operation and potential hazards.
- Before performing any cleaning or maintenance procedure, disconnect the appliance from the electrical supply by either unplugging it from the mains or turning off the main switch of the system.
- In the event of equipment breakage or malfunction, power it off immediately. For any necessary repairs, seek assistance only from an authorized service center and request the use of original spare parts. Failure to adhere to these regulations may jeopardize equipment safety and pose risks to the user.
- Do not reduce or disable the acoustic signal indicating generator activation. A functioning activation signal can help minimize or prevent injuries to patients or personnel in the event of accidental activation.
- Do not test the equipment's operation by generating power between the active and neutral electrode or between the active electrode and metal parts. Testing in this manner can be unsafe.
- If required, use appropriate fume extraction methods in the surgical field to manage the release of smoke or fumes generated during procedures.

ATTENTION: When using the equipment in an operating room, it is essential to utilize only immersion-tight foot switches (such as code 00304.00 for a single watertight pedal or code 00305.03 for a double watertight pedal). This ensures safety during surgical procedures.

PATIENT SAFETY

During high-frequency electrosurgery procedures, it's crucial to understand that the patient becomes a conductor of electrical voltage relative to the earth potential. Consequently, if there is contact between the patient and electrically conductive objects (such as metal objects, damp or wet sheets, cloths, etc.), it can result in the generation of electric current at the point of contact.

It is allowed to conduct thorough inspections of the device and its accessories before use, and to strictly adhere to all pertinent safety regulations.

CORRECT PATIENT POSITIONING

Prevent any deliberate or accidental contact between the patient and grounded metal components by ensuring the following:

- The patient does not come into contact with metal parts such as the operating table or supports.
- Ensure that respirator hoses do not rest on the patient's body.
- Always maintain coatings on the grounded operating table to dissipate electrostatic charges.
- Place the patient on a thick fabric with insulating properties, covered with an adequate number of layers.
- Ensure the patient does not touch damp sheets or mattresses.
- Prevent body secretions, cleaning fluids, or other liquids from soaking dry sheets.
- Keep the area beneath the patient free of liquid residue.
- Employ catheters to manage urinary excretions.
- Keep areas of the body prone to increased sweating or areas with skin-to-skin contact points dry using drapes (e.g., arm/body trunk, leg/leg, breasts, skin folds).
- Properly insulate all conductive and grounding supports and brackets.
- Adjust the anesthesia dosage to prevent excessive sweating.


CORRECT APPLICATION OF THE NEUTRAL ELECTRODE

In monopolar electrosurgery, the use of a neutral electrode, also known as a current leakage plate, is essential. It facilitates the safe return of the cutting or coagulation current to the electrosurgical unit. There are two types of neutral electrodes:

1. **Monopolar Neutral Electrode:** In this type, there is no control over the contact between the neutral electrode and the patient.
2. **Bipartite Neutral Electrode:** This type offers control over the contact between the neutral electrode and the patient.


Ensuring the correct placement of the neutral electrode is of utmost importance to prevent burns and minimize patient risks. Below are some valuable tips to achieve this:

1. Correct positioning

In the image alongside, the correct positioning of the split neutral electrode is illustrated. The patient-plate should be placed perpendicular to the surgical field. Avoid positioning it transversely and instead, prefer a vertical or diagonal orientation. This promotes a uniform distribution of the current over the surface of the neutral electrode, minimizing the risk of burns to the patient.

2. Incorrect Positioning

In the image alongside, the incorrect positioning of the split neutral electrode is illustrated. The parallel arrangement between the patient-plate and the surgical field causes a non-uniform distribution of current across the two surfaces of the neutral electrode, potentially leading to alarm notifications on the unit and preventing the correct activation of the device.

For both single-part and dual-part electrodes, before proceeding with the placement of the neutral electrode, clean and remove any residues of foreign substances from its surface.

Do not apply the neutral electrode on scars, bony prominences, or anatomical areas where prosthetic implants or monitoring electrodes are present. Instead, apply it on well-irrigated tissues, such as muscles and in proximity to the surgical site.

If a disposable neutral electrode is being used, adhere to the expiration dates. If a reusable neutral electrode is used, ensure that the fastening systems provide stability.

It is of paramount importance that the neutral electrode is firmly applied over its entire surface to prevent burns. When a neutral electrode partially detaches from the patient, the current density in the remaining electrode area increases. As the current density beneath the neutral electrode becomes uneven, non-uniform heating occurs, especially at the edges of the neutral electrode.

If the electrode were placed in an area subjected to pressure during the procedure, the compressive load would result in reduced skin perfusion. Consequently, the generated heat can only be partially dissipated, thereby increasing the risk of burns. Furthermore, there is an increased risk of pressure points (decubitus) formation due to the heating that occurs. This temperature rise leads to a higher demand for oxygen (O₂) and energy in the affected area, contributing to the potential development of pressure areas on the body.

HIGH-FREQUENCY ELECTROSURGERY IN LAPAROSCOPY

Laparoscopic surgery, or minimally invasive surgery, has revolutionized the field of surgical procedures, bringing significant benefits in terms of patient recovery times and healing. In this context, the use of high-frequency monopolar surgery (HF) is widely employed due to its flexibility in performing cuts, coagulations, and mixed cuts that combine both functions. However, this surgical approach carries risks for the patient, especially the risk of burns.

The risks of burns can be exacerbated by various factors, including limited visual field, inadequate maintenance of laparoscopic equipment, interference on the monitor, insufficient surgeon preparation or distractions, excessive smoke development, inadequate insulation, capacitive currents, and accidental contact of the active electrode tip with surrounding tissue. These factors can contribute to an increased risk of burns, internal injuries, tissue necrosis, and organ perforation.

Furthermore, within the surgical environment, where the active electrode closely interacts with conductive instruments and bodily tissue, there are factors that can promote the transfer of electrical currents to concealed areas. These factors include:

1. **Direct Coupling:** This occurs when the active electrode makes contact with another metal instrument, leading to the transmission of electrical current and an increased risk of burns to nearby tissue, such as the intestines or other organs.
2. **Insulation Failure:** In this case, the insulation surrounding the electrode can become compromised due to excessive voltage, improper use, or damage to the electrode shaft. This can happen during surgery or during the cleaning and sterilization of instruments. An invisible insulation breakdown, when the electrode is activated, poses a risk of unpredictable and potentially more insidious burns. Surprisingly, a minor insulation breakdown can be more hazardous than a major one, as it concentrates the current, making burns more likely.
3. **Capacitive Coupling:** This occurs when the active electrode induces electrical current in conductive materials, even if the insulation remains intact. During high frequency electrosurgery, the rapid changes in the electric field around the active electrode are only partially hindered by insulation, generating ionic currents that, upon contact with tissue, can cause sufficient heating to induce burns.

Addressing these risks with great care and implementing preventive measures is crucial to ensuring patient safety during high-frequency surgery in a laparoscopic setting.

To minimize the risk of burns during high frequency electrosurgery procedures in laparoscopy, consider the following preventive measures:

- **Comprehensive Staff Training:** Provide thorough and meticulous training for medical and healthcare personnel involved in electrosurgery procedures. It is essential to ensure they have a complete understanding of the procedures, associated risks, and preventive measures.
- **Detailed Inspection of Surgical Instruments:** Conduct a meticulous visual examination of surgical instruments, including the active electrode and laparoscope. This can help identify any defects or wear that could increase the risk of burns.
- **Use of Disposable Electrodes:** Although disposable electrodes may have thinner insulation that does not necessarily reduce the risk of insulation breakdown or capacitive coupling, they are free from wear and tear.
- **Avoidance of Hybrid Material Cannulas:** Steer clear of using cannulas made of hybrid materials, such as plastic and metal, as they can heighten the risk of direct coupling and capacitive coupling.
- **Adoption of Bipolar Technique:** While the bipolar technique may be less versatile than the monopolar one, it is considered safer because heat injuries are localized and occur only with prolonged current application.

In conclusion, burns are a genuine concern in high frequency electrosurgery procedures. However, with a deep understanding of potential causes and thorough preparation of the medical team, it is possible to limit their occurrence and effectively manage potentially risky situations.

PUTTING INTO SERVICE

- **Inspect the Equipment for Shipping Damage:** Check the equipment for any damage caused during shipping. If you find any damage, please report it immediately to the carrier and document the damage found. You should also inform LED SpA or your seller. When returning the equipment, use the original packaging or packaging that ensures safe transport.
- **Unpack and Review the Documentation:** Carefully remove the appliance from its packaging and study the provided documentation and operating instructions. Ensure that the mains voltage, as indicated on the nameplate data, matches the local mains voltage (mains frequency: 50-60Hz). If necessary, replace the fuses with the values specified on the data plate.
- **Connect to a Grounded Power Socket:** Connect the power cord to a mains socket with a reliable earth connection.


USING THE EQUIPMENT WITHOUT PROPER GROUNDING IS STRICTLY PROHIBITED.

- **Ensure Proper Installation:** The equipment must be placed on a flat surface with dimensions at least corresponding to the base of the equipment. Leave a minimum of 25 cm of space around the equipment.
- **Connect Mains Cable:** Attach the mains cable to the power socket located on the rear panel of the unit.
- **Establish Equipotential Connection:** Connect the equipotential connection point on the rear left side of the unit to the equipotential socket of the system.
- **Connect Pedal or Footswitch (Optional):** If you have a pedal or double footswitch (optional), connect it to the corresponding connector on the front panel of the unit.
- **Connect Handpiece:** If you are using a handpiece with buttons, connect it to the appropriate connection points. If you are using a handpiece without buttons, it must be connected to the "active" socket.
- **Operate in a Suitable Environment:** Ensure that the equipment is used only in a dry environment. Any condensation that forms must be allowed to evaporate before operating the equipment. Adhere to the ambient temperature and humidity levels, and do not exceed the specified limits.
- **Environmental conditions:**

	OPERATION	TRANSPORT / STORAGE
Temperature:	from 10 to 40 °C	from -10 to 50 °C
Relative Humidity:	from 30% to 75%	from 10 to 100%
Atmospheric Pressure:	from 70 to 106 kPa	from 50 to 106 kPa

- **Upon Activation:** When turning on the equipment using the switch located on the rear panel, the device will start with the function and power levels used during the last activation.
- **Monopolar Mode:** In monopolar mode, it is essential to connect the neutral electrode cable to the corresponding neutral electrode. If a split neutral electrode is used, ensure that the circuit is properly closed (connected correctly to the patient). This action will cause the red indicator light on the neutral electrode connector to stop flashing if the impedance value is acceptable.
- **Before Using the Equipment:** Before attempting to use the equipment, it is necessary to connect the patient plate cable. When using split electrodes, close the circuit by connecting the electrode to the patient. If the impedance value read by the equipment is acceptable, the OC indicator light will stop flashing, and you will receive an audible signal to indicate proper functionality.

CONNECTION AND USE OF ACCESSORIES

Depending on the unique requirements of each procedure, there may be a need for optional accessories. These accessories are designed to cater to specific demands, enhance precision, or optimize the outcomes of the procedure. Some examples of optional accessories may include:

- ***For the bipolar procedure:***

1 : Connection cable for bipolar accessories

2 : Bipolar Accessory (e.g., bipolar clamp)

- ***For vessel synthesis and coagulation (Vessel Sealing):***

3 : Connection cable for Artesy Sealer Clamp

4 : Artesy Sealer Clamp

CONNECTORS AND CONTROLS

LABEL ON REAR PANEL

Safety requirements for high-frequency (HF) surgical equipment necessitate the clear printing of data and symbols on the cabinet, or at the very least, on one of the panels of the power unit. These data and symbols serve a crucial role in defining the equipment's characteristics and monitoring its operational conditions to guarantee safe and effective utilization of the device.

TECHNICAL DATA


FREQUENCY	360kHz
CUT OUTPUT	200W - 300Ω
ENHANCED CUT OUTPUT	150W - 300Ω
BLEND OUTPUT	150W - 300Ω
FORCED OUTPUT	150W - 200Ω
SOFT OUTPUT	100W - 200Ω
FULGURATION OUTPUT	100W - 1000Ω
BIPOLAR PURE OUTPUT	120W - 50Ω
BIPOLAR TUR OUTPUT	200W - 50Ω
BIPOLAR COAG OUTPUT	120W - 50Ω
BIPOLAR SEALING OUTPUT	200W - 50Ω
POWER SUPPLY	100-240 Vac - 50/60 Hz selectable
ABSORPTION	750VA
FUSES	2xT 10AL, 250V
DUTY - CYCLE	intermittent 10 second emission / 30 second pause

MEANING OF GRAPHIC SYMBOLS

N°	SYMBOL	DESCRIPTION
1		Floating neutral electrode: not connected to ground at either high or low frequencies.
2		CF Class equipment protected against defibrillator-induced discharge.
3		Non-ionising radiation generator equipment.
4		Follow the instructions for use.
5		CE Mark (2017/745/EU) + Notified Body Number 0051 = IMQ Italy.
6		The product should not be disposed of in urban waste containers but must be disposed of through separate collection.
7		Manufacturer.
8		Serial Number.
9		Production date.
10		Unique Device Identification.
11		Medical Device.
12		Dealer.
13		No maintenance by the user.
14		Catalogue number (Code).
15		Temperature Limits.
16		Humidity Limits.
17		Atmospheric Pressure Limits.
18		This Way Up.
19		FRAGILE – Handle With Care.
20		Keep away from sunlight.
21		Protect against moisture.
22		Number of maximum stackable packages.
23		Weight.
24		Dimensions.
25		Number of Pieces.
26		Recycle.
27		Model/Trade Name.
28		Protection against harmful ingress of water or particulate matter.
29		Fuse.

BOX LABEL

With reference ISO15223-1 "Medical Devices-Symbols to be used with medical device, labels, labelling and information to be supplied" and ISO780 "Packaging – Pictorial marking for handling of goods" On box label of UNIT's carton box are present these indications:

- | | |
|--|---|
| 1. ISO15223-1 (5.1.1) MANUFACTURER | 12. ISO 7001: 2007 RECYCLING (Indicates the location of a bin or container) |
| 2. ISO15223-1 (5.7.10) UDI code = EAN code | 13. CE + Notified Body Number for MD Class (2017/745/UE) |
| 3. ISO15223-1 (5.1.6) CATALOGUE NUMBER | 14. ISO15223-1 (5.77) MD (Medical Device) |
| 4. ISO15223-1 (5.1.7) SERIAL NUMBER | 15. ISO15223-1 (5.1.3) DATE OF MANUFACTURER |
| 5. WEIGHT OF BOX | 16. ISO15223-1 (5.3.7) TEMPERATURE LIMITS (Indicates temperature limits within which the transport package shall be stored and handled) |
| 6. DIMENSIONS OF BOX | 17. ISO15223-1 (5.3.8) HUMIDITY LIMITS (Indicates humidity limits within which the transport package shall be stored and handled) |
| 7. ISO15223-1 (5.3.4) KEEP DRY (Transport package shall not be kept from moisture) | 18. ISO15223-1 (5.3.9) ATMOSPHERIC PRESSURE LIMITS (Indicates atmospheric pressure limits within which the transport package shall be stored and handled) |
| 8. ISO15223-1 (5.3.2) KEEP AWAY FROM SUNLIGHT (Transport package shall not be exposed to sunlight) | 19. WEEE PRODUCT (Directive 2012/19/EU) |
| 9. ISO15223-1 (5.3.1) FRAGILE (Contents of the package are fragile therefore it shall be handled with care) | |
| 10. STACKING LIMIT BY NUMBER (Indicates the maximum number of identical products that can be safely stacked on the bottom package) | |
| 11. ISO780 (3) THIS WAY UP (Indicates correct upright position of the transport package) | |

FRONT PANEL

1. USB PORT

To perform software updates, use the USB 2.0 port located on the front of the unit.

2. TOUCHSCREEN DISPLAY

The backlit LCD touchscreen display allows for the visualization and selection of all set and adjustable parameters within a specific procedure.

3. MONOPOLAR OUTPUT CONNECTOR

This is the connection point for the handpiece with dual buttons to activate cutting (CUT) and coagulation (COAG) functions. If handpieces without buttons or optional monopolar cables are used, they should be connected to the socket labeled "ACTIVE."

4. NEUTRAL ELECTRODE CONNECTOR

This is the connection point for the neutral electrode to be applied to the patient. Disposable or reusable, single-part or split neutral electrodes can be used.

5. BIPOLAR OUTPUT CONNECTOR

This is the connection point for bipolar accessories.

6. PEDAL CONNECTOR


On the front panel, there is a socket for connecting the pedal or optional dual pedal unit (pedalboard).

OPERATING MODES

CONTROL AND POWER ON

The unit is operated directly through the icons displayed on the touchscreen of the device. To confirm a selection, simply touch the icon.

Once the electrosurgical unit is powered on, the specific software begins to load on the screen. The progress of the software application process is indicated by the filling of the bar at the bottom of the screen. At the end of this process, the home screen will be displayed:

On the home screen, you can select and interact with the following options:


- Programs
- Surgery
- Settings
- Update (Software Update via USB Port).

Using the "*info*" option, you can view the currently installed software versions. By pressing the "*help*" button, you can access an informative summary that is helpful for correctly interpreting the indications displayed on the screen.

PROGRAMS


By selecting the *"Programs"* option on the home screen, two separate pages will open:

1. **Preset:** In this section, pre-configured programs that have already been set up and stored in the system will be displayed. These represent standard settings for a variety of common surgical procedures.

Simply click *"Start Program"* to activate the pre-set program and make it ready for use.

2. **Custom:** On this page, you will see programs that have been created and customized by the user or surgeon. These programs reflect specific settings that the user has configured to address particular needs or unique procedures.


At first, the screen will appear empty. For instructions on how to create a personalized program, please refer to the next section (*"Programs Section"*).

The distinction between these two pages allows operators to easily access a range of predefined or customized settings, depending on the specific requirements of the surgical procedure. This provides flexibility in selecting settings according to the context.

Furthermore, by pressing '*Help*' on the display, you can obtain information about the features of the selected program. To return to the home page, simply select the '*Back*' option.

SURGERY

When selecting the '*Surgery*' option on the home page, the following screen will be displayed:


Through this screen, you can access different sections:

- A. Monopolar Cutting Section
- B. Monopolar Coagulation Section
- C. Neutral Electrode Circuit Impedance Indication
- D. Timing Section
- E. Programs Section
- F. Home Icon
- G. Bipolar Cutting Section
- H. Bipolar Coagulation Section

Each of these options allows you to adjust various parameters, which can be saved as a new customized program.


MONOPOLAR SECTION

In the sections dedicated to adjusting cutting or coagulation, both in monopolar and bipolar modes, there are three icons:

1. *Icon for adjusting the output power*

By tapping the icon with the number, you can adjust the output power. Below is an example to illustrate the procedure:

Pressing the '+' and '-' buttons adjusts the power, while pressing the '✓' button confirms the selected value.

2. *Icon to select the deliverable current*

By tapping the icon with the function, you can select the desired one.

CUT CURRENT (CUT)

The ideal current for cutting is pure sinusoidal without modulation, i.e., with a 100% duty cycle.

BLENDED CURRENT (BLEND)

The blended current (BLEND) is suitable for coagulated cutting when you want to achieve deep coagulation combined with cutting. Its waveform has a lower modulation percentage compared to pure coagulation. This current consists of a sinusoidal current component suitable for cutting, combined with a low-voltage coagulation current (SOFT COAG). This creates a current that is suitable for coagulated cutting without the formation of eschar or charring and is particularly suitable for endoscopic procedures.

ENHANCED CUT CURRENT (ENHANCED CUT)

The enhanced cut current (ENHANCED CUT) is a sinusoidal current characterized by amplitude modulation and is suitable for cutting tissues, particularly adipose tissues.

FORCED COAG CURRENT (FORCED COAG)

The modulated current (FORCED COAG) has good surface coagulative properties but may also lead to the potential formation of eschar and slight tissue carbonization. The advantage of this type of coagulation lies in the speed at which the desired effect is achieved.

FORCED COAG is also referred to as Speedy.

SOFT COAG CURRENT (SOFT COAG)

The low-voltage, low modulation current (SOFT COAG) is suitable for coagulating deep tissue layers, where cellular albumin coagulation is achieved without carbonization or eschar formation. The coagulation process in this case is slower than in FORCED COAG coagulation.

SOFT Coag is also referred to as Pin Point, Dessicate, or Deep.

FULGURATION COAGULATION CURRENT (FULGURATE)

The high-voltage fulguration coagulation current (FULGURATE) is directed into the active electrode that is not in contact with the tissue being treated and is primarily responsible for coagulation. This method is ideal for treating large surfaces with widespread and superficial bleeding (e.g., liver resection) and/or for achieving coagulation at the open sternum level in cardiac surgery.

3. Icon for activating the footswitch

By tapping the icon with the footswitch, you can choose which one to activate, for example:

Footswitch Activated

Footswitch Not Activated

BIPOLAR SECTION

In the bipolar section, just like in the monopolar section, you can adjust the output power, select the deliverable current, and activate the footswitch. The difference lies in the mode-specific selectable functions for each mode. You will need to connect the bipolar accessories to the connectors for this function and use the pedal switch.

BIPOLAR CUT CURRENT (BIPOLAR CUT)

The current supplied by bipolar forceps is a high-voltage pure sinusoidal current suitable for cutting without coagulation.

TRANS URETHRAL BIPOLAR RESECTION CURRENT (BIPOLAR TUR)

The TUR current, using a specific bipolar accessory, is suitable for both cutting and coagulated cutting when forced coagulation along with cutting is desired. This current consists of a sinusoidal current component suitable for cutting, along with a high-voltage current component for coagulation.

BIPOLAR COAG CURRENT (BIPOLAR COAG)

This type of coagulation is achievable with bipolar forceps and allows for the delivery of radiofrequency output power at low impedances, such as those typically encountered in the tissue section that can be positioned between the tips of the forceps.

VESSEL SEALING

This type of function is suitable for arterial and venous vascular synthesis and coagulation using radiofrequency clamping. The vessel is gently clamped with slight pressure, radiofrequency is initiated, and pressure with the forceps is maintained until coagulation occurs.

For the operation of the bipolar modes, please refer to the "AUTOSTART AND AUTOSTOP" section.

AUTOSTART AND AUTOSTOP


In the Bipolar functions, you can access four different operating settings:

- **No preset output automation**

The output is only activated by pressing the pedal and stops when released.

- **START**


To activate the "*start bipolar*" or "*start sealing*" function, you need to press the pedal, establish contact between the active electrode and the tissue, and this will activate the output. If you wish to stop, simply release the pedal.

- **STOP**

To activate the "*stop bipolar*" or "*stop sealing*" function, you need to press the pedal; the output is activated (even if there is no contact between the tissue and the active electrode) and stops when the tissue is coagulated or if the pedal is released.


- **AUTOSTART/AUTOSTOP**

To activate the "*start/stop bipolar*" or "*start/stop sealing*" function, you need to press the pedal, establish contact between the active electrode and the tissue, and this will activate the output. The output stops when the tissue is coagulated or if the pedal is released.


DISPENSING SCREEN

In the dispensing state, the screen will display the function along with its respective level. This screen remains for a few seconds during which you can adjust the output level using the + and - buttons.

FLASH TIMER SECTION

Through the "FLASH TIMER" section, you can select the "*continuous*" option for continuous dispensing or "*timed*" to set the parameters for pulses.

For programming pulse timing, consider:

- The "*time*" section, where you set the duration of the pulses (from 10 milliseconds to 30 seconds).
- The "*repeat*" section, to choose the number of pulse repetitions (from 1 to infinite).

It's crucial to keep the pedal pressed throughout the time you want to use the timed function.

TIME		Supply Time		Step
		<i>From</i>	<i>To</i>	
10÷90	msec	10 msec	90 msec	10 msec
100÷900	msec	0.1 sec	0.9 sec	100 msec
1.0÷30	sec	1 sec	9.5 sec	0.5 sec

The "*timed*" function can be adjusted with all functions, but it does not interfere with the AUTOSTART/AUTOSTOP function in BIPOLAR mode.

PROGRAMS SECTION


In the *"Surgery"* screen, you can add a custom program or view the existing programs.

To create a new program, follow the following procedure:

- Set the desired values in the 'Surgery' screen.
- Press '+ save' and choose the name to assign, completing the process by clicking 'enter'.

The new customized programs can be viewed in the 'Programs' section.

Through the icon representing three lines, you can perform various actions on the customized program:

Rename it.

Modify it.

Delete it.

NEUTRAL ELECTRODE CONTROL

The neutral electrode circuit is continuously monitored by a special circuit. This circuit checks, especially when bipartite neutral electrodes are in use, for any loss of contact between the patient reference plate or changes in the conductivity characteristics of the neutral electrode. Such changes can lead to a reduction in circuit conductivity, increasing the risk of burns to the patient.

To reduce acoustic pollution, the audible alarm only triggers when the delivery pedal is held down.

It's important to note that when using single-section (single-part) neutral electrodes, the circuit only checks the connection of the neutral electrode to the unit. Therefore, it's crucial to ensure that the entire surface of the neutral electrode is correctly and securely applied to the patient.


LIGHTS

1. **Monopolar Output Connector**
2. **Monopolar Output Indicator**, which has three modes:
 - Off Indicator – Monopolar functions off (tap on the display "MONOPOLAR" to turn on and off the section)
 - Yellow Indicator – Cut current output
 - Blue Indicator – Coagulation current output
3. **Bipolar Output Connector**
4. **Bipolar Output Indicator**, which has three modes:
 - Off Indicator – Bipolar function not active (tap on the display "BIPOLAR" to activate or deactivate the function)
 - Yellow Indicator – Cut current output
 - Blue Indicator – Coagulation current output
5. **Connector For Connecting the Neutral Electrode**
6. **Neutral electrode connection indicator**, which has two modes:
 - White light - Correct positioning of the neutral electrode
 - Red light - Incorrect or missing neutral electrode positioning
7. **Connector For Connecting the Footswitch**
8. **Foot Pedal Indicator**, if off indicates the footswitch is not connected

SETTINGS

Starting from the home page and selecting the "*Settings*" option will bring up the following screen:

Through this screen, you will be able to access three sections:

- Language.
- Display, where you can adjust brightness, touchscreen sensitivity, volume, and switch between colour and black-and-white display.
- Surgery, which allows you to select the surgical intervention delay response and fast-medium-slow autostart/autostop.

With the start and autostart functions (see "AUTOSTART AND AUTOSTOP" section), you can select the delay time between contact between the active electrode and the tissue and the activation of the output (from 0.1 s to 2.0 s).

The sensitivity settings for autostart and autostop functions can be set to SLOW, MEDIUM, and FAST.

Pressing the "home" icon will take you back.

UPDATE

Starting from the home page and choosing the "Update" option, you can update:

- Software
- Images
- Protocols
- Firmware


To perform updates, you need to connect a compatible device to the USB connector, which contains the compatible file for the software, images, protocols, or firmware to be updated. Follow these steps:

1. Insert the compatible device into the equipment's USB connector.
2. On the screen, select the corresponding option among "Software," "Images," "Protocols," or "Firmware."
3. Confirm the selection through the popup that appears.
4. Follow the on-screen instructions to complete the update process.
5. Once the update is complete, you can select the "*Home*" or "*Back*" option to exit the procedure and return to the main screen.

It is crucial to carefully follow the on-screen instructions throughout the entire process to ensure a correct and safe update.

BACK PANEL

1. Power Socket
2. Power Switch
3. Fuse Holder / Voltage Selector
4. Equipotential Socket

POWER SUPPLY MODULE OF THE EQUIPMENT AND VOLTAGE SELECTOR

The equipment's power supply module serves as the connection point for supplying power to the internal electronics of the equipment. This power supply module includes the power connector and line fuses. The voltage selector is located inside the power supply module.

WARNING: Before powering on the equipment, the operator should ensure that the mains voltage indicated on the voltage selector matches the voltage it is connected to and that appropriate fuses for the selected voltage have been inserted.

POWER SWITCH

To turn on the power to the equipment, move the switch to position 1. When the power is on, the front panel will light up. Moving the switch to position 0 will turn off the power supply, serving as an emergency switch in case of any faults.

TECHNICAL CHARACTERISTICS

Tolerance	Description	DIATERMO MB 200 T
-	Electrosurgical unit code	GMA10100.T40
-	Detection system of tissue impedance (Bipolar coagulation – auto start/autostop)	●
-	Bipolar coag with automatic activation/disactivation	●
-	Minimum preselectable power	1
-	Selection of the power through touch screen	●
± 20%	Maximum output power CUT (W)	200W → 300Ω
± 20%	Maximum output power BLEND (W)	150W → 200Ω
± 20%	Maximum output power ENHANCED (W)	150W → 300Ω
± 20%	Maximum output power FORCED COAG (W)	150W → 200Ω
± 20%	Maximum output power SOFT COAG (W)	100W → 200Ω
± 20%	Maximum output power FULGURATION (W)	100W → 1000Ω
± 20%	Maximum output bipolar power BIPOLAR CUT (W)	120W → 50Ω
± 20%	Maximum output bipolar power BIPOLAR TUR (W)	200W → 50Ω
± 20%	Maximum output bipolar power BIPOLAR COAG (W)	120W → 50Ω
± 20%	Maximum output bipolar power BIPOLAR VESSEL SEALING (W)	200W → 50Ω
± 5%	Modulation factor BLEND (Hz)	50
± 5%	Modulation factor ENHANCED (Hz)	1.25
± 5%	Modulation factor FORCED COAG (kHz)	20
± 5%	Modulation factor FULGURATION (kHz)	20
± 5%	Modulation factor BIPOLAR TUR (Hz)	50
± 0.2	Crest Factor CUT	1.6
± 0.2	Crest Factor BLEND	2.0
± 0.2	Crest Factor ENHANCED CUT	2.5
± 0.2	Crest Factor FORCED COAG	3.0
± 0.2	Crest Factor SOFT COAG	1.6
± 0.2	Crest Factor FULGURATION	4.0
± 0.2	Crest Factor BIPOLAR CUT	1.6
± 0.2	Crest Factor BIPOLAR TUR	2.0
± 0.2	Crest Factor BIPOLAR COAG	1.6
± 0.2	Crest Factor BIPOLAR VESSEL SEALING	1.6
± 10%	Working frequency	360 kHz
± 15%	Maximum output voltage CUT (Vpp)	2000
± 15%	Maximum output voltage BLEND (Vpp)	1800
± 15%	Maximum output voltage ENHANCED CUT (Vpp)	1800
± 15%	Maximum output voltage FORCED COAG (Vpp)	2500
± 15%	Maximum output voltage SOFT COAG (Vpp)	1000
± 15%	Maximum output voltage FULGURATION (Vpp)	4000
± 15%	Maximum output voltage BIPOLAR CUT (Vpp)	650
± 15%	Maximum output voltage BIPOLAR TUR (Vpp)	1200
± 15%	Maximum output voltage BIPOLAR COAG (Vpp)	650
± 15%	Maximum output voltage BIPOLAR VESSEL SEALING (Vpp)	650
± 0.5	Size (LxHxP) mm	370x144x319
± 10	Weight (kg)	6
± 5%	Selectable mains power (Vac)	100 – 240
± 1%	Mains frequency (Hz)	50-60
± 0	Fuses (5x20) TIMED	2xT 10AL 250V
± 10%	Electrical input power (VA)	750
± 10%	Electrical input current (240Vac) (A)	3,15
± 10%	Electrical input current (100Vac) (A)	7,5
-	Five steps adjustable sound level	●
-	Self-check	●

Tolerance	Description	DIATERMO MB 200 T
-	Power accuracy output warning	●
-	Skin Plate Electronic Control ¹	●
-	Split or not split patient plate allowed	●
-	Working (CUSTOM) condition storing ²	More than 50
-	Electrical Class (EN60601-1)	Class I Applicated Part CF
-	MDR 2017/45/UE Class	II b
-	IP Protection Class (EN 60529)	IP32
-	EN55011 (CISPR 11) Class (Class/Group)	2 / A
-	Patient circuit	F
-	Duty Cycle (action / pause) in seconds	10 / 30
-	Output power control by footswitch or finger-switch	●
-	Defibrillation-proof	●
-	Equipotential binding	●
-	ABS cabinet	●

● = PRESENT

- = NOT PRESENT

HARDWARE REQUIREMENTS

Microcontroller	ARM Cortex M4
Clock Frequency	200 MHz
Flash	2048 KB
Ram	512 KB
SDRAM DDR2	512 MB
Nand Flash	1 Gb
Peripherals	UART, I2C, SPI, Watch-dog timer, USB2.0
Visual	Display touchscreen 7" 800x480 px

¹ Patient to plate contact monitoring system

² Continuous storing of the last settings

MAINTENANCE

GENERAL

There are no user-adjustable parts inside the equipment for calibration or service purposes. Opening the equipment casing voids the warranty. If repair or adjustment is necessary, the entire equipment should be sent to the LED SpA service center in APRILIA (LT), ITALY, along with a description of the issue. User maintenance primarily involves cleaning and sterilizing accessories and performing equipment checks before each use. Conducting functional and safety checks to verify parameters is the responsibility of specialized technical personnel.

CLEANING THE ENCLOSURE

Whenever possible, use only disposable accessories and dispose of them as special hospital waste. However, because some accessories may need to be reused, it is essential to thoroughly clean and sterilize them before each new use. Refer to the instructions for use (IFU) provided by the supplier of each reusable accessory for information on the maximum number of cycles and the type of sterilization required.

CLEANING AND STERILIZATION OF ACCESSORIES

If disposable non-sterile accessories are used, you must meticulously follow the Instructions for Use (IFU) provided by the manufacturer for the sterilization method and to dispose of them according to the currently applicable regulations.

In the case of reusable accessories, you must adhere to the maximum number of cycles and the sterilization method specified in the manufacturer's Instructions for Use for each accessory.

GUIDE TO TROUBLESHOOTING

If you encounter any problems, it's crucial to verify that the installation and accessory setup have been executed correctly. The error code shown on the seven-segment display, along with the description on the LCD display, can offer valuable insights for diagnosing and resolving the issue.

Problems	Probable Cause	Solution
The equipment does not power on.	There is an interruption or absence of the main power supply.	Please verify the connection of the main cable and replace any fuses with new ones of the appropriate type where necessary.
The Neutral Electrode alarm is consistently active.	There is an interruption or lack of contact in the neutral electrode circuit.	Check the connection of the cable to the neutral electrode. Replace the cable of connection of the neutral electrode.
The unit fails to respond to the activation command.	Possible causes for the issue could include a breakdown of the handpiece or pedal, or incorrect connections of these components.	Replace the handpiece or the pedal. Verify the connection of the handpiece or of the pedal.
Error Code 001.	Current delivery control is activated during the power-on process.	Disconnect the handpiece or the pedal and switch on the unit again.
Error Code 004.	Error in the data conversion circuit.	Call for Service.
Error Code 010.	Error in the output power activation circuit.	Call for Service.
Error Code 022.	Communication Error.	Call for Service.

REPAIRS

High-frequency cables or electrode handle cannot be repaired. Always replace a faulty part with a new one.

FUSE REPLACEMENT


Before replacing the fuses, disconnect the equipment from the power supply.

For fuse replacement, use 5x20 type fuses and follow these steps:

(A-B) Using a small screwdriver, remove the fuse holder cartridges from the power module.

(C) Insert the fuses, referring to this information:

Voltage: 100-240 V Fuse Type: Time-Delay T10 AL, 250V / 5 x 20 mm

PRE-USE UNIT CHECK

Whenever the unit is scheduled for use, it is essential to perform a safety check, considering the following points:

- Inspect the integrity of cables, connections, wire breakage, etc.
- Ensure that all electrical equipment is properly grounded.
- Verify that all accessories to be used are available and sterilized.
- Check the operation of the OC light by disconnecting the reference electrode cable. Activate the unit and verify the OC light and audible alarm.
- Test the warning lights and audible signals by activating the CUT and COAG power switches.


SAFETY AND MEASUREMENT FUNCTION CHECKS

Periodically (at least once a year), safety checks and measurements should be scheduled by the Department of Biomedical Engineering or other specialists. These include:


- Inspection of power cables and connectors.
- Visual inspection of mechanical protections.
- Evaluation of protections against hazards from liquid leakage, penetration, drips, humidity, hygiene products, and disinfectants.
- Review of data on the device panel.
- Verification of the availability of the instruction manual.
- Examination of high-frequency output.
- Measurement of grounding conductivity resistance.
- Measurement of high-frequency leakage current.
- Assessment of neuromuscular stimulation.
- Verification of output power correction.

Information about elimination of this product (Applicable in the countries with separate collection systems)	
	On the end of the life, the present product <u>mustn't</u> be eliminated as urban refusal, but it must be eliminated in a separated collection. If the product is eliminated in unsuitable way, it is possible that some parts of the product (for example some accumulators) could be negative for the environment and for the human health. The symbol on the side (barred dustbin on wheel) denotes that the products mustn't throw into urban refuses container, but it must be eliminated with separate collection. In case of abusive elimination of this product, could be foreseen sanctions.


DIAGRAMS


Diagrams of half and maximum output power versus impedance load 100-2000Ω CUT100%


Diagrams of half and maximum output power versus impedance load 100-2000Ω BLEND


Diagrams of half and maximum output power versus impedance load 100-2000Ω ENHANCED


Diagrams of half and maximum output power versus impedance load 100-2000Ω FORCED COAG


Diagrams of half and maximum output power versus impedance load 100-2000Ω SOFT COAG


Diagrams of half and maximum output power versus impedance load 100-2000Ω FULGURATION


Diagrams of half and maximum output power versus impedance load 10-1000Ω BIPOLAR CUT


Diagrams of half and maximum output power versus impedance load 10-1000Ω BIPOLAR TUR


Diagrams of half and maximum output power versus impedance
load 10-1000Ω BIPOLAR COAG


Diagrams of half and maximum output power versus impedance
load 10-1000Ω BIPOLAR VESSEL


Diagrams of output power versus nominal value
CUT100%


Diagrams of output power versus nominal value
BLEND


Diagrams of output power versus nominal value
ENHANCED


Diagrams of output power versus nominal value
FORCED COAG


Diagrams of output power versus nominal value
SOFT COAG


Diagrams of output power versus nominal value
FULGURATION


*Diagrams of output power versus nominal value
BIPOLAR CUT*


*Diagrams of output power versus nominal value
BIPOLAR TUR*


*Diagrams of output power versus nominal value
BIPOLAR COAG*


*Diagrams of output power versus nominal value
BIPOLAR VESSEL*


*Diagrams of maximum mains voltage output versus Vp
CUT*


*Diagrams of maximum mains voltage output versus Vp
BLEND*


*Diagrams of maximum mains voltage output versus Vp
ENHANCED*


*Diagrams of maximum mains voltage output versus Vp
FORCED COAG*


*Diagrams of maximum mains voltage output versus Vp
SOFT COAG*


*Diagrams of maximum mains voltage output versus Vp
FULGURATION*


*Diagrams of maximum mains voltage output versus Vp
BIPOLAR CUT*

*Diagrams of maximum mains voltage output versus Vp
BIPOLAR TUR*

*Diagrams of maximum mains voltage output versus Vp
BIPOLARCOAG*

*Diagrams of maximum mains voltage output versus Vp
BIPOLAR VESSEL*

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