Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

TRUSTEEL™

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ENGLISH

DESCRIPTION

Trusteel (Stainless Steel), is a sterile synthetic monofilament non absorbable surgical suture prepared from 316 LVM grade Stainless Steel. 316 LVM grade Stainless Steel alloy contains iron, copper, manganese. Molybdenum, Nickel and Chromium. Trusteel sutures comply with the requirements established by the European Pharmacopoeia (E.P) for Sterile Non-Absorbable sutures, and United States Pharmacopoeia (U.S.P) for Non-Absorbable Surgical Suture.

INTENDED USE

Trusteel suture is intended for use in abdominal wound closure, hernia repair, sternal closure and orthopaedic procedures including cerclage and tendon repair.

INDICATION

Indicated for surgical procedures which require sutures with high tensile strength, such as in conditions like abdominal wound closure, hernia repair, sternal closure and orthopaedic procedures including cerclage and tendon repair.

CONTRAINDICATIONS

The use of the suture is contra indicated in patients with known sensitivities or allergies to 316 LVM Stainless Steel, or its constituent metals such as chromium and nickel.

INTENDED USERS

Trained and registered health care professional only.

INTENDED PATIENT POPULATION

This suture can be used in patients irrespective of age and sex, in line with the intended use, indications and contraindications.

PERFORMANCE

Trusteel suture elicits a minimal initial inflammatory reaction in tissues and is not absorbed

APPLICATION

Suture should be selected and implanted depending on patient's condition, surgical experience, surgical technique and wound size.

WARNINGS

This suture is intended to be used by a trained and registered

health care professional.

Users should be familiar with surgical procedures and surgical suturing techniques involving non-absorbable stainless steel sutures before using Trusteel suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.

Users must apply their professional judgment when determining the appropriate sized suture based on the specific indication, wound size, preferred surgical technique, patient's condition and history.

As a non absorbable suture, it may act as a foreign body and as with any foreign body, presence of bacterial contamination may enhance bacterial infectivity. Thus, acceptable surgical practice should be followed with respect to drainage and closure of infected or contaminated wounds.

Wound dehiscence may occur when suture fails to provide adequate wound support inclosure of the sites where expansion, stretching or distension occur. Postoperatively and until healing is complete, fixation provided by this suture should be considered as temporary and may not withstand weight bearing or other unsupported stress.

The fixation provided by this suture should be protected. The postoperative regimen prescribed by the surgeon should be strictly followed to avoid adverse stresses applied to the suture.

Under certain circumstances, notably orthopaedic procedures, immobilization by external support may be employed at the discretion of the user.

Device removal when required, should be followed by adequate postoperative management.

PRECAUTIONS

Care should be taken to avoid damage when handling surgical needles. Grasp the needle in an area one third (1/3) to one half (1/2) of the distance from the attachment end to the point.. Grasping in the point area could impair the penetration performance & cause fracture of the needle. Grasping at the butt or attachment end could cause bending or breakage of the needle. Reshaping needles may cause them to lose strength and make less resistant to bending and breaking. Broken needles may result in extended or additional surgeries or residual foreign bodies.

Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury. Discard used needles in 'Sharps' containers.

When handling this suture, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. As with any suture material, adequate knot security requires the accepted surgical techniques of flat and square ties, with additional throws as warranted by surgical circumstances and the experience of the user.

Do not tie bands around pouches/packs of sutures as it may damage the pouches/pack.

Store the suture pouches/pack as supplied and avoid bending of the suture pouches/pack during transport or storage, as bending of suture before usage may lead to suture breakage.

Do not re-sterilize, do not reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or cross-contamination, which may lead to infection or transmission of blood-borne pathogens to patients and users.

SIDE EFFECTS

Side effects associated with the use of Trusteel suture include allergic response in patients with known sensitivities to 316 LVM Stainless Steel, or its constituent metals such as chromium and nickel; transient local irritation at the wound site, minimal initial inflammatory tissue reaction, pain, edema, and erythema at the wound site, and wound dehiscence. Like all foreign bodies, Trusteel suture may enhance an existing infection.

STERILITY

Trusteel suture is sterilized by ethylene oxide gas. Do not re-sterilize, do not reuse, Do not use if package is opened or damaged.

Discard opened, unused sutures.

SUPPLY

Trusteel suture is available in various U.S.P. sizes and E.P Metric sizes.

Suture, are supplied sterile in pre-cut lengths, both non-needled and attached to various needle type, shape and length. Sutures are packed in a printed box with quantity indicated on the box label

STORAGE

Recommended storage conditions: Below 30° C, away from moisture and direct sunlight. Do not use after expiry date.

DIPOSAL

Discard used sutures and needles contaminated with blood and tissue in the container meant for infectious waste. Unused expired pouches should be incinerated or disposal should be done as per local regulations.

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Gima S.p.A. Via Marconi, 1 - 20060 Gessate (MI) Italy gima@gimaitaly.com - export@gimaitaly.com www.gimaitaly.com