CANNULA TRACHEOSTOMICA NON CUFFIATA / CUFFIATA
UNCUFFED / CUFFED TRACHEOSTOMY TUBE
CANULE DE TRACHÉOTOMIE SANS BALLONNET / À BALLONNET
TRACHEOSTOMIEKANÜLE OHNE CUFF / MIT CUFF
CÁNULA DE TRAQUEOSTOMÍA SIN BALÓN / CON BALÓN
TUBO DE TRAQUEOSTOMIA SEM BRAÇADEIRA / COM BRAÇADEIRA
ΤΡΑΧΕΙΟΣΩΛΗΝΑΣ ΧΩΡΙΣ ΑΕΡΟΘΑΛΑΜΟ / ΜΕ ΑΕΡΟΘΑΛΑΜΟ
ΤΡΑΧΕΟΣΤΟΜΙΧΗ ΤΡΥΒΑ ΒΕΖ ΜΑΝΨΕΤ / Σ ΜΑΝΨΕΤ
TRACHEOSTOMICKÁ KANYLA BEZ MANŽETY / S MANŽETOU
TRAKEOSTOMIRØR UDEN MANCHET / MED MANCHET
MANSETITA / MANSETIGA TRAHHEOSTOOMIATORU
TRAKEOSTOMIAPUTKI ILMAN KUFFIA / JOSSA KUFFI
TRAHEOSTOMSKA CIJEV BEZ MANŠETE / S MANŠETOM
BALLONOS / BALLON NÉLKÜLI TRACHEOSZTÓMIÁS KANÜL
TRACHEOSTOMINIS VAMZDELIS BE MANŽETĚS / SU MANŽETE
TRAHEOSTOMIJAS CAURULE BEZ MANŠETES / AR MANŠETI
TRAKEOSTOMISLANGE UTEN MANSJETT / MED MANSJETT
MET CUFF / CUFFLOZE TRACHEACANULE
RURKA TRACHEOSTOMIJNA BEZ MANKIETU / Z MANKIETEM
TUB DE TRAHEOSTOMIE FĂRĂ MANȘETĂ / CU MANȘETĂ
TRACHEOSTOMICKÁ TRUBICA BEZ MANŽETY / S MANŽETOU
TRAHEOSTOMSKA CEVKA BREZ MANŠETE / Z MANŠETO
TRAKEALKANYL UTAN KUFF / MED KUFF

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Tracheostomy Tubes are available uncuffed/plain or cuffed. Tubes are supplied with compatible cuff, a one-way luer-lock valve, a responsive pilot balloon and radiopaque line.

INDICATIONS

Tracheostomy Tubes are indicated for airway management of tracheostomized patients.

CONTRAINDICATIONS

Care must be taken to avoid contact of a LASER beam or an electrosurgical active electrode with this and other tracheostomy tubes. Such contact, especially in the presence of oxygen-enriched mixtures, could result in rapid combustion of the tube with harmful thermal effects and with emission of corrosive and toxic combustion products, including hydrochloric acid (HCL).

WARNING

- Inflation of cuff by "feel" alone or by using a measured amount of air is not recommended since resistance is an unreliable guide during inflation.
- Deflate cuff prior to repositioning the tube. Movement of the tube with cuff inflated could result in patient injury or damage to the cuff, require a tube change.
- Diffusion of nitrous oxide mixture, oxygen or air may either increase or decrease cuff volume and pressure. Inflating the cuff with the gas mixture which will contact its external surface is recommended as a means to decrease the incidence of such diffusion.
- Do not overinflate cuff. Overinflation can result in tracheal damage, rupture of the cuff with subsequent deflation, or in cuff distortion which may lead to airway blockage.
- Do not use if package has been opened or damaged.
- Do not use needle to inflate cuff. Use Luer syringe.
- For single patient use only.
- Do not use after the expiry date.
- This is a disposable product for single use, any resterilization or reuse of the product can result in hazard to the patient.
- Precautions shall be taken when discarding device and disposal of the device shall be made in accordance with applicable national regulations for biological hazardous waste.

PRECAUTIONS

1. In the event abnormal anatomy is encountered or unusual positioning of the head or neck is required, care should be exercised to avoid the potential for kinking of the tracheostomy tube. The use of a reinforced tube should be considered for unusual cases.
2. Seat the connector firmly in both the tracheostomy tube and adapter on the ventilation equipment to prevent disconnection during use.
3. Non-standard dimensioning of some connector on ventilator or anesthesia equipment may make secure mating with the tracheostomy tube 15 mm connector difficult.
4. Expert clinical judgment should be exercised in the selection of the appropriate size tracheostomy tube for each individual patient.
5. Intubation and extubation should be performed following currently accepted medical techniques.
6. Each tube's cuff, pilot balloon and valve should be tested by inflation prior to use.
7. Care should be taken to avoid damaging the thin-walled cuffed during intubation.
8. Avoid using Lidocaine Topical Aerosol as it has been associated with the formation of pinholes in PVC cuffs (Jayasuriya, K.D., and Watson, W.F.: PVC Cuffs and Lignocaine-based Aerosol. Brit. J. Ann. 53:1368, 1981). The same authors report that lidocaine hydrochloride solution does not have this effect.
9. Tree-way stopcocks or other devices should not be left inserted in the inflation valve for extended periods of time. The resulting stress may break the valve housing causing the cuff to deflate.
10. Alcohol based solutions should not be in contact with the tube.

ADVERSE REACTIONS

Standard textbooks and scientific literatures should be consulted for any specific adverse reaction information.

DIRECTIONS FOR USE: Use aseptic technique

1. Remove the sterile tracheostomy tube from the package.
2. Test/inflate cuff using a luer-tip syringe.
3. Completely deflate cuff after test inflation.
4. With the obturator in place (if supplied), insert the tracheostomy tube through the stoma and position in the tracheal in accordance with currently accepted medical technique.
5. Following insertion, remove the obturator (if used) and slowly inflate the cuff with the minimum amount of air required to effect a seal.
6. Remove the syringe from the one-way valve.
7. Attach the swivel adapter to the 15 mm termination of the tracheostomy tube and connect the ventilator circuit to the swivel adapter. The anti-disconnect lugs on the swivel should be utilized to help guard against accidental disconnect.
8. Secure the tracheostomy tube.

STORAGE CONDITIONS

Protect product from moisture and excessive heat.

Avoid prolonged exposure to ultraviolet and fluorescent lighting.

Store in a manner to prevent crushing. Stock rotation by First-In-First-Out basis.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

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