36 rue Brunel
75017 PARIS - France
Phone: +33 (0)1 44 95 14 95

FRA	Dispositifs médicaux pour hystérosalpingographie et chromopertubation	2
ENG	Medical devices for hysterosonography and chromopertubation	3
SPA	Productos sanitarios y accesorios para histerosalpingografía y cromopertubación	4
BUL	Медицински изделия и принадлежности за хистеросалпингография и хромопертубация	5
CES	Zdravotnické prostředky a příslušenství pro hysterosalpingografii a chromopertubaci	6
DAN	Medicinsk udstyr og tilbehør til Hysterosalpingografi og kromopertubation	7
ELL	Ιατροτεχνολογικά προϊόντα και παρελκόμενο για υστεροσαλπιγγογραφία και Χρωματογραφία	8
EST	Meditsiiniseadmed ja tarvikud hüsterosalpingograafia ja kromopertubatsiooni läbiviimiseks	9
FIN	Lääketieteelliset laitteet ja tarvikkeet hysterosalpingografiaa ja kromopertubaatiota varten	10
DEU	Medizinprodukte und Zubehör für Hysterosalpingographie und Chromopertubation	11
HRV	Medicinski proizvodi i pribor za histerosalpingografiju i kromopertubaciju	12
HUN	Orvosi eszközök és tartozékok hiszteroszalpingográfiához és kromopertubációhoz	13
ITA	Dispositivi medici e accessorio per isterosalpingografia e cromoperturbazione	14
LAV	Medicīniskās ierīces un piederumi histerosalpingogrāfijai un hromopertubācijai	15
LTU	Histerosalpingografijos ir laparoskopijos naudojant dažus medicinos priemonės ir priedai	16
NLD	Medische hulpmiddelen en toebehoren voor hysterosalpingografie en chromopertubatie	17
POL	Wyroby medyczne i akcesoria do histerosalpingografii i chromopertubacji	18
POR	Dispositivos médicos e acessórios para histerossalpingografia e cromotubação	19
RON	Dispozitive medicale și accesorii pentru histerosalpingografie și cromopertubare	20
SLK	dravotnicke pomôcky a príslušenstvo pre hysterosalpingografiu a chromopertubáciu	21
SLV	Medicinski pripomočki in dodatki za histerosalpingografijo in kromopertubacijo	22
SWE	Medicinsk utrustning och tillbehör för hysterosalpingografi och kromopertubation	23

Commercial Code	Brand name
1404020	Sliding Rigister Ø 20mm
1404025	Sliding Rigister Ø 25mm
1404030	Sliding Rigister Ø 30mm
1404035	Sliding Rigister Ø 35mm
Accessoire/ Accessory	
1507000	Bloc Connector

FRA

Description

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ENG

Description

The Sliding Register is used for HSG and chromopertubation. It includes:

- A tapered and transparent cup with a central cone-shaped cannula, fitted onto a cylindrical handle, external diameters 20, 25, 30 or 35 mm
- A lateral catheter attached to the cup to be connected to a vacuum system,
- A central injection tubing slides into the cup and has 4 markings spaced 2cm apart.

The cup is available in 4 different diameters: (20, 25, 30, 35 mm) allowing to adapt to different morphologies of the cervix.

It is a sterile, single-use device.

*The Bloc Connector is an accessory of the Sliding Register & the Register that allows their connection to a vacuum source.

Sterilization

Sterilisation with ethylene oxide

Indication

The Sliding Register is a medical device designed to be used during hysterosalpingography which is a radiographic examination. The indications of that using are the following case:

- Infertility with or without the presence of repeated miscarriages
- Congenital or anatomic abnormalities
- Anomalies of the menstrual cycle.
- Pain in the pelvis tract
- Abnormal menses
- Preoperative control for uterine or tubal surgery

The Sliding Register is also used for chromopertubation during laparoscopy, and allows to verify tubal permeability

Precautions of use

- The product is intended for use by health professionals with knowledge of gynaecology/obstetrics, in compliance with aseptic measures.
- Laboratoire CCD strongly recommends that the health professional take note of the information given in the IFU.
- The product is delivered in a sterile state. Before use, make sure that the packaging has not been opened, and/or damaged, check the expiry date, and the integrity of the product.
- The product is intended for single use only, and must not be cleaned, or resterilised. In case of cleaning, and/or re-sterilisation, its characteristics can no longer be guaranteed, and may cause cross-infection regardless of the method used.
- Use only under the conditions indicated.

Instructions for use

1. Place the patient on the fluoroscopic machine in a gynecologic examination position.
2. Clean the external genital area with antiseptic solution
3. Insert a sterile speculum into the vagina for good visualization of the cervix.
4. Disinfect the cervix with an antiseptic solution
5. Before introducing the catheter:
 - Allow at least 10 cm of the sliding injection catheter to stick out from the cup.
 - Connect the injection catheter to a syringe containing the solution to inject.
 - Purge the whole set. The solution's progression can be viewed easily within the sliding catheter.
 - Connect the lateral catheter to the vacuum system (Bloc connector or syringe luer lock).
6. Introduce the extremity of the injection catheter 1 or 2 cm into the cervical orifice.
7. Push gently the cup towards the cervix by sliding it along the injection catheter maintained immobile. The cup is easily guided along the catheter and fits automatically over the centre of the cervix providing a perfect seal. The guide-marks every 2 cm on the proximal section of the injection catheter indicate the length of the catheter inside the cervix. This length can be changed during the examination, as the catheter can slide in or out from the cup.
8. Create the negative pressure progressively while controlling the vacuum to seal the cervix without any risk of bruising. Use the clamp on the lateral catheter to avoid air backflow.
9. To avoid the discomfort of the patient, the opening of the speculum can be reduced during the procedure.
10. Inject slowly the radiopaque solution or sterile methylene blue inside the uterine cavity. Use the valve on the central catheter to adjust the amount of liquid injected.
11. The cup and cannula set is radiolucent and allows control of the proximal zone of the endocervix; the radiologist can thus operate at a distance from the radiological field while conveniently mobilizing the uterus.
12. Unclamp the side tubing and carefully remove the entire device at the end of the procedure.

Complications

- Discomfort or moderate pain
- Infection.
- Mild to moderate uterine cramps.
- Allergic reaction to the contrast agent.
- Uterine perforation and tubal rupture: very rare

Contraindications

Do not use these products:

- In case of pregnancy, or suspicion of pregnancy,
- In case of acute pelvic inflammation or infection
- In case of unexplained pelvic tenderness
- In case of presence of an intrauterine contraceptive device
- In case of allergy to contrast product

Storage conditions

- Single packaging.
- The product must not be used in the event of damaged, or open primary packaging before the medical procedure.
- Store away from direct sunlight.
- Store in a dry place.

Validity date

Use before the expiry date stated on the packaging.

Elimination

After use, this product presents a potential biological hazard. It should be disposed of in accordance with current guidelines for the handling, and disposal of potentially infectious materials.

Contact for medical information and materiovigilance

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