GIMA

PROFESSIONAL MEDICAL PRODUCTS

Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

MISURATORE DI PRESSIONE LEO CON SOFTWARE
LEO BLOOD PRESSURE MONITOR
WITH SOFTWARE
TENSIOMÈTRE LEO AVEC LOGICIEL
DE LEO MEDIDOR DE PRESIÓN
CON SOFTWARE
MEDIDOR DA TENSÃO LEO COM SOFTWARE

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Contents

Chapter 1	Safety Precautions	49
	1.1 Operation for AC Adapter (Separate Sale)	50
	1.2 Operation for Battery	51
Chapter 2	Main Unit	53
Chapter 3	Button Functions	56
Chapter 4	External Interfaces	57
Chapter 5	Dry Battery/AC Adapter Installation	58
	5.1 Dry Battery Installation	58
	5.2 Using the AC Adapter.....	59
Chapter 6	Setting the Date and Time	60
Chapter 7	Unit	61
Chapter 8	User Switch	62
Chapter 9	Applying the Arm Cuff	63
Chapter 10	Usage Method of Sphygmomanometer	64
	10.1 Accurate Measurement Way	64
	10.2 BP Measurement	65
Chapter 11	Memory Function	66
	11.1 Review Memory Values	66
	11.2 Delete Memory Values.....	67
Chapter 12	Alarm Function	68
	12.1 Physiological Alarm	68
	12.2 Technical Alarm	69

Chapter 13 SpO₂ Measurement Function	69
Chapter 14 SpO₂ Measurement Method	72
Chapter 15 Installation of the Software.....	73
15.1 Demand of Editor	73
15.2 Installation of Software	74
Chapter 16 Keys and Symbols	74
Chapter 17 Error Message.....	76
Chapter 18 Troubleshooting	77
Chapter 19 Maintenance, Cleaning and Keeping.....	79
Chapter 20 NIBP Specification.....	82
Chapter 21 SpO₂ Specification	84
Appendix	85

Chapter 1

SAFETY PRECAUTIONS

- Before use, carefully read “Safety Precautions” for a correct use.
- To prevent users suffered hurt or damnification due to improper use, refer to “Safety Precautions”, and use this product properly.

For safety reasons, be sure to comply with safety precautions.

Note

If not to use correctly, it exists that a potentially hazardous situation which may result in injury to the user or patient or damage to the equipment or other property.

Note

Human and the surrounding environment are in the range of electromagnetic radiation which generated by this device, there is no energy hazard to human or the surrounding environment, also not affected by the electromagnetic interference of the external device. This device has slight energy radiation, should be avoided too close to the high precision requirements equipment.

Note

Self-diagnosis and treatment using measured results may be dangerous. Follow the instructions of your physician.

Please hand measurement results to the doctor who know your health to accept diagnosis.

For severe blood circulation disorder or arrhythmia patients, please use the device under the guidance of a doctor.

Otherwise it may lead to acute hemorrhage, or measurement error as a result of squeezed arm.

For Infant and the person who can't express oneself, please use the device under the guidance of a doctor.

Otherwise it may cause accident or dissension.

Please do not use for any other purpose.

Otherwise it may cause accident or holdback.

Please use special cuff.

Otherwise it is possible that measurement result is incorrect.

Please do not keep the cuff in the over-inflated state for a long time.

Otherwise it may cause risk.

Do not disassemble, repair and change the device.

Otherwise it can not measure correctly.

Do not use the device in the case of there are flammable anesthetic gasses mixing with the air or nitrous oxide.

1.1 Operation for AC Adapter (Separate Sale)

Note

Please use the dedicated AC adapter of this device.

Otherwise it may cause trouble.

The dedicated AC adapter must use AC 100V-240V.

Otherwise it may cause fire or electric shock.

When there is breakage of dedicated AC adapter plug or wire, please do not use it.

Otherwise it may cause fire or electric shock.

Please do not plug or unplug the adapter on the socket with wet hands.

Otherwise it may cause electric shock or injury.

1.2 Operation for Battery

Note

Please use 4 “AA” manganese or alkaline batteries, do not use batteries of other types.

Otherwise it may cause fire.

New and old batteries, different kinds batteries can not be used together.

Otherwise it may cause battery leakage, heat, rupture, and damage to the device.

Please don't put wrong positive and negative of battery. When the batteries power exhausts, replace with four new batteries at the same time.

Please take out the dry batteries when you do not use the device for a long time.

Otherwise it may cause battery leakage, heat, rupture, and damage to the device.

If electrolyte of the batteries immodestly get in your eyes, immediately rinse with plenty of clean water.

It will cause blindness or other hazards, Should immediately go to the nearest hospital for treatment.

If electrolyte of the batteries immodestly glues on the skin or the clothes, immediately rinse with plenty of clean water.

Otherwise it may hurt the skin.

Advice

Do not strike or drop the device.

Do not inflate before the cuff wraps around the arm.

Do not inflect the cuff and the air tube forcibly.

Description of functions

The Sphygmomanometer apply to measure the non-invasive blood pressure and SpO₂ of human (adult, children, neonate), uses three-user mode, each user can store 100 items records of measurement results at most. Each record includes detailed measuring time, systolic pressure, diastolic pressure, average

pressure, pulse rate and record number, etc. With 2.8 inch color LCD screen, clear interface, the function of data review is complete. User can implement ON/OFF, manual measuring, system setup, parameters change and other operations with seven buttons which are located on the front panel of the device. The Sphygmomanometer uses audible and visual alarm, when the battery power is low, the buzzer will intermittent buzzing and LCD screen displays “Low Power” to prompt user replace batteries. When the measurement data exceeds the set alarm limit, the font color of measurement results will change to red and the audible alarm will occur, user can open and close the alarm sound according to needs. With timing shutdown function, if there is no operation and SpO₂ measurement, the device will automatically turn off after 2 minutes. With USB interface, Users can send measurement results to PC. Refer to the help or explication of the related software for specific operation.

Purpose

The device apply to measure the non-invasive blood pressure and SpO₂ of human. Record parameter value of blood pressure to provide the reference for the health care professional.

Warning


Please use the device on the adult object who can read the user manual and the error message shown on the screen. Read the user manual before use the device in order to take actions according to the manual when something wrong with the device. For pediatric and neonatal population, measurement should be only performed by qualified personnel. And please make sure to select the right user mode and cuff before use .

Chapter 2

MAIN UNIT

The production is in the package. Open the package and confirm whether the production is whole.



Accessories

(Specification: limb circumference 22-32cm (middle part of upper arm), please choose suited cuff when measuring pediatric or other).


USB Data Line

Software CD

User Manual

Optional Accessories

AC adapter

Input: voltage: AC 100V~240V
frequency: 50Hz/60HZ
rated current: AC 150mA
Output: DC 6.0V \pm 0.2V 1.0A
or DC5.0V \pm 0.2V 1.0A


Integrated SpO₂ probe: Y10UCH150

(This part is only suitable for European Union market)

A. SpO₂ measurement

Measurement range: 0%~100%

Error: 70~100%: \pm 2%; Below 70%: unspecified

- B. PR (pulse rate) measurement
Measurement range: 30bpm~250bpm
Error: ± 2 bpm or $\pm 2\%$ (select the larger)
- C. Optical sensor: red light (wavelength: 660nm, 6.65mW)
infrared light (wavelength: 880nm, 6.75mW)

Note

- The optional probe of the Sphygmomanometer is an integrated SpO₂ probe, the measuring part is integrated with the probe.
- The service life of the integrated SpO₂ probe is three years.

**Cuff**

There are several suitable cuffs (range of limb circumference, middle part of upper arm)

the range of limb circumference is 6-11cm

the range of limb circumference is 10-19cm

the range of limb circumference is 18-26cm

the range of limb circumference is 32-43cm

**Note**

- The cuff is a consumable. Calculate by measuring 6 times a day(3 times each morning and evening), the service life of the cuff is about 1 year.(using our experimental conditions).
- In order to correctly measure blood pressure, please replace the cuff in time.
- If the cuff leaks, please contact our company to buy a new one. The cuff purchased separately does not include the airway tube plug. When replacing, please do not throw the airway tube plug away, install it on the new cuff.

**Note**

When the product and accessories described in this manual are about to exceed the period of use, they must be disposed according to relevant product handling specification. If you want to know more information, please contact our company or representative organization.

Chapter 3

BUTTON FUNCTIONS

All the operations to the Electronic Sphygmomanometer are through the buttons. The names of the buttons are above them. They are:

ON/OFF Press this button to turn on/off the device.

START/STOP Press it to inflate the cuff and start a blood pressure measurement. During measuring, press it to cancel the measurement and deflate the cuff.

At all levels interface, the three buttons correspond respectively with the text prompts below the LCD screen, pressing any button will carry on the corresponding function, eg: **MENU - ENTER - LIST - USER** etc.

Up and down buttons, respectively carry on the functions of moving the cursor up and down, changing the parameters and switching the status.


Chapter 4

EXTERNAL INTERFACES

Note

Please hold the airway tube plug to remove the NIBP cuff.


**Note**

The cuff socket and AC adapter socket of the sphygmomanometer only connect to the stated equipment which accord with the requirements of IEC 60601-1: 1988 “Medical electrical equipment-Part 1: General requirements for safety”.

Chapter 5

DRY BATTERY/AC ADAPTER INSTALLATION

The product can use dry battery or AC adapter as power source.

5.1 Dry Battery Installation

1. Demount the battery cover in the direction of the arrow.
2. Install “AA” dry batteries according to \oplus \ominus polarities.
3. Slide to close the battery cover.
Icon “

**Note**

Dispose of the exhausted batteries according to applicable local regulations about environmental.

5.2 Using the AC Adapter

1. Connect the device and the AC adapter. Insert the AC adapter plug into the AC adapter socket on the right side of the device.
2. Please insert the power plug of AC adapter into the 100V- 240V socket.

**Note**

- Hold and pull the Housing to remove the AC Adapter from the electrical outlet. Do not remove by pulling on the cord.
- Remove the AC Adapter plug from the unit.
Please be sure to use dedicated AC adapter.

**Note**

You'd better take off batteries when use the AC adapter as power source. If there is any damage to the AC adapter, you should use batteries to run the device.

When adapter and batteries are both used at the same time, the battery power will not be consumed. Switch adapter and battery as power supply when the device is off, otherwise, the device may shutdown due to power failure.

Chapter 6

SETTING THE DATE AND TIME

It is necessary to set date and time after turning on the device.

The Electronic Sphygmomanometer can automatically stores measurement results, date and time.

If dry battery power exhausts or removed, the time to stop.

At the moment, please reset date and time.


The Electronic Sphygmomanometer stores measure results of three users automatically, and up to 100 items for every user. If the date and time are set correctly, the date and time when measuring will be correct in the memory, otherwise it may not be correct. The results can be uploaded to PC via USB and processed with the PC software.


Note

Correctly use data upload function:

1. First, make sure the time is right, there are two modes of time setting:
 - (1) When using the Sphygmomanometer for the first time or after the Sphygmomanometer has been placed without power supply for a certain time (more than 5 minutes), after power on, there is a prompt of time error on the main interface, set date and time with **UP**, **DOWN** and **ENTER** button.
 - (2) Press **MENU** button on the main interface to enter system menu, then enter **SYSTEM TIME** item, the current time will be displayed on the screen. Set date and time with **UP**, **DOWN** and **ENTER** button.

2. After completing all time setting, select **CONFIRM** option and press **ENTER** button to confirm the setting value. If you do not want to change the time, select **EXIT** option and press **ENTER** button to return to the previous menu.

Note


Please choose the computer which should be ensured compliance with the requirements of IEC60950, or else it may damage the device.

Chapter 7 UNIT

There are two units: “mmHg” and “kPa”.

The default is: “mmHg”.

Enter the **SYSTEM SETUP** submenu in **SYSTEM MENU**, then select **UNIT** option to switch units between “mmHg” and “kPa”.

Chapter 8 USER SWITCH

The Electronic Sphygmomanometer stores the measure results of three users automatically, and up to 100 items for every user.

Press **USER** button in main interface to switch users.
Or press **USER PURVIEW** item in **SYSTEM SETUP** menu to switch users.



Note

When the **USER PURVIEW** is set to **ALL**, current user can be switched under main interface; when set to a certain user, it will not be able to switch under main interface.
User type can be set to adult, pediatric and neonatal three different kinds, setting method is as follows:

Chapter 9 APPLYING THE ARM CUFF

Both left and right arm can be measured.
Bare your arm or cloth close-fitting clothing during measurement.
Carry out the operation in a room with comfortable temperature.

**When measuring, take the thick clothes off instead of rolling up the sleeves.
In order to measure accurately, pay attention to applying the cuff properly (left arm).**


1. Insert the arm cuff air plug
in the cuff socket of sphygmomanometer.


2. Stretch cuff into a barrel for the arm
can conformable enter into the barrel.


3. Left arm penetrate through the cuff,
the air tube of the cuff will pass the top of your palm.


4. Wrap the cuff to your upper arm. Make the air tube
inside the forearm and aligned with your middle finger.

5. The bottom of the cuff should be
approximately 2cm~3cm above your elbow.

6. Be fixed with cloths, and wrapped tight cuff,
the arm and the cuff should not have gaps.


Chapter 10

USAGE METHOD OF SPHYGMOMANOMETER

10.1 Accurate Measurement Way

Measurement in quiet and relaxing state.

1. Place your elbow on a table.
2. The cuff is level with your heart.
3. The palm of the hand is up, and the body relax.

Advice

Try to measure your blood pressure at the same time each day with the same arm and the same pose for consistency.

The high and low location of cuff will cause changes in measurement results.

Do not touch Electronic Sphygmomanometer, cuff and windpipe during measurement.

Measurements should be taken in a quiet place and keep the body relaxed.

Remain still for 4~5 minutes before measurement.

Relax the body, do not let the muscle activity.

Do not talk and movement during the measurement.

Wait more than 4~5 minutes between measurements.


Do not use precision instruments near the Sphygmomanometer.

Please use the sphygmomanometer at an environment of suitable temperature and humidity (refer to Chapter 20), otherwise it will cause measurement error.

10.2 BP Measurement

1. Press **START/STOP** button to take a measurement.

During measurement, please keep correct pose and quiet state, do not move.

Wish to stop measurement

During measurement, press **START/STOP** button, device will stop inflating, and release the air from the cuff.

2. Confirm measurement value

Measurement value can be automatically stored. ([Memory Function] refer to Chapter 11)

*Self-diagnosis and treatment using measured results may be dangerous. Follow the instructions of your physician.

Note

- Wait at least 4-5 minutes between two measurements.
When repeating measurements, because the arm appears congestion, it may not get correct blood pressure measurement. After the blood flow, take a measurement again.
 - When some factors affect the measurement results in measurement process, error messages will appear on the screen, you can obviate the malfunction and restart a measurement.
3. In no physiological alarm state, press any button to carry on the corresponding button function; in audio alarm state, press any button (except **ON/OFF** button) to clear up the audio alarm.
 4. Take off the cuff, press **ON/OFF** button to turn the device off.


Chapter 11

Memory Function

The sphygmomanometer is designed to store the blood pressure, pulse rate values and the date and time when measured, which are up to 100 groups. If there have been stored 100 groups, the earliest results will be deleted when saving the 101 group of measurement results.

11.1 Review Memory Values


1. In the main interface (interface when boot-strap), press **MEMORY** button to review the most recent measurement values in large-print with the serial number from 1 to 100.


2. Press **UP / DOWN** button to circularly switch the former measurement values.


*The right figure shows that there is no measurement result can be displayed.

3. Press **LIST** button to display the data list interface.

4. Press **TREND** button to display trend interface.


End to display the measurement values.

Press **EXIT** button to return to the main interface or hold **ON/OFF** button to turn the power off.

11.2 Delete Memory Values

Users can delete all memory values instead of separately delete the special item.

1. Press **MENU** button to enter the system menu, select **DELETE DATA** option to enter its interface, select the user whose data will be deleted, after confirming again, all measurement results of the selected user will be deleted.

2. Finish Operation
Select **CONFIRM** or **EXIT** to return to the previous menu, or hold **ON/OFF** button to turn the power off.

Chapter 12

ALARM FUNCTION


Alarms are classified into two categories: technical alarm and physiological alarm.

12.1 Physiological Alarm

The sphygmomanometer has the function of overrun alarm, user can press **MENU** button to enter system menu, select **ALARM SETUP** option to enter its interface, then set the limit of blood pressure value, When the BP measurement result is higher than the high limit or lower than the low limit and the

alarm is ON, the physiological alarm will occur; in **ALARM SETUP** interface, select **SpO₂ ALARM SETUP** option to enter its interface, when the SpO₂ measurement result is higher than the high limit or lower than the low limit and the alarm is ON, the physiological alarm will occur.

In the state of physiological alarm, press any button to cancel the alarm and it does not affect the next alarm; the alarm can be closed permanently with alarm switch of the alarm setup menu until the alarm switch be opened again.

12.2 Technical Alarm

When power is about to exhaust and alarm is on, the alarm will occur. This alarm can not be cancelled unless being closed or the power replaced.

Chapter 13 SpO₂ MEASUREMENT FUNCTION

(This chapter is only suitable for European Union market)

Precaution for SpO₂ measurement

Note

- Make sure the nail covers the light. The probe wire should be on the backside of the hand.
- SpO₂ value always displays in the fixed place.
- Testee can not use enamel or other makeup.
- Testee's fingernail can not be too long.
- The fingers which are too cold or too thin may affect the measure accuracy, please insert the thicker finger such as thumb or middle finger deeply enough into the probe.
- The SpO₂ probe is suitable for children above four years old and adults(Weight should be between 15Kg to 110Kg). The device may not work for all patients. If you are unable to achieve stable readings, discontinue use it. Do not use the device on infant or neonatal patients.
- The update period of data is less than 5 seconds, which is changeable according to different pulse rate values.
- Pulse rate waveform through the normalization, when pulse rate waveform becomes smooth and stable, the read value is optimal value, and the waveform at the moment is the most standard one.
- The maximum temperature for the contact surface of the device with the body is less than 41°C, and the temperature is measured by a temperature measuring device.
- The SpO₂ probe has been calibrated before leaving factory.

Warning

- Pulse oximeter can overestimate the SpO₂ value in the presence of Hb-CO, Met-Hb or dye dilution chemicals.
- Cables of ES (Electrosurgery) equipment and SpO₂ probe must not be tangled up.
- Do not put the SpO₂ probe on extremities with arterial catheter or venous syringe.


- Do not perform SpO₂ measuring and NIBP measuring on same arm at one time, because obstruction of blood flow during NIBP measuring may adversely affect the reading of SpO₂ value.
- Please don't measure this device with function tester for the device's related information.
- The person who is allergic to rubber can not use this device.
- For the special patients, there should be a more prudent inspecting in the placing process. The device can not be clipped on the edema and tender tissue.
- The light (the infrared is invisible) emitted from the device is harmful to the eyes, so the user and the maintenance man should not stare at the light.
- Check if the cable of SpO₂ probe is in normal condition before measuring. After unplugging the SpO₂ probe cable from the socket, "SpO₂%" and "bmp" on the screen will disappear.
- Do not use the SpO₂ probe once the package or the probe is found damaged. Instead, you shall return it to the vendor.
- The uncomfortable or painful feeling may appear if using the device ceaselessly, especially for the microcirculation barrier patients. It is recommended that the sensor should not be applied to the same finger for over 2 hours. Prolonged and continuous monitoring may increase jeopardy of unexpected change of dermal condition such as abnormal sensitivity, erubescence, vesicle, repressive putrescence, Particularly in newborns or in a Perfusion disorders and changes or immature skin form of the patient. According to skin quality change, correct optical path alignment and attachment methods to regularly check the place of SpO₂ probe, and change the attachment position when the quality of skin decline. More frequent examinations may be required for different patients.
- Please refer to the correlative literature about the clinical restrictions and caution.
- This device is not intended for treatment.
- Please don't use the SpO₂ probe during MRI or CT scanning.
- When the device is ON, if power interrupt for more than 30s, the SpO₂ probe needs no operation after the power is restored, after the device is turned on, ensure that the SpO₂ probe can be used normally.

Chapter 14

SpO₂ MEASUREMENT METHOD

(This chapter is only suitable for European Union market)

1. Attach the SpO₂ probe to the appropriate site of the patient finger as the following figure.

Place SpO₂ probe

2. Plug the connector of the SpO₂ probe cable into the USB socket in the lower right of the device.
The main interface will switch to SpO₂ interface. The update period of data is less than 5 seconds, which is changeable according to different individual pulse rate. This operation brings no affection to other functions.

Warning

Uncomfortable or painful feeling may appear if using the device ceaselessly, especially, for the microcirculation barrier patients. It is recommended that the sensor should not be applied to the same finger for over 2 hours.


Measurement Limitation

During operation, the accuracy of oximeter readings can be affected by:

- High-frequency electromagnetic interference such as electrosurgical apparatus connected to the system.
- Intravenous dyestuff.

- Excessive patient movement.
- Outside ray radiation.
- Improper SpO₂ probe installation or incorrect contact position of the patient.
- Temperature of SpO₂ probe(optimum temperature between 28°C and 40°C).
- Place the SpO₂ probe on an limb that has a blood pressure cuff, arterial catheter, or intravascular line.
- Concentrations of dysfunctional hemoglobin, such as carboxyhemoglobin (COHb) and methemoglobin (MetHb).
- SpO₂ is too low.
- Bad circular injection of the part being measured.
- It is required to use SpO₂ probe which is provided by our company, contact with our sale department when changes SpO₂ probe.

Chapter 15

Installation of the Software

15.1 Demand of Editor

Pentium IV 1.8G or more

Operation System: Windows XP

EMS memory: 256M or more

Hard Disk: 40G or more

Display: 17 inch or more

CD-ROM

USB: 2 or more

Resolution of printer: 600 DPI

15.2 Installation of Software

1. Place the CD-ROM in the CD-ROM compartment located on your computer.
2. If Auto Play for CDs is enabled, place CD in reader and follow instructions when they appear in the screen; otherwise follow install instructions below:
Open Windows Explorer.
Click on the root CD-ROM directory.
Double click file setup software.
Follow the instructions in the screen.
Refer to “Software Help” for details about the operation method of the PC software.

Chapter 16 KEYS AND SYMBOLS

Symbol	Description	Symbol	Description
	Instruction manual/booklet	PED	Pediatric
		NEO	Neonatal
		INFO	Information
SYS	Systolic pressure		Open the alarm sound indication
MAP	Mean pressure		
DIA	Diastolic pressure		Close the alarm sound indication
PR	Pulse rate (bpm)		
ADU	Adult		Low-power

	Full-power	SN	Serial number
	<ol style="list-style-type: none"> 1. no NIBP data to review 2. no finger inserted to SpO₂ probe (This item is only suitable for European Union market) 3. An indicator of signal inadequacy 		<p>This item is compliant with Medical Device Directive 93/42/EEC of June 14, 1993, a directive of the European Economic Community</p>
			European Representative
	Class II equipment		<p>USB socket connect SpO₂ probe (This item only suitable for European Union market)</p>
	WEEE (2002/96/EC)		
	Type BF Applied Part		

Chapter 17

ERROR MESSAGE

Error message will be displayed in the screen if there is something wrong when measuring. The causes and solutions are shown as follows:

Error Message	Causes	Solutions
Self-test failure System failure	Function abnormal	Please contact us
Loose cuff	Cuff is not connected correctly.	Correctly connect cuff (refer to Chapter 9)
Air leakage	Cuff plug fall off	Make sure the cuff plug is securely inserted in the windpipe (refer to Chapter 9)
Air pressure error	Air pressure error	Refer to the troubleshooting
Weak signal	The pulse signal is too weak or the cuff is loose.	Correctly connect cuff (refer to Chapter 9)
Overpressure	Cuff is blocked or squeezed	Correctly connect cuff (refer to Chapter 9)
Excessive movement Over range Saturated signal	The signal extent is too big owing to the arm or body moving or other reasons when measuring	Keep arm, body still, measure again
Time out	It takes too much time	

Chapter 18

TROUBLESHOOTING

Abnormal Phenomenons	Causes	Solutions
BP measurement values too high or too low	Cuff is not connected correctly	Correctly connect cuff (refer to Chapter 9)
	Talk or move arms when measuring	Keep quiet and restart a measurement
	The turnup clothing presses the arm	Take off the clothing which presses the arm, and restart a measurement
No pressure	Cuff leakage	Buy a new cuff
	The cuff windpipe is not correctly connected with cuff	Correctly connect
	Cuff is not inflated	Stop using the device and contact us
Cuff deflates in short time	Loose cuff	Correctly apply cuff
It can not carry on measurement when press the measurement button		Switch on the power once again and restart a measurement
Power off suddenly when inflating	No use for a long time, the power of batteries can be exhausted owing to the changed temperature	Replace all four batteries with new ones

Hold the on/off button but can not start the device	Power of batteries can be exhausted	Replace all four batteries with new ones
	The battery polarities is reversed	Check the battery installation for proper placement of the battery polarities
Cuff inflation start before press the measurement button or never stop inflating when measuring		Pull out the cuff to deflate. Stop using the device and contact us
Cuff never deflation		Pull out the cuff to deflate. Stop using the device and contact us
Air pressure error	No deflation or deflation error or inflation without stop	Pull out the cuff to deflate. Stop using the device and contact us
	Others	Keep arm, body still, measure again
No press value displayed or the value unchanged or change erratically when cuff inflated		Pull out the cuff to deflate. Stop using the device and contact us
Other phenomenon		Switch on the power once again and restart an operation. Replace the batteries. If no, please contact us

Chapter 19

MAINTENANCE, CLEANING AND KEEPING

*Please follow the instructions in the user manual. If you do not comply, our company will not assume responsibilities of the quality.

Warning

Must remove batteries and cut off the alternating power supply before cleaning the sphygmomanometer, accessories must be clean separately.

Caution

- Do not high pressure disinfect the sphygmomanometer and accessories.
- Do not immerse the sphygmomanometer and accessories into any liquid.
- Do not use any sphygmomanometer and accessories that may be damaged or deteriorated.

Cleaning

- If there is something dirty on the sphygmomanometer and accessories, please wipe it with a dry and soft cloth.
- If the sphygmomanometer and accessories are very dirty, you can clean them with a soft cloth dampened with water or neutral detergent after the cloth is full twisted.

Maintenance

*Please follow the instructions in the user manual. If you do not comply, our company will not assume responsibilities of the quality.

- Frequently clean the Electronic Sphygmomanometer.
- The device should be inspected regularly (or according to the inspection

procedure of hospital). The device can be inspected by inspection organizations which designated by the state, or inspected by professionals, also, users can contact with our company.

Warning

Do not submerge the device in water.

Advice

- Do not use any naphtha, thinner or gas to clean.

- Do not attempt to clean or wash the cuff.

Keeping

Advice

Do not place the machine in the following areas:

- Easy to splash water areas.
- Direct sunlight, extreme hot, humidity, dust, causticity gas areas.
- Lean or the area which can cause vibration, impact.
- Chemicals or corrosive gas storage areas.
- Remove the batteries if the unit will not be used for long time.

**Warning**

It is recommended that you check if there is any damage on the Sphygmomanometer or the accessories regularly, if you find any damage, stop using it, and contact our Customer Service immediately.

**Warning**

The disposal of scrap instrument and its accessories and packing (including batteries, plastic bags, foams and paper boxes) should follow the local laws and regulations.

EMC declaration

-  When this device is installed or putted into service, EMC should be paid more attention, as the portable and mobile RF communications equipment with higher EM interference can affect this device.
-  The internal components and cables should not be changed, as this may decrease the IMMUNITY of the device.
-  The Electronic Sphygmomanometer should not be used adjacent to or stacked with other equipment.

Classification

EMC: Group I Class B.

According to the MDD 93/42, the classification of this medical device: II a.

The type of protection against electroshock: Class II equipment.

The degree of protection against electroshock: type BF applied part.

The degree of protection against ingress of water: IPX0.

According to the mode of operation: Continuous.

According to the degree of safety of application in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

Chapter 20

NIBP SPECIFICATION

Name	Electronic Sphygmomanometer	
Display mode	2.8" color LCD Display	
NIBP specification		
Measurement method	Oscillometric method	
Working mode	Automatic	
Operation mode	Continuous operation	
Measurement Range	Pressure	adult 0~290mmHg (0~38.6kPa)
		pediatric 0~235mmHg (0~31.3kPa)
		neonatal 0~140mmHg(0~18.6kPa)
	Pulse: 40~240/min	
Overpressure protection	adult mode	295±5mmHg (39.33±0.67kPa)
	pediatric mode	240±5mmHg (32±0.67kPa)
	neonatal mode	145±5mmHg (19.33±0.67kPa)
Inflation	adult	160±5mmHg (21.33±0.67kPa)
	pediatric	120±5mmHg (16±0.67kPa)
	neonatal	70±5mmHg (9.33±0.67kPa)
Alarm Range	adult mode	SYS ALM: 40~270mmHg - DIA ALM: 10~215mmHg
	pediatric mode	SYS ALM: 40~200mmHg - DIA ALM: 10~150mmHg
	neonatal mode	SYS ALM: 40~135mmHg - DIA ALM: 10~100mmHg

Resolution	Pressure: 1mmHg(0.133kPa)
Accuracy	Static pressure: ± 3 mmHg(± 0.4 kPa)
Error	The BP Value of the device is equivalence with the measurement value of Stethoscopy. The error meets the requests in the ANSI/AAMI SP-10:2002+A1:2003 +A2:2006
Operating temperature/ humidity	+5°C~40 °C 15%RH~80%RH
Transport	Transport by general vehicle or according to the order contract, avoid pounded, shake and splash by rain and snow in transportation
Storage	Temperature: -20°C~+55°C; Relative humidity: $\leq 95\%$; No corrosive gas and drafty
Atmospheric pressure	70hPa~1060hPa
Power supply	4 "AA" alkaline batteries, AC Adapter (AC, 100V-240V, optional)
Battery life	When the temperature is 23°C, limb circumference is 270mm, the measured blood pressure is normal, 4 "AA" alkaline batteries cab be used about 300 times
Dimensions	130(L)*110(W)*80mm(H)
Unit Weight	300gram
Safety classification	class II, type BF applied equipment
Service life	The service life of the device is five years or 10000 times of BP measurement
Accessories	Standard Configure: Adult Cuff: limb circumference 22-32cm (middle of upper arm) Software CD, User Manual, USB data line, four "AA" alkaline batteries Separate Sale: Pediatric Cuff: limb circumference 10-19cm (middle of upper arm) Neonatal Cuff: limb circumference 6-11cm (middle of upper arm) AC Adapter Input: AC 100-240V 50-60Hz - AC 500mA Output: DC 6.0V \pm 0.2V 1.0A - DC 5.0V \pm 0.2V 1.0A

Chapter 21

SpO₂ SPECIFICATION

(This chapter is only suitable for European Union market)

Name	SpO ₂ Probe (Accessory Separate Sale)	
Model	Y10UCH150	
Measurement Range	SpO ₂ Measuring Range: 0%~100%; Pulse Rate Measuring Range: 30bpm~250bpm	
Resolution		
SpO ₂	1%	
PR	1bpm	
Measurement Accuracy		
SpO ₂	70%~100% ±2% 0%~69% undefined	
PR	±2bpm or ±2% (select larger)	
Measurement Performance In Weak Filling Condition		
pulse-filling ratio: 0.4%	SpO ₂ error	±4%
	pulse rate error	±2bpm or ±2% (select larger)
Optical Sensor		
Red light	wavelength is 660nm, 6.65mW	
Infrared	wavelength is 880nm, 6.75mW	

Appendix

Guidance and manufacturer's declaration - electromagnetic emissions for all EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration - electromagnetic emission		
The device is intended for use in the electromagnetic environment specified below. The customer of the user of the device should assure that it is used in such and environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic immunity for all EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration - electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.

<p>Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11</p>	<p><5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec</p>	<p><5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec</p>	<p>Mains power quality should be that of a typical commercial or hospital environment. The device can continue the operation during power mains interruptions due to the usage of battery.</p>
<p>Power frequency (50/60Hz) magnetic field IEC61000-4-8</p>	<p>3A/m</p>	<p>3A/m</p>	<p>Mains power quality should be that of a typical commercial or hospital environment.</p>
<p>NOTE: U_T is the a.c. mains voltage prior to application of the test level.</p>			

Guidance and manufacturer's declaration - electromagnetic immunity for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration - electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	<p>Recommended separation distance</p> $d = \left[\frac{3.5}{V^1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E^1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E^1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

A Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

B Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM - for EQUIPMENT or SYSTEM that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the device

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{V^1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{E^1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E^1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Disposal: *The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment. For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.*

GIMA WARRANTY CONDITIONS

Congratulations for purchasing a GIMA product. This product meets high qualitative standards both as regards the material and the production. The warranty is valid for 12 months from the date of supply of GIMA. During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel traveling expenses and packaging not included. All components subject to wear are not included in the warranty. The repair or replacement performed during the warranty period shall not extend the warranty. The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use. GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc. The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed. The defected products must be returned only to the dealer the product was purchased from. Products sent to GIMA will be rejected.

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