5 Operating Guide

- 5.1 Insert the two batteries properly to the direction, and then replace the cover.
- 5.2 Open the clip as shown in Figure 5.

Figure 5. Put finger in position

- 5.3 Let the patient's finger put into the rubber cushions of the clip (make sure the finger is in the right position), and then clip the finger.
- 5.4 Press the button once on front panel.

- 5.5 Do not shake the finger and keep the patient at ease during the

process. Meanwhile, human body is not recommended in movement status.

- 5.6 Get the information directly from screen display.
- 5.7 In boot-strap state, press button, and the device is reset.
- 5.8 Under non-measurement state, it will enter standby mode automatically when there is no operation within 5s.

Fingernails and the luminescent tube should be on the same side.

6 Maintain, Transport and Storage

6.1 Cleaning

- A. Before cleaning, ensure good room ventilation, remove the inflammable and combustible materials around it, and avoid open flame and electrical power source, etc.
- B. Remove the internal battery, then wipe the enclosure and finger clip cushions respectively with a medical gauze dipped in isopropanol (70%) (note: isopropanol produced by any formal company can be used.) for at least 1 minute. Dry them naturally or clean them with a clean and dry cloth.
- C. Determine whether the cleaning is adequate by visual inspection. If stains are found, please repeat the above steps.

6.2 Maintenance

- A. Check the main unit and all accessories periodically to make sure that there is no visible damage that may affect patient's safety and monitoring performance. It is recommended that the device should be inspected weekly at least. When there is obvious damage, stop using it.
- B. Please clean the device before/after using it according to the User Manual (6.1).
- C. Please replace the batteries in time when low-battery appears.
- D. Please take out the batteries if the device is not used for a long time.
- E. The device need not to be calibrated during maintenance.

6.3 Transport and Storage

- A. The packed device can be transported by ordinary conveyance or according to transport contract. During transportation, avoid strong shock, vibration and splashing with rain or snow, and it can not be transported mixed with toxic, harmful, corrosive material.
- B. The packed device should be stored in room with no corrosive gases and good ventilation. Temperature: -40°C~+60°C; Relative humidity: ≤95%.

7 Troubleshooting

Trouble	Possible Reason	Solution
The values can not be displayed normally or stably.	The finger is not properly inserted. The finger is shaking or the patient is moving. The device is not used in environment required by the manual. The device works abnormally.	Please insert the finger properly and measure again. Let the patient keep calm. Please use the device in normal environment. Please contact the after-sales.
The device can not be turned on	1) The battery is drained away or almost drained away. 2) The battery is assembled incorrectly. 3) The device's malfunction.	1) Please change batteries. 2) Please assemble the battery again. 3) Please contact the local service center.
The display disappears suddenly.	The device enters into the energy saving mode. Low battery. The device works abnormally.	Normal. Please change batteries. Please contact the after-sales.

8 Key of Symbols

Symbol	Description
	Type BF applied part
	Follow instructions for use
%SpO ₂	The pulse oxygen saturation(%)
PRbpm	Pulse rate (bpm)

	The battery voltage indication is deficient (change the battery in time avoiding the inexact measure)
	1. No finger inserted 2. Signal inadequacy indicator
	Battery positive electrode
	Battery cathode
	1. Exit standby mode. 2. Reset
	Serial number
	Alarm inhibit
	WEEE disposal
IP22	Covering Protection rate
	Temperature limit
	Humidity limit
	Atmospheric pressure limit
	This way up
	Fragile, handle with care
	Keep away from rain
	Recyclable
	Manufacturer
	Manufacture Date
	Authorized representative in the European Union
	Medical Device compliant with Directive 93/42/EEC
P/N	Material code
	Lot number
	Medical device
	Keep away from sunlight
	Imported by
	Unique device identifier

Note : Your device may not contain all the following symbols.

9 Function Specification

SpO ₂ [see note 1]	
Display range	0% ~ 99%
Measured range	0% ~ 100%
Accuracy [see note 2]	70%~100%: ±2%; 0%~69%: unspecified.
Resolution	1%
PR	
Display range	30 bpm ~ 250 bpm
Measured range	30 bpm ~ 250 bpm
Accuracy [see note 3]	±2 bpm during the pulse rate range of 30 bpm ~ 99 bpm and ±2% during the pulse rate range of 100 bpm ~ 250 bpm.
Resolution	1 bpm
Accuracy under low perfusion [see note 4]	Low perfusion 0.4%: SpO ₂ : ±4%; PR: ±2 bpm during the pulse rate range of 30 bpm ~ 99 bpm and ±2% during the pulse rate range of 100 bpm ~ 250 bpm.
Light interference	Under normal and ambient light conditions, the SpO ₂ deviation ≤ 1%
Pulse intensity	Continuous bar graph display, the higher display indicates the stronger pulse.

Optical sensor [see note 5]	
Red light	Wavelength: about 660 nm, optical output power: < 6.65 mW
Infrared light	Wavelength: about 905 nm, optical output power: < 6.75 mW
Safety class	Internally powered equipment, type BF applied part
International Protection	IP22
Working voltage	DC 2.6 V—3.6 V
Working current	≤ 25 mA
Power supply	1.5V (AAA size) alkaline batteries × 2 or rechargeable battery
Battery life	The device can continuously work for 20 hours when it was powered by two new batteries within the warranty period.
Dimension and Weight	
Dimension	57 mm(L) × 31 mm(W) × 32 mm(H)
Weight	About 50 g (with the batteries)

Note 1: the claims of SpO₂ accuracy shall be supported by clinical study measurements taken over the full range. By artificial inducing, get the stable oxygen level to the range of 70 % to 100 % SpO₂, compare the SpO₂ values collected by the secondary standard pulse oximeter equipment and the tested equipment at the same time, to form paired data, which are used for the accuracy analysis.

There are 12 healthy volunteers (male: 6, female: 6; age: 18~50; skin color: dark black: 3, medium dark: 1, light: 7, white: 1) data in the clinical report.

Note 2: because pulse oximeter equipment measurements are statistically distributed, only about two-thirds of pulse oximeter equipment measurements can be expected to fall within ±Arms of the value measured by a CO-OXIMETER.

Note 3: Patient simulator has been used to verify the pulse rate accuracy, it is stated as the root-mean-square difference between the PR measurement value and the value set by simulator.

Note 4: percentage modulation of infrared signal as the indication of pulsating signal strength, patient simulator has been used to verify its accuracy under conditions of low perfusion. SpO₂ and PR values are different due to low signal conditions, compare them with the known SpO₂ and PR values of input signal.

Note 5: optical sensors as the light-emitting components, will affect other medical devices applied the wavelength range. The information may be useful for the clinicians who carry out the optical treatment. For example, photodynamic therapy operated by clinician.

EMC

This equipmen is suitable for professional healthcare facility environments and home healthcare environments

Warning :

- Don't near active HF SURGICAL EQUIPMENT and the RF shielded room of an ME SYSTEM for magnetic resonance imaging, where the intensity of EM DISTURBANCES is high.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the this equipment, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Note :

- This equipment needs special precautions regarding EMC and needs to put into service according to the EMC information provided below.
- The basic performance: SpO₂ measured range: 70% ~ 100%, absolute error: ±2%; PR measured range: 30 bpm ~ 250 bpm, accuracy: ±2 bpm or ±2%, whichever is greater.
- When the device is disturbed, the data measured may fluctuate, please measure repeatedly or in another environment to ensure its accuracy.
- Other devices may affect this device even though they meet the requirements of CISPR.

Table 1:

Guidance and Declaration - Electromagnetic Emissions	
Emission test	Compliance
Radiated RF emissions CISPR 11	Group 1
Radiated RF emissions CISPR 11	Class B
Harmonic distortion IEC 61000-3-2	Not applicable
Voltage fluctuations and flicker IEC 61000-3-3	Not applicable

Table 2:

Guidance and Declaration - Electromagnetic Immunity		
Immunity test	IEC60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/	Not Applicable

output lines		
Surge IEC 61000-4-5:	±1 kV line(s) to line(s) ±2 kV line(s) to earth 0 % UT;	Not Applicable
Voltage dips and Voltage interruptions IEC 61000-4-11	0,5 .cycle .At0°,45°,90°,135°,180°, 225°,270°and315°. 0 % UT; 1 cycle and 70 % UT ; 25/30 cycles /Single phase:at 0°. 0 % UT; 250/300 cycle	Not Applicable
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	3 V 0,15MHz - 80 MHz 6 V in ISM and amateur radio bands between 0,15MHz to 80 MHz 80%AM at 1kHz	Not Applicable
Radiated RF IEC61000-4-3	10V/m 80 MHz-2,7GHz 80%AM at 1kHz	10V/m 80 MHz-2,7GHz 80%AM at 1kHz
OTE UT is the a.c.mains voltage prior to application of the test level		

Table 3:

Guidance and manufacturer's declaration - electromagnetic Immunity						
	Test Frequency (MHz)	Band (MHz)	Service	Modulation	IEC60601 1-1-2 Test level (V/m)	Compliance level (V/m)
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380 – 390	TETRA 400	Pulse modulation on b) 18 Hz	27	27
	450	430 – 470	GMRS 460, FRS 460	FM c) ±5 kHz deviation 1 kHz sine	28	28
	710	704 – 745	LTE Band 13,17	Pulse modulation on b) 217 Hz	9	9
	810	800 – 870	GSM 800/900, CDMA 900, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation on b) 18 Hz	28	28
	1720	1700 – 1845	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation on b) 217 Hz	28	28
	2450	2400 – 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation on b) 217 Hz	28	28
	5240	5100 – 5500	WLAN 802.11 a/n	Pulse modulation on b) 217 Hz	9	9
	5785	5800				

Disposal: The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies

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