CATETERE URETRALE STERILE MONOUSO
STERILE URETHRAL CATHETER FOR SINGLE USE
CATHÉTERA URÉTRAL STÉRILE À USAGE UNIQUE
CATÉTER URETRAL ESTÉRIL DE UN SOLO USO
CATETER URETRAL ESTÉRIL DE UTILIZAÇÃO ÚNICA
STERILER HARNRÖHRENKATHETER FÜR DEN EINMALGEBRAUCH
ΑΠΟΣΤΕΙΡΩΜΕΝΟΣ ΟΥΡΗΘΡΙΚΟΣ ΚΑΘΗΤΗΡΑΣ ΜΙΑΣ ΧΡΗΣΗΣ
STERYLNY CEWNIK UROLOGICZNY DO JEDNORAZOWEGO UŻYTKU
STERILNÍ URETRÁLNÍ KATÉTR NA JEDNO POUŽITÍ
STERIL URINRÖRSKATETER FÖR ENGÅNGSBRUK
STERIILI JA KERTAKÄYTTÖINEN VIRTSAPUTKIKATETRI
STERILNI URETRALNI KATETER ZA ENKRATNO UPORABO
STERILNÝ MOČOVÝ KATÉTER NA JEDNO POUŽITIE
CATETER URETRAL STERIL DE UNICĂ FOLOSINȚĂ
STERIELE URETHRALE KATHETER VOOR EENMALIG GEBRUIK
STERILNI KATETER ZA URETRU ZA JEDNOKRATNU UPORABU
STERIL URINRÖRSKATETER FÖR ENGÅNGSBRUK
STERILT URINRØRSKATETER TIL ENGANGSBRUG
СТЕРИЛЕН УРЕТРАЛЕН КАТЕТЪР ЗА ЕДНОКРАТНА УПОТРЕБА
VIENKARTINIO NAUDOJIMO STERILUS ŠLAPLĖS KATETERIS
STERILS URĪNIZVADKANĀLA KATETRS VIENREIZĒJAI LIETOŠANAI
STERIILNE UREETRA KATEETER ÜHEKORDSEKS KASUTAMISEKS

Gima 22003 - 22004 - 22005 - 22006 - 22007 - 22022 - 22023 - 22024 - 22025 - 22026

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ENGLISH

Product name:

Sterile Urethral Catheter for Single Use (Nelaton catheter)

Type:

Paediatric male, Paediatric female, Male, Female

Sizes:

Straight Tip/ Standard tip

Intended purpose:

For temporary emission by inserting it to bladder through urethra when patients can not micturate by themselves.

Intended users:

It can be used by trained doctors, nursing staff, patients or accompanying personnel.

Structural Composition:

The product is composed of drainage funnel and tube

Indications:

Urinary obstruction and urine retention, etc.

Contraindications:

Those who are allergic to PVC materials should not use this product.

Patient group:

Patients with urinary retention caused by urinary tract obstruction. Rescue critically ill patients.

Patients in need of diagnosis and treatment of bladder disease. People undergoing surgery.

Patients with spinal cord injury.

Performance characteristics:

Single-use and sterile

Method:

a. Wash your hands and put on gloves to clean around the urethral opening. Those who have foreskin should pull the foreskin down to expose the urethral opening and clean it.

b. After preparing for catheterization, open the individual package from the opening.

- c. Lubricate the catheter with medical liquid paraffin or medical silicone oil.
- d. Insert the urinary catheter into the bladder slowly through the urethra for normal catheterization.

Any residual risk and any undesirable side-effects:
none

Installation information:
none

Clinical benefits:

The product help patients to temporary emission.

Caution:

a. When using the product, it should strictly implement the requirements of aseptic operation specifications and relevant regulations, and can be used by trained doctors, nursing staff, patients or accompanying personnel.

However, when patients or accompanying staff conduct urinary catheterization for the first time, please use it under the guidance of professional medical staff. In the course of use, the operator or patient should be treated in a timely manner according to medical regulations.

b. If you find any possible quality doubt during use, please stop using it immediately and notify our company quickly.

c. The product should be sterilized with EO and the product should be sterile.

d. The product is disposable and should be destroyed after use.

e. Please check the individual package before using it.

f. Don't use if package is opened or damaged.

g. Sterile and validity for 5 years. Don't use if overdue.

h. This product needs to be used in conjunction with a container for urine. When using it, it should avoid splashing urine and polluting the environment.

i. Discard after use. After use, the products shall be collected, transported, stored and disposed to ensure that the products are harmless to human health and the environment. Or the product can be disposed in accordance with the local medical regulations for waste disposal.

j. If the product is used in a household environment, after use, the product shall be disposed as household hazardous waste.

The product shall be collected, transported, stored and disposed by a dedicated establishment or an undertaking.

Or the product can be disposed in accordance with the local medical regulations for waste disposal.

k. All serious incidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

Storage:

The product should be stored in a well-ventilated room temperature.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

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