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**ENGLISH****OXYGEN THERAPY AND AEROSOL THERAPY MASKS****INSTRUCTIONS FOR USE****INDICATIONS**

Therapeutic administration of oxygen.

**NB:** Disposable product - do not reuse.

The oxygen source must be in accordance to the regulations in force.

Devices must be connected and activated by qualified personnel.

**PREPARATION AND CONNECTION**

**MODEL: OS/100 - OS/100P - OS/100N**

**MEDIUM CONCENTRATION MASK**

Attach tubing connector to oxygen source. Attach other end to plastic output of mask. The mask has been designed to operate with a flow of between 5 and 10 LPM and to deliver oxygen flows of between 40% and 60%.

**MODEL: OS6K - OS/60K - OS/62K - OS/70K - OS/72K**

**VARIABLE CONCENTRATION MASK WITH CUP**

Connect prescribed diluter to corrugated tubing (for OS/62K - OS/72K select the concentration and check that the mobile part of the regulator is fully entered in its seat). Connect cup to diluter. Attach sure flow tubing to diluter. Attach other end to plastic output of mask. The recommended flows and relative percentages of the oxygen delivered are printed on the diluter or on the label.

**MODEL: OS/50 - OS/50E - OS/50P**

**NON REBREATHING MASK WITH PLASTIC BAG AND CHECK VALVE**

Flatten plastic bag. Attach tubing connector to oxygen source. Attach other end to plastic output of mask. Adjust the capacity in order to prevent the bag from ever deflating more than

half during inspiration. The oxygen flow delivered is contained between 90% and 100%.

**MODEL: OS/80 - OS/80P**

**AEROSOL MASK WITH NEBULIZER**

Unscrew cap from base and add prescribed medication. Reassemble. Secure mask to top of nebulizer of opening provided. Attach one end of tubing to nebulizer output and the other end to oxygen or compressed air source. Verify that the connection is correct (as reported below).

For the duration of the administration follow the medical prescription. For an optimal nebulisation the recommended flow is between 6 and 9 LPM with a pressure of 3.44 bar (344 kPa, equivalent to 50 PSI). Efficient nebulisation is guaranteed up to a pressure of 1.73 bar (173 kPa, equivalent to 25 PSI).

**ALL ABOVE MENTIONED MODELS: OXYGEN FLOW SETTING AND MASK DRESSING**

Make sure tubing is safely fixed to oxygen source. Deliver oxygen and regulate the flow as prescribed by physician. Place mask on patient's face covering both mouth and nose. Pass elastic strap over patient's head and ears, on the neck. Adjust strap tension to hold mask in position by pulling. Model metal nose piece to patient's face.

**MODEL: OS/110 - OS/110P - OS/110K**

**FOR TRACHEOTOMY PATIENTS**

OS/110 - OS/110P: connect aerosol tubing (not supplied with this package) between mask and gas source. Select correct filter flow at gas source and check for gas flow through mask device.

OS/110K: Connect prescribed diluter to corrugated tubing. Connect cup to diluter. Attach sure flow tubing to diluter. Attach other end to plastic output of mask. The recommended flows and relative percentages of the oxygen delivered are printed on the diluter.

OS/110 - OS/110P - OS/110K: position the elastic strap behind the neck and gently pull the end of the strap until the mask is secure. The mask inlet swivel is 360° to allow tubing to position itself for supine or upright patients. When using suction, loosen mask (by pulling strap in opposite direction or unsnapping one side of strap) and remove mask from suction area. Replace mask as previously noted. Warning: Be sure all connections are secure.

**CONTRA-INDICATIONS**

No contra-indications are described for oxygen therapy. For medicine administration refer to instructions supplied with medicine.

**WARNINGS**

If device is reused, contrarily to the present instructions, this may compromise following:

1. Not cleanliness of the device and possible presence of biological residues that might cause cross-infections.
2. Alteration of materials.
3. Loss of initial functional features of product.

- Before connecting to oxygen source verify its functioning.

- The accidental crushing of kinking of the connection tubing do not allow oxygen to flow properly and may cause the disjunction of the connector. It is therefore recommended that the users control carefully.

- An oxygen flow of more than 10 litres per minute may cause a disturbing noise in the system. Also possibly liquid, coming from the humidifier, can filtrate into the mixture that is sent to the patient.

- All above models are disposable and cannot be used for a continuous period of more than 30 days.

**WASTE DISPOSAL**

Waste coming from hospitals must be disposed of at thermal destruction only in authorised places. It is suggested to treat the waste coming from private houses with disinfectant solutions (Sodium hypo-chlorite, Chloro-oxidant, etc.) before disposal even though the regulations in force might not require this.

**GUARANTEE-NOTICE**

FIAB guarantees that the product complies with Directive 93/42/EEC and has been manufactured according to the procedures of FIAB Quality System certified ISO 13485. No responsibility may be ascribed to the producer who shall not be held liable for medical costs, direct or indirect damage due to lacking function or mis-function of above models, when used differently from the instructions for use. We recommend to report opportunely any

malfunction or defect of the product to FIAB Quality Assurance Service.

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**PRÉPARATION ET RACCORDEMENTS**

**MODELES: OS/100 - OS/100P - OS/100N  
A CONCENTRATION MOYENNE**

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