
GIMA

**PALLONE RIANIMATORE MONOUSO
DISPOSABLE RESUSCITATOR BAG
INSUFFLATEUR DE REANIMATION MONOUSAGE
BALON DE REANIMACION MONOUSO**

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Disposable resuscitator bag

The resuscitator bag Gima is conceived for use in addition to artificial respiration and cardiopulmonary reanimation.

Instructions for use

- Connect the oxygen tube to an adjustable oxygen source.
- In case it is connected to a reservoir, adjust the oxygen emission so that the reservoir completely expands during inhalation and almost completely empties during expiration thus recharging the bag.
- Before using the product, check for its efficiency preferably connecting it to a test lung: check that the connections, the reservoir, if present, and the valves complete the whole ventilation cycle correctly.
- In case of need connect the face mask to the apposite connector.
- Follow the necessary ventilation and cardiac resuscitation.
- Press the bag and observe the thorax expand during inhalation.
- Release the pressure on the bag and observe the thorax reduce during expiration.
- While ventilating check for:
 - Cyanosis symptoms
 - The correctness of the applied ventilation method
 - The airway pressure
 - The correct operation of all valves
 - The correct operation of the reserve and the oxygen tube

During ventilation the non-return valve could contaminate with vomit, blood or other secretions. In this case, stop ventilation and proceed as follows:

- Rapidly compress the bag several times so that to expel the substance through strong jets of air.
- If the non return valve is not clean yet: dip the valve in the water and rapidly compress the bag several times with strong jets of air. If the valve is still obstructed, throw out the bag.

Prescriptions

- Do not use the bag in presence of toxic substances.
- Remove the oxygen reserve and the supplementary valve if oxygen is not being given. The non-removal reduces the recharge speed and consequently the ventilation frequency.
- Do not give supplementary oxygen in presence of free flames.
- Do not use oil, grease or other hydrocarbon containing substances on the bag and its components. The supplementary oxygen, given under pressure, can cause explosions if mixed with hydrocarbons.
- The use of this device is reserved to qualified medical personnel trained in lung ventilation and cardiac resuscitation advanced techniques.
- Before using the device on patients, make sure you are able to assemble, disassemble and use it.
- During ventilation, always check the airways pressure with a pressure gauge
- Qualified personnel trained in PEEP (Positive End Expiratory Pressure) only can use PEEP technique with this kind of device.

Technical specifications

Bag capacity: 1600 ml adult mod. (34277); 500 ml child mod. (34248); 280 ml infant mod. (34249)
Ejection capacity: 700 ml (one hand) – 900 ml (two hands) adult mod. - 300 ml - 350 ml child mod.
- 150 ml - 225 ml infant mod.

Reserve capacity (optional): 2500 ml for adult and child models - 500 ml infant mod.

Ventilation cycle rate: 20 cycles/min. adult and child models - 40 cycles/min. infant mod.

Oxygen concentration: with reserve 99% - without reserve 45% for adult and child; 90% infant

Exerted pressure intervals: 60 ± 10 cmH₂O max. for adult mod.
 40 ± 5 cmH₂O for child and infant mod.


The results have been obtained in the following conditions:

Adult: Vt-600ml, conformity 0.02L/cmH₂O, resistance 20cm H₂O/L/s

Child: Vt- 300, conformity 0.01L/cmH₂O, resistance 20cm H₂O/L/s

Infant: Vt-70ml, conformity 0.01L/cmH₂O, resistance 20cm H₂O/L/s

The bag performance can vary according to the operator and different factors such as room temperature, patient's lung capacity, ventilation rate and operator's hands size.

Insufflateur de reanimation monousage

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