Italian
Medical
Touch

MANUALE D'ISTRUZIONI USER MANUAL

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The medical equipment is designed liner is used to collect organic fluids by aspiration. If may connected to medical suction equipment which must comply with UNI EN ISO 10079 standards. The medical device may be sterilised and is equipped with a floating overflow valve, located inside the lid with the function of blocking the exit of material from the jar.

GENERAL WARNING

READ THE INSTRUCTION MANUAL CAREFULLY BEFORE USING THE DEVICE

ONLY QUALIFIED PERSONNEL MAY USE THIS DEVICE

1. When opening the pack, carefully check that the equipment is intact. It should be free from damage which makes it difficult to use.
2. Pay particular attention to:
 - Avoid use of medical equipment by children and / or incapacitated persons without correct supervision;
 - Store and use the equipment in areas which are protected from atmospheric agents and away from any heat sources
 - Solvents or abrasive products may not be used for cleaning, which may damage the equipment's plastic components
 - Only use original spare parts.
3. **The medical device equipment must be used for which it has been designed.** Any other uses are considered improper and dangerous. The manufacturer accepts no liability for any consequences which may be caused by improper use.
4. **Disposal of waste and of the medical equipment must be carried out in accordance with the specific legislation in force in each country.**

Manufactured by:

CA-MI s.n.c. di Attolini Mario & C.
Via Ugo La Malfa n°31
43010 Pilaastro (PR) Italia
Tel. 0521 / 637133 – 631138
Fax. 0521 / 639041
E-mail: vendite@ca-mi.it
Web: www.ca-mi.it
Rev. 0 del 06/06/07

CARATTERISTECHNICAL CHARACTERISTICS

Jar / Lid Material	Polycarbonate
Dimension of entry fittings	Ø 9.5 (external diameter) Ø 6.5 (internal diameter)
Dimension of exit fittings	Ø 9.5 (external diameter) Ø 6.5 (internal diameter)
Dimension of suction tube	Min. Ø 6x10 – Max. Ø 8x14 (recommended dimensions)
Dimension of patient's tube	Min. Ø 6x10 – Max. Ø 8x14 (recommended dimensions)
Storage conditions	Room temperature: - 10°C ÷ + 40°C
Maximum interval on scale	100ml

SIMBOLOGY

	Device conform to the MDD 93/42/CEE
	Warning, read the instruction manuale

CLEANING ACCESSORIES

Wearing disposable latex gloves, proceed to clean the accessories as follows:

1. Grip the jar with your hand and turn the cap in a counter-clockwise direction.
2. Remove the red cage covering the float located in the jar cap (conical part with 3 openings)
3. Once the cage is removed, the float with the red seal comes out of the jar cap.
4. To disinfect the accessories, use either denatured alcohol or hypochlorate-based solutions, easily purchasable from chemist shops.
5. Once cleaning has been completed, insert the float in the cage, with the seal facing towards the cage opening.
6. Insert the opening of the cage into the seat in the jar cap.

It is possible autoclave the accessories cover and jar: insert the parts in autoclave and to effect a cycle of sterilization with vapor to the temperature of 121°C (pressure relative 1 bar) having care to position turned upside-down the graduated jar (with turned fund toward the tall one). After the sterilization and the cooling to the temperature the components' environment, to verify that this last don't result damaged; disassembled therefore the container for inhaled liquids following the inverse operations to the dismantlement.

The product's mechanical resistance is guarantee up to **30 cleaning and sterilising cycles**.

Over this limit, deterioration in the physical and mechanical properties of the plastic material may be seen. For this reason it's recommended to replace the components collection jar and lid.

ACCESSORIES SUPPLIED

DESCRIPTION
Jar
Lid with floating overflow valve
Seal for jar