Laringoscopi in fibra ottica
Fiber Optic Laryngoscopes
Laryngoscope à fibre optique
LWL-Laryngoskope
Laringoscopios de fibra óptica
Laringoscópios em fibra ótica
Λαρυγγοσκόπιο οπτικής ίνας

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Operating Instruction

Attention

Thank you for purchasing our product. This product meets the most stringent requirements regarding the selection of manufacturing materials and also the final control. Operating instructions manual relates to our brand Lunalite & Maxlite of Laryngoscope. These instructions should be followed to ensure durability of this product. The operator must carefully read and understand this manual thoroughly to keep the product performance durable and reliable for longer period.

After opening the packages, first of all it is necessary to check all the components against the standard configuration. Check that they are all present and in perfect conditions.

Features

The M.D. Lunalite & Maxlite Fiber Optic laryngoscope blades are manufactured from antimagnetic 18/8- quality stainless steel type AISI 303/304 which is highly resistant to corrosion and conform to the ISO 7376 standard.

The Maxlite F. O. Blades are built with an integrated F. O. bundle with no cavities to trap or body fluids, thus allowing the blade to be easily and decontaminates. This contributes largely to the elimination of cross infection.

High quality and bigger Fiber bundles in Maxlite F. O. blade ensure excellent light transmission. 6,000 lux with 2.5V Xenon lamp and 13,000 lux with 3.5V Xenon lamp. The LED Handles (2.5V, 3.5V) provide 3x brighter illumination than Xenon Handles (2.5V, 3.5V).

Maxlite Blades are maintenance free and autoclavable up to 134 C / 5 min approximately for 4,000 times. Gas Plasma sterilisation is permitted for Maxlite F. O. Laryngoscope.

Using the Maxlite Blades

Do not grip the lever when removing the blade. Do not apply any pressure, in the directions shown in the illustration Fig 1, which could force the blade and lever apart.

Damage to the linkage may occur, resulting in incorrect action, or stiffness of the adjustable tip. The Maxlite blade is assembled to the handle in the normal manner. The lever to operate the tip section will then extend behind the handle.

The lever should not be touched during the initial stage of use, until the tip of the laryngoscope has been inserted into the vallecula.

Once this stage has been reached, movement of the lever towards the handle will elevate the tip of the blade, and therefore lift the epiglottis, without the need to increase the force exerted by the main part of the blade. Release the lever before withdrawing the blade.

1. Engage The blade by aligning The slot of the blade on to the hook pin of the handle and apply a sufficient force 1N-5N to make engagement as show in Fig 2.
2. Apply force upward to bring the blade in operating position as in Fig 3.
3. To bring the blade in stand by position apply force downward as in Fig 3.

Light Guide Replacement Procedure

1. Remove locking screw with a screw driver as show in Fig 4.
2. Pull out the green block and slide out the light guide.
3. Fix new light guide of similar size and replace The locking screw.
4. Be sure that the screw is properly screwed in Maxlite Blades.
5. To remove the light guide, first remove screw A, to allow removal of lever B as shown in Fig 5. The Maxlite blade levering mechanism may require articulation, before The light guide can be detached from The blade.

Great care should be taken while performing this procedure to avoid structural damage to the fiber blade.

Fiber Optic Handle Lamp Replacement Procedure

1. Unscrew head from the barrel counter-clock wise.
2. Remove the lamp from head.
3. Replace the new MD lamp.
4. Screw the Head clockwise to the barrel

Always try to keep the lens of lamp clean for better performance.

Battery Replacement Procedure

1. Unscrew bottom cap of handle and remove batteries.
 2. For greater longevity alkaline batteries are recommended as replacement. Ordinary carbonic batteries may also be used.
 3. Replace with appropriate size batteries making sure the +ve and -ve terminals are placed as described below
- * 2.5 V handles can be used with dry batteries as well as rechargeable batteries.
 - * 3.5 V handles can be used with rechargeable battery only.

Care & Maintenance

Remove batteries before cleaning, high level disinfecting or sterilization of the laryngoscope system.

Blades

1. Cleaning Procedure:

Immediately after use, the laryngoscope system should be rinsed under cool running tap water until all visible soil is removed. Ensure that all hard-to-reach areas are flushed with the running tap water. Immerse sealed laryngoscope system in a presoak enzymatic cleaner solution, which was prepared in accordance to manufacturers recommendations for a minimum of two minutes. Remove device from enzymatic cleaner solution and rinse with lukewarm running tap water for a minimum of one minute to remove all residues and visible soils. Then, immerse device in enzymatic detergent. Remove bottom cap and brush thoroughly using a soft-bristle brush, while ensuring that all hard to reach are reached visible soils / residue removed. Dry with lint free, clean cloth or filtered pressurized air.

Follow with HIGH-LEVEL DISINFECTION or STEAM STERILIZATION PROCEDURE.

Ultrasonic cleaning is strictly prohibited

2. Disinfection:

Soaking in solutions or Thermo chemically in a washer Sterilizer up to 65°C maximum may perform disinfection. Manufacturers instruction regarding duration and concentration of solutions should be strictly adhered. After Disinfection, rinse throughly in sterile water and dry with a lint free clean cloth.

3. Cold Soak Solution:

To achieve a high-level disinfection, Cidex® OPA or 2.4% Glutaraldehyde solution may be used according to manufacturers instructions. Dry with lint free, clean cloth or filtered pressurized

air. Reassemble all parts, load handle with batteries and test the system for proper function. if not functional, review battery/ lamp testing instructions below

Do not immerse Blades in Bleach, Betadine or A Potassium Hydroxide solutions. Doing so will severl damage instruments also avoid metal to metal contact after soaking, the blades should be rinsed under sterile water to remove chemical residues and dry with lint free clean cloth or filtered pressurized air.

4. Sterilization:

Before performing any of the procedures described below, The blade should be cleaned as described in the cleaning Procedure.

Note: Its is recommended to remove the fiber optic light guide from the blade before sterilization it effect the polishing of fiber and decreases the light output.

5. Gas Sterilization:

Gas sterilization by Ethylene oxide up to a maximum temperature of 65°C and 8 p.s.i, may be performed, which is preferred especially if sterilization is to be performed regularly

6. Steam Sterilization:

Steam Sterilization can also be performed. Insert device in appropriate autoclave pouch.

(A) GRAVITY DISPLACEMENT STEAM

- * **Temperature:** 121°C (250°F) 132°C (270°F)
- * **Cycle Time:** 30 Minutes 4 Minutes
- * **Dry Time:** 15 Minutes 20 Minutes

NOTE: Do not exceed temperature of 135°C and pressure of be avoided as these processes will damage the instrument.

Flash autoclaving and hot air sterilization should be avoided as these processes will damage the instrument.

Handles

Cleaning/Sterilization:

Battery Handle with stand the same cold soak solution and autoclave ranges outlined in the blade section.

However, the batteries and lamp must be removed prior to disinfection / sterilization. Battery Handles can with stand exposure to the ethylene oxide. Lamp may be cleaned with cotton ball dampened in alcohol (IPA).

The main handle and a cap will also withstand and exposure to ethylene oxide.

Do not allow excess fluid to seep into electrical contact, batteries must be removed before cleaning & sterilization.

Disposal: *The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.*

For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.

GIMA WARRANTY CONDITIONS

Congratulations for purchasing a GIMA product.

This product meets high qualitative standards both as regards the material and the production.

The warranty is valid for 12 months from the date of supply of GIMA.

During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel traveling expenses and packaging not included.

All components subject to wear are not included in the warranty.

The repair or replacement performed during the warranty period shall not extend the warranty.

The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use.

GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc.

The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed.

The defected products must be returned only to the dealer the product was purchased from.

Products sent to GIMA will be rejected.

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