PULSOSSIMETRO OXY-4
OXIMETER OXY-4
OXYMÈTRE DE POULS OXY-4
PULOXIMETER OXY-4
PULSIOXÍMETRO OXY-4
MEDIDOR DE OXI-PULSAÇÕES OXY-4
OXIMETER OXY-4
OKSYMETR OXY-4

MANUALE D'USO E MANUTENZIONE
USE AND MAINTENANCE BOOK
INSTRUCTIONS DE FONCTIONNEMENT ET ENTRETIEN
BETRIEBS UND WARTUNGS ANWEISUNGEN
MANUAL DE USO Y MANTENIMIENTO
MANUAL DE USO E MANUTENÇÃO
ΕΓΧΕΙΡΙΔΙΟ ΧΡΗΣΗΣ ΚΑΙ ΣΥΝΤΗΡΗΣΗΣ
PODRĘCZNIK EKSPLOATACJI I KONSERWACJI

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Instructions to User

Read these instructions carefully before using this equipment. These instructions describe the operating procedures to be followed strictly. Failure to follow these instructions can cause measuring abnormality, equipment damage and personal injury. The manufacturer is NOT responsible for the safety, reliability and performance issues and any monitoring abnormality, personal injury and equipment damage due to user's negligence of the operation instructions. The manufacturer's warranty service does not cover such faults.

- The uncomfortable or painful feeling may appear if using the device ceaselessly, especially for the microcirculation barrier patients. It is recommended that the sensor should not be applied to the same finger for over 2 hours.
- For the individual patients, there should be a more prudent inspecting in the placing process. The device can not be clipped on the edema and tender tissue.
- The light (the infrared is invisible) emitted from the device is harmful to the eyes, so the user and the maintenance man, can not stare at the light.
- Testee can not use enamel or other makeup on the finger.
- Testee's fingernail can not be too long.
- Please peruse the relative content about the clinical restrictions and caution.

1. Safety

1.1 Instructions for Safe Operations

- Check the main unit and all accessories periodically to make sure that there is no visible damage that may affect patient's safety and monitoring performance. It is recommended that the device should be inspected once a week at least. When there is obvious damage, stop using the oximeter.
- Necessary maintenance must be performed by qualified service engineers ONLY. Users are not permitted to maintain it by themselves.
- The oximeter cannot be used together with the devices not specified in User's Manual.

1.2 Warnings

- *Explosive hazard—DO NOT use the oximeter in environment with inflammable gas such as some ignitable anesthetic agents.*
- *DO NOT use the oximeter while the testee measured by MRI and CT*
- *To dispose the device, the local law must be followed.*

1.3 Attentions

- *Keep the oximeter away from dust, vibration, corrosive substances, explosive materials, high temperature and moisture.*
- *If the oximeter gets wet, please stop operating it. When it is carried*

from cold environment to warm and humid environment, please do not use it immediately.

- DO NOT press the keys on front panel with sharp materials.

- High temperature or high pressure steam disinfection to the oximeter is not permitted. Refer to User's Manual for instructions of cleaning and disinfection.

- DO NOT have the oximeter immersed in liquid. When it needs cleaning, please wipe its surface with disinfect solution by soft material. Do not spray any liquid on the device directly.

- When cleaning the device with water, the temperature of water should be lower than 60°C.

2. Overview

The pulse oxygen saturation is the percentage of HbO₂ in the total Hb of the blood, so-called the O₂ concentration in the blood. It is an important bio-parameter to the respiration. Many of the respiration disease will cause hypoxemia, even endanger the patient's health. As a result, monitoring the SpO₂ is indispensable in the clinical rescuing. The traditional method to measure SpO₂ is to analyze the sample of the patient's blood, so can get the partial pressure of oxygen and calculate the SpO₂ by use the blood-gas analyzer. This method is inconvenient and can not be used to monitor continuously. For the purpose of measuring the SpO₂ more easily and accurately, Creative developed the Fingertip Oximeter. The device can measure the pulse rate simultaneously.

The Fingertip Oximeter is compact, convenient to use and carry and with low power consumption. You just need to put the fingertip into the sensor of the device, the SpO₂ value will appear on the screen immediately.

2.1 Features

- Large true color OLED display of SpO₂, PR Pulse Bar, PI & Plethysmogram;
- Innovative 4 directions display;
- Automatic power on/off;
- Audible & visible over-limit indication;
- Shift parameter display between PR and PI;
- 2AAA alkaline batteries with low power consumption;
- Low battery voltage indication.

2.2 Major Applications and Scope

The Fingertip Oximeter can detect SpO₂ and pulse rate through patient's finger, and indicate the pulse intensity by the bar-display. This device is applicable to home, hospital (including internal medicine, surgery, anesthesia, pediatrics, emergency

room etc.), oxygen bar, the community medical center, alpine area and it also can be used before or after sports, and the like.

This device is not appropriate to be used for continuous monitoring.

2.3 Principle of Measurement

The measurement of pulse oximeter is that it uses a multi-functional oxyhemoglobinometer to transmit some narrow spectrum light bands through blood samples and to measure attenuation of spectrum with different wavelengths according to the characteristic that RHb, O₂Hb, Met Hb and COHb absorb the light of different wavelength, thereby determining O₂Hb saturation of different fractions. O₂Hb saturation is called "fractional" O₂Hb saturation.

Fractional O₂Hb saturation = $[\text{O}_2\text{Hb} / (\text{RHb} + \text{O}_2\text{Hb} + \text{Met Hb} + \text{COHb})] \times 100$

Oppositely, pulse oxygen oximeter measures functional O₂Hb saturation:

Functional O₂Hb saturation = $[\text{O}_2\text{Hb} / (\text{RHb} + \text{O}_2\text{Hb})] \times 100$

Present SpO₂ oximeter transmits light of two wavelengths only, red light (wavelength 660 nm) and infrared (wavelength 940 nm), to differentiate HbO₂ from HbR. One side of the sensor contains two LEDs, and the other side contains a photoelectric detector. SpO₂ oximeter measures HbO₂ saturation in the blood by the light plethysmograph when the pulse beats. The result is quite precise when HbO₂ saturation is over 70% ~ 95%.

2.4 Caution

- A. The finger should be placed properly (see the attached illustration of this manual), or else it may cause inaccurate measurement.
- B. The SpO₂ sensor and photoelectric receiving tube should be arranged in a way with the subject's arteriole in a position there between.
- C. The SpO₂ sensor should not be used at a location or limb tied with arterial canal or blood pressure cuff or receiving intravenous injection.
- D. Do not fix the SpO₂ sensor with adhesive or else it may result in venous pulsation and inaccurate measure of SpO₂.
- E. Make sure the optical path is free from any optical obstacles like rubberized fabric.
- F. Excessive ambient light may affect the measuring result. It includes fluorescent lamp, dual ruby light, infrared heater, direct sunlight and etc.
- G. Strenuous action of the subject or extreme electrosurgical interference may also affect the accuracy.
- H. Testee can not use enamel or other makeup on the finger.

3. Battery Installation

1. Press the lock button of the battery cover, meanwhile, pull the cover back and take it out.

2. Refer to Figure 2, insert two AAA size batteries into the battery compartment properly.

3. Replace the cover. Please make sure that the batteries are correctly installed since the incorrect installation may cause the device inoperable.

4. Operation

4.1 Start Measurement

1. Open the clip as shown in Figure 3.
2. Put finger into the rubber cushions of the clip (make sure the finger is in the correct position), and then clip the finger.
3. The device will power on automatically in 2 seconds, and start to display software version number.
4. Next enter into data display screen (as shown in Figure 4). The user can read the values and view the waveform from display screen.


Figure 1


Figure 2 Battery Installation

Figure 3: Put finger into the Oximeter

Figure 4 A1


Figure 4 A2

Figure 4 B1


Figure 4 B2

Figure 4 C1

Figure 4 C2

Figure 4 D1

Figure 4 D2

Screen Description:

"%SpO2": The title of SpO2; "99": SpO2 value, unit:%;

"PR": The title of Pulse Rate; "65": Pulse Rate value, unit: bpm (beat per minute);

"❤️": Pulse beat icon;

"▬": Pulse bar-graph;

"PI%": The title of Perfusion Index; "1.4": Perfusion Index value, unit: %;

"🔋": Battery power indicator

5. Change display direction

Four directions display alternately. Short time press "Display Key" to flip the screen 90° each time in a cyclical manner as shown in Figure 4. When the screen displays towards the left side, the plethysmogram will be viewed.

6. Shift parameter display between PR and PI during measurement
Long time press the "Display Key", shift the parameter display between PR and PI. But when the PR is shifted to PI display and no button operation is performed after 20 seconds, the PI will change to PR display automatically.

5. Technical Specifications

A. SpO₂ measurement:

dual-wavelength LED sensor with wavelength:

Red light: 663 nm, Infrared light: 890 nm.

Maximal average optical output power: $\leq 1.5\text{mW}$

Measuring range: 35%~100%

Measuring accuracy:

Arms $\leq 3\%$ for SpO₂ range from 70% to 100%

Note: Arms is defined as root-mean-square value of deviation according to ISO 9919.

SpO₂ low over-limit: 90%

B. Pulse Rate measurement:

Measuring range: 30bpm~240bpm

Measuring accuracy: $\pm 2\text{bpm}$ or $\pm 2\%$ (whichever is greater)

Pulse Rate over-limit: high over-limit: 120bpm; low over-limit: 50bpm

C. Perfusion Index (PI) Display

Range: 0.2%~20%

D. Audible & visual over-limit indication

When measuring, if SpO₂ value or pulse rate value exceeds the limit, the device will beep automatically and the value which exceeds its limit will flash on the screen. The Oximeter will shut down automatically in 8 seconds with no signal.

E. Display: Color OLED Display

F. Power supply requirement:

2 x LR03 (AAA) alkaline batteries

Working voltage: 2.2V~3.3VDC

Operating current: $\leq 40\text{mA}$

G. Environment requirement

Operating Temperature: 5 ~40°C

Operating Humidity: 30~80%

Atmospheric pressure: 70~106kPa

H. The performance under low perfusion condition

The accuracy of SpO₂ and PR measurement still meets the specification described above when the pulse modulation amplitude is as low as 0.6%.

I. Resistance to interference of surrounding light:

The difference between the SpO₂ value measured in the condition of indoor natural light and that of darkroom is less than $\pm 1\%$.

J. Dimensions: 60 mm (L) × 33 mm (W) × 30 mm (H)

Net Weight: 35g (including battery)

K. Classification:

The type of protection against electric shock: Internally powered equipment.

The degree of protection against electric shock: Type BF applied parts.

The degree of protection against harmful ingress of liquids: Ordinary equipment without protection against ingress of water.

Electro-Magnetic Compatibility: Group I, Class B

6. Accessories

- A. A lanyard
- B. Two batteries
- C. A pouch
- D. A User Manual
- E. Quality Certificate

Note: The accessories are subject to change. Detailed items and quantity see the Packing List.

7. Repair and Maintenance

- A. Please change the batteries when the low-voltage indicator lightens.
- B. Please clean the surface of the device before using. Wipe the device with alcohol first, and then let it dry in air or wipe it dry.
- C. Please take out the batteries if the oximeter will not be used for a long time.
- D. The best storage environment of the device is -20°C to 55°C ambient temperature and 10% to 100% relative humidity.

High-pressure sterilization cannot be used on the device.

Do not immerse the device in liquid.

It is recommended that the device should be kept in a dry environment.

Humidity may reduce the useful life of the device, or even damage it.

Disposal: *The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.*

For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.

8. Troubleshooting

Trouble	Possible Reason	Solution
The SpO2 and Pulse Rate can not be displayed normally	<ol style="list-style-type: none">1. The finger is not properly positioned.2. The patient's SpO2 is too low to be detected.	<ol style="list-style-type: none">1. Place the finger properly and try again.2. Try again; Go to a hospital for a diagnosis if you are sure the device works all right.
The SpO2 and Pulse Rate display instable	<ol style="list-style-type: none">1. The finger is not placed inside enough.2. The finger is shaking or the patient is moving.	<ol style="list-style-type: none">1. Place the finger properly and try again.2. Let the patient keep calm.
The device can not turn on	<ol style="list-style-type: none">1. The batteries are drained or almost drained2. The batteries are not inserted properly.3. The device's malfunction.	<ol style="list-style-type: none">1. Change batteries.2. Reinstall batteries.3. Please contact the local service center.
The indicator light is off suddenly	<ol style="list-style-type: none">1. The device will power off automatically when it gets no signal for 8 seconds.	<ol style="list-style-type: none">1. Normal.2. Change batteries.

GIMA WARRANTY CONDITIONS

Congratulations for purchasing a GIMA product.

This product meets high qualitative standards both as regards the material and the production. The warranty is valid for 12 months from the date of supply of GIMA. During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel traveling expenses and packaging not included. All components subject to wear are not included in the warranty. The repair or replacement performed during the warranty period shall not extend the warranty. The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use. GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc. The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed.

The defected products must be returned only to the dealer the product was purchased from. Products sent to GIMA will be rejected.

Table 1: Guidance and manufacturer's declaration - electromagnetic emissions

The device is intended for use in the electromagnetic environment specified below. The user of this device should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment
RF emissions CISPR 11	Group 1	This device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	This device is suitable for use in all establishments. It uses internal power supply and has no connection with power supply network.
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	N/A	

Table 2: Guidance and manufacturer's declaration - electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The user of this device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	N/A	N/A

Surge IEC 61000-4-5	±1 kV line (s) to line(s) ±2 kV line(s) to earth	N/A	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	N/A	N/A
Power frequency (50/60Hz) magnetic field IEC61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: UT is the a.c. mains voltage prior to application of the test level.

Table 3: Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The user of this device should assure that it is used in such an environment.

NOTE 1 : At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 : These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

A: Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, and electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating this device.

B: Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Immunity test	IEC 60601 test level	Compliance Level	Electromagnetic environment
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms da 150 kHz 80 MHz</p> <p>3 V/m da 80 MHz 2,5 GHz</p>	<p>N/A</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3 \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ,a should be less than the compliance level in each frequency range .b Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Table 4: Recommended separation distances between portable and mobile RF communications equipment and the device

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2,5 GHz
0,01	N/A	0,12	0,23
0,1	N/A	0,38	0,73
1	N/A	1,2	2,3
10	N/A	3,8	7,3
100	N/A	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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