: HEINE GAMMA 3.1®, 3.2®, 3.3®, C3®

HEINE Optotechnik GmbH & Co. KG
Kientalstr. 7 · 82211 Herrsching · Germany
Tel. +49 (0) 81 52 / 38 - 0
Fax +49 (0) 81 52 / 38 - 202
E-Mail: info@heine.com · www.heine.com
med 0711 1/10.11

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Instructions for use

Before using the HEINE GAMMA® stethoscope, read these instructions carefully and keep them in a safe place for future reference.

Range of applications

HEINE GAMMA® stethoscopes are used for auscultation of the activity of human organs and in connection with non-invasive blood pressure measurements.

Areas of application / indications

Stethoscopes are used by physicians and trained personnel for auscultation of the circulatory system, respiratory system as well as the gastrointestinal system in order to interpret the respective sounds. During blood pressure measurement a stethoscope is used for the detection of Korotkow sounds in order to determine the systolic and diastolic blood pressure of adults and children. The stethoscope consists of chest piece, stethoscope tube and binaurals with separate ear olives. HEINE GAMMA® Stethoscopes are available with various chest pieces and tubes. Please see the following description of the different models:

GAMMA 3.1 Pulse Stethoscope

Especially useful for blood pressure measurement. The chest piece fits comfortably under the cuff of the sphygmomanometer.

GAMMA 3.2, 3.3 Acoustic Stethoscope

For the auscultation of high and low frequency tones. Switchable.

GAMMA C3 Cardio Stethoscope

For the auscultation of high and low frequency tones. Switchable. Special features: dual membrane, bell and double channel tube. The dual membrane permits auscultation of high and low frequency sounds and the bell facilitates evaluation of lower frequencies. (Cardiology)

Preparation for use

For the best auditory quality, please pay attention to the following:

To prevent injury

The ear pieces should be firmly attached to the binaurals.

Adjustment of the tension of the binaurals

Adjust to fit correctly.

Adjustments:

• Binaurals are too tight for comfort

To **reduce** the spring tension, firmly grasp each section of the binaurals between thumb and fore finger and widen carefully.

• Binaurals are too loose

To **increase** the spring tension, grasp both sides of spring with one hand and squeeze them together.

This process can be repeated as often as needed.

Switching the double chest piece

Depending on the needed frequency range, select either the membrane or the open bell side. In order to change hold the connecting tube in one hand and the chest piece in the other, now turn it 180° until there is a noticeable click.

If the connecting tube points towards the membrane

Auskultation with the bell

If the connecting tube points towards the bell

Auskultation with the membrane

Important hints

General information

- Measurements should only be interpreted by qualified persons.
- Avoid unnecessary impact on and vibration of the instrument.
- Avoid damaging membrane and soft components with sharp or pointed objects.
- Do not insert the ear pieces too far into ear canal, check for cleanliness regularly and clean as needed.
- Do not sterilize thermally.
- Do not dip in cleaning fluid and ensure that no fluids enter the tubing system.
- Never use a stethoscope without ear pieces.
- Never use a stethoscope without a membrane.
- Use only HEINE original spare parts.

Storage and transport conditions

- Store the instrument at a temperature between -10°C and 40°C
- Humidity air humidity below 85%
- Do not expose the instrument, in particular the soft componets, to direct ultra-violet radiation

Side effects

No side effects are known.

Maintenance and care

For cleaning we recommend a clean, damp cloth with a mild soap solution. If needed, use a **non-aggressive** disinfectant for the membrane. To extend the flexibility of the stethoscope tube we recommend the use of a vinyl lubricant.

All HEINE GAMMA® stethoscopes are latex-free.

Symbols used:

Please follow the instructions attached

The Product complies with the requirements of the Council Directive 93/42/EEC of 14 June 1993 (Medical Device Directive)

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