Abbott

NycoCard™

HbA1c

EN	5-8	ES	25-28	NO	45-48
CS	9-12	FR	29-32	PL	49-52
DA	13-16	HU	33-36	PT	52-56
DE	17-20	IT	37-40	RO	57-60
EL	21-24	NL	41-44	RU	61-64

REF 1116813

1116927 Rev. A 2019/02

SYMBOLS/ SYMBOLY/ SYMBOLER/ SYMBOLE/ ΣΥΜΒΟΛΑ/ SÍMBOLOS/ SYMBOLES/
JELMAGYARÁZAT/ SIMBOLI/ SYMBOLEN/ SYMBOLER/ SYMBOLE/ SÍMBOLO/ SIMBOLURI/
СИМВОЛЫ

Conformity to the European directive 98/79/EC
on *in vitro* diagnostic medical devices
Splnění evropské směrnice 98/79/ES o
zdravotnických prostředcích pro diagnostiku
in vitro
Overensstemmelse med EU-direktivet 98/79/EC
om *in vitro* diagnostisk udstyr til medicinsk brug
Konformität mit der europäischen Richtlinie
98/79/EC über medizinische Geräte zur *In-vitro*-
Diagnostik
Συμμορφώνεται με την ευρωπαϊκή οδηγία
98/79/ΕΚ περί ιατροτεχνολογικών βοηθημάτων
που χρησιμοποιούνται στη διάγνωση *in vitro*
Conformidad con la Directiva europea 98/79/CE
sobre dispositivos médicos para el diagnóstico
in vitro
Conformité avec la directive européenne 98/79/
CE sur les appareils médicaux de diagnostic
in vitro
Az Európai Parlament és a Tanács 98/79/EK
irányelve az *in vitro* diagnosztikai orvostechnikai
eszközökről
Conformità con la direttiva europea 98/79/CE sui
dispositivi medicali per la diagnostica *in vitro*
Voldoet aan de Europese richtlijn 98/79/EG
betreffende medische hulpmiddelen voor *in-
vitro*diagnostiek
Overensstemmelse med EU-direktiv 98/79/EC for
in vitro diagnostisk medisinsk udstyr
Zgodność z dyrektywą 98/79/WE w sprawie
wyrobów medycznych używanych do diagnostyki
in vitro
Conformidade com a Diretiva Europeia 98/79/CE
de dispositivos médicos de diagnóstico *in vitro*
Conformitate cu directiva europeană 98/79/CE
privind dispozitivele medicale pentru diagnostic
in vitro
Знак соответствия европейской директиве
98/79/EC по медицинским устройствам,
предназначенным для диагностики *in vitro*

Lot number
Číslo šarže
Lot-nummer
Chargen-Nr.
Αριθμός παρτίδας
Número de lote
Numéro de lot
Gyártási tételszám
Numero di lotto
Partijnummer
Lotnummer
Numer partii
Número de lote
Număr de lot
номер партии

Caution
Varování
Forsigtig
Achtung
Προσοχή
Advertencia
Faire attention
Figyelmeztetés
Attenzione
Opgelet
Advarsel
Ostrzeżenie
Atenção
Avertizare
Предупреждение

In Vitro diagnostic medical device
Diagnostické lékařské zařízení pro testy *in vitro*
In vitro diagnostisk medicinsk udstyr
Medizinisches Gerät zur *In-vitro*-Diagnostik
ιατροτεχνολογικό βοήθημα που χρησιμοποιείται
στη διάγνωση *in vitro*
Dispositivo médico para el diagnóstico *in vitro*
Dispositif medical de diagnostic *in vitro*
In vitro diagnosztikai orvosi eszköz
Dispositivo medico per la diagnostica *in vitro*
In-vitro-diagnostisch medisch apparaat
Medisinsk udstyr til *in vitro* diagnostikk
Wyrob medyczny do diagnostyki *in vitro*
Dispositivo médico de diagnóstico *in vitro*
Dispozitiv medical pentru diagnostic *in vitro*
Медицинское устройство для диагностики *in vitro*

Catalogue number
Katalogové číslo
Katalognummer
Bestellnummer
Αριθμός καταλόγου
Número de catálogo
Numéro de référence du catalogue
Katalógusszám
Numero di catalogo
Catalogusnummer
Katalognummer
Numer katalogowy
Número do catálogo
Număr de catalog
Номер по каталогу

2°C 8°C

Storage temperature 2-8°C (36-46°F)
 Skladovací teplota 2-8 °C (36-46 °F)
 Opbevaringstemperatur 2-8°C (36-46°F)
 Lagertemperatur 2-8°C (36-46°F)
 Θερμοκρασία φύλαξης 2-8°C (36-46°F)
 Temperatura de almacenamiento 2-8°C (36-46°F)
 Température de conservation 2-8°C (36-46°F)
 Tárolási hőmérséklet 2-8 °C (36-46 °F)
 Temperatura di conservazione 2-8°C (36-46°F)
 Bewaarstemperatuur 2-8°C (36-46°F)
 Lagringstemperatur 2-8°C (36-46°F)
 Temperatura przechowywania 2-8°C (36-46°F)
 Temperatura de conservação 2-8°C (36-46°F)
 Temperatură de depozitare 2-8°C (36-46°F)
 Температура хранения 2-8 °C (36-46 °F)

Expiry date (year-month-day)
 Datum exspirace (rok-měsíc-den)
 Udløbsdato (år-måned-dag)
 Verfallsdatum (Jahr-Monat-Tag)
 Ημερομηνία λήξης (έτος-μήνας-ημέρα)
 Fecha de caducidad (año-mes-día)
 Date de péremption (année-mois-jour)
 Lejárati dátum (év-hónap-nap)
 Data di scadenza (anno-mese-giorno)
 Uiterste gebruiksdatum (jaar-maand-dag)
 Utløpsdato (år-måned-dag)
 Data przydatności (rok-miesiąc-dzień)
 Prazo de validade (ano-mês-dia)
 Data expirării (an-lună-ziuă)
 Срок годности (год-месяц-день)

Consult the instruction for use
 Viz návod k použití.
 Se brugervejledningen
 Bitte die Gebrauchsanleitung beachten
 Συμβουλευτείτε τις οδηγίες χρήσης
 Consulte las instrucciones de uso
 Consulter le mode d'emploi
 További információt a használati útmutatóban talál
 Consultare le istruzioni per l'uso
 Raadpleeg de gebruiksaanwijzing
 Se bruksanvisningen
 Sprawdzić w instrukcji użytkowania
 Consulte as instruções de utilização
 Consultați instrucțiunile de utilizare
 См. инструкцию по использованию

Contents sufficient for 24 tests
 Obsah postačuje pro 24 testů
 Indholdet er tilstrækkeligt til 24 test
 Inhalt ausreichend für 24 Tests
 Περιεχόμενα επαρκή για 24 δοκιμές
 Contenido suficiente para 24 pruebas
 Contenu suffisant pour 24 dosages
 24 vizsgálatra elegendő tartalom
 Contenuto sufficiente per 24 test
 Inhoud voldoende voor 24 testen
 Innholdet rekker til 24 tester
 Zawiera ilość wystarczającą na 24 badań
 Conteúdo suficiente para 24 testes
 Conține o cantitate suficientă pentru 24 teste
 Содержимого достаточно для проведения 24 анализов

Manufacturer
 Výrobce
 Producent
 Hersteller
 Κατασκευαστής
 Fabricante
 Fabricant
 Gyártó
 Produttore
 Fabrikant
 Produzent
 Producent
 Fabricante
 Producător
 Изготовитель

Date of Manufacture (year-month-day)
 Datum výroby (rok-měsíc-den)
 Produktionsdato (år-måned-dag)
 Herstellungsdatum (Jahr-Monat-Tag)
 Ημερομηνία κατασκευής (έτος-μήνας-ημέρα)
 Fecha de fabricación (año-mes-día)
 Date de fabrication (année-mois-jour)
 Gyártás dátuma (év-hónap-nap)
 Data di fabbricazione (anno-mese-giorno)
 Productiedatum (jaar-maand-dag)
 Produksjonsdato (år-måned-dag)
 Data produkcji (rok-miesiąc-dzień)
 Data de fabrico (ano-mês-dia)
 Data fabricației (an-lună-ziuă)
 Дата изготовления (год-месяц-день)

For professional near-patient testing and laboratory use. For use with the NycoCard™ Reader II.

⚠ The NycoCard Reader II is calibrated for use with the NycoCard Tests only.

PRODUCT DESCRIPTION

Intended use

NycoCard™ HbA1c is an *in vitro* diagnostic medical device for quantitative determination of glycated hemoglobin (HbA1c) in human whole blood. For monitoring of long-term metabolic control in persons with diabetes mellitus^{1,2,3}.

Principle of the assay

NycoCard HbA1c is a boronate affinity assay. The kit contains test devices with a porous membrane filter, test tubes prefilled with reagent and a washing solution. The reagent contains agents that lyse erythrocytes and precipitate hemoglobin specifically, as well as a blue boronic acid conjugate that binds *cis*-diols of glycated hemoglobin. When blood is added to the reagent, the erythrocytes immediately lyse. All hemoglobin precipitates. The boronic acid conjugate binds to the *cis*-diol configuration of glycated hemoglobin. An aliquot of the reaction mixture is added to the test device, and all the precipitated hemoglobin, conjugate-bound and unbound, remains on top of the filter. Any excess of coloured conjugate is removed with the washing solution. The precipitate is evaluated by measuring the blue (glycated hemoglobin) and the red (total hemoglobin) colour intensity with the NycoCard Reader II, the ratio between them being proportional to the percentage of HbA1c in the sample.

Kit contents, 24 Tests

TD/Test Device	1 x 24 units
Plastic device containing a membrane filter.	
R1/Reagent	1 x 24 x 0.2 mL
Glycinamide buffer containing dye-bound boronic acid and detergents.	
R2/Washing Solution	1 x 2.0 mL
Morpholine buffered NaCl solution and detergents.	

Material required but not provided with the kit

- Capillary tube or pipette (5 µL) for sample collection
- Capillary tube holder
- Pipette (25 µL) and pipette tips
- NycoCard Reader II

Warnings and precautions

- For *in vitro* diagnostic use.
- Do not interchange components from different test kit lots.
- Do not use test kits after the expiry date.
- R1/reagent and R2/washing solution contain sodium azide, which is a toxic agent.
- R2/washing solution contains morpholine, which is a toxic agent.
- Use gloves.

TEST CHARACTERISTICS

Analytical specificity

NycoCard HbA1c measures the total glycated hemoglobin and reports the HbA1c value.

The following hemoglobin (Hb) variants have been analysed and found not to affect the NycoCard HbA1c Test result: HbAC, HbAE, HbAF, HbAJ, HbAS and Hb Yamagata^{4,5}. Carbamylated and pre-glycated hemoglobin does not affect the test.

Standardization

NycoCard HbA1c is traceable to the IFCC Reference Method for Measurement of HbA1c⁶.

Reference range

For methods reporting DCCT traceable values, the upper limit of non-diabetic reference range is 6.4% HbA1c^{11,12}.

Measuring range

Hemoglobin measuring range: 7-18 g/dL.

The NycoCard Reader II displays the HbA1c concentration in mmol/mol, percentage (%) or in combination, depending on the serial number of the instrument.

SN $\geq$ 78450: Consult the NycoCard Reader II User Manual for instructions on how to change the HbA1c measuring unit.

HbA1c measuring units

Two different measuring units are in use for reporting HbA1c test results⁷⁻¹⁰:

- mmol/mol – the HbA1c values are aligned to the IFCC Reference Method for Measurements of HbA1c.
- percentage (%) – the HbA1c values are aligned to the assay used in the DCCT study, also known as NGSP-HbA1c.

A linear relationship has been established between the two units⁹:

$$\text{DCCT-HbA1c (\%)} = 0.092 \times \text{IFCC-HbA1c (mmol/mol)} + 2.15$$

NycoCard™ READER II	HbA1c (%)		HbA1c (mmol/mol)	
	Measuring range	Unit	Measuring range	Unit
SN < 59681	3-18	0.1	---	---
SN 59681-78449	4-15	0.1	---	---
SN $\geq$ 78450	4-15	0.1	20-140	1.0

Precision

In professional use, a coefficient of variation (CV) < 5% is usually obtained.

Limitations of the test

- Elevated amounts of glucose, bilirubin, lipids and fructosamine were added to blood samples with normal and elevated HbA1c values. No interference was obtained.
- Hemolysed samples with plasma Hb >3 g/100 mL will interfere with the test system.

STORAGE

Unopened kits

The expiry date of the kit applies to storage at 2-8°C in the original container. Avoid direct sunlight or exposure to temperatures above 25°C. Do not freeze.

Opened kits

R1/Reagent must be stored dark at 2-8°C. Equilibrate the R1/reagent to room temperature (20-25°C) before use. Equilibration to room temperature can be achieved by holding the tube in a closed hand for 30 seconds. The reagent tube can be stored for maximum 6 hours at room temperature before use. Avoid direct sunlight.

TD/Test device can be stored at room temperature (15-25°C). Store the test devices in the original bag and avoid humidity below 20% and above 70%. The TD/test device should have room temperature when used.

R2/Washing solution can be stored at room temperature (15-25°C). The R2/washing solution should have room temperature when used.

Sample material

Blood samples with anticoagulant (EDTA, heparin and NaF) can be stored up to 10 days at 2-8°C before analysis. Do not freeze.

TEST PROCEDURE

Important procedural notes!

- Do not interchange components from different kits or kit lots.
- Bring the R1/reagent to room temperature (20-25°C) before use.
- Do not touch the membrane with the pipette tip.
- Change the pipette tip between each pipetting step.
- The procedure steps should be performed successively without any interruptions.

Sample material

Capillary and venous whole blood. Avoid measuring hemolysed samples (see “Limitations of the test”). Mix venous samples well by inverting the tube 8-10 times before collecting a sample.

Test procedure (illustrations on page 65)

1 Prepare sample

Add 5 µL whole blood to the test tube with R1/reagent. Mix well. Leave the tube for minimum 2 minutes, maximum 3 minutes. Use a timer.

Note! Equilibrate the R1/reagent to room temperature (20-25°C) before use.

2 Apply sample

Remix to obtain a homogenous suspension. Apply 25 µL of the mixture to a TD/test device. Hold the pipette approx. 0.5 cm above the test well and empty the pipette quickly in the middle of the test well. Allow the mixture to soak completely into the membrane. Wait for minimum 10 seconds.

Note! Avoid air bubbles.

3 Apply R2/washing solution

Apply 25 µL R2/washing solution to the TD/test device. Allow the reagent to soak completely into the membrane. Wait for minimum 10 seconds.

Note! Avoid air bubbles.

4 Read the test result

Read the test result within 5 minutes using the NycoCard Reader II. Follow the NycoCard Reader II Instruction Manual.

INTERPRETATION OF RESULTS

Interpret the NycoCard HbA1c Test results with careful consideration to the patient's medical history, clinical examinations and other laboratory results. If the test result is questionable or if clinical signs and symptoms

appear inconsistent with the test result, retest the sample or confirm the result using another method. Analyse control materials frequently to verify the performance of the NycoCard Reader II Test System.

QUALITY CONTROL

A quality control material with NycoCard HbA1c specific target values should be used to confirm the efficacy of the reagents and the correct performance of the test. The measured value should be within the acceptance limits stated for the control material. Do not analyse patient samples when control values are outside range.

Do not use lyophilised control materials with the NycoCard HbA1c Assay.

Keeping a permanent record of all quality control results is recommended.

TROUBLE SHOOTING

Possible cause	Corrective action
1 – Unexpected low test results	
R1/reagent was too cold when used.	Equilibrate R1 to 20-25°C before use.
R1 has been exposed to light or room temperature > 8 hours.	Discard the test kit. Retest the sample with a new kit.
A hemolysed sample material was used.	Collect a new sample.
2 – Unexpected high test results	
Too little R2 was added to the TD/test device.	Make sure that 25 µL of R2 is applied to the test device.
There are white spots on the TD membrane.	Apply the sample quickly, in the middle of the test well. Avoid air bubbles.
3 – The R2/washing solution does not soak into the TD/test device membrane	
A hemolysed sample material was used.	Collect a new sample.
4 – The NycoCard™ READER II displays “Hb conc. too low”	
The Hb concentration of the sample is below the measuring range.	Retest the sample. Fill 2 capillaries with sample (= 10 µL) and drop both capillaries into one tube of R1. The NycoCard Reader II will display the correct test result; no manual calculation is required.
5 – The NycoCard™ READER II displays “Reduce Hb conc.”	
The Hb concentration of the sample is above the measuring range.	Retest the sample. Mix R1 from two tubes (= 400 µL) and apply 1 capillary (5 µL) with sample. The NycoCard Reader II will display the correct test result; no manual calculation is required.

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1. Lenzi S *et al.* The Clinical Usefulness of Glycated Hemoglobin in Monitoring Diabetes Mellitus: A Long-Term Study Clin Chem 1987; 33:55-56.
2. The Diabetes Control and Complications Trial Research Group. The Effect of Intensive Treatment of Diabetes on the Development and Progression of Long-Term Complications in Insulin-Dependent Diabetes Mellitus. N Eng J Med 1993; 329:977-986.
3. Sacks BD *et al.* Guidelines and Recommendations for Laboratory Analysis in the Diagnosis and Management of Diabetes Mellitus. Clin Chem 2011; 57(6):e1-e47.
4. European Reference Laboratory for Glycohemoglobin. Manufacturer Check-Up Certification of NycoCard HbA1c, 1998.
5. Chae Lim Jung *et al.* Hemoglobin Yamagata: Hemoglobin Variant Detected by HbA1c Test. Korean J Lab Med 2009; 29:536-540.
6. Jeppsson JO *et al.* Approved IFCC Reference Method for the Measurement of HbA1c in Human Blood. Clin Chem Lab Med 2002; 40(1):78-89.
7. American Diabetes Association, European Association for the Study of Diabetes, International Federation of Clinical Chemistry and Laboratory Medicine. Consensus Statement on the Worldwide Standardization of the Hemoglobin A1C Measurement. Diabetes Care 2007; 30(9):2399-2400.
8. Nordin G *et al.* Recommendation for the term and measurement unit for "HbA1c". Clin Chem Lab Med 2007; 45:1081-1082.
9. Hoelzl W *et al.* IFCC reference system for the measurement of haemoglobin A1c in human blood and the national standardization schemes in the United States, Japan and Sweden: a method-comparison study. Clin Chem 2004; 50(1):166-174.
10. Weykamp C *et al.* The IFCC Reference Measurement System for HbA1c: a 6-year progress report. Clin Chem 2008; 54(2):240-248.
11. The International Expert Committee. International Expert Committee Report on the Role of the A1c Assay in the Diagnosis of Diabetes. Diabetes Care 2009; 32(7):1327-1334
12. American Diabetes Association. Diagnosis and Classification of Diabetes Mellitus. Diabetes Care 2011; 34 (Suppl 1):S62-S69.

1

2

3

4

Abbott

Abbott Diagnostics Technologies AS

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P.O. Box 6863 Rodeløkka

NO-0504 Oslo, Norway

www.abbott.com/poct

ISO 13485 certified company

© 2019 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

1116927 Rev. A 2019/02