∴ **HEINE® EN 50 UNPLUGGED**
HEINE® mPack UNPLUGGED

HEINE Optotechnik GmbH & Co. KG
Kientalstr. 7 · 82211 Herrsching · Germany
Tel. +49 (0) 81 52 / 38 - 0
Fax +49 (0) 81 52 / 38 - 202
E-Mail: info@heine.com · www.heine.com
med 6507 1/04.2013

HEINE® EN 50 UNPLUGGED HEINE® mPack UNPLUGGED

DEUTSCH 4

ENGLISH 8

FRANÇAIS 12

ESPAÑOL 16

ITALIANO 20

SVENSKA 24

NEDERLANDS 28

SUOMI 32

PORTUGUÊS 36

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HEINE® EN 50 UNPLUGGED HEINE® mPack UNPLUGGED

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Instructions for use


HEINE® EN 50 UNPLUGGED HEINE® mPack UNPLUGGED

 Read these instructions carefully before using the HEINE® EN50 UNPLUGGED or the HEINE® mPack UNPLUGGED and keep them in a safe place for future reference.

Product overview


HEINE® mPack UNPLUGGED:

Connector for the instrument or an UNPLUGGED power unit for mobile charging

Capacity indicator

Serial-No.


Charging contacts for charging via the HEINE® EN50 UNPLUGGED

HEINE® EN 50 UNPLUGGED:


Charging cradle II

Contact pins

Charging cradle I

HEINE UNPLUGGED power unit (plug-in transformer):

Mains adapters

Mains adapters:

EU

USA / Japan

UK

Australia

Safety

Intended use: The HEINE® mPack UNPLUGGED is intended exclusively for supplying power to HEINE instruments on headbands specially designed for the purpose. There are two chargers available for charging the Li-polymer battery inside the HEINE® mPack UNPLUGGED:

- The HEINE® EN50 UNPLUGGED, with two charging cradles, is intended for wall mounting.
- The HEINE UNPLUGGED plug-in transformer is intended as a charger for mobile operation.

Safety information: Never open the housing – except when changing the battery of the mPack UNPLUGGED. Do not make any alterations to the units. All repairs to the unit or the connecting leads should only be carried out by qualified personnel.

Safety in use:

- Use only an original HEINE plug-in transformer to connect the unit to the local mains supply.
- Make sure that the local mains voltage is the same as the nominal voltage shown on the identification plate of the HEINE plug-in transformer.
- Unplug the mains connector if you will not be using the unit for a long period if you are cleaning it or if a fault occurs. It should not be assumed there is no voltage present just because the Power LED is not lit.
- Do not use the plug-in transformer if it shows visible signs of damage.
- Take care that the battery housing does not get damaged. Damaged batteries must not be used.
- Pay attention to the mounting instructions when wall-mounting the HEINE® EN50 UNPLUGGED charging station. Take care to mount the charging station outside the patient environment.
- Ensure that the units are beyond the reach of small children.
- Operation in areas where an explosion hazard is present is not permitted.
- Protect your unit from very high humidity, direct sunlight and high or extremely low temperatures.
- Do not install the unit near radiators or other heat sources.
- We can only guarantee the safety and performance of this unit if it is fitted with original HEINE accessories and spare parts.
- The unit is not waterproof! Take care that no liquids get into the unit.
- Condensation can form when moving from a cold to a warm environment. Wait until the unit has warmed up and is absolutely dry before using it.
- The capacity of the mPack UNPLUGGED Battery can reduce with increasing use. This is a normal procedure with time and inevitable. Careful use of the battery can reduce this effect. We recommend that the instrument is use should be switched of after use and the battery recharged regardless of its state of charge.

Symbols

Observe the safety information.

The battery in the mPack UNPLUGGED should be recycled in accordance with the applicable legislation.

In accordance with local regulations this product should be disposed of as an electronic device separately.

Locking/unlocking of the mPack UNPLUGGED to the headband.

Power supply connector (9VDC) for charging the mPack UNPLUGGED or for direct operation.

Instrument connector

Protection from electric shock (Type BF equipment)

Double insulated equipment of Protection Class II

DC

The CE marking indicates that the product conforms to the provisions of Directive 93/42/EEG for medical products.


WITH RESPECT TO ELECTRIC SHOCK FIRE, AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH UL60601-1/CAN/CSA C22.2 No. 601.1



Setting up

Fitting the HEINE® mPack UNPLUGGED to the instrument headband:

First remove the rear segment cover of the headband. To do this, hold the headband in the region of the rear segment cover with your left hand and pull the segment cover to the right with your right thumb.

Then fit the mPack UNPLUGGED to the rear part of the headband by supporting the headband with your left hand and placing the mPack UNPLUGGED on the headband so that the LEDs are arranged to the left. The mPack UNPLUGGED will engage securely with a little light pressure from your left thumb. Connect the brightness control to the mPack UNPLUGGED.

The batteries must be charged before using the mPack UNPLUGGED for the first time. You can charge the mPack UNPLUGGED either from the UNPLUGGED plug-in transformer or from the EN50 UNPLUGGED charging station.


When supplied, new batteries are only partially-charged. The mPack UNPLUGGED must be charged before the instrument is used for the first time. The LED indicator corresponds to the actual battery capacity only after a full charge and discharge cycle. (Also applies when changing the battery.)

Note:

- To enable the LED display of the mPack UNPLUGGED to show the true battery capacity, the unit must first be completely discharged and then recharged. Switch on the instrument connected to the mPack UNPLUGGED until the battery is completely discharged, then recharge the mPack UNPLUGGED until charging is completed. (Also applies when changing the battery.)

Wall mounting the HEINE® EN50 UNPLUGGED:

Use the four provided screws to fix the mounting bracket to the wall. Place the EN50 UNPLUGGED on top of the mounting bracket and fasten with the two top screws.

Operation

Charging via the HEINE® EN50 UNPLUGGED (wall-mounted charger):

Connect the plug-in transformer to a suitable mains outlet (pay attention to the mains supply voltage). The green LED on the plug-in transformer will light. The mPack UNPLUGGED can be charged in charging cradle I or charging cradle II. Using gentle pressure, insert the mPack UNPLUGGED into the charging cradle until it engages securely.

Charging cradle I is intended for parking the complete instrument and for charging the mPack UNPLUGGED. The instrument is automatically switched off as soon as the complete instrument is parked

in charging cradle I. The mPack UNPLUGGED will engage in the charging cradle with light pressure, thus ensuring safe charging. After approx. 1 second, the LEDs start to light in sequence. During charging, the LEDs blink.

(See the "LED display" chapter for further information.) When you remove the instrument from the charging station and switch it on with the brightness control, all five LEDs will light if the battery is fully charged. As the level of charge falls, the LEDs go out in turn. If the orange LED starts to flash, you should recharge the battery.

Note: You can park the instrument in the charging station at any time, even after short periods of operation. To increase the working life of the battery, recharging only occurs if the residual charge is less than 75%.

Charging by means of the UNPLUGGED power unit:

Please note that charging with the UNPLUGGED mains transformer is only possible when the brightness control is switched off.

Connect the plug-in transformer to a suitable mains outlet (pay attention to the mains supply voltage). The green LED on the plug-in transformer will light. Next connect the cable to the right or left connector on the mPack UNPLUGGED.

The charging process begins. All LEDs will go out as soon as the charging process is complete. Disconnect the charger. You can now switch the instrument on by means of the brightness control.

The following applies to either charging method:

- The charging process is dependent on the residual charge of the battery and requires a maximum of approximately 2 hours.
- The charging process switches off automatically when the battery is full, all LEDs go out.
- Overcharging of the battery is not possible.

Direct operation:

Direct operation is possible if there is no capacity in the mPack UNPLUGGED battery but you want to continue to work with the instrument. To do this, simply connect the UNPLUGGED plug-in transformer to the connector on the left or right hand side and switch the instrument on with the brightness control.

Note: The battery will not be charged at the same time during "direct operation".

LED display

Charging

- During charging the LED's operate sequentially, in use they remain switched on.
- If all LEDs are out, the charging process is complete and the battery is fully charged.
- mPack UNPLUGGED up to 2010 behave differently. Up to 2010 the LED are switched on and off in the same way. However, the number of LED's illuminated shows the state of charge of the battery.

Operating the mPack UNPLUGGED

- The number of LEDs illuminated indicates the level of charge in the battery. The full operating period is achieved when at least 4 LEDs are lit.
- If the orange LED is flashing, you should charge the mPack UNPLUGGED.
- If all the LEDs are out, the battery is fully discharged.

Servicing and cleaning

- Always disconnect from the mains before cleaning.
- It is possible to wipe with a damp cloth. If a cleanser is used, it must be approved for use on plastics.
- Disinfectant cleaning can be done using an alcohol-based surface disinfectant.
- Spray disinfection, soaking or disinfection in a machine is not permitted.
- Check the function of each unit regularly. Use the units only when they are functioning correctly.

- Do not damage the contact surfaces on the mPack UNPLUGGED or the contact pins in the EN50 UNPLUGGED charger with pointed or sharp-edged objects.
- Clean the contacts of EN 50 UNPLUGGED and mPack UNPLUGGED with a cloth moistened with alcohol.

Changing the battery


In the course of time, batteries gradually wear out and require longer charging times; this is normal. If you charge the batteries regularly and notice a shortening of the operating time or an extension of the charging time, it is probably time to replace them. The new Li-Pol battery can be inserted simply by opening the lid of the mPack UNPLUGGED:

- Isolate the mPack by removing it from the charger cradle or by unplugging the UNPLUGGED mains transformer.
- Unscrew the eight screws in the lid of the mPack UNPLUGGED with the screwdriver supplied and remove the lid, which contains the Li-Pol battery.
- Make sure that no metal objects such as screws fall into the unit.
- Protect the electronics from electrostatic discharge (ESD), by avoiding contact with electronic components.
- Before replacing the lid, remove the protective cap from the three-pole battery contacts.
- Re-attach the lid with the eight screws, making sure that no wires are damaged when doing so.

Note: The battery should be disposed of correctly and should be completely discharged before disposal.

Note:

- Do not short-circuit the Li-Pol battery. This can occur when a metal object such as a coin or paper clip connects the two battery contacts together. This can damage the battery.
- Do not expose to fire (Explosion hazard).

Product specification

Technical Data

HEINE® EN50 UNPLUGGED:

Mains voltage:	100–240 V~ / 50–60 Hz
Current:	700 mA
Output:	9 VDC; 3.0 A
Fuse:	Integral overload protection
Duration of use:	Continuous
Classification:	Type BF

HEINE® mPack UNPLUGGED:

Input:	9 VDC; 1,5 A
Output:	Variable, max. 6.5 V / 10 W
Fuse:	Integral overload protection
Weight:	approx. 190 g
Dimensions:	162 mm x 52 mm x 24 mm
Duration of use:	depends on bulb type and brightness, e.g. 5W bulb / HEINE OMEGA® 500, approx. 120 min. at maximum brightness)
Charging time:	approx. 2 hours

HEINE Li-Pol Rechargeable battery:

Cell type:	Li-polymer
Nominal voltage:	7,4 V
Capacity:	approx. 1850 mAh

Environmental conditions:

Operation:	+10 °C to +40 °C
Storage:	+5 °C to +45 °C
Transport:	-20 °C to +45 °C
Rel. Humidity:	10% to 75%
Atmospheric pressure:	700 hPa to 1060 hPa

Applicable Standards

ISO 7000: 2004	(Graphical symbols for use on equipment)
EN ISO 15223-1	(Grafische Symbole zur Kennzeichnung)
IEC 60601-1: 1996 EN 60601	(Medical Electrical Equipment – PART 1)
UL 60601-1: 2003	(General requirements for Safety for Medical Electrical Equipment)
IEC 60601-1-2	(Medical Electrical Equipment – COLLATERAL STANDARD: ELECTROMAGNETIC COMPATIBILITY – REQUIREMENTS AND TESTS)

What if...?

In the event of a fault, it is generally recommended that the mPack UNPLUGGED be connected briefly to the charger.

Fault description	LED display on mPack UNPLUGGED	Cause of problem	Remedy
Green Power LED on power unit does not light	-	- no mains supply voltage present - mains supply voltage outside permissible range	-
Instrument does not function	No LEDs on	Battery is discharged	Charge the mPack UNPLUGGED
		Electronics defect in the mPack UNPLUGGED or brightness controller	Refer to your service center. Check whether direct operation is still possible.
		Ambient temperature too high	Do not expose the mPack UNPLUGGED to ambient temperatures greater than 40 °C.
	LEDs on	Bulb in instrument is defective	Change the bulb
		Fault on brightness controller	- Turn the rotary knob on/off briefly. - Plug in the connecting cable again
mPack UNPLUGGED does not charge	fourth green LED does not light	remaining capacity in the battery is greater than 75%; frequent, brief recharging is to be avoided, to preserve the battery.	-
	even after a long charging time, only the orange LED is flashing	Due to a disturbance in the instrument, e.g. when changing the bulb, the LED display on the mPack UNPLUGGED no longer corresponds to the actual state of charge.	If the instrument is functioning, fully discharge and recharge the battery once; the LED display will then be resynchronized.
	- no mains voltage present - mains voltage does not correspond to the nominal mains voltage		-
Instrument does not function	LEDs blink	Defective bulb or bulb not correctly installed.	Change the bulb or instal corretly. Switch HC 50 off and then on again.

General Conditions of Warranty

Instead of the statutory warranty time period of 2 years, HEINE will grant a guarantee of 5 years from the date of the consignment of the goods ex works, concerning its equipment (excluding disposables, e.g. bulbs, single-use articles, and rechargeable batteries).

The guarantee covers irrefragable workmanship, on condition of the proper use of the equipment and the observation of the operating instructions. During the warranty- and guarantee time period, errors and deficiencies arising on the equipment will be rectified free of charge, in so far as such are evidenced by defective materials, processing and/or constructional errors. Should buyer complain of a material deficiency during the warranty time period, then the onus of proof is always to be on the orderer, that the product was defective already upon receipt of the goods. The statutory warranty and the guarantee do not apply to loss or damage caused by wear and tear, negligent use, the non-employment of original HEINE components and/or spares (in particular bulbs, as these have been especially developed for HEINE instruments in accordance with the following criteria: colour temperature, useful service life, safety, optical quality and performance. The statutory warranty and the guarantee do not apply to interventions by persons not authorised by HEINE or when the operating instructions are not observed by the customer. Any modification of a HEINE product with parts or additional parts which do not conform to the original HEINE specification will invalidate the warranty for the correct function of the product and further invalidate any warranty claims which result from such a change or modification. Further claims, in particular claims for replacement of loss or damage, which are experienced otherwise than directly on the HEINE product itself, are hereby excluded.

For **repairs**, please contact your supplier, who will return the goods to us or our authorized agent.

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