TotalShield® II Surgical Helmet

English

GENERAL INFORMATION INDICATIONS FOR USE

The *TotalShield II* Surgical Helmet or *TotalShield II* Surgical Helmet with LED Light is for use with the *TotalShield II* Zippered Surgical Toga or *TotalShield II* Surgical Hood as the *TotalShield II* Surgical Helmet System that is intended to be worn by surgical personnel to provide a barrier between the operating environment and the surgical personnel in order to protect against contamination and/or exposure of infectious body fluids and harmful microorganisms.

DESCRIPTION

The *TotalShield II* Surgical Helmets are components of *TotalShield II* Surgical Helmet System which includes the *TotalShield II* Zippered Surgical Toga and Surgical Hood. The *TotalShield II* Helmet draws external ambient air in via an integrated Battery-operated fan and circulates it throughout the helmet to aid in user comfort. The *TotalShield II* Surgical Helmet includes optional Battery powered LED lighting.

SPECIFICATIONS

REFERENCE NUMBER	00992000100	00992000200	00992000500	00992000600
Product Description	Surgical Helmet	Surgical Helmet with LED Light	Surgical Helmet w/ Lock	Surgical Helmet LED w/ Lock
Height	10.2 in. (25.9 cm)	10.2 in. (25.9 cm)	10.2 in. (25.9 cm)	10.2 in. (25.9 cm)
Width	9.1 in. (23.1 cm)	9.1 in. (23.1 cm)	9.1 in. (23.1 cm)	9.1 in. (23.1 cm)
Length	12.9 in. (32.8 cm)	12.9 in. (32.8 cm)	12.9 in. (32.8 cm)	12.9 in. (32.8 cm)
Beam Intensity	Not Applicable	≥ 40,000 lux (Working Distance 40 cm)	Not Applicable	≥ 40,000 lux (Working Distance 40 cm)
Fan Speed	Four Speeds 10.4 cfm, 12.8 cfm, 16.0 cfm, 19.5 cfm (nominal)			
Supply/Input	14.4 V, 3100 mAh and 6200 mAh			
Useful Life	Approximately 3 years or 900 uses			
Note	Transport or store in an environment that limits exposure to dust, moisture and temperature extremes.			
Warning: Inspect the device prior to each use. If you see the following, it is recommended to replace your device: • Damaged control buttons or LED adjustment slide. • Damaged LED lens or shroud • Helmet exhibits evidence of blood or bodily fluid soiling. (Follow cleaning instructions before using the helmet for future procedures.)				
Approval		Medical – General Medical Equipment as to Electric Shock, Fire and Mechanical Hazards Only, in Accordance with: Conforms To: AAMI STD ES60601-1:2005+A1 (2012) IEC STD 60601-1-6: 2010 Ed.3+A1 Certified To: CSA STD C22.2 No. 60601-1:2014 Ed.3		

ANSI/AAMI Std ES 60601-1 CLASSIFICATION

Type of protection against electric shock	Internally Powered Equipment
Degree of protection against electric shock	No applied parts
Mode of operation	Continuous operation
IPX0 Classification according to the degree of protection against ingress of water	IPX0 (Ordinary Equipment)

Affixed CE mark indicates compliance with the following European Directives:

- Directive 93/42/EEC – Medical Device
- Regulation 2016/425/EEC – Personal Protective Equipment (PPE);

Helmets were developed in accordance with EN 166:2001

PPE has been EC-type examined by the following Notified Body

BSI Group The Netherlands B.V.
Say Building, John M. Keynesplein 9,
1066 EP Amsterdam, The Netherlands

MARKING DEFINITIONS – Z166F

Code	Definition
Z	Manufacture Code, Zimmer
166	European Standard – Eye Protection Specification
F	Mechanical Strength, Low Energy Impact

WARNINGS

- Always wear a surgical head covering before using the *TotalShield II* Surgical Helmet.
- Use only Zimmer Biomet *TotalShield II* approved components and accessories unless otherwise specified.
- Not intended to be used as a respiratory protection system.
- The Surgical Helmet System is not made with natural rubber latex.
- Use only by qualified and trained medical professionals. The medical professional performing the procedure is responsible for the final decision on whether to use Surgical Helmet System.
- Follow normal institutional protocol for maintaining sterile technique at all times.
- DO NOT cut or puncture the helmet power cable.
- Always hold the connector at the end of the helmet power cable while connecting to the Holster with a Battery.
- DO NOT use the equipment in the presence of flammable anesthetics or gases.
- Users of a *TotalShield II* Surgical Helmet should be aware that precautions should be taken in regards to Electromagnetic Compatibility (EMC). The device should be installed and used according to the EMC information provided in the instructions for use. See EMC Guidance Tables included in this manual.

CONTRAINDICATIONS

None known.

RESIDUAL RISK

Failure to adhere to the manufactures instructions for use could lead to equipment malfunction and/or user injury.

ACCESSORY INFORMATION

Use only *TotalShield II* Surgical Helmet System accessories.

REFERENCE NUMBER	DESCRIPTION
00992020000	Six-Bay Smart Charger
00992020002	Two-Bay Smart Charger
00992010200	Rechargeable Li-Ion Battery
00992010400	Battery Holster
00992000306	Rear Comfort Support
0099200307	Crown Comfort Support


REFERENCE NUMBER	DESCRIPTION
00992000400	Helmet Pads
00992000401	Fan Cover Repair Kit
00992000402	Headband Repair Kit
00992000403	Hook and Loop Repair Kit
00992000404	Headband Repair Kit w/ Lock

CONTACT INFORMATION

To obtain accessories or additional information regarding your unit, contact your Zimmer distributor or write or phone. Visit www.ZimmerBiomet.com for the Declaration of Conformity.

MAIL	TELEPHONE
Zimmer Surgical, Inc. 200 West Ohio Avenue Dover, Ohio 44622 U.S.A.	1-330-343-8801 (local) 1-800-830-0970 (toll-free)

FEATURES

Feature		What it does
A	 LED Power ON/OFF Button	Turns the LED light module ON or OFF. The LED Light automatically starts when the power cable is connected to the Holster with a Battery. Available only on the <i>TotalShield II</i> Surgical Helmet with LED Light.
B	 Fan Speed Control Button	Press this button to cycle through the four available fan speeds. The fan speeds cycle from speed 1 to 2 to 3 to 4 and then back to speed 1. The fan automatically starts when the power cable is connected to the Holster with a Battery.
C	LED Light Adjustment Slide	Allows adjustment of LED light. Available only on the <i>TotalShield II</i> Surgical Helmet with LED Light.
D	LED Light Module	Provides a large area of illumination for deep cavity visualization and enhanced tissue differentiation. Available only on the <i>TotalShield II</i> Surgical Helmet with LED Light.
E	LED Light Shroud	The shroud prevents glare from the shield. Available only on the <i>TotalShield II</i> Surgical Helmet with LED Light.
F	Front Air Outlets	The outlets are designed to reduce heat and direct air away from the User's eyes bringing comfort that supports the User's ability to focus during long procedures.
G	Rear Air Outlets	The outlets are designed to direct air along the User's neck bringing comfort that supports the User's ability to focus during long procedures.
H	Foam Pads	Seven removeable foam pads provide cushion and a comfort fit.
I	Attachment Points	Secures the Toga or Hood to the Surgical Helmet.
J	Fan Inlet	Draws external ambient air in and circulates it throughout the helmet.
K	Rear Adjustment Knob	The Boa® Fit System makes it easy to single-handedly dial in a personalized fit by turning the knob clockwise to tighten and counterclockwise to loosen. Turn the headband rear adjustment dial to secure the helmet for a comfortable fit.
L	Power Cable	Provides power from the Battery Holster to the Helmet.
M	Power Cable Plug	Connects Helmet to the Battery Holster to operate the system.
N	Rear Strap	Provides up and down movement for improved fit. <i>TotalShield II</i> Surgical Helmet w/Lock and <i>TotalShield II</i> Surgical Helmet LED w/Lock have rear strap locked in down position.

DONNING INSTRUCTIONS

CAUTION Always wear a surgical head covering before using the *TotalShield II* Surgical Helmet

1. Apply the foam pads inside the helmet as necessary. Hook and loop dots on the headband are used to secure the foam pads to the headband.
2. Place the helmet on the head. Ensure clear field of vision.
3. Turn the headband rear adjustment knobs to secure the helmet for a comfortable fit.
4. Connect the helmet power cable to the Battery Holster. Insert a fully charged Battery into the Battery Holster and fasten the Battery Holster on the waist. Helmet fan operates automatically.

WARNING DO NOT place fingers on the fan inlet while the fan is in operation

5. See Battery Charge Status Window to check the status of the Battery charge. Beeps indicate that the Battery is low on charge. Press the fan speed control buttons to adjust the fan speed for personalized air flow. An audible sound will provide indication of speed.

NOTE Users are now ready to scrub prior to surgery and donning a Surgical Toga or Hood. Refer to the instructions for use included in the *TotalShield II* Zippered Surgical Toga or Surgical Hood.

NOTE When doffing the Surgical Toga or Hood, take care to avoid soiling the *TotalShield II* Surgical Helmet with bodily fluids through contact with the exterior of the toga or hood.

Additional Instructions for Surgical Helmet with LED Light Only

6. Press the LED button to turn on and off the light.

CAUTION Battery Life will decrease with the use of the LED light. Turn OFF the light when not in use.

7. Adjust the LED light adjustment lever to direct the light to the desired location.

NOTE Helmet power can be disconnected by removing the power cable from the Battery Holster.

MAINTENANCE

WARNING Do not attempt to dismantle, modify or repair this equipment or its internal components

1. Upon receipt and prior to each use, inspect the device for any damage. Do Not use if any damage is detected. If the initial inspection results are satisfactory operate the device.
2. After each use, it is recommended to inspect the helmet for evidence of blood or bodily fluid soiling. If present, or a known exposure to bodily fluids has occurred, it is recommended the helmet be discarded. When used as directed, soiling of this type is unlikely and therefore cleaning methods have not been validated to address blood and other infectious bodily fluids.

3. When used as directed, the *TotalShield II* Surgical Helmet has a low likelihood of being soiled. Therefore, cleaning is recommended on an as-needed/desired basis.

CLEANING INSTRUCTIONS

- Prior to cleaning, disconnect the Holster and Battery from the helmet.
- Foam pads are not compatible with cleaning agents. Remove the foam pads when cleaning the helmet. Replace the foam pads as desired or if soiled, damaged, or missing.
- If the helmet exhibits evidence of blood or bodily fluid soiling, or a known exposure to bodily fluids has occurred, it is recommended the helmet be discarded. When used as directed, soiling of this type is unlikely and therefore cleaning methods have not been validated to address blood and other infectious bodily fluids.
- The helmet is compatible with mild detergent solution with a neutral PH or isopropyl alcohol.
- To achieve low level disinfection:
 - Wipe the exterior of the device with a damp cloth using 70% isopropanol to remove visible soil.
 - Keep the surfaces wet with 70% isopropanol for a minimum of 1 minute.

CLEANING WARNINGS

- The LED light lens is to be cleaned only with a lens tissue.
- Do Not allow liquid and/or moisture into the electrical connections or fan opening of the helmet.
- Do Not sterilize or immerse.

BATTERY SPECIFICATIONS

CAUTION

- Only use *TotalShield II* Rechargeable Li-ion Batteries REF 00992010200 with *TotalShield II* Surgical Helmets.
- Refer to the instructions for use included in the *TotalShield II* Rechargeable Li-ion Batteries for warnings and handling recommendations.

REFERENCE NUMBER	REF 00992010200
Product Description	<i>TotalShield II</i> Rechargeable Li-ion Battery
Electrical Output	14.4 V — — —
Capacity	3100 mAh
Charging Time	2.5 hours
Battery Life Time	
Surgical Helmet	~10 hours
Surgical Helmet w/ LED	~4 hours


BATTERY Charging Instructions


CAUTION

- Only use the *TotalShield II* Battery Charger REF 00992020000 and 00992020002 to charge the *TotalShield II* Rechargeable Li-ion Batteries REF 00992010200.
- Refer to the instructions for use included in the *TotalShield II* Battery Charger for warnings and recommendations.


NOTES

- For optimal performance, charge the Batteries completely before initial use and after each use. Check the status window before each use to confirm the charging status (See image below). LED Status window will show Battery status.

1. Place the Charger on a flat, level surface away from sources of heat and moisture.
2. Connect the Power Supply to the mains AC supply using the Cord provided.
NOTE For the Two-bay Battery Charger (REF 00992020002) only, plug the DC connector from the Power Supply into the back of the Charger.
3. Place the Battery in the Charger bay located on the *TotalShield II* Battery Charger
4. Align the arrows of the Battery and the Charger to ensure that the Battery makes contact with the Charger
5. The LED status window on the Battery Charger will indicate the status of the Battery (refer to the LED indication section for details) and the Charger will automatically begin charging.
6. Remove the Battery when the Charger LED status window indicator is green solid.

LED INDICATION

NOTE The Battery status is determined by the LEDs in the status window.

Green Flashing	Battery Charging
Green Solid	Battery Fully Charged
Blue Flashing	Battery in calibration mode
Blue Solid	Battery fuel gauge calibrated
Red Flashing	Battery fuel gauge in need of recalibration
Red Solid	Error

BATTERY RECALIBRATION INSTRUCTIONS

NOTES

- During recalibration the Charger may become warm.
 - The calibration cycle will perform best if the Battery has a full charge in the beginning.
 - The Battery can either be calibrated and charged or just charged. Calibration takes longer than just charging.
 - It is recommended that the recalibration button is not pressed after the calibration has started.
1. When calibration is required, the red LED will flash when the Battery is placed into the Charger bay.
NOTE Recalibration is used to reset the Battery gauge to full accuracy.
 2. Press the recalibration button on the front of the Charger (see Figure 9) to recalibrate the Battery. If the button is not pressed the Charger will automatically charge the Battery.
 3. When the blue LED flashes, the Battery is recalibrating.
 4. When the calibration is completed, a solid blue LED will be present.

TROUBLESHOOTING GUIDE


SYMPTOM	POSSIBLE CAUSES/SOLUTIONS
Insufficient and/or no airflow	• Increase the fan speed using the fan speed control button. • Ensure the helmet power cable is securely connected to the Battery Holster. • Replace the Battery with a fully charged Battery. • Replace the Battery Holster. • Ensure the Surgical Hood or Toga filter is properly placed over the fan inlet. • Faulty Helmet, replace the Helmet if necessary. • Electromagnetic disturbance causes fan to stop working properly. Restart helmet by removing and reinserting battery. Replace the helmet if necessary.
Face Shield is covered with fog	• Increase the fan speed using the fan speed control button. • Ensure the helmet power cable is securely connected to the Battery Holster. • Replace the Battery with a fully charged Battery. • Replace the Battery Holster. • Ensure the Surgical Hood or Toga vent is properly placed over the fan inlet. • Faulty Helmet, replace the Helmet if necessary.
Insufficient and/or no LED light	• Press the LED power button to turn on the light. • Ensure the helmet power cable is securely connected to the Battery Holster. • Replace the Battery with a fully charged Battery. • Replace the Battery Holster. • Faulty Helmet, replace the Helmet if necessary. • Electromagnetic disturbance causes fan to stop working properly. Restart helmet by removing and reinserting battery. Replace the helmet if necessary.
Helmet beeps	• Remove the Battery and re-insert the Battery. • Replace the Battery. • Replace the Holster and re-insert the Battery.

REPORTING PROBLEMS

The user and/or patient should report any suspected serious incident related to the device by informing the manufacturer and the competent authority of the member state in which the serious incident has occurred

DISPOSAL/RECYCLE

Follow the hospital and/or local regulations to recycle or dispose of electrical equipment at the end of its useful life.

Li-ion

NOTICE Battery Pack Disposal

The Battery is of a Lithium-Ion type. When the Battery reaches the end of its useful life, the Battery should be disposed of by a qualified recycler or hazardous materials handler. Do not mix this Battery with the solid waste stream. Contact your local recycler or hazardous material handler for recycling or disposal information.

ENVIRONMENTAL CONDITIONS

Surgical Helmet

The following are the recommended environmental conditions for operation, storage and transportation of the *Total/Shield II* Surgical Helmets

CONDITION	OPERATION	TRANSPORTATION AND STORAGE
Ambient temperature	68°F (20°C) to 73°F (23°C) Typical 50°F (10°C) to 110°F (38°C) Maximum	0°F (-18°C) and 113°F (45°C) Typical
Relative humidity	30% to 60%	up to 80%
Altitude	≤2,000 m – Not intended for high altitude environments	

NOTE Transport or store in an environment that limits exposure to dust, moisture and temperature extremes.

BATTERY AND CHARGER TRANSPORTATION AND STORAGE CONDITION

Recommended Storage and Transportation Conditions 0°F (-18°C) to 113° F (45°C), up to 80% RH

CAUTION For optimal Battery life the Battery packs should be stored in an environment with low humidity, free from corrosive gas at a recommended temperature <70°F (21°C). Extended exposure to temperatures above 113°F (45°C) could degrade Battery life and performance.

NOTE Transport or store the Batteries and Chargers in an environment that limits exposure to dust, moisture and temperature extremes.

Electromagnetic Compatibility (EMC) GUIDANCE TABLES

WARNING The *TotalShield II* Surgical Helmets are suitable for use in Professional Healthcare Facility Environments except for near active HF SURGICAL EQUIPMENT and the RF shielded room of a ME SYSTEM for magnetic resonance imaging.

WARNING Use of the *TotalShield II* Surgical Helmets adjacent to or stacked with other equipment should be avoided because it could result in improper operation.

WARNING Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of *TotalShield II* Surgical Helmets. Otherwise, degradation of the performance of this equipment could result.

WARNING Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

The following tables provide guidance on electromagnetic compatibility for the *TotalShield II* Surgical Helmets

EMC GUIDANCE AND DECLARATION - EM EMISSIONS

Guidance and manufacturer's declaration – electromagnetic emissions		
The <i>TotalShield II</i> Surgical Helmets are intended for use in the electromagnetic environments specified below. The customer or the user should assure that they are used in such an environment. The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or reorienting the equipment.		
EMISSIONS TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
RF emissions CISPR 11	Group 1	The <i>TotalShield II</i> Surgical Helmets use RF energy only for their internal function. Therefore, their RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations/Flicker emissions IEC 61000-3-3	N/A	The <i>TotalShield II</i> Surgical Helmets are suitable for use in all establishments other than domestic and those directly connected to the public low-voltage Power Supply network that supplies buildings used for domestic purposes.

EMC GUIDANCE AND DECLARATION - EM IMMUNITY/DISTURBANCES

Guidance and manufacturer's declaration – electromagnetic immunity The <i>TotalShield II</i> Surgical Helmets are intended for use in the electromagnetic environments specified below. The customer or the user should assure that they are used in such an environment.			
IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT – GUIDANCE
Electrostatic discharge (ESD) IEC 61000-4-2	± 2, ± 4, ± 8 kV contact	± 2, ± 4, ± 8 kV contact	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic materials, the relative humidity should be at least 30%.
	± 2, ± 4, ± 15 kV air	± 2, ± 4, ± 15 kV air	
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for Power Supply lines	N/A	~ AC Power Mains power quality should be that of a typical commercial or hospital environment.
	± 1 kV for input/output lines	N/A	
Surge IEC 61000-4-5	± 1 kV line(s) to line(s)	N/A	~ AC Power Mains power quality should be that of a typical commercial or hospital environment.
	± 2 kV line(s) to earth	N/A	
Voltage dips, short interruptions and voltage variations on Power Supply input lines IEC 61000-4-11	100% drop, 0.5 periods, 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°	N/A	~ AC Power Mains quality should be that of a typical commercial or hospital environment. If the user requires continued operation during ~ Power Mains interruptions, it is recommended to be powered from an uninterruptible Power Supply or a Battery.
	30% dip, 25/30 periods	N/A	
	100% dip, 1 period	N/A	
	Voltage Interruptions (all input current)	N/A	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

EMC GUIDANCE AND DECLARATION - EM EMISSIONS/RF

Guidance and manufacturer's declaration – electromagnetic immunity			
The <i>TotalShield II</i> Surgical Helmets are intended for use in the electromagnetic environments specified below. The customer or the user should assure that they are used in such an environment.			
IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the <i>TotalShield II</i> Surgical Helmets, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.7 GHz where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <i>d</i> is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance levels in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m	
NOTES - At 80 MHz and 800 MHz, the higher frequency range applies. - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephone and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the <i>TotalShield II</i> Surgical Helmets are used exceeds the application RF compliance levels above, the <i>TotalShield II</i> Surgical Helmets should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the <i>TotalShield II</i> Surgical Helmets.			
b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

EMC GUIDANCE AND DECLARATION – IMMUNITY/SEPARATION DISTANCES

Recommended separation distances between portable and mobile RF communication equipment and <i>TotalShield II</i> Surgical Helmets			
The <i>TotalShield II</i> Surgical Helmets are intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the <i>TotalShield II</i> Surgical Helmets can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the <i>TotalShield II</i> Surgical Helmets as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter in Watts (W)	Separation distance according to frequency transmitter in Meters (m)		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.7 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTES - At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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© 2019, 2020 Zimmer Surgical, Inc.