# GIMA

REF 28183

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**SUPERVEGA EVO** is a surgical aspirator with the following electrical features: 14V  $\overline{\text{---}}$  4A with AC/DC adapter mod. UE60-140429SPA1 of FUHUA (input: 100-240V~ - 50/60Hz - 100VA) or Internally powered equipment (Li-Ion Battery 14.8V  $\overline{\text{---}}$  5.2A) or with cigarette lighter adapter (12V  $\overline{\text{---}}$  4A). **SUPERVEGA EVO** is a desk-type electric suction unit for the aspiration of body liquids, oral, nasal and tracheal aspiration in adults or children. This device has been designed to offer ease of transport and continuous use, thanks to an electronic system that manages the power supply. The large LCD display facilitates the use of the device and increases suction by adjusting the control buttons.

The signal on the LCD screen, placed on the front panel, allows you to monitor the charge status of the Lithium battery. The lithium battery and the innovative Feedback system, which guarantees a smart use by automatically controlling and adjusting the suction power, allow the battery to increase its autonomy and decrease the noise produced.

The "PROXIMITY" function allows the device to be activated/deactivated via an infra-red proximity sensor (detecting the presence of the hand from a distance of tenths of centimetres without touching the suction unit) and it prevents and avoids possible cross-contamination between patients as they are treated in turns.

Made of highly heat resistant, electrically insulated plastic material in conformity with the latest European safety standard, the product is supplied with a complete polycarbonate autoclavable jar with overflow valve.

## GENERAL WARNING

**READ INSTRUCTION MANUAL CAREFULLY BEFORE USE.**

**THE DEVICE IS FOR USE BY QUALIFIED PERSONNEL (SURGEON / PROFESSIONAL NURSE / ASSISTANT)**

**THE USE OF THE DEVICE AT HOME IS RESTRICTED TO AN ADULT IN FULL POSSESSION OF MENTAL FACULTIES AND / OR HOME CARERS**

**THE INSTRUMENT MUST NOT DISASSEMBLED. FOR TECHNICAL SERVICE ALWAYS CONTACT GIMA S.P.A.**

## IMPORTANT SAFETY RULES

1. Check the condition of the unit before each use. The surface of the unit should be carefully inspected for visual damage. Check the mains cable and **do not connect to power** if damage is apparent;
2. Before connecting the appliance always check that the electric data indicated on the data label and the type of plug used, correspond to those of the mains electricity to which it's to be connected;
3. Respect the safety regulations indicated for electrical appliances and particularly:
  - Use original components and accessories provided by the manufacturer to guarantee the highest efficiency and safety of the device;
  - The device can be used only with the bacteriological filter;
  - Never immerse the appliance into water;
  - Do not place or store the aspirator in places where it may fall or be pulled into the bathtub or washbasin. In the event it is accidentally dropped, do not attempt to remove the device from the water whilst the plug is still connected: disconnect the mains switch, remove the plug from the power supply and contact the GIMA technical service department. Do not attempt to make the device work before it has been thoroughly checked by qualified personnel and/or the GIMA technical service department.
  - Position the device on stable and flat surfaces in a way that the air inlets on the back aren't obstructed;
  - To avoid incidents, do not place the aspirator on unstable surfaces, which may cause it to accidentally fall and lead to a malfunction and/or breakage. Should there be signs of damage to the plastic parts, which may expose inner parts of the energised device, **do not connect the plug to the electrical socket**. Do not attempt to make the device work before it has been thoroughly checked by qualified personnel and/or the GIMA technical service department.
  - Don't use in the presence of inflammable substances such as anaesthetic, oxygen or nitrous oxide;
  - Don't touch the device with wet hands and always prevent the appliance coming into contact with liquids;
  - Don't leave the appliance connected to the power supply socket when not in use;
  - Don't pull the power supply cable to disconnect the plug remove the plug from the mains socket correctly;
  - Store and use the device in places protected against the weather and far from any sources of heat. After each use, it is recommended to store the device in its own box away from dust and sunlight.
  - Don't use the device thoracic drainage.
  - In general, it is inadvisable to use single or multiple adapters and/or extensions. Should their use be necessary, you must use ones that are in compliance with safety regulations, however, taking care not to exceed the maximum power supply tolerated, which is indicated on the adapters and extensions.
4. For repairs, exclusively contact technical service and request the use of original spare parts. Failure to comply with the above can jeopardise the safety of the device;
5. **Use only for the purpose intended.** Don't use for anything other than the use defined by the manufacturer. The manufacturer will not be responsible for damage due to improper use or connection to an electrical system not complying with current regulation.
6. The medical device requires special precautions regarding electromagnetic compatibility and must be installed and used in accordance with the information provided with the accompanying documents: the SUPERVEGA EVO device must be installed and used away from mobile and portable RF communication devices (mobile phones, transceivers, etc.) that may interfere with the said device.
7. Instrument and accessory discharging must be done according to current regulations in the country of use.

**SYMBOLS**

	Insulation Class II (when used with AC / DC adapter or when connected to a car cigarette lighter cable)	
	CE marking in conformity with EC directive 93/42/CEE and subsequent changes	
	General warnings and/or specifications	
	Consult the instruction manual	
	Manufacturer: GIMA S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) - Italy	
	Date of manufacture	
	Applied Part type BF (suction probe)	
	Keep away from sunlight	
	Keep in a cool, dry place	
	Atmospheric Pressure	
	Battery (Li-Ion Battery 14,8V --- 5.2A)	
	Direct Current	
	ON / OFF	
	Batch production	
	Serial number	
	Model / Ref number	
 10R-041439	Device approved according to European Directive 95/54/EEC and international standards ECE-R10	
	Degree of protection an electrical device provides in the case of accidental or intentional contact with the human body or with objects, and protection in the case of contact with water.	
	<b>1st DIGIT PENETRATION OF SOLIDS</b>	<b>2nd DIGIT PENETRATION OF LIQUIDS</b>
	Protected against solids having a dimension greater than Ø 12mm	Protected against the vertical flow of drops of water

**DISPOSAL OF WASTE BATTERIES - (Directive 2006/66/EC)**

The battery provided with this product shall not be treated as household waste. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. The recycling of the materials will help to conserve natural resources.

At the end of their life hand the batteries over to the applicable collection points for the recycling of waste batteries. For more detailed information about recycling of this product or battery, please contact your local Civic Office, your household waste disposal service or the shop where you purchased the product.

**CLEANING DEVICE**

Use a soft dry cloth with not – abrasive and not – solvent detergents. To clean the device external parts always use a cotton cloth dampened with detergent. Don't use abrasive or solvent detergents.

**GB**

**PARTICULAR CARE SHOULD BE TAKEN TO ENSURE THAT THE INTERNAL PARTS OF THE EQUIPMENT DO NOT GET IN TOUCH WITH LIQUIDS. NEVER CLEAN THE EQUIPMENT WITH WATER.**

During all clearing operations use protection gloves and apron (if need be, also wear a face mask and glasses) to avoid getting in contact with contaminating substances (after each utilization cycle of the machine).

## ACCESSORIES SUPPLIES

DESCRIPTION	CODE
COMPLETE ASPIRATION JAR 1000ml	RE 210001/02
CONICAL FITTING	RE 210420
TUBES SET 6 mm x 10 mm	51100/01
HYDROPHOBIC AND ANTIBACTERIAL FILTER	SP 0121
SUCTION PROBE CH20	25723
AC/DC ADAPTER (UE60-140429SPA1)	SP 0208/01
POWER SUPPLY CORD FOR AC/DC ADAPTER	SP 0020/03
CIGARETTE LIGHTER CABLE	SP 0007/02

The filter is produced with (PTFE) hydrophobic material to prevent fluids entering the pneumatic circuit. It should be changed immediately if it becomes wet or if there is any sign of contamination or discolouration. It should also be changed if the unit is used with a patient whose risk of contamination is unknown. **Don't use the suction unit without the protection filter.** If the suction unit is used in an emergency or in a patient where the risk of contamination is not known the filter must be changed after each use. Available under request with different versions with complete jar 2000ml.

**Suction catheter:** Single-use device to be used on a single patient. Do not wash or re-sterilize after use. Reuse may cause cross-infections. Don't use after lapse of the sell-by date

Check the expiry date on the original packaging of the suction catheter and check the integrity of the sterile packaging. GIMA declines any liability for injury to the patient correlated to the deterioration of the above-mentioned sterile packaging due to handling of the original packaging by third parties.

**WARNING:** Suction tubes for insertion in the human body purchased separately from the machine should comply with ISO 10993-1 standards on material biocompatibility.

**Aspiration jar:** The mechanical resistance of the component is guaranteed up to 30 cycles of cleaning and sterilization. Beyond this limit, the physical-chemical characteristics of the plastic material may show signs of decay. Therefore, we recommend that you to change it.

**Silicone tubes:** the number of cycles of sterilization and/or cleaning is strictly linked to the employment of the said tube. Therefore, after each cleaning cycle, it is up to the final user to verify whether the tube is suitable for reuse. The component must be replaced if there are visible signs of decay of the material constituting the said component.

**Conical fitting:** the number of cycles of sterilization and the number of cleaning cycles is strictly linked to the employment of the said component. Therefore, after each cleaning cycle, it is up to the final user to verify whether the fitting is suitable for reuse. The component must be replaced if there are visible signs of decay in the material constituting the said component.

**Service life of the device:** More than 30000 hours of operation in accordance with the standard conditions of testing and operation. Shelflife: maximum 5 years from the date of manufactured.

## CLEANING OF ACCESSORIES

Before using the device, the manufacturer advises you to clean and/or sterilize the accessories.

Washing and / or cleaning the autoclavable jar as to be carried out as follows:

- Wear protection gloves and apron (glasses and face mask if necessary) to avoid contact with contaminating substances;
- Disconnect the tank from the device and remove the said container from the support of the device.
- Separate all the parts of the cover (overflow device, washer).
- Disconnect all tubes from the jar and the protection filter
- Wash each part of the container from secretions under cold running water and then clean every single part in hot water (temperature not exceeding 60°C)
- Once again, carefully wash each single part using, if necessary, a non-abrasive brush to remove any deposits. Rinse with hot running water and dry all parts with a soft cloth (non-abrasive). It is possible to wash with commercial disinfectants by carefully following the instructions and dilution values supplied by the manufacturer. After cleaning, leave the parts to dry in an open, clean environment.
- Dispose of the aspiration catheter according to that provided by local laws and regulations.

The silicone aspiration tubes and the conical fitting may be carefully washed in hot water (temperature must not exceed 60°C). After cleaning, leave the parts to dry in an open, clean environment.

When cleaning is complete, reassemble the container for liquid aspirations according to the following procedure:

- Place the overflow valve into its seat in the cover (under VACUUM connector)
- Insert floating valve keeping the o-ring towards the opening of the cage
- Place the o-ring into its seat around the cover
- After completing assembling operations always make sure that cover seals perfectly to avoid vacuum leakages or liquid exit

After disposing of disposable parts and disassembling the jar wash in running cold water and rinse thoroughly. Then soak in warm water (temperature shall not exceed 60°C). Wash thoroughly and if necessary use a non-abrasive brush to remove incrustations. Rinse in running warm water and dry all parts with a soft cloth (non-abrasive). The jar and the cover can be autoclaved by placing the parts into the autoclave and running one sterilization stem cycle at 121°C (1 bar relative pressure – 15 min) making sure that the jar is positioned upsidedown. Mechanical resistance of the jar is guaranteed up to 30 cycles of sterilization and cleaning at the indicated conditions (EN ISO 10079-1). Beyond this limit the physical-mechanical characteristics of the plastic may decrease and replacement of the part is therefore recommended. After sterilization and cooling at environment temperature of the parts make sure that these are not damaged. The aspiration tubes can be sterilized on autoclave using a sterilization cycle at 121°C (1 bar relative pressure – 15 min). The conical connector can be sterilized on autoclave using a sterilization cycle at 121°C (1 bar relative pressure – 15 min).

**DO NOT WASH, STERILIZE OR PUT IN AUTOCLAVE THE ANTIBACTERIAL FILTER**

**PERIODICAL MAINTENANCE CHECKS**

The **SUPERVEGA EVO** suction equipment does not need maintenance or lubrication. It is, however to inspect the unit before each use. With regard to training, given the information contained in the user manual and since it is easy to understand the said device, it doesn't appear to be necessary. Unpack the instrument and **always check** integrity of plastic parts, they might have been damaged during previous use. Turn switch on and verify that loud noises are not present, these can indicate wrong functioning. Make sure you do not hear excessively annoying noises that may highlight a malfunction. The device is protected by a safety fuse (**F 10A L 250V**) situated in the cigarette lighter cable. When replacing, always check the type and value as indicated. Internally, the device is protected (see electrical specifications) by two fuses (**T 15A L 125V**) that cannot be reached from the outside. Therefore, contact the manufacturer to request the assistance of an authorized and qualified technician when they need to be replaced. If it's replaced make sure that its replacement is always the same type and value, as indicated. The device is equipped with a Lithium-Ion Battery, which cannot be reached from the outside. Only refer to GIMA technical service personnel to replace it. The device should be checked at least once every 12 months by technical service. Every 24 months it is compulsory to have a safety inspection and technical maintenance performed on the device.

**ONLY USE BATTERIES RECOMMENDED BY GIMA. USE OF BATTERIES OTHER THAN THOSE RECOMMENDED WILL MAKE THE WARRANTY VOID AND NULL. ONLY QUALIFIED PERSONNEL ARE AUTHORISED TO REPLACE THE INTERNAL BATTERY. ANY OPERATION PERFORMED BY NON-TRAINED PERSONNEL MAY CAUSE DANGER (E.G. EXCESSIVE TEMPERATURE)**

<b>Fault type</b>	<b>Cause</b>	<b>Solution</b>
1. Fixed Red Back-Light	Low battery	If it is intended to be used in ambulances, position the device on the support bracket and leave the battery charge until the <b>Fixed White Back-Light</b> gives the signal (TAB. I).
2. No LED is On and the Back-Light failed to switch on	Locking device	Internal technical problem. Contact technical assistance.
3. No aspiration	Jar Cap badly screwed down	Unscrewed the cap, then rescrew it correctly
4. No aspiration	Lid seal not in its seat	Unscrew the cap and insert the seal properly in its seat
5. The float doesn't close	If the cap has been washed, ensure that the float is not partially detached	Insert the float into it's place
6. The float doesn't close	The float it's covered by dirty material	Unscrewed the cap, leave the and put in on autoclave
7. Low suction	Foam inside the jar	Fill the jar to 1/3 full of ordinary water
8. No aspiration due to flow leakage of mucus	Filter blocked	Replace filter
9. The Vacuum power on the patient side is either very low or absent	<ul style="list-style-type: none"> <li>• Protection filter blocked or damaged</li> <li>• Connection tubes blocked, kinked or disconnected</li> <li>• Shut-off valve blocked or damaged</li> <li>• Pump motor damaged</li> </ul>	<ul style="list-style-type: none"> <li>• Replace the filter</li> <li>• Replace or reconnect the tubes, check the jar connections</li> <li>• Empty the jar, or disconnect the tube from the jar and unblock the shut-off valve. The unit will only work in the upright position</li> <li>• Contact the technical service</li> </ul>
10. Noisy	Technical internal problem	Contact the technical service
<b>Faults 1 - 2 - 3 - 4 - 5 - 6 - 7 - 8 - 9 - 10</b>	<b>None of the remedies has achieved the desired results</b>	<b>Contact the seller or GIMA After-sales Assistance Service</b>

**Lithium-Ion Battery Charging Cycles:** The Lithium-Ion battery contained inside the device is guaranteed for a number greater than 300 charging cycles. In the vicinity of 300 charging cycles you can ask for verification of operating status to the manufacturer or to require replacement of the battery pack in such a way that you always have the component in perfect condition.

If the overflow security system it's activated, don't proceed with the liquid aspiration. If the overflow security system doesn't work there are two cases:

1° case - If the overflow security system doesn't work the aspiration will be stopped by the bacteriological filter who avoid the liquid penetration inside the device.

2° case - If both the security system doesn't work, there is the possibility that liquid comes inside the device, in this case return the device to GIMA technical service.

*GIMA S.p.A. will provide upon request electric diagrams, components list, descriptions, setting instructions and any other information that can help the technical assistance staff for product repair.*

**BEFORE EVERY CHECKING OPERATION, IN CASE OF ANOMALIES OR BAD FUNCTIONING, PLEASE CONTACT GIMA TECHNICAL SERVICE. GIMA DOES NOT GIVE GUARANTEE IF INSTRUMENT, AFTER THE TECHNICAL SERVICE CHECKING, APPEARS TO BE TAMPERED.**

## INSTRUCTION FOR USE

- The device must be checked before each use in order to detect malfunctions and / or damage caused by transport and / or storage.
- The working position must be such as to allow one to reach the control panel and to have a good view of the LCD screen, the suction tank and the antibacterial filter.

**WARNING:** For proper use, place the aspirator on a flat, stable surface in order to have the full volume of use of the jar and better efficiency of the overflow device.

### Operation using with AC/DC switching adapter

- Connect the short silicon tube, with antibacterial filter, to the suction connector (VACUUM). The other tube, with one end connected to the antibacterial filter with the other end to jar's lid connector where has been fixed the red float. When the 90% of the volume of the jar is reached there is the activation of the security float (the float close the aspiration connector on the jar) to the avoid liquid penetration inside the device.
- Connect the long silicon tube to the other jar's lid connector (PATIENT PORT);
- Connect the switching adapter to the device with the appropriate connector and insert the power cable plug to the power socket. To start the treatment press the ON/OFF button and check its green back-light. When the device is switched on the LCD screen activates displaying a bar-graph, graduated scale, and negative pressure value, which must be set by the operator;
- The negative pressure is set with keys - and +, which are found below the LCD screen: Press one of the two keys with the motor on. The level of the negative pressure will flash at the centre of the display. Release the key and, after 3 seconds, the value will be operating (fixed digit on the display) and saved. The desired value remains fixed until the operator varies the negative pressure, whilst the bar-graph moves along the semi-circumference stopping in correspondence (graph) of the value set.

The maximum negative pressure value that can be set is: -75kPa (-0.75 bar)

The minimum negative pressure value that can be set is: -10kPa (-0.10 bar)

- To suspend and / or terminate the treatment press the ON / OFF button and and pull the plug out from the power socket.
- Unscrew the jar's lid and fill the jar 1/3 full or ordinary water (this for an easy cleaning operations and an rapid reaching of the functionally vacuum) then rescrew the lid on the jar correctly.
- To extract the accessories and start with cleaning.
- The battery is fully charged in about 70 minutes with continuous operation.


**WARNING:** Before using the device, check the battery power status. Before each use proceed with charging the battery. To maintain the device in good conditions, recharge the battery every 3 months (when not in use).

**WARNING:** Ensure that the FLUID SIDE or IN marker on the filter is on the side facing the collection jar lid and fitted into the "VACUUM". A wrong connection causes immediate destruction in case of contact with sucked liquids.

**WARNING:** The power supply cable plug is the element of separation from the electrical mains system: even if the units equipped with a special on / off switch button, the power supply plug must be kept accessible once the device is in use so as to allow a further method of disconnection from the mains supply system.


**Filter assembling**

The device features an energy saving function which allows the device to increase its autonomy and decrease the number of revs of the motor. In order to do so, it controls the negative pressure which is generated only when the operator suctions body fluids. If the operator does not carry out the suction process the device lowers the number of revs of the motor, thus protecting the duration of the internal battery.

**NEVER USE THE DEVICE WITHOUT JAR AND / OR PROTECTION FILTER**

**Operation using cigarette lighter DC 12V**

- Connect the device's external plug 12V to the lighter plug with the cigarette lighter cable. Check the battery power status of the vehicle before the cigarette lighter cable. Press the switch to start suction
- Press the switch to the I position to turn it on.

**WARNING: Only use the originally supplied or recommended replacement cigarette lighter cables (view the chapter "Important Safety Rules")**

**WARNING:** Before using the device, check the battery power status. Before each use proceed with charging the battery. To maintain the device in good conditions, recharge the battery every 3 months (when not in use).

**Recharging operations:** to be able to charge the internal battery it is necessary to connect the universal switching adapter (mod. UE 60140429SPA1 of FUHUA) to the electric network for approx. 360 minutes with the main switch to position 0. The battery's autonomy when fully charged is approx. 70 minutes with continued operations.

**LIGHT INDICATORS:**

The device is provided with a light indicator (directly on the LCD screen) which allows you to view the operation of the device, autonomy of the battery, and the charge phase in progress.

The light indications, which appear during operation, are indicated in table I.

The charge cycle starts automatically it is necessary to connect the universal switching adapter (mod. UE 60140429SPA1 of FUHUA) to the electric network. The charge occurs only when the motor of the suction unit is off.

DISPLAY IN CHARGE: this is on with white back-light and it remains on for the entire time. It switches off if there is no external voltage. The charge phase lasts until the device is powered and the battery is kept at its maximum charge.

The large digits display "CHA" .

**TAB. I – INDICATOR LIGHTS DURING OPERATIONS**

Back-light signals	Phase	Problem / Cause	Solution
Fixed white back-light and progressive activation of the bar-graph	During the Charge phase	Battery charge in progress	Wait
Fixed white back-light and fixed activation of the bar-graph	During the Charge phase	Charging cycle complete	Remove the power supply
Fixed green back-light	During battery operation	Primary status / Battery fully charged	Battery operation guaranteed
Fixed yellow back-light	During battery operation	Intermediate status / Battery not fully charged	Battery operation guaranteed / start the charge cycle upon activation of the red back-light
Fixed red back-light	During battery operation	Low battery	Start charge phase <b>ATTENTION:</b> A long and continuous beep is produced during this signal to warn the operator about the low charge level of the battery.
Flashing red back-light	Automatic switch-off of the device due to discharged battery	Battery fully discharged	The flashing back-light switches on upon reactivation of the device. Immediately start the battery charge cycle.

**TAB. II – WRITTEN SIGNALS / BUTTON LED SIGNAL**

Button LED signal / written signal	Function	Colour	Position
ON/OFF button	Power on	Green	Near the front panel key
Proximity Button	Switch the ON/OFF key	Blue	LED located above the "Proximity" key
CHA	Displays the battery charge status	Fixed white back-light and progressive activation of the bar-graph	Central in digits of large dimensions
Small three-digit number (000)	Displays the voltage of the device if in charge phase or it displays the full scale of the bar-graph.	Fixed white or green back-light depending on whether external power supply or battery and progressive activation of the bar-graph are used.	Sideways in small digits
Bar-graph (arc bar-graph)	It displays the trend of the negative pressure or of the charge status	Progressive activation of the bar-graph (black)	Arc-shaped in the middle of the screen
Large three-digit number (000)	It displays the negative pressure set by keys "up" and "dw"	Fixed white or green back-light depending on whether external power supply or battery is used.	Central in digits of large dimensions

**PROXIMITY FUNCTIONS:**

- Pressing the "Proximity" (Lights BLUE LED) is activated on / off the motor by an infrared proximity sensor that detects the presence of the operator's hand from a few inches away. This allows the operator to use the device without touching or focus on pressing the button.
- The on / off button is also active with the function of proximity inserted and can be used as an alternative.
- To remove the "Proximity" you have to re-press the button.

The "Proximity" is retained, or if it was active before power off, then on again to return to this, but if it was not active will remain disabled. The function set to power down the card after 2 minutes after turning off the engine if it is not then on again.

The unwanted approach of the hand to the Proximity sensors causes the device to switch off. To reactivate the function, place the hand close to the sensors.

**NOTE:** It is up to the end user to activate or deactivate the Proximity sensor. With the Function off, the device is activated/deactivated by pressing the ON/OFF key.

**RISK OF ELECTROMAGNETIC INTERFERENCE AND POSSIBLE REMEDIES**

This section contains information regarding the conformity of the compliance with the EN 60601-1-2 Standard. The SUPERVEGA EVO surgical aspirator is an electro-medical device that requires particular precautions regarding electromagnetic compatibility and which must be installed and commissioned according to the electro-magnetic compatibility information supplied. Portable and mobile radio communication devices (mobile phones, transceivers, etc.) may interfere with the medical device and should not be used in close proximity with, adjacent to or on top of the medical device. If such use is necessary and unavoidable, special precautions should be taken so that the electro-medical device functions properly in its intended operating configuration (for example, constantly and visually checking for the absence of anomalies or malfunctions). The use of accessories, transducers and cables different to those specified, with the exception of transducers and cables sold by the appliance and system manufacturer as spare parts, can lead to an increase in emissions or in a decrease of the immunity of the device or system. The following tables supply information regarding the EMC (Electromagnetic Compatibility) characteristics of the electro-medical device.

<b>Guidance and manufacturer's declaration - electromagnetic Emissions</b>		
The surgical aspirator SUPERVEGA EVO is intended for use in the electromagnetic environment specified below. The customers or the user of the surgical aspirator SUPERVEGA EVO should assure that it's used in such an environment.		
<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>
Irradiated / Conducted emissions CISPR11	Group 1	The surgical aspirator SUPERVEGA EVO only used RF energy only for its internal functioning. Therefore, its RF emissions are very low and are not cause interference in proximity of any Electronic appliances.
Irradiated / Conducted emissions CISPR11	Class [B]	The surgical aspirator SUPERVEGA EVO can be used in all environments, including domestic and those connected directly to the public mains distribution that supplies power to environments used for domestic scopes.
Harmonic emissions EN 61000-3-2	Class [A]	
Voltage fluctuations / flicker emissions EN 61000-3-3	Complies	

<b>Guidance and manufacturer's declaration - Immunity Emissions</b>			
The surgical aspirator SUPERVEGA EVO is intended for use in the electromagnetic environment specified below. The customers or the user of the surgical aspirator SUPERVEGA EVO should assure that it's used in such an environment.			
<b>Immunity Test</b>	<b>Level indicated by the EN 60601-1-2</b>	<b>Compliance Level</b>	<b>Electromagnetic environment - guidance</b>
Electrostatic discharge (ESD) EN 61000-4-2	± 6kV on contact ± 8kV in air	The device doesn't change its state	Floors should be wood, concret or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst EN 61000-4-4	± 2kV power supply lines ± 1kV for input / output lines	The device doesn't change its state	Mains power quality should be that of a typical commercial environment or hospital.
Surge EN 61000-4-5	± 1kV differential mode	The device doesn't change its state	Mains power quality should be that of a typical commercial environment or hospital.
Loss of voltage, brief voltage interruptions and variations EN 61000-4-11	5%U <sub>T</sub> (>95% dip U <sub>T</sub> ) for 0.5 cycle 40%U <sub>T</sub> (>60% dip U <sub>T</sub> ) for 5 cycle 70%U <sub>T</sub> (>30% dip U <sub>T</sub> ) for 25 cycle <5%U <sub>T</sub> (>95% dip U <sub>T</sub> ) for 5 sec	-	Mains power quality should be that of a typical commercial environment or hospital If the user of the surgical aspirator SUPERVEGA EVO request that the appliance operates continuously, the use of a continuity unit is recommended.
Magnetic field EN 61000-4-8	3A/m	The device doesn't change its state	The power frequency magnetic field should be measured in the intended installation location to assure that it's sufficiently low.
Nota U <sub>T</sub> is the value of the power supply voltage			

**Guidance and manufacturer's declaration – Immunity Emissions**

The surgical aspirator SUPERVEGA EVO is intended for use in the electromagnetic environment specified below. The customers or the user of the surgical aspirator SUPERVEGA EVO should assure that it's used in such an environment.

<b>Immunity Test</b>	<b>Level indicated by the EN 60601-1-2</b>	<b>Compliance Level</b>	<b>Electromagnetic environment - guidance</b>
Conducted Immunity EN 61000-4-6	3Vrms 150kHz to 80MHz (for non life-supporting devices)	$V_1 = 3 \text{ V rms}$	The portable and mobile RF communication devices, including cables, must not be used closer to the SUPERVEGA EVO device, than the separation distance calculated by the equation applicable to the transmitter frequency. Recommended separation distance $d = [3.5 / V_1] \sqrt{P}$ $d = [12 / E_1] \sqrt{P}$ from 80 MHz to 800MHz $d = [23 / E_1] \sqrt{P}$ from 800 MHz to 2.5 GHz  Where P is the maximum nominal output voltage of the transmitter in Watt (W) depending on the manufacturer of the transmitter and the recommended separation distance in metres (m). The intensity of the field from the fixed RF transmitters, as determined by an electro-magnetic study of the site <sup>1)</sup> , could be lower than the level of conformity of each frequency interval <sup>2)</sup> . It is possible to check for interference in proximity to devices identified by the following symbol: 
Radiated Immunity EN 61000-4-3	3V/m 80MHz to 2.5GHz (for non life-supporting devices)	$E_i = 3 \text{ V / m}$	

Note 1: At 80 MHz and 800 MHz the interval with the highest frequency is applied

Note 2: These guide lines may not be applicable in all situations. The electro-magnetic propagation is influenced by the absorption and by reflection from buildings, objects and people.

a) The field intensity for fixed transmitters such as the base stations for radiotelephones (mobile and cordless) and terrestrial mobile radio, amateur radio devices, radio AM and FM transmitters and TV transmitters can not be theoretically and accurately foreseen. To establish an electro-magnetic environment generated by fixed RF transmitters, an electro-magnetic study of the site should be considered. If the field intensity measured in the place where the device will be used surpasses the above mentioned applicable level of conformity, the normal functioning of the device should be monitored. If abnormal performance arises, additional measures such as changing the device's direction or positioning may be necessary.

b) The field intensity on an interval frequency of 150 kHz to 80 MHz should be less than 3 V/m.

**Recommended separation distance between portable and mobile radio-communication devices and the monitor**

The SUPERVEGA EVO surgical aspirator is intended to operate in an electro-magnetic environment where RF irradiated interferences are under control. The client or operator of the SUPERVEGA EVO device can help prevent electro-magnetic interference by keeping a minimum distance between the portable and mobile RF communication devices (transmitters) and the SUPERVEGA EVO device, as recommended below, in relation to the radio-communication maximum output power.

<b>Maximum nominal output power of the Transmitter W</b>	<b>Separation distance from the frequency transmitter (m)</b>		
	150 kHz to 80 MHz $d = [3.5 / V_1] \sqrt{P}$	80 MHz to 800 MHz $d = [12/E_1] \sqrt{P}$	800 MHz to 2.5 GHz $d = [23/E_1] \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters with a maximum nominal output power not shown above, the recommended separation distance in metres (m) can be calculated using the equation applicable to the transmitter frequency, where P is the maximum nominal output power of the transmitter in Watt (W) depending on the transmitter's manufacturer.

Note 1: At 80 MHz and 800 MHz the interval with the highest frequency is applied

Note 2: These guide lines may not be applicable in all situations. The electro-magnetic propagation is influenced by the absorption and by the reflection from buildings, objects and people.

**WARRANTY CONDITIONS**

This product is guaranteed for a period of **24 months** from the date of purchase. The warranty includes the repair or replacement of defect spare parts free of charge, if the defect has been clearly described by the customer and determined by technical service. Inspections on the part of the seller, performed at the request of the customer and intended to determine whether the device is fully functional, are not covered by the free-of-charge warranty service. This service will be charged to the customer depending on the effort required. The consumables components are not subject to warranty. Consumable components are silicon tubes, filters, seals, conical adaptor and suction catheter. Also excluded from warranty is all damage resulting from improper handling, wilful damage or improper care of the device.

The warranty shall expire if repairs and servicing are not carried out by technical service.

**RULES FOR RETURNING AND REPAIRING**

**UNDER NEW EUROPEAN RULES, GIMA REQUIRES THE FOLLOWING PROCEDURES TO BE CARRIED OUT TO PROTECT THE INSTRUMENT AND THE SAFETY OF ALL WHO COME IN CONTACT WITH IT.**

Before returning an instrument for repair, the external surfaces and all accessories **MUST** be carefully disinfected with a cloth soaked in methylated spirits or hypochlorite-based solution. The instrument and accessories should then be placed in a bag with a note outlining the disinfection undertaken.

Failure to follow this procedure will result in the instrument being returned to the purchaser unrepaired.

Instruments returned for repair **MUST** be accompanied by a description of the problem. GIMA will not be responsible for damage caused through improper use. To avoid such damage, please read the instruction carefully.

Where GIMA determines that an instrument is faulty, a replacement will be provided only if a SALES RECEIPT and STAMPED GUARANTEE are provided. GIMA will not be responsible for damage accessories. These may be replaced at the customer's expense.

**GIMA S.p.A.**

Via Marconi, 1 – 20060 Gessate (MI) – Italia

**ITALIA:** Tel. 02 953854.1 – Fax. 02 95381167

e-mail: [gima@gimaitaly.com](mailto:gima@gimaitaly.com) – [www.gimaitaly.com](http://www.gimaitaly.com)

**INTERNATIONAL:** Tel. +39 02 953854209/221/225 – Fax. +39 02 95380056

e-mail: [export@gimaitaly.com](mailto:export@gimaitaly.com) – [www.gimaitaly.com](http://www.gimaitaly.com)

GIMA S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) - Italy