



## 01. General device description Urilyzer 100 (Pro)

# 01. General device description – Urilyzer 100 (Pro)

### Content:

1. General Description
2. Specifications
  - 2.1. Technical Specifications
  - 2.2. Product Presentations
  - 2.3 Pack List



## 01. General device description Urilyzer 100 (Pro)

### 1. General description

The Urilyzer 100 analyzer is a lightweight instrument for reading CombiScreen Plus urinalysis strips (CombiScreen 11SYS Plus, CombiScreen 7SYS Plus and CombiScreen 5SYS Plus). It is a reflectance photometer that analyzes the intensity and colour of light reflected from the reagent areas of the strip. It is designed for In Vitro Diagnostic (IVD) use by qualified physicians and laboratory staff; however no special training is required to operate the instrument.

The analyzer can be set up to be as simple or sophisticated as you prefer. You may simply insert a dipped urinalysis strip into the analyzer and the result will be reported. By modifying the user options the measurements can be automatically started, printed and transferred.

Alternatively, you have the option to enter the Sample ID, Patient ID, and colour and clarity of the specimen manually. This added information will be reported along with the test results. The system also allows full customization to determine which fields appear on the printouts. You can also enable the user management features, thus the Operator ID is also recorded with the test results. The advanced security functions may also prevent unauthorized use through configurable settings. These features can be set independently.

The touch screen displays instructions and prompts you through the operation of the analyzer. In addition, you enter information through the touch screen. A barcode reader or an external keyboard can also be connected for accurate data entry to reduce transcription errors.



## 01. General device description Urilyzer 100 (Pro)

## 2. Specifications

The Urilyzer 100 and the Urilyzer 100 Pro have got several features, which facilitate the daily routine work and ensure a reliable and highly standardized urine analysis:

### 2.1 Technical Specifications

|                      |  |
|----------------------|--|
| <b>Type:</b>         | reflectance photometer with 4 discrete wavelengths 505, 530, 620, 660 nm   |
| <b>Throughput:</b>   | maximum 50 strips/hour (in normal mode)  |
| <b>Display:</b>      | 3.5" QVGA touch-screen LCD (resolution: 240x320)   |
| <b>Memory:</b>       | Urilyzer 100: 1000 test results / 500 QC results ()<br>Urilyzer 100 Pro: 3000 test results / 1000 QC results   |
| <b>Printer:</b>      | internal thermo printer (roll diameter max. 60mm)  |
| <b>Dimensions:</b>   | Width: 208 mm (8.2 inches)<br>Depth: 280 mm (11 inches)<br>Height: 69 mm (2.7 inches)  |
| <b>Weight:</b>       | 1.14 kg (unpacked, without power supply)   |
| <b>Power supply:</b> | 100...240V AC $\pm$ +10% -15%, 50/60Hz $\pm$ 5% external mains adapter   |
| <b>Interfaces:</b>   | PS2 (external keyboard, barcode reader)<br>serial RS232<br>USB Type B<br>USB Type A<br>Ethernet (Urilyzer 100 Pro only)<br>microSD card holder (Urilyzer 100 Pro only) |

#### Operating conditions:

|  |                  |                |                  |
|--|------------------|----------------|------------------|
| <b>Altitude:</b>                               | 3000 m (9842 ft) |                |                  |
|  | Optimum          | Operating      | In storage       |
| <b>Temperature:</b>                            | 20 °C to 26 °C   | 15 °C to 32 °C | -10 °C to +60 °C |
| <b>Relative humidity:<br/>(non-condensing)</b> | 35% to 55%       | 20% to 80%     | 20% to 85%       |



## 01. General device description Urilyzer 100 (Pro)

### 2.2 Product Presentations

#### Urilyzer 100



#### Urilyzer 100 Pro





## 01. General device description Urilyzer 100 (Pro)

### 2.3 Pack List

|             |   |
|-------------|---|
| Analyzer    | 1 x Urilyzer 100 <b>or</b><br>1 x Urilyzer 100 Pro  |
| Accessories | 2 x Test strip holder<br>1 x Check strip<br>1 x Power supply (AC Adapter 100–240 V, 50/60 Hz)<br>1 x Power supply cord<br>1 x Roll Thermal printer paper<br>1 x User Manual on CD<br>Quick reference guide(s) |