# GIMA

PROFESSIONAL MEDICAL PRODUCTS

Gima S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy  
Italia: tel. 199 400 401 - fax 199 400 403  
Export: tel. +39 02 953854209/221/225 fax +39 02 95380056  
gima@gimaitaly.com - export@gimaitaly.com  
[www.gimaitaly.com](http://www.gimaitaly.com)

**PINZE BIPOLARI**  
**BIPOLAR FORCEPS**  
**PINCES BIPOLAIRES**  
**PINZAS BIPOLARES**  
**PINÇAS BIPOLARES**  
**ΔΙΠΟΛΙΚΕΣ ΛΑΒΙΔΕΣ**

**MANUALE D'USO E MANUTENZIONE**  
**USE AND MAINTENANCE BOOK**  
**INSTRUCTIONS DE FONCTIONNEMENT ET ENTRETIEN**  
**MANUAL DE USO Y MANTENIMIENTO**  
**MANUAL DE USO E MANUTENÇÃO**  
**ΕΓΧΕΙΡΙΔΙΟ ΧΡΗΣΗΣ ΚΑΙ ΣΥΝΤΗΡΗΣΗΣ**

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## FEATURES

The bipolar forceps GIMA are surgical instruments to be used with equipment for bipolar electrosurgery. All models are covered with insulating material, can be sterilized in autoclave and are equipped with universal coupling for use with Diatermo GIMA and the most common models on the market.

## PRESCRIPTIONS

The product must be used by qualified medical personnel only, trained in electrosurgery techniques and maintenance operations on this kind of products.

Particular attention must be paid as regards the instruments care. The incorrect use or the lack of proper maintenance can make the instruments quickly deteriorate.

*These instruments can cause shocks or burns if used improperly!  
Avoid putting these instruments in contact with others, which are not insulated!*

*Do not use electro-surgery instruments on patients with Pacemakers, or in presence of flammable or anaesthetic liquids!*

Diathermic knife electrodes used with these instruments cause the destruction of tissues and are extremely dangerous if they are not correctly used. It is therefore recommended to carefully read the operating instructions supplied by the diathermic knife electrode manufacturer.

The product must be sterilized before using it.

## UNPACKING

*Always remember that packing elements (paper, cellophane, stitches, adhesive tape, etc.) can cut and/or hurt if they are not carefully handled.*

*They shall be removed with adequate means and shall not be left at the mercy of irresponsible persons; the same is valid for tools used to remove packages (scissors, knives, etc.).*

After opening the packages, first of all it is necessary to check all pieces and parts composing the product. Check that they are all present and in perfect conditions.

## MAINTENANCE

It is important to carefully examine the product before using it to monitor the presence of breaks, defects or wear and tear. Particular attention must be paid to moving parts, tips and connectors. The instruments must be kept in suitable boxes or protected trays, with their tips covered with gauzes or apposite caps in dry and clean place.

In case the product is defective or excessively worn out, do not use it and do not try to repair it. For repair apply to your own retailer, GIMA or qualified personnel.

## CLEANING

Remove all residues that have formed during the use of the product with a specific brush for cleaning instruments. The brush must be soft and not metallic. Then use a delicate detergent for cleaning. Thoroughly rinse the device with sterile purified water in order to eliminate all remaining detergent residues and dry well using a sterile cloth. Do not completely immerse the device in liquids.

## STERILIZATION

Separately wrap each clamp or put them in a container so they do not come into contact with each other or with other instruments.

**Steam:** steam sterilization in autoclave at 121°C (250°F) for 20 minutes.

**Flash:** steam sterilization in autoclave at 134°C (275°F) for 10 minutes.

**Gas (ETO):** carefully follow the instructions of the manufacturer of the sterilization unit.

Be careful that the product does not get in contact with acids or aggressive disinfectants, which can corrode it.

The product is reusable and in conformity with ANSI/AAMI standards on sterilizing.

We guarantee the product for a minimum of 20 sterilizing cycles, if they are sterilized in compliance with the following criteria. Most instruments fully pass this number of cycles in presence of a care and maintenance protocol.

It is recommended to follow the instructions of the autoclave manufacturer and arrange the instruments in such a way that their surface is directly and fully exposed to steam. Make them cool before using it.

## GIMA WARRANTY CONDITIONS

Congratulations for purchasing a GIMA product.

This product meets high qualitative standards both as regards the material and the production.

The warranty is valid for 12 months from the date of supply of GIMA.

During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel traveling expenses and packaging not included.

All components subject to wear are not included in the warranty.

The repair or replacement performed during the warranty period shall not extend the warranty.

The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use.

GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc.

The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed.

The defected products must be returned only to the dealer the product was purchased from.

Products sent to GIMA will be rejected.

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Tecno Instruments Pvt. Ltd.