**DIAPASON PER OTORINO GIMA - TUNING FORK GIMA FOR
 OTORHINOLARYNGOLOGY - DIAPASON POUR ORL GIMA - DIAPASÓN PARA
 OTORRINO GIMA - DIAPASÃO PARA MÉDICOS ESPECIALISTAS EM OUVIDOS GIMA -
 ΔΙΑΠΑΣΩΝ ΟΤΟΡΙΝΟΛΑΡΥΓΓΟΛΟΓΙΑΣ GIMA**

**Manuale d'uso e manutenzione
 Use and maintenance book
 Instructions de fonctionnement et entretien
 Manual de uso y mantenimiento
 Manual de uso e manutenção
 ΕΓΧΕΙΡΙΔΙΟ ΧΡΗΣΗΣ ΚΑΙ ΣΥΝΤΗΡΗΣΗΣ**

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Thank you for choosing the tuning fork GIMA for otorhinolaryngology. The absolute precision of the sounds is guaranteed by the special alloy with which the instrument has been made.

PRESCRIPTIONS

The product shall be used by qualified medical staff only, trained to the used of this kind of product.

The tuning fork GIMA for otorhinolaryngology is very simple to use and with some expedients its service life can be prolonged.

 *To make the tuning fork vibrate, do not use stiff objects, made of steel or other metals, make use of the digital pulp of your hand or alternatively a piece of wood.*

The tuning forks cannot be sterilized. To clean them, make use of a soft cloth.

In case of breaks, it is advisable not to use it. Apply to your dealer.

FEATURES AND OPERATION

The tuning fork GIMA for otorhinolaryngology is used for different audiometric tests (Rinnie test, Weber test, etc.), which allow testing the sensitivity to the sound, which is transmitted both through the air (COA) and through the bones (CO). To use it, it is sufficient to make it vibrate with a soft surface, such as digital pulp of the hand.

Conduction	Description	Parts of the apparatus where it is possible to detect possible auditive defects
COA	Through the air an acoustic stimulus is sent to the ear.	External auditive channel, medium or internal ear, acoustic nerve or central auditive pathways.
CO	It is carried out by placing the stem of the tuning fork in contact with the patient's head; the vibration caused by the instrument will directly stimulate the internal ear.	Orecchio interno, nervo acustico, vie centrali.

For a correct auditory test, it is necessary an active cooperation of the patient. In fact, the measurements carried out with the tuning forks are based on the answers given by the same patient during the test.

Table of vibration absolute numbers of major scale double vibrations c1 = 256

Code GIMA	Destination DEU - GB	FRA - ITA	C Ut (Do)	D Ré	E Mi	F Fa	G Sol	A La	H Si B
31230	c - h	ut ² - si ²	128	144	160	170 2/3	192	213 1/3	240
31231	c ¹ - h ¹	ut ³ - si ³	256	288	320	341 1/3	384	426 2/3	480
31232	c ² - h ²	ut ⁴ - si ⁴	512	576	640	682 2/3	768	853 1/3	960
31233	c ³ - h ³	ut ⁵ - si ⁵	1024	1152	1280	1365 1/3	1536	1706 2/3	1920
31234	c ⁴ - h ⁴	ut ⁶ - si ⁶	2048	2304	2560	2730 2/3	3072	3413 1/3	3840
31235	c ⁵	ut ⁷	4096						

TRANSPORT AND STORAGE

The equipment packed for transport and storage can be exposed to the following environment conditions:

- 1) Room temperature: from -40°C to +70°C
- 2) Relative Humidity: from 10% to 100% (condensate included).
- 3) Atmospheric pressure: from 500 to 1060 hPa.

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