
GIMA

OTOSCOPI PARKER
OTOSCOPES PARKER
OTOSCOPES PARKER
OTOSCOPIOS PARKER

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
FEATURES

The otoscope PARKER is medical equipment suitable to directly examine the external acoustic duct and the tympanic membrane. The otoscope is available in two different models: the first one is called PARKER otoscope and the second one PARKER COLOUR.

The two models differ in the colour and type of connection between the head and the handle. The otoscope PARKER (cod. 31440) is provided with a screw cap, while the model PARKER COLOUR has got a bayonet cap and it is available in different colours: green (cod. 31435), blue (cod. 31436) and yellow (cod. 31437).

Both models are equipped with a 2,5 V bulb and are supplied by two 1,5 V alkaline batteries type “C”. The otoscopes PARKER are provided with a rheostat to adjust the intensity of light, a removable glass allowing a 3X magnification and an output for air insufflation to check the mobility of the tympanic membrane.

The otoscopes PARKER are delivered in an imitation leather case and are equipped with three autoclavable specula Ø 2,5 - 3,5 - 4,5 mm.

For the otoscopes PARKER following accessories/spare parts are available:

Accessories codes for PARKER otoscopes	Minimum quantity
31447 Set 3 veterinary autoclavable specula	1 set
31448 Set 3 autoclavable specula	1 set
31449 2,5 V Bulb for PARKER otoscopes	1 pc.
31491 Disposable mini specula Ø 2,5 mm – black colour	500 pcs.
31492 Disposable mini specula Ø 4 mm – black colour	500 pcs.
31494 Adapter for disposable mini alligator speculum	1 pc.
31493 Adapter for disposable mini speculum	1 pc.
31495 Coloured disposable mini specula (various diameter)	1000 pcs.
31499 Small pump for air insufflation	1 pc.

The product or its components cannot be used for purposes different from the one specified in the present manual.

PRESCRIPTIONS


The product shall be used exclusively with alkaline batteries type “C” correctly installed. It shall not be dipped into water! When cleaning, make sure that the lamp socket or the battery compartment is not soaked with water.

The product shall be used by qualified medical personnel only, in out-patients’ department. It is not conceived for the employment in operation room. Due to its powerful luminous intensity it is inadvisable to directly point it at eyes.

Do not use the equipment in case it is damaged. Apply to your retailer. Avoid precarious repairs. Repairs shall be carried out with original spare parts only, which shall be installed according to the intended use.

Since the product is made of corrosion-proof materials suitable for the environmental conditions foreseen for its normal use, does not require special care, however it is necessary to store it in a closed place making sure that is protected from dust and dirt to assure its hygienic conditions. Moreover, it is recommended to store the product in a place which can be reached easily by the personnel in case of necessity.

UNPACKING

Always remember that packing elements (paper, cellophane, stitches, adhesive tape, etc.) can cut and/or hurt if they are not carefully handled. They shall be removed with adequate means and shall not be left at the mercy of irresponsible persons; the same is valid for tools used to remove packages (scissors, knives, etc.).

After opening the packages, first of all it is necessary to check all pieces and parts composing the product. Check that they are all present and in perfect conditions.

INSTALLATION

Before using it, it is necessary to insert the batteries into the apposite compartment. Proceed as follows:

- Open the battery compartment. To do this counterclockwise turn the apposite ring nut located on the handle.
- Once opened, insert the batteries paying attention to the direction of poles (see picture).
- To close the battery compartment, clockwise screw down the ring nut and check that the batteries are in contact with poles.
- In case the otoscope head unit separates from the handle, proceed as follows:
- For PARKER versions with screw cap, it is necessary to slightly unscrew the knob located on the side and pull the head while tightly holding the otoscope handle.

To reinstall the head unit on the handle, it is sufficient to insert it again into the apposite seat, taking care that the notch present on the head bottom matches with the knob. At this point it is sufficient to exert a light pressure towards the bottom and screw the knob.

- For PARKER COLOUR versions with bayonet cap, it is necessary to hold the otoscope handle and head group, counter-clockwise turn this last one until a click is heard, at this point it is possible to take the head off.

To insert the head, make the bayonet cap grooves match, push the head towards the bottom and clockwise rotate it.

- The operation necessary to remove and replace the speculum is very simple:

Counterclockwise rotate the speculum and draw it out of its seat. To reposition the speculum, insert it into the head group, make the groove match with the small pivot and clockwise rotate.

OPERATION

First of all it is necessary to choose the size and the type of speculum suitable to the acoustic duct to examine.

Besides the ones provided as standard, on request it is possible to deliver disposable specula with different sizes and colours (cod. 31491 - 31492 - 31495), which can be assembled on the head by means of an adapter (cod. 31493). An adapter can be also supplied to install a disposable alligator speculum (cod. 31494) (see picture), suitable to examine children.

The examiner shall turn on the light by means of the black knurled ring nut located on the handle upper side. Both PARKER models are equipped with a rheostat. Thus pressing the white pushbutton and slowly clockwise rotating the ring nut it is possible to increase the intensity of light. Counter-clockwise rotate the ring nut to turn down the light until it is off.

After the use, make sure that the light is off.

MAINTENANCE

Periodically check the battery conditions, making sure that no sign of corrosion or oxidation is present, in case of necessity replace them with new alkaline ones. Carefully handle the batteries as the liquids they contain can irritate skin and eyes.

Before being used, thoroughly check the product. The same operation shall be carried out after cleaning.

- Check that the connection between the head and the handle is perfect and that the rheostat correctly works.
- If the light is intermittent or in case it does not turn on, check the bulb, the batteries and the electric contacts.

In case the bulb stops working, it is possible to replace it with original spare parts (cod. 31449). To replace the bulb, it is necessary to take the speculum out, counter-clockwise turn the bulb and screw the new one down (see picture). It is advisable to use only original bulbs.

Before removing the bulb, make sure that the instrument has been off for some minutes; otherwise you run the risk of being burnt

- Check that the glass is clean; on the contrary it shall be removed. To do this counter-clockwise rotate it and at the same time pull it towards you.

After cleaning, reintroduce the glass into its seat making the groove match with the small pivot and clockwise rotating it (see picture).

- It is advisable to always have new bulbs and batteries at disposal for the maximum light.
- The otoscope head and handle can be cleaned with a cloth soaked with disinfectant, however they cannot be dipped into liquids.
- The speculum can be sterilized in boiling water or in autoclave, in this case the temperature shall not be over 121°C; alternatively cold disinfectants can be used with the same concentration used to sterilize medical instruments.

The product has been conceived to last forever; it is resistant to most part of chemicals except for methyl ethyl ketones (M.E.K.), diluents and nail polish removers. All internal parts are brass and stainless steel plated preventing rusting. If the lamp is exposed to humidity remove the batteries, thoroughly dry to avoid short-circuits.

GUARANTEE

Thank you for having acquired our product.

This product meets the most stringent requirements regarding the selection of quality manufacturing materials and also the final control. The product has a 12-month warranty, valid from the date it is delivered by GIMA.


During the guarantee period free repair and/or replacement of any defective parts due to faulty manufacture will be given, labour, postal, transport and packaging charges etc. are not included. The guarantee therefore excludes components subject to wear and tear such as parts in rubber or PVC, etc. No compensation will be given for loss of use of the product.

Furthermore replacement or repairs effected during the guarantee period do not lengthen the duration of the guarantee.

This guarantee is void in the event of: repairs effected by unauthorised persons or with spare parts not approved by GIMA, damage or defects caused by negligence, blows, abnormal use of the appliance or faulty installation.

The guarantee will be invalidated if the registration number has been removed, cancelled or altered. Faulty appliances must be returned only to the retailer where the item was bought. Any item dispatched directly to us will be refused.

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Disposal: The product must not be disposed of along with other domestic waste.

The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.

For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.

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ATTENTION: The operators must carefully read and completely understand the present manual before using the product.

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
FEATURES

The ophthalmoscope PARKER is a medical equipment suitable to examine the interior of the eye and especially the retina. The ophthalmoscope is equipped with a 2.7 V halogen light, for a better quality of image and powered by two alkaline batteries type "C", 1,5 Volt.

The ophthalmoscopes PARKER are provided with a rheostat for the adjustment of the light intensity and a wheel, which allows choosing among the 24 available lenses, from -25 to +40 diopters.

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PRESCRIPTIONS

The product shall be used exclusively with alkaline batteries type "AA" correctly installed. It shall not be dipped into water! When cleaning, make sure that the lamp socket or the battery compartment is not soaked with water.

The product shall be used by qualified medical personnel only, in out-patients' department. It is not conceived for the employment in operation room.

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
UNPACKING

Always remember that packing elements (paper, cellophane, stitches, adhesive tape, etc.) can cut and/or hurt if they are not carefully handled.

They shall be removed with adequate means and shall not be left at the mercy of irresponsible persons; the same is valid for tools used to remove packages (scissors, knives, etc.).

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- Open the battery compartment. To do this counterclockwise turn the apposite ring nut located on the handle.

- Once opened, insert the batteries paying attention to the direction of poles.

- To close the battery compartment, clockwise screw down the ring nut and check that the batteries are in contact with poles.

- In case the ophthalmoscope head unit separates from the handle, proceed as follows:

- Slightly unscrew the bigger knob located on the side and pull the head while tightly holding the ophthalmoscope handle.

To reinstall the head unit on the handle, it is sufficient to insert it again into the apposite seat, in connection with the notch present on the head bottom. Clockwise screw down the knob.

OPERATION

First of all, the examiner shall turn on the light by means of the black knurled ring nut located on the handle upper side. The ophthalmoscope is equipped with a rheostat. Thus pressing the white pushbutton and slowly clockwise rotating the ring nut it is possible to increase the intensity of light. Counterclockwise rotate the ring nut to turn down the light until it is off.

Then the examiner shall select the power of the corrective lens, by rotating the apposite knurled ring nut; the available values range between +25 and -40 diopters. The power of the selected lens can be read inside the hole placed on the bottom of the head.

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MAINTENANCE

Periodically check the battery conditions, making sure that no sign of corrosion or oxidation is present, in case of necessity replace them with new alkaline ones.

Carefully handle the batteries as the liquids they contain can irritate skin and eyes.

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- Check that the connection between the head and the handle is perfect and that the rheostat correctly works.

- If the light is intermittent or in case it does not turn on, check the bulb, the batteries and the electric contacts.

In case the bulb stops working, it is possible to replace it with original spare parts (cod. 31479).

To replace the bulb, it is necessary to counterclockwise unscrew the smaller knob placed on the bottom of the head. It is possible to perform this operation both with the head connected to the handle and separated from it. Once the upper section of the head has been removed, counterclockwise unscrew the bulb and replace it with a new one (see picture). When fixing the upper section of the head again, before screwing down, be sure that the screw is in connection with the apposite notch. It is advisable to use only original bulbs.

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
- It is advisable to always have new bulbs and batteries at disposal for the maximum light.

- The ophthalmoscope head and handle can be cleaned with a cloth soaked with disinfectant, however they cannot be dipped into liquids.

The product has been conceived to last forever; it is resistant to most part of chemicals except for methyl ethyl ketones (M.E.K.), diluents and nail polish removers. All internal parts are brass and stainless steel plated preventing rusting. If the lamp is exposed to humidity remove the batteries,

thoroughly dry to avoid short-circuits.

GUARANTEE

Thank you for having acquired our product. This product meets the most stringent requirements regarding the selection of quality manufacturing materials and also the final control. The product has a 12-month warranty, valid from the date it is delivered by GIMA.

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