DEUTSCH
ENGLISH
FRANÇAIS
ESPAÑOL
ITALIANO

HEINE Otoscoptes

HEINE Optotechnik GmbH & Co. KG
Kientalstr. 7 · 82211 Herrsching · Germany
Tel. +49 (0) 81 52 / 38 - 0
Fax +49 (0) 81 52 / 38 - 202
E-Mail: info@heine.com · www.heine.com
med 0514 2014-04-09

Bedienung

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Please read and follow these instructions for use of and keep them for future reference.

General Conditions of Warranty

Instead of the statutory warranty time period of 2 years, HEINE will grant a guarantee of 5 years from the date of the consignment of the goods ex works, concerning its equipment (excluding disposable, e.g. bulbs, single-use articles, and rechargeable batteries). Date of consignment means that HEINE hands the goods over to the transport carrier, freight forwarder or any other person designated by the Customer for the transport of the goods without loading the collecting vehicle. The guarantee covers irreplaceable workmanship, on condition of the proper use of the equipment and the observation of the operating instructions. During the warranty- and guarantee time period, errors and deficiencies arising on the equipment will be rectified free of charge, in so far as such are evidenced by defective materials, processing and/or constructional errors. Should buyer complain of a material deficiency during the warranty time period, then the onus of proof is always to be on the orderer, that the product was defective already upon receipt of the goods. The statutory warranty and the guarantee do not apply to loss or damage caused by wear and tear, negligent use, the non-employment of original HEINE components and/or spares (in particular bulbs, as these have been especially developed for HEINE instruments in accordance with the following criteria: colour temperature, useful service life, safety, optical quality and performance. The statutory warranty and the guarantee do not apply to interventions by persons not authorised by HEINE or when the operating instructions are not observed by the customer.

Any modification of a HEINE product with parts or additional parts which do not conform to the original HEINE specification will invalidate the warranty for the correct function of the product and further invalidate any warranty claims which result from such a change or modification. Further claims, in particular claims for replacement of loss or damage, which are experienced otherwise than directly on the HEINE product itself, are hereby excluded.

Intended Use


HEINE Otoscoptes are used for the examination of the ear, the external auditory canal and the tympanum. The instruments have an optical magnifying equipment and a battery- or rechargeable battery-powered illumination equipment. The instruments may only be used by authorized and qualified medical personnel.

For U.S. only:
Federal law restricts this device to sale by or on the order of a Physician or Practitioner.

Warnings and Safety Information

CAUTION! Indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent. (Background color yellow; foreground color black).

NOTE! Indicates valuable advice in terms of installation, operation, maintenance or repair. Notes are important, but not related to hazardous situations.

Product overview

- | | | |
|---------------------|---------------------|-------------------------------|
| 1 Slit | 4 Connector plug | 7 Insufflation port connector |
| 2 Viewing window | 5 HEINE Handle | 8 SANALON S tips (re-use) |
| 3 Insufflation port | 6 Insufflation bulb | 9 AllSpec tips (single use) |

Setting up

To set up the HEINE Ootoscoptes, screw the instrument head into the HEINE battery handle or plug it on the HEINE rechargeable handle.

Verify that the lamp voltage complies with the supply voltage of the handle. The coloured marking on the bottom of the lamp shows you the lamp voltage:

White ring = HEINE XHL 2.5V bulb only for use with the HEINE Battery handle
Red ring = HEINE XHL 3.5V bulb only for use with the HEINE Rechargeable Handle

This does not apply to the otoscoptes of the mini3000 series. Their lamp voltage is always 2.5V independently of the handle.

Operation

The HEINE otoscoptes may only be inserted into the auditory canal when a tip or a speculum is fitted. Place the speculum on the otoscope so that the projection inside the speculum fits into the bayonet slit (1). Tighten the speculum by twisting it slightly to the right.

The magnifying lens is inside the viewing window (2). To facilitate instrumentation, the viewing window can be swivelled to both sides or folded up.

The otoscoptes have a connection port (3) for the optionally available insufflation bulb (6, 7). The test of tympanic mobility can be carried out when the viewing window is closed using the insufflation bulb. Dose the pressure carefully with the insufflation bulb.

The setting up and operation of the HEINE handles are described in a separate instruction of use.

Hygienic Reprocessing

Instructions on hygienic reprocessing must be adhered to, based on national standards, laws and guidelines.

Classification according to KRINKO: non-critical

Spaulding Classification USA: non-critical

Allow the device to cool down before reprocessing.

When contamination is suspected, a hygienic reprocessing has to be carried out immediately or at least once a working day.

The described cleaning and disinfection measures do not replace the specific rules applicable for the establishment.

HEINE Optotechnik only approves the agents and procedures mentioned below. Other procedures/ reprocessing agents have to be validated.

Cleaning and disinfection may only be carried out by personnel with sufficient hygienic knowledge.

Observe the disinfectant manufacturer's instructions, especially the specifications regarding concentration and material compatibility and the safety guidelines.

Do not use spray or immersion disinfection, dripping wet or heavily foaming tissues. Do not use ultrasonic reprocessing.

Otoscope head and insufflation bulb

The otoscope and the insufflation bulb shall be cleaned outside manually with a soft cloth and inside with a cotton-wool bud (wipe cleaning and wipe disinfection).

The otoscope and the insufflation bulb are not suitable for mechanical reprocessing (cleaning, disinfection or sterilisation) or for immersion in liquids.

Recommended agents:

Cleaning agent: Neodisher® MediClean

Disinfectant agent: alcoholic (e.g. Incides® Tissus)

Reusable SANALON S Tips

The reusable tips shall be cleaned and disinfected manually (in an immersion bath) or by machine.

Recommended agents:

Cleaning agent: Neodisher® MediClean

Disinfection agent EU: rotasept®, USA: Cidex®OPA Solution

The reusable tips can be reprocessed up to 360 cycles.

Disposable tips AllSpec/UniSpec

AllSpec-Tips are for single use only. Do not attempt to clean them as the tip may be damaged and could injure the patient.

Changing the bulb

HEINE mini3000® LED F.O.

With the HEINE mini3000® LED F.O. Otokospe the LED cannot be changed.

HEINE mini3000®

In the case of the HEINE mini3000® Otokospe, swivel the viewing window (2) to the side and pull out the bulb backwards. Insert the new bulb as far as possible into the socket.

HEINE BETA® 200 F.O. / HEINE BETA® 400 F.O. / HEINE K180® F.O. / HEINE mini3000® F.O.

• Detach the otoscope from the handle and pull the bulb out of the guide tube of the connector plug (4) taking it at the narrow collar.

• Push in the new bulb as far as it will go.

HEINE BETA® 100 / HEINE® K100

• Remove the specula and pull the shield off the bulb.
• Unscrew the bulb.
• Replace the shield. This reduces stray light and reflections from the inside of the specula and should always be fitted.

Verify that the lamp voltage complies with the supply voltage of the handle.

ENGLISH

Maintenance and Service

The instruments do not require maintenance or service.

General Notes and Warnings

Check the correct operation of the device before use! Do not use the device if there are visible signs of damage.

Do not use the device in explosive or oxygen loaded environments!

Do not use the device near strong magnetic fields like MRI scanners!

Do not modify the instrument. Repairs shall only be carried out by qualified persons.

The speculum / tip must be examined for sharp edges before use.

The maximum on-time of the instruments shall be 1 minute in combination with an off-time of 10 minutes.

Do not shine light of the instrument directly into the eyes.

Avoid direct contact with mucous membranes or injured skin.

Store and use the device in dry and dust-free environments only! If you don't use the device for a longer period of time, please remove the batteries in advance.

The appendix contains the tables „Guidance and Manufacturer's declaration – Electromagnetic immunity“, technical data and an explanation of the used symbols.

Disposal

The product must be recycled as separated electrical and electronic devices. Please observe the relevant state-specific disposal regulations.

Electromagnetic Compatibility

Medical electric devices are subject to special precautionary measures with regard to electromagnetic compatibility (EMC). Portable and mobile high frequency communication equipment can affect medical electric devices.

This device is intended for use by medical professionals in the electromagnetic environment specified below. The user of this device should assure that it is used in such an environment. The use of accessories, converters or cables other than the ones specified by HEINE might lead to increased emission reduced electrical immunity of the medical equipment.

The device may not be stacked directly near or used directly beside other devices. If the device is to be operated in a stack or with other devices, the device should be watched to ensure it operates properly in this location.

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Guidance and manufacturer’s declaration – electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environment.		
Emission test	Compliance	Electromagnetic environment – Guidelines
RF emissions CISPR11	Group 1	The device uses RF energy only for its internal function. Therefore, RF-emission is very low and it is unlikely that any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. Warning: This device is intended only for use by medical professionals. This is a device of class A CISPR 11 in the domestic environment, this device may cause radio interference, so that it may be necessary in this case, to take appropriate remedial measures, as e.g. orientation, new arrangement or shielding of the device or restrict the connection to the site.
Harmonic Emissions IEC 61000-3-2	not applicable	
Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3	Passed	

Guidance and manufacturer declaration - Electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – Guideline
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	Floors should be wood, concrete or covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC 61000-4-4	± 2 kV for mains cables ± 1 kV for input and output lines	± 2 kV for mains cables ± 1 kV for input and output lines	The supply voltage quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV voltage phase – phase ± 2 kV voltage phase – earth	± 1 kV voltage phase – phase ± 2 kV voltage phase – earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U _r , (>95% dip in UT) for 1/2 period 40% U _r , (60% dip in UT) for 5 periods 70% U _r , (30% dip in UT) for 25 periods <5% U _r , (>95% dip in U _r) for 5 seconds	<5% U _r , (>95% dip in UT) for 1/2 period 40% U _r , (60% dip in UT) for 5 periods 70% U _r , (30% dip in UT) for 25 periods <5% U _r , (>95% dip in U _r) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered by a UPS (uninterruptible power supply) or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Comment: U _r is the a.c. supply voltage prior to application of the test level.			

Technical specification

Nominal voltage	2,5V or 3,5V
Nominal current	Typ. 440-760mA
Protection class	Internal power supply
Application part	Type BF
Environmental conditions for operation	10°C to 35°C 10% to 75% rel. humidity 700hPa to 1060hPa
Environmental conditions for storage and transport	-20°C to 50°C 40% to 95% rel. humidity 500hPa to 1060hPa
Device classification according to IEC 62471:2006	Exempt group

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Guidance and manufacturer’s declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment Guidelines
Conducted RF IEC 61000-4-6	3 Veff 150 kHz to 80 MHz	3 V eff	Portable and mobile RF communication equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated HF IEC 61000-4-3	3 V/m 80MHz to 2,5GHz	3 V/m	Recommended separation distance: d = 3,5/3 · SQRT (P/W) d = 3,5/3 · SQRT (P/W) 80 MHz to 800 MHz d = 7/3 · SQRT (P/W) 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey*, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:

Note 1: At 80Hz and 800MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted