**TERMOMETRO DIGITALE BASALE**

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**BASAL DIGITAL THERMOMETER**

Carefully read the following instructions before use. Basal temperature (BBT: Basal Body Temperature) is the temperature of the human body registered upon awakening before carrying out any form of activity. Measuring the basal temperature is one of the most simple and economic instruments to assess period trends and, in particular, to understand if ovulation has occurred. This is because basal temperature varies according to the period phase and according to the quantity of progesterone in circulation: during the follicular phase, when progesterone is practically absent, the basal temperature is much lower, while during the luteous phase, when the level of progesterone increases, the basal temperature is much higher compared with the previous phase. Across a period of approximately 28 days, around the 14th day a drastic temperature fall is registered (already in its lower phase) followed, the next day, by a drastic thermal increase that is completed within 48 hours. The temperature then remains high for the remaining part of the period. The falling peak informs the woman that ovulation is about to begin. The temperature then increases when the luteous body starts to produce progesterone. It is important to register several periods in order to understand how the trend of the basal temperature works. The best days for conception vary from the day in which the downward trend is registered to approximately 2 days after. The basal digital thermometer is a device for basal body temperature measurement suitable both for hospital and home use. Compared to the traditional glass and mercury models, it offers the advantage of simple reading, quick and accurate measurement, memory and alarms, etc. In particular it is safer for people and environment because of the absence of mercury. The thermometer has been calibrated by the manufacturer and further calibrations are not necessary. For any doubt on the device accuracy, please get in contact with your retailer.

**Technical features**

Measuring range: 32-43 °C  
Accuracy: ± 0,01 °C  
Battery life: approx. 100 hours of uninterrupted use

**Use**

- Before the use, disinfect the probe with alcohol or place a cover on the probe.
- Press the power button; a "beep" sound is issued.
- After 2 seconds the letter °C or °F starts blinking, now it is possible to start measuring.
- During the measurement, remain in the same place/position until the end.
- The measurement lasts about 2 minutes and when it is complete, the thermometer will sound an acoustic signal. The detected temperature is shown on the display. The thermometer will automatically turn off after about 10 minutes. However to prolong the life of the battery, it is best to press the power button soon after the use.
- To display the memorized value of the last measurement, press the ON/OFF button for 2 seconds with thermometer off.

The detected values vary according to the method of use.  
**Rectum measurement: From a medical standpoint, this is the most accurate and recommended method.** Delicately insert the probe in the rectum for Approx. 2 cm. Do not move until the measurement is completed. (40-60 seconds).  
**Armpit measurement: Compared to measurement in the rectum, armpit measurement generally varies by about 0.5 °C ~ 1.5 °C.** Before starting, keep your arm close to your body so the armpit remains close for about 5 minutes. Place the probe in the centre of the armpit in a stable position; do not move your arm until the measurement end. (80-120 seconds).  
**Tongue measurement: Compared to measurement in the rectum, the value measured is generally 0.3 °C ~ 0.8 °C lower.** Place the probe under the tongue in a stable position in contact with the mucous, breathe with nose and do not open your mouth until measuring has been concluded (50-70 seconds).  
To make comparison of the values measured in different periods easier it is **necessary to always use the same measurement method.**  
To obtain a more accurate measurement we suggest continuing measurement for 3 minutes after the beep.

**Disinfect the probe only with alcohol or similar, DO NOT boil it! DO NOT sterilize in autoclave!**  
The product is water-resistant therefore it can be washed under running water. In any case DO NOT dip the product for too long!  
Avoid exposing the unit to direct sunlight or atmospheric agents.  
Do not measure except body temperature!

**Battery replacement**  
- The battery must be replaced when the symbol or blinks on the display.  
- Delicately withdraw the back cover and remove the old battery by means of a toothpick or a needle tool.  
- Insert the new 1,5 V battery (LR/SR-41).  
- Replace the back cover

**Do not disassemble the unit except for battery replacement!**  
Do not leave the batteries or the back cover within children's reach!  
The product or battery disposal must be in compliance with the standards in force in your country.

**Disposal: The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.**  
For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.

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