

## USE INSTRUCTION FOR MINI CYCLE

- 1. Take mini cycle body out of box, then handle it with care. Look at picture 1.
- Lay out mini cycle body on the ground in the same level, then put knob screw A onto it as picture shows. Lookat picture 2.
- 3. Turn Knob screw A so as to make the axis B pressure properly. (toward left to put on pressure, otherwise decrease pressure) Turn on the calculator button C and rotate axis B, immediately it is running in use. Look atpicture 3.(see the calculator instruction for detail.)
- 4. Pull the pin 0 toward inside after use and fold it welllike picture 5 does.
- (1) Inspect it and make sure it is placed steadily and in same level before use;
- (2) Don't forget to let the pressure be minimum when start to use, then increase it to reach what it doesneed;
- (3) Applicable people :adult,children,old people;

WARNING: Please don't try to stand up on the mini cycle, that will damage the machine.



## **EXERCISE COMPUTER INSTRUCTION MANUAL** SPECIFICATIONS:

-00·00-99·59MIN TIME COUNT-REPS/MIN-----0-9999 CA LORIE

KEY FUNCTION:

MODE: This key lets you to select and lock on to a particular function you want, and hold the key for 4 secondsto have all function values reset.

OPERATION PROCEDURES:

AUTO ON/OFF: The system turns on when any key is depressed or when it sensor an input from the speedsensor. The monitor will be automatically shut off if there is no any signal coming in 4 minutes

NOTE 1.Battery Spec:1.5V, LR44 (1pc) - 2.If improper display on monitor, please re-install the batteries. - 3. Battery install procedure









GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

Fabbricante / Manufacturer / Fabricant / Fabricante / Fabricante Hersteller / Producent / Producător / Gyártó / Παραγωγός: 
ZHEJIANG TODO HARDWARE MANUFACTURE CO., LTD No.31th Floor, Jinmao Building, CBD Centre, Yongkang, Zhejiang, China Made in China

