# GIMA

**BORSA ACQUA CALDA**  
**HOT WATER BAG**  
**BOUILLOTTE**  
**WARMWASSERTASCHE**  
**BOLSA AGUA CALIENTE**  
**BOLSA DE ÁGUA QUENTE**  
**ΘΕΡΜΟΦΟΡΑ**

**MANUALE D'USO E MANUTENZIONE**  
**USE AND MAINTENANCE BOOK**  
**INSTRUCTIONS DE FONCTIONNEMENT ET ENTRETIEN**  
**BETRIEBS UND WARTUNGS ANWEISUNGEN**  
**MANUAL DE USO Y MANTENIMIENTO**  
**MANUAL DE USO E MANUTENÇÃO**  
**ΕΓΧΕΙΡΙΔΙΟ ΧΡΗΣΗΣ ΚΑΙ ΣΥΝΤΗΡΗΣΗΣ**

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## INTENDED USE

Before using it, ALWAYS ask for your doctor's opinion.

BEFORE EACH USE and from time to time apply to your doctor, who will show you possible collateral effects due to the specific use of the Device.

The present device aims at easing symptoms and/or dysfunctions generated on the human body by diseases; the bag can only contain water at the temperature prescribed in the specifications.

It can only be applied to external healthy parts of the body.

Other uses are considered improper and therefore FORBIDDEN:

- Do not use it in contact with injured parts or on parts inside the human body.
- Do not use it as invasive tool or together with invasive tools.
- Do not use together with or as part of other medical instruments, or together with drugs or substances considered improper.
- Do not modify or alter its components and/or features.
- Do not blow into the device and do not inflate it.

## PRECAUTIONS AND WARNINGS

- STRICTLY OBSERVE THE INTENDED USE OF THE DEVICE.
- THE DEVICE CANNOT BE USED BY CHILDREN, DISABLE OR IMPAIRED PEOPLE WITHOUT SUPERVISION.
- THE WATER TEMPERATURES SHOWN IN THE SPECIFICATIONS ARE THE MINIMUM AND MAXIMUM TEMPERATURES ALLOWED IN ORDER TO AVOID PHYSICAL ALTERATION. WARM UP WATER AT A TEMPERATURE SUITABLE FOR YOUR BODY SENSITIVENESS.
- IN CASE OF CONTEMPORARY USE OF SIMILAR DEVICES DO NOT INVERT DETACHABLE PARTS AND/OR ACCESSORIES.
- WHILE POURING WATER INSIDE THE BAG, MAKE USE OF HEAT PROTECTION EQUIPMENT (GLOVES) TO AVOID THE RISK OF BURNS IN CASE OF ACCIDENTAL SPRAYS.
- AFTER THE BAG IS FILLED, MAKE THE EXCEEDING AIR LEAK BEFORE PLUGGING IT. CLOSE THE BAG AND CHECK THE LOADING INLET TIGHTNESS.
- DO NOT USE OR LAY THE DEVICE ON/WITH: HOUSEHOLD ELECTRIC APPLIANCES, HOT OR POINTED SURFACES.
- BEFORE USING THE DEVICE ALWAYS MAKE SURE THAT ALL COMPONENTS ARE INTEGRAL AND NOT DAMAGED (FOR EXAMPLE CRACKING, DENTING, TEARING), THEIR COLOUR IS UNCHANGED AND THEY ARE FREE FROM SURFACE MOULDS, ETC.
- IT IS ADVISABLE TO BIND THE CLOSING PLUG AT THE LOADING INLET WITH SUITABLE ROPE.

## FEATURES

<b>Device class</b> I (93/42/CEE Directive)	<b>Device service life (normal use)</b> See label	<b>Capacity</b> 2 litres
<b>Disposable device</b> NO	<b>Custom-sized device</b> NO	<b>Device delivered STERILE</b> NO
<b>Water temperature</b> From + 15 °C to + 50 °C	<b>Overall dimensions (cm)</b> 37 x 20 x 44	<b>No-load mass (N/Kg.)</b> 4,12 / 0,42

For a longer-lasting service life, it is advisable to store it in the original packing and in dry environment with temperatures ranging from +1 °C to + 40 °C, far from heat sources (sun or adverse weather conditions), sheltered from dust and dangerous materials and far from the children's reach.

### MAINTENANCE AND DISPOSAL

Before using the device for the **FIRST TIME** and each time it is necessary, wash the interior of the bag with warm drinkable water and the external surface with warm water and neutral soap; do not use alcohol, solvents or acids. Do not use brushes, but only soft sponges so that not to scratch the rubber surface. After the use, completely empty the bag and make it drain/dry. For long periods of non-use, it is advisable to spread the device with talcum powder.

In case of replacement or breakages, dispose it according to the local legislation on the subject of waste disposal. No toxic substances are present. The device is composed of rubber and plastic material.

### GUARANTEE

Thank you for having acquired our product.

This product meets the most stringent requirements regarding the selection of quality manufacturing materials and also the final control. The product has a 12-month warranty, valid from the date it is delivered by GIMA.

During the guarantee period free repair and/or replacement of any defective parts due to faulty manufacture will be given, labour, postal, transport and packaging charges etc. are not included. The guarantee therefore excludes components subject to wear and tear such as parts in rubber or PVC, etc.

No compensation will be given for loss of use of the product.


Furthermore replacement or repairs effected during the guarantee period do not lengthen the duration of the guarantee.

*This guarantee is void in the event of: repairs effected by unauthorised persons or with spare parts not approved by GIMA, damage or defects caused by negligence, blows, abnormal use of the appliance or faulty installation.*

The guarantee will be invalidated if the registration number has been removed, cancelled or altered.

Faulty appliances must be returned only to the retailer where the item was bought. Any item dispatched directly to us will be refused.

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