**Please** read these safety instructions carefully before use, and keep them for future question.

### 1. Safety indications

- ♥ ALWAYS consult with your physician or supplier to determine proper adjustment and use of the device.
- ♥ Don't store outside and don't let under the sun (Beware: don't let the product in a car in case of high outside temperature !!).
- ♥ Store in a dry room in a flat position to prevent deforming, stretching and tearing.
- ♥ Store away from sharp objects and pointed edges to prevent tears and punctures.
- ♥ If stored below 15°C, unpack the heel protector and expose it for 30 minutes in a room with a temperature > 15°C.
- ♥ Transport flat pads to avoid tears and punctures.
- ♥ Look for cracks in foam base before use.
- ♥ Before using the heel protector, make sure that the surface where you will place it, doesn't have sharp or pointed edges which could destroy the foam base.
- ♥ Don't use in extreme temperature conditions (above 38°C or 100°F or below 0°C or 32°F).
- ♥ No modification must be made to this device.
- ♥ ALWAYS check if the heel protector is properly and securely positioned in place before using.
- ♥ In case you didn't use your heel protector for a long time, let it checked by the technical staff, before using.
- ♥ ALWAYS observe the weight limit on the labelling of your product.
- ♥ Use only as a heel protector and in a proper way: don't knee or step on the product.
- ♥ If your skin develops redness, discontinue the use of the heel protector immediately and consult your doctor or therapist.
- ♥ Follow the advice of the professional who prescribed this medical device, respect duration and frequency of use.
- ♥ In case of discomfort consult your professional.

### 2. Indication

It is easier to prevent pressure ulcers than to treat them. Bedsores can be avoided by changing frequently lying or seating positions to avoid stress on the skin and to minimize the risk of pressure on vulnerable bony areas. Other strategies include skin care, regular skin inspections and good nutrition.

Using a gel or foam support can help relieve pressure and make sure that the body is correctly positioned in the chair.

### 3. Description and Fitting instructions

The heel protector is a contoured moulded visco-elastic foam support. The foam increases the contact surface and maximizes pressure relief. The heel protector's shape ensures a good posture of the leg and completely offloads heels.

This heel protector is delivered with a polyurethane cover. Hook and loop closure.

 Viscoelastic foam is a «slow recovery» foam. During compression, the foam progressively conforms to the shape of the body, and after the weight is removed, the foam slowly reassumes its initial shape.

It also reacts to body temperature and ambient temperature, softening with heat and more easily adjusting to body contours.

Foam has following characteristics:

- Micro perforated cells bounce slower, better viscoelasticity.
- Air circulates through the foam.
- Viscoelasticity remains constant over a wide temperature range. Pleasant touch due to very fine cell structure.
- Low emission and low odor foams.

Ask your technical staff how to install the cushion.

### 4. Washing instructions

To clean Cover:

Remove the cover from the cushion. Cover is machine washable at 30°C. Dry flat and away from any source of heat. Do not use detergent or bleach, do not autoclave.

To clean Cushion:

Remove the cover, then wipe lightly with a damp cloth. Do not use soap and do not submerge in water. After cleaning make sure to wipe dry and air dry the heel protector before using. Re-install the cover over the heel protector.

### 5. Composition

Polyurethane CosyPUr viscoelastic foam developed by BASF. Cover is made of polyurethane.

### 6. Disposal

Dispose of this product in accordance with your regional waste disposal regulations.

### 7. Technical specifications

This product has a 2-year warranty, from the date of purchase, against any manufacturing defect.

This product is biocompatible. This product is CE certified.

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