# LARINGOSCOPI A FIBRE OTTICHE
# *FIBRE-OPTIC LARYNGOSCOPES*
# LARYNGOSCOPES À FIBRES OPTIQUES
# *LARYNGOSKOP MIT FIBER OPTIK LICHTLEITER*
# LARINGOSCOPIOS DE FIBRAS ÓPTICAS
# *LARINGOSCÓPIO A FIBRAS ÓPTICAS*
# ΛΑΡΥΓΓΟΣΚΟΠΙΑ ΜΕ ΟΠΤΙΚΕΣ ΙΝΕΣ
# منظار الحنجرة العامل بالألياف البصريّة

## MANUALE D'USO E MANUTENZIONE
## *USE AND MAINTENANCE BOOK*
## INSTRUCTIONS DE FONCIONNEMENT ET ENTRETIEN
## *BETRIEBS UND WARTUNGS ANWEISUNGEN*
## MANUAL DE USO Y MANTENIMIENTO
## *MANUAL DE USO E MANUTENÇÃO*
## ΕΓΧΕΙΡΙΔΙΟ ΧΡΗΣΗΣ ΚΑΙ ΣΥΝΤΗΡΗΣΗΣ
## *دليل الإستعمال والرعاية*

**ATTENZIONE:** Gli operatori devono leggere e capire completamente questo manuale prima di utilizzare il prodotto.
***ATTENTION:*** *The operators must carefully read and completely understand the present manual before using the product.*
**AVIS:** Les opérateurs doivent lire et bien comprendre ce manuel avant d'utiliser le produit.
***ACHTUNG:*** *Die Bediener müssen vorher dieses Handbuch gelesen und verstanden haben, bevor sie das Produkt benutzen.*
**ATENCIÓN:** Los operadores tienen que leer y entender completamente este manual antes de utilizar el producto.
***ATENÇÃO:*** *Os operadores devem ler e entender completamente este manual antes de usar o produto.*
**ΠΡΟΣΟΧΗ:** Οι χειριστές αυτού του προϊόντος πρέπει να διαβάσουν και να καταλάβουν πλήρως τις οδηγίες του εγχειριδίου πριν από την χρήση του.
***الحذر:*** *على العمال قراءة وفهم هذا الدليل بكامله قبل البدء باستعمال المنتج.*


CE

**GIMA Spa - Via Marconi, 1 - 20060 Gessate (MI) - Italia**

**ITALIA**: Tel. 199 400 401 (8 linee r.a.) - Fax 199 400 403
E-mail: gima@gimaitaly.com - www.gimaitaly.com

**INTERNATIONAL**: Tel. ++39 02 953854209 - Fax ++39 02 95380056
E-mail: export@gimaitaly.com - www.gimaitaly.com

*M34487-M-Rev.2-04.10*

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
## FEATURES

The GIMA laryngoscopes are useful medical appliances to examine the larynx and to execute tracheal intubations.
The Green GIMA line is made by anti-magnetic stainless steel which is highly resistant to corrosion and uses fibre-optic blades and handles compatible with all major producers (certified ISO 7376-3). Light is transmitted from the handle to the tip of the blade through a fibre-optic bundle that can be taken apart to be sterilized independently or to be replaced in case of malfunction.
The handle uses two alkaline batteries "C" type 1,5 volt or rechargeable batteries, not supplied with the item.
The product or its components cannot be used for purposes different from the one specified in the present manual.

## PRESCRIPTIONS

*The product shall be used exclusively with alkaline batteries type "C" correctly installed. It shall not be dipped into water! When cleaning, make sure that the lamp socket or the battery compartment is not soaked with water.*

*Do not use the equipment in case it is damaged. Apply to your retailer.*
*Avoid precarious repairs. Repairs shall be carried out with original spare parts only, which shall be installed according to the intended use.*

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**Disposal:** The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.
For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.

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