

HEIGHT ROD

THE SOEHNLE PROFESSIONAL HEIGHT ROD – WHEN IT'S A QUESTION OF HEIGHT.



Height measuring appliances

Folding height rod 5003



- > Ultrasound measuring technology, temperature-compensated
- > Tilt sensor to correct if measuring instrument is out of horizontal and error message
- > Small space requirement due to fold function
- > Powered by two 1.5V AA batteries
- > Plastic housing
- > LCD, digit height: 9mm
- > Acoustic acknowledgement after measurement is completed
- > Unit switchover to "inches"

| Order Number | Measuring Range | Divisions | Dimensions | Weight |
|--------------|-----------------|-----------|---------------|--------|
| | cm | cm | mm | g |
| 5003.01.001 | 50-300 | 1 | 345 x 92 x 28 | 250 |

Body height tape measure 5002.01



- > Body height tape measure for wall mounting
- > Small and practical

| Order Number | Measuring Range | Divisions | Dimensions | Weight |
|--------------|-----------------|-----------|----------------|--------|
| | cm | mm | mm | g |
| 5002.01.001 | 0-220 | 1 | 150 x 30 x 140 | 130 |

Body height rod 5002.02



- > For mounting on scales 7800, 7820 and 7700
- > Wall mounting possible
- > Folding head stop

| Order Number | Measuring Range | Divisions | Dimensions | Weight |
|--------------|-----------------|-----------|-----------------|--------|
| | cm | cm | mm | g |
| 5002.02.001 | 105-210 | 0.5 | 50 x 1 175 x 22 | 800 |

Baby height rod 5001.10



- > Easy to use, quick and precise
- > Slide made of white plastic
- > Rod made of metal
- > Easy-to-clean

| Order Number | Measuring Range | Divisions | Dimensions | Weight |
|--------------|-----------------|-----------|------------|--------|
| | cm | cm | mm | g |
| 5001.10.001 | 40-80 | 0.5 | 500 x 340 | 300 |